

Waived Combination Survey

Table of Contents

Kit Contents.....	1
Important: Before You Begin.....	1
Detailed Testing Instructions.....	2
Reporting Your Results.....	2
Biohazard Warning.....	5
For Assistance.....	6
Master Lists.....	7

Kit Contents

- Two whole blood specimens
 - HE-01 and HE-02
- Two urine specimens
 - CM-01 and CM-02

Important: Before You Begin

New for this Mailing

Beginning with this mailing, the result form is available in an HTML format via e-LAB Solutions Suite. Depending on the program's reporting requirements, sections of the result form may include logic that flags missing required information (eg, units of measure, methods). Your laboratory must still verify and approve all data prior to the due date.

Reporting Code Changes

The following manufacturers have either deleted or updated codes for this mailing:

- Arkray

Storage and Stability Instructions

1. Upon receipt store specimens at 2 - 8°C.
2. Specimens are stable for 14 days at this temperature.

Critical Reporting Information

1. For any testing that you do not routinely perform in your laboratory, leave all reporting areas for that test blank unless otherwise noted.
2. A unit of measure (UOM) must be filled in for the hemoglobin (HGB) challenge for proper evaluation.
3. If you report hematocrit (HCT) using microhematocrit centrifugation, spin these specimens for exactly 5 minutes. This instruction does **not** apply to CritSpin/StatSpin or HemataSTAT II/C-70 users
4. **Urinalysis: Failure to select an inappropriate instrument, method, or manufacturer may result in comparison with the wrong peer group and unacceptable performance.** To accommodate various reporting styles, multiple reporting sections are provided for protein, glucose, ketones, bilirubin, blood/hemoglobin, and leukocyte esterase. **Choose ONLY 1 reporting method to enter your results.** Report results as you would a patient result.
 - If you use the Siemens Reagent Strip on the Siemens Clinitek 500, report your



Important: See the Biohazard Warning at the end of these instructions.
Customer Contact Center 800-323-4040 option 1 (domestic) or 001-847-832-7000 option 1 (international)

ATTACHMENT C

method as 2936 Siemens Reagent Strips (All) and your instrument as 2292 Siemens Clinitek 500.

- If you use the Roche Chemstrip and read this visually, report your method as 2937 Roche Chemstrip/combur and your instrument as 2044 Visual.

Detailed Testing Instructions

Handling Instructions

Whole blood specimens for hemoglobin and hematocrit testing:

1. Allow the vial to stand at room temperature for 20 minutes before mixing.
2. Briskly roll the vial back and forth between the palms of the hands for 60 seconds, then gently invert 10 times.
3. Mix until the cells are completely and uniformly suspended, ie, no cell button is seen on the bottom of the tube. Do not shake the vial or use a mechanical mixer.
4. Allow the vial to rest undisturbed about 15 seconds for dispersion of small air bubbles.
5. Immediately prior to use, gently invert the vial 10 times.
6. **Hemocytometer Users:** Preparation time must be extended to ensure proper lysing of fixed cells. The initial dilution should sit for 20 - 30 minutes in the Unopette. During this time period, prepare a Petri dish moisture chamber by moistening a folded paper towel or other absorbent material with water and placing it in the bottom of the Petri dish. The hemocytometer should be placed on top of the moist towel and allowed to sit in the covered moisture chamber for 10 - 15 minutes before reading. The moist chamber will prevent specimen evaporation.

Urine specimen testing:

Urinalysis

1. Allow specimen to come to room temperature for 20 minutes before mixing.
2. Mix specimen well before beginning tests. **Complete testing within 8 hours.** Perform only those tests that are ordinarily performed in your laboratory.
3. Do not perform confirmatory tests for positive dipstick results, even though your laboratory protocol requires them for patient specimens.
4. The analyte used in urine proficiency specimens for bilirubin and urobilinogen testing may give slightly atypical (pink) color development by some methods when examined visually. If this occurs, the specimen should be graded based upon the intensity of the color.

Urine hCG

1. Allow specimen to come to room temperature before mixing.
2. Mix specimen well before beginning tests.
3. **Complete testing within 8 hours.**

Reporting Your Results

General Reporting Instructions

1. Each mailing, verify the accuracy of your reporting codes (eg, manufacturer, method, instrument, reagent) by reviewing the online result form or the Method Summary Page attached to the front of your result form.
2. The inclusion of reporting codes on the result form does not imply US FDA approval.
3. **Exception Codes:** If you must report an analytical problem for a test or individual analyte, leave the result area for that section blank and select one of the following

bubbles on the result form within that section. The exception code bubble that you select will apply **only** to the result area(s) left blank. Documentation on the use of these codes is the responsibility of the laboratory and should be kept internally.

- **11 Unable to analyze**

Use code 11 to indicate why specimens were not analyzed (eg, instrument not functioning, reagents not available).

- **22 Result is outside the method/instrument reportable range**

Use code 22 if you obtain a high or low result outside the reportable range of your method or instrument. Do not use this code if there is an option to select **greater than or less than**.

- **33 Specimen unsatisfactory**

To use code 33, you **must** contact the CAP.

If you select an exception code bubble **and** enter data on the result form, the data will be graded.

4. An option for **greater than** values is provided in certain sections of the result form. Select the option if appropriate and enter your highest detectable limit in the boxes provided or refer to the use of exception codes.
5. Quantitative results for this Survey can now be automatically transmitted using e-LAB Solutions Connect. To learn more, visit cap.org and search for e-LAB Solutions Connect.
6. For reporting Urinalysis and Specific Gravity method and instrument codes, refer to list below for examples:

If you use:	Then for method use:	Then for instrument use:
Refractometer	2533 Refractometer	2044 Visual
Siemens Reagent Strip on Siemens Clinitek 500	2936 Siemens Reagent Strip	2292 Siemens Clinitek 500
Siemens Reagent Strip and read visually	2936 Siemens Reagent Strip	2044 Visual
Roche Chemstrip on Roche Chemstrip 101	2937 Roche Chemstrip	3241 Roche Chemstrip 101
Roche Chemstrip and read visually	2937 Roche Chemstrip	2044 Visual

7. **Bayer Acetest users:** Results should be read at 30 seconds as instructed for urine specimens in the Bayer package insert. Positive results may change if incubated longer due to the use of surrogate material for the Ketone analyte.
8. **Reducing Substances:** If your laboratory performs reducing substances, you must perform both challenges, otherwise leave blank if your laboratory does not perform this test.
9. **Urobilinogen:** If you select the reporting bubble for code **418**, note that $< 35 \mu\text{mol/L}$ is to be used by ARKRAY Aution MAX users only.
10. **Roche users:**
 - **Chemstrip, Specific Gravity (Visual):** For increased accuracy, add 0.005 to the reading if the pH is ≥ 7.0 . (If the strip is read on an instrument, the correction is applied automatically.)
 - **pH:**
 - **pH Visual:** Results should be reported in whole numbers; 0.5 increments are not allowed.
 - **Cassette and all other pH Instrumentation:** Results should be reported in whole numbers from 5.0 to 9.0, with the exception of 6.5.
 - **Glucose:**
 - **Roche Glucose Visual Method only:** Results may be reported as 500 mg/dL.

- **All Other Roche Glucose Methods:** 500 mg/dL is not an appropriate response. Analyzers are not able to discern results of 500 vs. 1000 mg/dL.
 - **Ketones:** All Roche Urinalysis instruments should not report 80 - 100 mg/dL. Analyzers are not able to discern ketones between 60 and 150.
 - **Protein, Urisys 1100:** 75 mg/dL is not an appropriate response.
11. **Siemens users:**
- **Multistix, Specific Gravity (Visual):** According to the Siemens package insert, for increased accuracy add 0.005 to the visual specific gravity reading when the urine pH is equal to or greater than 6.5. Report the corrected specific gravity on your result form. (Results from strips read by Clinitek are automatically corrected.)
 - **Clinitek Atlas users, Protein:** For reporting qualitative protein values of ≥ 300 , select the bubble for 300 - 600 mg/dL or 3+.
12. When reporting a high qualitative value, **do not** choose a response greater than your method allows.

Conversion Table (Protein, Glucose, and Ketones)

Protein Method 2 reporting options in SI units	
Conventional (mg/dL)	SI (g/L)
Negative	Negative
10 - 20	0.1 - 0.2
30 - 70	0.3 - 0.7
75 - Roche Users only	0.75
100 - 200	1.0 - 2.0
300 - 600	3.0 - 6.0
> 600 or ≥ 1000	> 6.0 or ≥ 1000
Glucose Method 2 reporting options in SI units	
Conventional (mg/dL)	SI (mmol/L)
Negative	Negative
30 - 100	1.6 - 6.0
150 - 300	8.3 - 16.6
500	28
> 500 or ≥ 1000 or ≥ 2000	> 28 or ≥ 55 or ≥ 111
Ketones Method 3 reporting options in SI units	
Conventional (mg/dL)	SI mmol/L
Negative	Negative
5 - 10	0.5 - 1.0
15 - 20	1.5 - 2.0
40 - 60	3.9 - 6.0
80 - 100	7.8 - 10.0
≥ 150	≥ 15.0

Per CLIA, as published by the United States Federal Register

- Proficiency Testing (PT) specimens must be tested with the laboratory's regular workload, using routine methods, and testing the PT specimens the same number of times it routinely tests patient specimens.
- If referral for testing is routinely performed for patient specimens, the practice cannot be followed for PT specimens. Referral is considered to be movement of the specimen from a laboratory with a CLIA identification number to another laboratory that has a different CLIA identification number.
- Laboratories must ensure that personnel do not share results or refer PT specimens for any reflex or testing outside their CLIA identification number.

Disclaimer

Survey specimens, their progeny, unmodified derivatives, or modifications thereof may not be transferred or incorporated into a program intended for sale. Survey specimens, their progeny, unmodified derivatives, or modifications thereof, reagents, and disposable equipment used in PT, when disposed of, should be autoclaved or incinerated and disposed of as hazardous waste. Disposal must follow local regulations, if more stringent than regulations enforced by the CDC or the FDA.

Submitting Results

1. Results **must** be received at the CAP no later than midnight, Central Time by the due date on the result form.
2. Your site administrator must establish a Web account via cap.org and opt in your laboratory. Information about opting in and a unique PIN was emailed to all laboratory directors. If your laboratory director does not have this information, contact the CAP.
3. Laboratory staff who enter results online must first create a personal Web account, affiliate themselves with a laboratory, then log in.

Biohazard Warning

All Survey specimens should be treated as if potentially infectious and should be handled as if they are capable of transmitting disease.

Survey specimens are prepared from blood or other source material obtained from human donors or animals.

When working with Survey specimens, precautions should be taken to protect yourself and others from accidental exposure to infectious agents such as HIV, HBV, and HCV.

HIV can be transmitted through accidental parenteral inoculation, mucous membranes, or non-intact skin contact with HIV infected blood or body fluids. HBV and HCV can be transmitted through accidental parenteral inoculation, mucous membranes, non-intact skin contact, aerosolization or ingestion.

Precautions described in CDC and FDA recommendations and OSHA blood borne pathogen rules should be followed at all times when handling Survey specimens and reagents.

Such precautions include the following:

- Gloves should be put on **before opening the container** and should be kept on throughout the period specimens are handled. Replace gloves if contaminated, or if their ability to function as a barrier is compromised.
- At high altitudes, specimens should be opened in a hood or biologic safety cabinet.
- There should be no eating, drinking, or smoking in the laboratory.
- Hands should be washed after removing gloves and before leaving the testing area.
- Survey specimens and reagents should be kept in separate refrigerators from those containing blood or blood components for transfusion.
- Survey specimens, reagents, and disposable equipment used in testing should be autoclaved or incinerated and disposed of as hazardous waste. Disposal must follow local regulations, if more stringent than regulations enforced by the CDC or the FDA.

Warning: This mailing may contain chemicals known to the State of California to cause cancer, and to cause birth defects or other reproductive harm.

If there has been an accident in which you have been exposed to the Survey's materials, please call the CAP Hot Line at 800-443-3244 (domestic) or 001-847-470-2812 (international) at any time. You can access Safety Data Sheets (SDS/MSDS) by logging on to cap.org, clicking on the Laboratory Improvement tab, then Catalog and Ordering Information.

For Assistance

For replacement materials, please contact the CAP within **10** calendar days of the ship date for information. **Provide your CAP number and contact information with all correspondence.** Participants in countries serviced by a designated CAP distributor should contact their distributor's customer service department.

Telephone: 800-323-4040 option 1
(Monday - Friday, 7:00 AM – 5:30 PM US Central Time)
International Participants: 001-847-832-7000 option 1

Email: contactcenter@cap.org

Website: cap.org

Address: CAP Surveys Program
325 Waukegan Road
Northfield, IL 60093-2750
USA

Hematology Instrument Master List

Deleted codes

None

New/Updated codes

None

1939 Alere/Stanbio HemoPoint H2	1272 Critspin/Statspin (glass tubes)	1410 Microhematocrit (all PCV)
1208 Clarity HbCheck	1273 Critspin/Statspin (plastic tubes)	1745 STI HemataSTAT Easy Read
1523 Clay Adams Accu-stat	2020 HemoCue B-Hb Analyzer	2030 STI HemataSTAT II/C-70
1393 Clay Adams Readacrit	1990 HemoCue Hb 201+ Systems	3430 STI UltraCrit
1392 Clay Adams Triac	3435 HemoCue Hb 301 Systems	0010 Other, specify on result form

Urinalysis and Specific Gravity Method Master List

Deleted codes
~~3238 DiaScreen Reagent Strips~~
New/Updated codes

None

1733 77 Elektronika Reagent Strips	3425 Diagnostic Test Group Clarity Urocheck Reagent Strips	2492 Roche cobas u Pack Cassette
3367 ACON Mission Reagent Strips	1238 Eiken Uropaper	1221 Roche Urisys Cassette
1534 ARKRAY Aution Sticks	2831 IRIS vChem Reagent Strips	1313 Schein Urispec
2547 Bayer Acetest	1114 Jant Accustrip URS 10 Reagent Strips	2673 Siemens Clinitek Atlas Reagent
2544 Bayer Clinitest, 2 drop	1381 Jant Accutest URS 10 Reagent Strips	1622 Siemens Clinitek Novus/Pro Cassette
2543 Bayer Clinitest, 5 drop	3422 McKesson/Teco Reagent Strips	2538 Siemens Ictotest
2922 Biorex DIAZO-CHEK	2468 Medline Urinalysis Reagent Strips	2936 Siemens Reagent Strips (All)
1798 Biorex K-Check	1095 M-N (Macherey-Nagel) Medi-Test Uryxxon test strips	3427 SMC Direct RefuAH Reagent Strips
1626 BioSys/Consult Diag/PSS Select/YD Diag Reagent Strips	2546 pH Meter	1285 Sulfosalicylic acid
2574 BTNX Rapid Response Reagent Strips	2533 Refractometer	2552 Sysmex Reagent Strips
1030 CTMI/Germaine Labs/Thermo Fisher Reagent Strips	3071 Roche Chemstrip LN	2534 Urinometer
	2937 Roche Chemstrip/Combur	0010 Other, specify on result form

Urinalysis and Specific Gravity Instrument Master List

Deleted codes
~~4106 DiaScreen 50~~
New/Updated codes

None

1935 77 Elektronika urine chemistry analyzers	1062 Germaine Labs AimStrip/CTMI CT-120 Urine Analyzer	1579 Siemens Clinitek 10 or 100
3368 ACON Mission U120, U500	1864 Germaine Labs AimStrip Urine Auto Analyzer	1845 Siemens Clinitek 50
2126 ARKRAY Aution Hybrid	1689 Jant Accustrip URS 10 Reader	1580 Siemens Clinitek 200 or 200+
2675 ARKRAY Aution JET	1096 M-N (Macherey-Nagel) Uryxxon series	2292 Siemens Clinitek 500, Advantus
2676 ARKRAY Aution MAX	3423 McKesson 107-101 UA	1578 Siemens Clinitek 2000
1509 ARKRAY/Thermo BioStar PocketChem UA	3241 Roche Chemstrip 101	1465 Siemens Clinitek Atlas
3444 BioSys/Consult Diag/PSS Select/YD Diag Urine Analyzer	1945 Roche Chemstrip Criterion, Criterion II	2681 Siemens Clinitek Novus
2573 BTNX Rapid Response series	2097 Roche Chemstrip Urine Analyzer	2532 Siemens Clinitek Status
3424 Diagnostic Test Group Clarity Urocheck 120	1171 Roche Chemstrip Mini UA	3426 SMC Direct RefuAH U120
3371 Eiken US-2100R/2200	1790 Roche cobas u411, Urisys 1800	2186 Sysmex UX-2000
3372 Eiken US-3100R/3100Rplus	2490 Roche cobas u601	2603 Teco Uritek TC-101
	1789 Roche Urisys 1100	2044 Visual
	1821 Roche Urisys 2400	0010 Other, specify on result form

Urine hCG Manufacturer Master List

Deleted codes

None

New/Updated codes

None

1735 Abbott TestPack Combo	1205 Henry Schein OneStep+ Combo	2872 Polymedco Poly stat (urine only)
2027 Abbott TestPack Plus Combo w/ or w/out OBC	1382 Henry Schein OneStep+ Urine Cassette	1372 Pulse Scientific
1736 Abbott TestPack Plus - Urine	1383 Henry Schein OneStep+ Urine Strip	1414 Quidel QuickVue+One-Step Combo
1784 ABI SureStep	1117 Immunostics Immuno hCG Detector	1475 Quidel QuickVue One-Step Combo
1782 ABI SureStep Combo	1118 Immunostics Immuno hCG Detector Combi (serum/urine)	1129 Quidel QuickVue One-Step Urine
2938 Acceava Basic II (urine)	1119 Immunostics Immuno hCG Detector Stix	1338 Quidel RapidVue (urine only)
3168 Acceava Combo II	2474 Instant Technologies iPregnancy Cassette (serum/urine)	2479 Ramco Quik-Trak (urine only)
2648 ACON One Step (urine only)	2475 Instant Technologies iPregnancy Cassette (urine only)	2480 Ramco Quik-Trak Combi Test (serum/ urine)
2055 Alere hCG Cassette (20) (Acceava)	2476 Instant Technologies iPregnancy Strip (urine only)	3137 Randox Direct
1959 Alere hCG Cassette (25) (Clearview)	1370 Jant Accustrip Value+	1800 Roche AccuStat Combo
2403 Alere hCG Combo Cassette (20/10)	3409 Jant Accutest Value+ (serum/urine)	2108 SA Scientific (serum/urine)
1299 Alere hCG Combo Cassette (25) (Clearview)	3410 Jant Accutest Value+ (urine only)	2107 SA Scientific (urine only)
2404 Alere hCG Dipstick	2981 LifeSign Status (serum/urine)	2483 SA Scientific (SAS) Ultra
1161 Beckman Coulter ICON 20	2980 LifeSign Status (urine only)	1836 Sekisui OSOM
2618 Beckman Coulter ICON 25	1180 LSC (LABSCO) PEP hCG Strip	1837 Sekisui OSOM Combo
2010 BioTron Diagnostics Foremost Combo	1178 LSC (LABSCO) PEP hCG Test	1270 Sekisui OSOM Ultra Combo
2354 Cardinal Health SP Brand Cassette	1060 Mainline Confirms (combo)	2532 Siemens Clinitek Status
2355 Cardinal Health SP Brand Combo	1468 Mainline Confirms (urine only)	2853 Signify hCG Card
2356 Cardinal Health SP Brand Dipstick	1477 McKesson Medi-Lab Performance Cassette	2854 Signify hCG Urine
2647 Cen-Med Elite Plus One Step hCG	1478 McKesson Medi-Lab Performance Combo	2855 Signify hCG Urine/Serum
2812 Clarity Test Cassette	1476 McKesson Medi-Lab Performance Dipstick	2442 SMC Direct RefuAH Combo
2813 Clarity Test Strip	2276 MediChoice Combi	2441 SMC Direct RefuAH Urine
2120 Chembio STAT-PAK	2277 MediChoice Urine Only	1438 Stanbio QuPID
1961 Clearview hCG Combo II	1091 Medline hCG Combo Cassette	1439 Stanbio QuPID Plus
1430 Consuli Diagnostics Cassette	2465 Medline hCG Urine Cassette	3425 Stanbio QuPID Plus E.R.
1432 Consult Diagnostics Combo	2466 Medline hCG Urine Dipstick	1135 Stanbio QuStick
1433 Consult Diagnostics Dipstick	3209 Meridian ImmunoCard STAT	3424 Stanbio True 20 Plus
1006 DE Healthcare Products TruView Cassette	1491 NDC Pro Advantage Urine Cassette	1467 Sure-Vue (serum/urine)
3140 DE Healthcare Products TruView One-Step	1288 NDC Pro Advantage Urine/Serum Cassette	1466 Sure-Vue (urine only)
3517 Ekla Novaplus hCG Combo	1492 NDC Pro Advantage Urine Strip	2052 Sure-Vue STAT (serum/urine)
3516 Ekla Novaplus hCG urine	3061 Polymedco Poly stat (serum/urine)	3428 Technologist Choice (serum/urine)
2860 Formosa One Sure Pregnancy Kit		1844 Trinity Biotech UniGold PT
1586 Germaine Labs AimStep		2559 Wampole UCG-Slide Test
1587 Germaine Labs AimStep Combo		3158 Wondfo One Step hCG Urine
1679 Germaine Labs AimStrip		3159 Wondfo One Step hCG Urine/Serum
		0010 Other, specify on result form




HCC2-A

2017

Results must be received at the CAP no later than midnight, Central Time by the due date below:

Specific Gravity			Alert! When reporting a high quantitative value, do not choose a response greater than your method allows.		
Method Code	010 <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/>	Instrument Code	020 <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/>	Exception Code	030 <input type="radio"/> 11 <input type="radio"/> 33
CM-01 040 <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/>			CM-02 060 <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/>		
050 <input type="radio"/> (4) Fill if less than <input type="radio"/> (5) Fill if greater than			070 <input type="radio"/> (4) Fill if less than <input type="radio"/> (5) Fill if greater than		

pH					
Alert! Roche pH Visual users: Results should be reported in whole numbers; 0.5 increments are not allowed. Alert! Roche Cassette and all other Roche pH users: Results should be reported in whole numbers from 5.0 to 9.0 with the exception of 6.5.					
Method Code	080 <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/>	Instrument Code	090 <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/>	Exception Code	100 <input type="radio"/> 11 <input type="radio"/> 33
CM-01 110 <input type="radio"/> 130 ≤ 3.5 <input type="radio"/> 136 6.5 <input type="radio"/> 131 4.0 <input type="radio"/> 137 7.0 <input type="radio"/> 132 4.5 <input type="radio"/> 138 7.5 <input type="radio"/> 133 5.0 <input type="radio"/> 144 8.0 <input type="radio"/> 134 5.5 <input type="radio"/> 177 8.5 <input type="radio"/> 135 6.0 <input type="radio"/> 145 ≥ 9.0			CM-02 120 <input type="radio"/> 130 ≤ 3.5 <input type="radio"/> 136 6.5 <input type="radio"/> 131 4.0 <input type="radio"/> 137 7.0 <input type="radio"/> 132 4.5 <input type="radio"/> 138 7.5 <input type="radio"/> 133 5.0 <input type="radio"/> 144 8.0 <input type="radio"/> 134 5.5 <input type="radio"/> 177 8.5 <input type="radio"/> 135 6.0 <input type="radio"/> 145 ≥ 9.0		





HCC2-A

2017

Results must be received at the CAP no later than midnight, Central Time by the due date below:

Handwritten: JNK

Protein			Alert! When reporting a high quantitative value, do not choose a response greater than your method allows.		
Alert! Clinitek Atlas users: For Protein values of ≥ 300, including values of 600 or ≥ 1000 mg/dL, select the bubble for 300 - 600 mg/dL or 3+.					
Reporting Method 1					
Method Code	<input type="text" value="010"/>	Instrument Code	<input type="text" value="020"/>	Exception Code	<input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 040 110 Negative <input type="radio"/> 207 (2+) <input type="radio"/> 533 Trace <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+)			CM-02 <input type="radio"/> 050 110 Negative <input type="radio"/> 207 (2+) <input type="radio"/> 533 Trace <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+)		
Reporting Method 2 (Unit of Measure - mg/dL)					
Alert! See conversion table on the kit instructions for information about equivalent conventional and SI values.					
Method Code	<input type="text" value="060"/>	Instrument Code	<input type="text" value="070"/>	Exception Code	<input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 090 110 Negative <input type="radio"/> 257 100 - 200 <input type="radio"/> 256 10 - 20 <input type="radio"/> 258 300 - 600** <input type="radio"/> 255 30 - 70 <input type="radio"/> 259 > 600 or ≥ 1000 <input type="radio"/> 109 75 - Roche users only*			CM-02 <input type="radio"/> 100 110 Negative <input type="radio"/> 257 100 - 200 <input type="radio"/> 256 10 - 20 <input type="radio"/> 258 300 - 600** <input type="radio"/> 255 30 - 70 <input type="radio"/> 259 > 600 or ≥ 1000 <input type="radio"/> 109 75 - Roche users only*		

Handwritten: WGA

* Not appropriate for Roche Urisys 1100 users.

** Clinitek Atlas Users: report ≥ 300 mg/dL, including values of 600 or ≥ 1000 mg/dL, as 300 - 600 mg/dL.

59124





HCC2-A

2017

Results must be received at the CAP no later than midnight, Central Time by the due date below.

Glucose			Alert! Select only <u>one</u> reporting method. Do not select more than one reporting area.		
Reporting Method 1					
Method Code	010	<input type="text"/>	Instrument Code	020	<input type="text"/>
			Exception Code	030	<input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 040 110 Negative <input type="radio"/> 207 (2+) <input type="radio"/> 533 Trace <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+)			CM-02 <input type="radio"/> 050 110 Negative <input type="radio"/> 207 (2+) <input type="radio"/> 533 Trace <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+)		
Reporting Method 2 (Unit of Measure = mg/dL)					
Alert! See conversion table on the kit instructions for information about equivalent conventional and SI values.					
Method Code	060	<input type="text"/>	Instrument Code	070	<input type="text"/>
			Exception Code	080	<input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 090 110 Negative <input type="radio"/> 260 30 - 100 <input type="radio"/> 296 150 - 300 <input type="radio"/> 297 500* <input type="radio"/> 261 > 500 or \geq 1000 or \geq 2000			CM-02 <input type="radio"/> 100 110 Negative <input type="radio"/> 260 30 - 100 <input type="radio"/> 296 150 - 300 <input type="radio"/> 297 500* <input type="radio"/> 261 > 500 or \geq 1000 or \geq 2000		

*Roche Glucose Visual Method: Results may be reported as 500 mg/dL.

All Other Roche Glucose Methods: Should not report 500 mg/dL. Analyzers are not able to discern results of 500 vs. 1000 mg/dL.

Reducing Substance		
Alert! If your laboratory performs reducing substances, you must perform 2 challenges, otherwise leave blank if your laboratory does not perform this test.		
Method Code	110	<input type="text"/>
		Exception Code
		120 <input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 130 120 Negative <input type="radio"/> 215 3/4%** <input type="radio"/> 124 1/4% (trace) <input type="radio"/> 218 1% <input type="radio"/> 212 1/2% <input type="radio"/> 221 \geq 2%	CM-02 <input type="radio"/> 140 120 Negative <input type="radio"/> 215 3/4%** <input type="radio"/> 124 1/4% (trace) <input type="radio"/> 218 1% <input type="radio"/> 212 1/2% <input type="radio"/> 221 \geq 2%	

** Please note 3/4% is not an appropriate response for the Bayer Clinitest, 2-drop method.





HCC2-A

2017


Results must be received at the CAP no later than midnight, Central Time by the due date below:

Ketones			Alert! Select only <u>one</u> reporting method. Do not select more than one reporting area.		
Reporting Method 1					
Method Code	010 <input type="text"/>	Instrument Code	020 <input type="text"/>	Exception Code	030 <input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 120 Negative <input type="radio"/> 366 Moderate <input type="radio"/> 533 Trace <input type="radio"/> 367 Large <input type="radio"/> 365 Small			CM-02 <input type="radio"/> 120 Negative <input type="radio"/> 366 Moderate <input type="radio"/> 533 Trace <input type="radio"/> 367 Large <input type="radio"/> 365 Small		
Reporting Method 2					
Method Code	060 <input type="text"/>	Instrument Code	070 <input type="text"/>	Exception Code	080 <input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 110 Negative <input type="radio"/> 207 (2+) <input type="radio"/> 533 Trace <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+)			CM-02 <input type="radio"/> 110 Negative <input type="radio"/> 207 (2+) <input type="radio"/> 533 Trace <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+)		
Reporting Method 3 (Unit of Measure = mg/dL)					
Alert! See conversion table on the kit instructions for information about equivalent conventional and SI values.					
Method Code	110 <input type="text"/>	Instrument Code	120 <input type="text"/>	Exception Code	130 <input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 110 Negative <input type="radio"/> 349 40 - 60 <input type="radio"/> 348 5 - 10 <input type="radio"/> 350 80 - 100* <input type="radio"/> 347 15 - 25 <input type="radio"/> 351 ≥ 150			CM-02 <input type="radio"/> 110 Negative <input type="radio"/> 349 40 - 60 <input type="radio"/> 348 5 - 10 <input type="radio"/> 350 80 - 100* <input type="radio"/> 347 15 - 25 <input type="radio"/> 351 ≥ 150		

Blank

Blank

* Roche Urinalysis instruments should not report 80 - 100 mg/dL. Analyzers are not able to discern ketones between 60 and 150.





HCC2-A

2017


Results must be received at the CAP no later than midnight, Central Time by the due date below:

Bilirubin*		Alert! Select only <u>one</u> reporting method. Do not select more than one reporting area.	
Reporting Method 1			
Method Code ⁰¹⁰	<input type="text"/>	Instrument Code ⁰²⁰	<input type="text"/>
		Exception Code ⁰³⁰ <input type="radio"/> 11 <input type="radio"/> 33	
CM-01 <input type="radio"/> 120 Negative <input type="radio"/> 369 Trace/Small <input type="radio"/> 368 Positive/Moderate <input type="radio"/> 367 Large		CM-02 <input type="radio"/> 120 Negative <input type="radio"/> 369 Trace/Small <input type="radio"/> 368 Positive/Moderate <input type="radio"/> 367 Large	
Reporting Method 2			
Method Code ⁰⁶⁰	<input type="text"/>	Instrument Code ⁰⁷⁰	<input type="text"/>
		Exception Code ⁰⁸⁰ <input type="radio"/> 11 <input type="radio"/> 33	
CM-01 <input type="radio"/> 110 Negative <input type="radio"/> 207 (2+) <input type="radio"/> 533 Trace <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+)		CM-02 <input type="radio"/> 110 Negative <input type="radio"/> 207 (2+) <input type="radio"/> 533 Trace <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+)	
Reporting Method 3 (Unit of Measure = mg/dL)			
Method Code ¹¹⁰	<input type="text"/>	Instrument Code ¹²⁰	<input type="text"/>
		Exception Code ¹³⁰ <input type="radio"/> 11 <input type="radio"/> 33	
CM-01 <input type="radio"/> 110 Negative <input type="radio"/> 376 6.0 - 10.0 <input type="radio"/> 288 0.5 - 1.0 <input type="radio"/> 377 >10.0 <input type="radio"/> 290 2.0 - 4.0		CM-02 <input type="radio"/> 110 Negative <input type="radio"/> 376 6.0 - 10.0 <input type="radio"/> 288 0.5 - 1.0 <input type="radio"/> 377 >10.0 <input type="radio"/> 290 2.0 - 4.0	

BLANK

GREEN

*Bilirubin and Urobilinogen Testing (visual read methods): Proficiency testing specimens may give a slightly atypical color development when examined visually. If this occurs, the specimen should be reported based on the intensity of the color.





HCC2-A

2017

Results must be received at the CAP no later than midnight, Central Time by the due date below:

Blood or Hemoglobin			Alert! Select only <u>one</u> reporting method. Do not select more than one reporting area.		
Reporting Method 1					
Method Code ⁰¹⁰	Instrument Code ⁰²⁰	Exception Code ⁰³⁰			
<input type="text"/>	<input type="text"/>	<input type="radio"/> 11 <input type="radio"/> 33			
CM-01 <input type="radio"/> 120 Negative <input type="radio"/> 369 Trace/Small <input type="radio"/> 368 Positive/Moderate <input type="radio"/> 361 Marked Positive/Large			CM-02 <input type="radio"/> 120 Negative <input type="radio"/> 369 Trace/Small <input type="radio"/> 368 Positive/Moderate <input type="radio"/> 361 Marked Positive/Large		
Reporting Method 2					
Method Code ⁰⁶⁰	Instrument Code ⁰⁷⁰	Exception Code ⁰⁸⁰			
<input type="text"/>	<input type="text"/>	<input type="radio"/> 11 <input type="radio"/> 33			
CM-01 <input type="radio"/> 110 Negative <input type="radio"/> 206 (3+) <input type="radio"/> 533 Trace <input type="radio"/> 205 (4+) <input type="radio"/> 208 (1+) <input type="radio"/> 203 (5+) <input type="radio"/> 207 (2+)			CM-02 <input type="radio"/> 110 Negative <input type="radio"/> 206 (3+) <input type="radio"/> 533 Trace <input type="radio"/> 205 (4+) <input type="radio"/> 208 (1+) <input type="radio"/> 203 (5+) <input type="radio"/> 207 (2+)		
Reporting Method 3 (Unit of Measure = Ery/μL)					
Method Code ¹¹⁰	Instrument Code ¹²⁰	Exception Code ¹³⁰			
<input type="text"/>	<input type="text"/>	<input type="radio"/> 11 <input type="radio"/> 33			
CM-01 <input type="radio"/> 110 Negative <input type="radio"/> 193 50 - 150 <input type="radio"/> 192 5 - 25 <input type="radio"/> 194 250			CM-02 <input type="radio"/> 110 Negative <input type="radio"/> 193 50 - 150 <input type="radio"/> 192 5 - 25 <input type="radio"/> 194 250		
Reporting Method 4 (Unit of Measure = mg/dL)					
Method Code ¹⁶⁰	Instrument Code ¹⁷⁰	Exception Code ¹⁸⁰			
<input type="text"/>	<input type="text"/>	<input type="radio"/> 11 <input type="radio"/> 33			
CM-01 <input type="radio"/> 110 Negative <input type="radio"/> 287 0.20 - 0.50 <input type="radio"/> 362 \pm 0.03 <input type="radio"/> 364 \geq 1.0 <input type="radio"/> 363 0.06 - 0.10			CM-02 <input type="radio"/> 110 Negative <input type="radio"/> 287 0.20 - 0.50 <input type="radio"/> 362 \pm 0.03 <input type="radio"/> 364 \geq 1.0 <input type="radio"/> 363 0.06 - 0.10		

BLANK

SPN

29258





HCC2-A

2017

Results must be received at the CAP no later than midnight, Central Time by the due date below:

Leukocyte Esterase			Alert! Select only <u>one</u> reporting method. Do not select more than one reporting area.		
Reporting Method 1					
Method Code	010 <input type="text"/>	Instrument Code	020 <input type="text"/>	Exception Code	030 <input type="radio"/> 11 <input type="radio"/> 33
CM-01 040 <input type="radio"/> 120 Negative <input type="radio"/> 366 Moderate <input type="radio"/> 533 Trace <input type="radio"/> 367 Large <input type="radio"/> 365 Small			CM-02 050 <input type="radio"/> 120 Negative <input type="radio"/> 366 Moderate <input checked="" type="radio"/> 533 Trace <input type="radio"/> 367 Large <input type="radio"/> 365 Small		
Reporting Method 2					
Method Code	060 <input type="text"/>	Instrument Code	070 <input type="text"/>	Exception Code	080 <input type="radio"/> 11 <input type="radio"/> 33
CM-01 090 <input type="radio"/> 110 Negative <input type="radio"/> 207 (2+) <input type="radio"/> 533 Trace <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+)			CM-02 100 <input type="radio"/> 110 Negative <input type="radio"/> 207 (2+) <input type="radio"/> 533 Trace <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+)		
Reporting Method 3 (Unit of Measure = μL)					
Method Code	110 <input type="text"/>	Instrument Code	120 <input type="text"/>	Exception Code	130 <input type="radio"/> 11 <input type="radio"/> 33
CM-01 140 <input type="radio"/> 110 Negative <input type="radio"/> 414 25 <input type="radio"/> 415 75 or 100 <input type="radio"/> 416 250 or 500			CM-02 150 <input type="radio"/> 110 Negative <input type="radio"/> 414 25 <input type="radio"/> 415 75 or 100 <input type="radio"/> 416 250 or 500		

10970





HCC2-A

2017

Results must be received at the CAP no later than midnight, Central Time by the due date below:

Nitrite		
Method Code ⁰¹⁰ <input type="text"/>	Instrument Code ⁰²⁰ <input type="text"/>	Exception Code ⁰³⁰ <input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 209 Negative <input type="radio"/> 210 Positive	CM-02 <input type="radio"/> 209 Negative <input type="radio"/> 210 Positive	

Urobilinogen* Alert! Select only <u>one</u> reporting method. Do not select more than one reporting area.		
Reporting Method 1		
Method Code ⁰⁶⁰ <input type="text"/>	Instrument Code ⁰⁷⁰ <input type="text"/>	Exception Code ⁰⁸⁰ <input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 110 Negative <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+) <input type="radio"/> 207 (2+)	CM-02 <input type="radio"/> 110 Negative <input type="radio"/> 206 (3+) <input type="radio"/> 208 (1+) <input type="radio"/> 205 (4+) <input type="radio"/> 207 (2+)	
Reporting Method 2 (Unit of Measure = mg/dL)		
Method Code ¹¹⁰ <input type="text"/>	Instrument Code ¹²⁰ <input type="text"/>	Exception Code ¹³⁰ <input type="radio"/> 11 <input type="radio"/> 33
CM-01 <input type="radio"/> 417 Negative or 0.0 - 0.2 mg/dL or < 3.2 µmol/L <input type="radio"/> 418 1.0 or < 2.0 mg/dL or 16 µmol/L or < 35 µmol/L** <input type="radio"/> 419 2.0/3.0 mg/dL or 34 µmol/L or ≥ 35 µmol/L or < 70 µmol/L <input type="radio"/> 420 4.0 or 4.0/6.0 mg/dL or ≥ 70 µmol/L or < 140 µmol/L <input type="radio"/> 421 ≥ 8.0 or 12.0 mg/dL or ≥ 140 µmol/L or 200 µmol/L	CM-02 <input type="radio"/> 417 Negative or 0.0 - 0.2 mg/dL or < 3.2 µmol/L <input type="radio"/> 418 1.0 or < 2.0 mg/dL or 16 µmol/L or < 35 µmol/L** <input type="radio"/> 419 2.0/3.0 mg/dL or 34 µmol/L or ≥ 35 µmol/L or < 70 µmol/L <input type="radio"/> 420 4.0 or 4.0/6.0 mg/dL or ≥ 70 µmol/L or < 140 µmol/L <input type="radio"/> 421 ≥ 8.0 or 12.0 mg/dL or ≥ 140 µmol/L or 200 µmol/L	

*Bilirubin and Urobilinogen Testing (visual read methods): Proficiency testing specimens may give a slightly atypical color development when examined visually. If this occurs, the specimen should be reported based on the intensity of the color.
 **Note: < 35 µmol/L for ARKRAY Aution MAX users only.





HCC2-A

2017

Results must be received at the CAP no later than midnight, Central Time by the due date below:

Urine, hCG	
Method Code ⁰¹⁰ 2052	Exception Code ⁰²⁰ <input type="radio"/> 11 <input type="radio"/> 33
CM-01 ⁰³⁰ <input type="radio"/> 120 Negative <input type="radio"/> 210 Positive	CM-02 ⁰⁴⁰ <input type="radio"/> 120 Negative <input type="radio"/> 210 Positive





Attestation/Use of Other Form

Attestation Statement

- As stated in the February 28, 1992 United States *Federal Register* under Subpart H 493-801 (b) (1), "the individual testing or examining the samples and the laboratory director must attest to the routine integration of the samples into the patient work load using the laboratory's routine methods." The laboratory director or designee and the testing personnel must sign on the result form.
- Retain a signed copy of this page in your laboratory for your records and inspection purposes. Attestation is a requirement for all CAP-accredited laboratories.
- If your laboratory requires additional space for signatures, copy this form as needed.

We, the undersigned, recognizing that some special handling may be required due to the nature of proficiency testing (PT) materials, have as closely as is practical, performed the analyses on these specimens in the same manner as regular patient specimens. We confirm that results were not shared or PT specimens referred or tested outside our CLIA identification number.

Director (or Designee) (signature required)

Survey Mailing Information (eg, CA2017)

010

030

Testing Personnel (signature required)

Testing Personnel (signature required)

Testing Personnel (signature required)

040

070

100

Use of Other

This section is provided to list methodology information not found on master lists. For online entry, you can enter only 255 characters. CAP Accreditation Program Participants: Do not use this section to make changes to your Test/Activity Menu. You can access a Test/Activity Menu Change Form through e-LAB Solutions Suite.

130

Signatures will not display when viewed online.

Customer Contact Center 800-323-4040 option 1 (domestic)
or 001-847-832-7000 option 1 (international)

AO

21022

