









IRIS iQ[®] 200 Auto-Particle Recognition (APR) Guide					
Particle	Size	Shape	Contrast	Texture	Image
Non-Squamous Epithelial Cells (NSE)	20um - 40um	Round, cuboidal, columnar-like	White for cytoplasm	Rough with cytoplasmic inclusions, vacuoles, multiple nuclei	
	Transitional epithelial cells (TREP)		Black for multiple nuclei or granules		
	12um - 20um	Polygonal, cuboidal, flat cell edge	Big contrast	Rough with cytoplasmic inclusions, vacuoles, big nuclei	
	Renal epithelial cells (REEP)		White for two-thirds cytoplasm Black for large round nuclei		
	20um - 60um	Oblong or cigar shape	White for cytoplasm	Rough, grainy cytoplasm	
	Proximal tubules cells		Black for small eccentric nucleus		
	14um - 25um	Oval to round	Less contrast	Rough, grainy cytoplasm	
Distal tubules cells	White for cytoplasm Grey for small dense centered or slightly eccentric nucleus				

IRIS iQ[®]200 Auto-Particle Recognition (APR) Guide

Particle	Size	Shape	Contrast	Texture	Image
Red Blood Cell (RBC)	6um – 8 um	Round, biconcave	Black White at center	Smooth	
White Blood Cell (WBC)	10um – 14um Neutrophils - Monocytes	Spherical rounded shape	Shimmer or glitter	Rough with granular cytoplasm	
White Blood Cell Clump (WBCC)	10um – 100um Neutrophils Clump	Spherical rounded shape	Shimmer or glitter	Rough with granular cytoplasm	
Squamous Epithelial Cells (SQEP)	40um – 60um	Thin, flagstone- shaped, flat sides	Big contrast White for cytoplasm Black for nucleus or granules	Less rough with finely granulated	



300-5129 B

Specifications subject to change without notice

Clinical Support 818-709-1244
Option 2

9172 Ejon Avenue
Chatsworth
California 91311

P 818.709.1244
F 818.700.9661
www.prolifis.com

IRIS iQ[®] 200 Auto-Particle Recognition (APR) Guide					
Particle	Size	Shape	Contrast	Texture	Image
Mucus (MUCS)	Multiple Sizes	Polygonal shape, with or without fine edge	Less contrast Depend on thread present	Smooth	
Sperm (SPRM)	Oval head 2.5um – 3.5um width, 4um – 5um length Tail <45um	Oval head Slender tail, uniform, straight	Big contrast Black for whole body	Smooth	

300-5129 B






Specifications subject to change without notice

Clinical Support 818-709-1244
Option 2

9172 Eton Avenue
Chatsworth
California 91311

P 818.709.1244
F 818.700.9661
www.proiris.com

IRIS iQ[®] 200 Auto-Particle Recognition (APR) Guide

Particle	Size	Shape	Contrast	Texture	Image
Hyaline Cast (HYAL)	Multiple Sizes	Cylindrical with parallel sides, round ends. Occasionally, tapered or tail-like end	Less contrast with few, if any, granules in matrix	Smooth or finely wrinkled	
Unclassified Cast (UNCC)	Multiple Sizes	Cylindrical with parallel sides, round ends. Occasionally, tapered or tail-like end	Large contrast with more granules in matrix	Rough with granulated matrix	
Unclassified Crystals (UNCX)	6um – 20um	Variety shape, diamond, lemon, barrel or cubes	Less contrast, majority black in colour	Smooth	
Bacteria (BACT)	8um – 100um	Rod or cocci	Less contrast with in black colour	Smooth	
Yeast	Single yeast cell: 5um - 7um	Single yeast cell: ovoid, round	Good contrast with the sharp edge	Smooth	
Budding (BYST) Pseudohyphae (HYST)	Pseudohyphae yeast cell: varying lengths	Pseudohyphae yeast cell: budding form			

9172 Eton Avenue
Chaisworth
California 91311

P 818.709.1244
F 818.700.9661
www.proiris.com


*Suggested Bacteria Grading Chart




RARE (1+)



FEW (2+)



MOD (3+)



MANY (4+)


Reminder: Verify bacteria in WBC, WBCC, Mucus and Squamous Epithelial categories.

*Suggested grading for use when searching for non-auto classified Bacteria.


Iris

9172 Eton Avenue
Chatsworth, CA 91311
A Division of IRIS International, Inc.
© Copyright 2008 Iris Diagnostics.
All rights reserved.
700-3064 Rev. A
(contains 300-4486 verbiage)

Analysis Flowchart: iChemVELOCITY





Strip Reader Module (SRM) Diagram




Analysis Process Summary

- Urine reacts with imbedded reagents on the test strip pads.
- Color change: measured by the instrument **[SRM]**.
- The test strip: transported to the waste container.
- Results: consolidated with Microscopy (as user defined) and reported.



Strip Conveyor System (SCS)

- Primary function of SCS:
 - Presents strips to the Strip Reader Module
 - Used strips are discarded into Strip Waste Container.




Specific Gravity, Color Clarity Module (CGM)

Specific Gravity

- Measured by Refractometer


Color and Clarity

- Color determined by transmitted light.
- Clarity determined by scattered light.



Strip Provider Module

- Primary functions of the Strip Provider Module (SPM):
 - Provides storage area to bulk load chemistry strips.
 - Controls humidity-ensures onboard strip stability.
 - Provides means for separating a single strip from the bulk of strips.
 - Provides a means for orienting the strips pads-face up.
 - Delivers single, properly oriented strips to the Strip Conveyor System (SCS).



Sample Probe


Primary Function of the Sample Probe:

- Mix the specimen inside the sample tube.
- Aspirates approximately 1mL of sample.
- Sample
 - Fills the specific gravity and color-clarity flowcells.
 - Approximately 10 μ L-dispensed onto each strip pad.
 - Short sample detection


Strip Loader [View of Inside Label]

Primary Function:
Module used to load
strips into the
iChemVELOCITY
analyzer


Strip Loader redesigned!
Image for strip
orientation only!



Label, Urine Chemistry
Strip Orientation



iChemVELOCITY: Internal View



iChemVELOCITY: Internal View

