

MANUAL DILUTION LOCAL PROCEDURE

Safety Message

Use universal precautions when handling patient samples.

Purpose

To provide instructions on the use of the Manual Dilution Worksheet LCS 7060-C.

Policy

- The Manual Dilution Worksheet LCS 7060-C is to be used to document the manual dilution of results for tests listed in the Regional Dilution Policy results exceeding instrument analytic measuring range (AMR).
- The lowest possible dilution will be used to obtain a result within the instrument's AMR.
- All instrument generated reports will be attached to the Manual Dilution Worksheet LCS 7060-C
- Submit completed worksheet and print-outs to manager for review.

Procedure

Step	Action																																										
1	Attach a reprint of the accession # onto the worksheet.																																										
2	Select the lowest possible dilution to obtain a result within the instrument's AMR.																																										
3	Perform the dilution with the appropriate verified diluent from the instrument's Chemistry Information Sheet manual.																																										
4.	Program the dilution into the instrument. <ul style="list-style-type: none"> • With or without ORDAC • Enter dilution factor 																																										
5.	Front-load diluted sample.																																										
6.	Print results from instrument.																																										
7.	Enter results in Remisol manually. <ul style="list-style-type: none"> • If results are greater than procedural limits in Cerner: <table border="1" data-bbox="836 1186 1169 1753"> <thead> <tr> <th>Procedure</th> <th>Result</th> </tr> </thead> <tbody> <tr><td>RBS</td><td>> 1200</td></tr> <tr><td>Acetaminophen</td><td>> 3000</td></tr> <tr><td>ALT</td><td>> 5000</td></tr> <tr><td>AST</td><td>> 5000</td></tr> <tr><td>CBZ Level</td><td>> 200.0</td></tr> <tr><td>CK Total</td><td>> 41000</td></tr> <tr><td>LDH</td><td>> 2700</td></tr> <tr><td>Magnesium</td><td>99.9</td></tr> <tr><td>PHB</td><td>> 800.0</td></tr> <tr><td>PTN</td><td>> 400.0</td></tr> <tr><td>Salicylate</td><td>> 1000.0</td></tr> <tr><td>Theo</td><td>> 400.0</td></tr> <tr><td>Valproic</td><td>> 1500.0</td></tr> <tr><td>R Vanco</td><td>> 400.0</td></tr> <tr><td>Xa Unifrac</td><td>> 1.10</td></tr> <tr><td>Digoxin NIS</td><td>> 45.0</td></tr> <tr><td>Pk Gent</td><td>> 120.0</td></tr> <tr><td>Lact Blood</td><td>> 110.0</td></tr> <tr><td>U Pro Pro_Cr</td><td>> 1500</td></tr> <tr><td>U Creat</td><td>> 800</td></tr> </tbody> </table> <ul style="list-style-type: none"> • Add in results comments section: <ul style="list-style-type: none"> ○ Verified by dilution comment e.g., CK 56,000 IU/L, Verified by dilution • Validate results 	Procedure	Result	RBS	> 1200	Acetaminophen	> 3000	ALT	> 5000	AST	> 5000	CBZ Level	> 200.0	CK Total	> 41000	LDH	> 2700	Magnesium	99.9	PHB	> 800.0	PTN	> 400.0	Salicylate	> 1000.0	Theo	> 400.0	Valproic	> 1500.0	R Vanco	> 400.0	Xa Unifrac	> 1.10	Digoxin NIS	> 45.0	Pk Gent	> 120.0	Lact Blood	> 110.0	U Pro Pro_Cr	> 1500	U Creat	> 800
Procedure	Result																																										
RBS	> 1200																																										
Acetaminophen	> 3000																																										
ALT	> 5000																																										
AST	> 5000																																										
CBZ Level	> 200.0																																										
CK Total	> 41000																																										
LDH	> 2700																																										
Magnesium	99.9																																										
PHB	> 800.0																																										
PTN	> 400.0																																										
Salicylate	> 1000.0																																										
Theo	> 400.0																																										
Valproic	> 1500.0																																										
R Vanco	> 400.0																																										
Xa Unifrac	> 1.10																																										
Digoxin NIS	> 45.0																																										
Pk Gent	> 120.0																																										
Lact Blood	> 110.0																																										
U Pro Pro_Cr	> 1500																																										
U Creat	> 800																																										

MANUAL DILUTION LOCAL PROCEDURE

Procedure

Step	Action
8.	Complete the Manual Dilution Worksheet LCS 7060-B and submit for manager review with instrument print-outs attached.

References

- Manual Dilution Policy, SCPMG QMS-0024, Rev 3
- Chemistry Information Sheets, Synchron DxC analyzer, Beckman Coulter.

Attachment

- Attachment A: Manual Dilution Worksheet LCS 7060-C

Authors

- Emmi Raymundo, CLS, MT(ASCP)
- Mina Acosta, CLS, MT(ASCP)
- Jay Raymund Castaneto, CLS, BSMT

MANUAL DILUTION LOCAL PROCEDURE

Document History Page

Change type: New, Major, Minor etc.	Changes Made to SOP – describe	Signature responsible person/date	Laboratory Operations Director Reviewed/ date	Laboratory Medical Director Reviewed/ Date	Date change implemented
New		M.Acosta 03-15-16	J.Wolf 03-21-16	S.Wirio MD 03-25-16	03-25-16
Major	<ul style="list-style-type: none"> Updated Manual Dilution Worksheet to reflect ALT and AST linear limits in Cerner. 	R.Castaneto 08-08-18	J.Wolf 08-09-18	S.Wirio MD 08-09-18	08-13-18
Major	<ul style="list-style-type: none"> Updated manual Dilution Worksheet to reflect Urine CREm linear limits in Cerner. 	<i>R. Castaneto</i> 10/17/18	<i>J. Wolf</i> 10/18/18	<i>S. Wirio</i> 11/6/18	11-6-18

MANUAL DILUTION LOCAL PROCEDURE ATTACHMENT A



Tests listed below have been approved for dilution by SCPMG Lab Systems:

ALT**	TDMs:
AST**	Acetaminophen
CK*	Digoxin
Glucose (Auto-ORDAC)	Salicylate
hCG	Tegretol/Carbamazepine
Lactic Acid	Theophylline
LDH	Valproic Acid
Magnesium	Vancomycin
Urine CREM***	
Urine Total Protein (mTP)	Anti-Xa
	STA-R Evolution

CLS Initials: _____
 Date: _____
 Analyte being Diluted: _____
 Results Reported: _____

1. Indicate dilution used by drawing a circle around dilution factor.
2. Answer question whether ORDAC was used in programming dilution on DXC.
3. Identify the diluent used.
4. Attach all instrument forms to this worksheet and place in QC review bin for manager review.

DILUTION FACTOR	SAMPLE	DILUENT
2	100 µL	100 µL
3	100 µL	200 µL
5	100 µL	400 µL
10	100 µL	500 µL
15	100 µL	1400 µL
20	50 µL	950 µL
50	10 µL	490 µL
100	10 µL	990 µL

Comments: _____

Manual ORDAC requested?

YES _____
 NO _____

Diluent used: Saline _____
 DI Water _____
 QC Material _____
 Normal Pooled plasma _____
 other: _____

NOTE: Verify diluent in Chemistry Information Sheet Manual

Cemer Procedural Limit[†]:

Procedure	Result
RBS	< 1.300
Acetaminophen	< 3000
ALT	> 5000
AST	> 5000
CBZ Level	< 260.0
CK Total	< 41000
LDH	< 2700
Magnesium	59.9
PHB	< 800.0
PTN	< 400.0
Salicylate	> 1000.0
Theo	< 460.0
Valproic	> 1500.0
R Vanco	< 400.0
Xa Unfrac	< 1.10
Oxgovan NIS	< 45.0
PK Gent	< 120.0
Lact Blood	< 110.0
U Plo Pio Cr	< 1500
U Creat	> 800

**Urine CREM max dilution is 1:2
 **ALT and AST max dilution is 1:15
 * CK may be diluted greater than 1:10
 Manual Dilution Worksheet LCS 7060-C