BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Tube

For collection, storage, and transport of urine specimens.

Not made with natural or dry rubber latex.

For *In Vitro* Diagnostic Use.

INTENDED USE

BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Tubes are single use tubes designed for collection, storage and transport of urine specimens for culture and sensitivity testing of bacteria for *in vitro* diagnostic use. The tubes are used in settings where a urine specimen is collected by a trained healthcare worker.

SUMMARY AND EXPLANATION

Bacteria quantification in urine is widely used as an aid in evaluating a patient for urinary tract infections (UTI).^{1,2,3,4} Colony forming units of 100,000 microorganisms or greater per milliliter of urine are generally considered indicative of infection.⁴

Urine frequently supports the proliferation of bacteria, which may multiply at the same rate as in the nutrient broth.⁵ Therefore, a urine sample delayed in transit or left at room temperature for an extended period of time may give an erroneous result.^{6,7}

As a means of preventing growth of the microorganisms from sources exogenous to the bladder, refrigeration or culturing within two hours of micturition is recommended.^{4,6,7} It is not always within the control of the laboratory to maintain the parameters necessary for accurate results.

BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Tubes contain a lyophilized urine maintenance formula and are evacuated to draw approximately 4.0 mL of urine. The lyophilized urine maintenance formula can maintain the bacterial population in the urine specimen for a period of up to 48 hours at room temperature at levels comparable to those urine specimens without additive, held under refrigeration for the same period of time.

PRODUCT DESCRIPTION

BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Urine Tubes are 4.0 mL draw volume, 13 x 75 mm, with a lyophilized maintenance formula and BD Hemogard™ Closure with olive shield color. A minimum fill line of 3.0 mL is indicated on the label. The tubes have sterile interiors.

Preservative

The mean concentration of the preservative in the urine sample in BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Tube is:

Boric Acid 2.63 mg/mL

Sodium Borate 3.95 mg/mL

Sodium Formate 1.65 mg/mL

PRECAUTIONS:

1. Caution should be used in handling the screw cap lid of the BD Vacutainer® Urine Collection Cup, which contains a needle under the label. Treat the screw cap of the cup as a contaminated sharp and discard in biohazard container approved for sharps disposal per your facility's recommended procedure.
2. Caution should be used in handling the transfer straw that contains a needle and should be disposed of as a contaminated sharp after use.
3. All biologic specimens and devices used to collect or store clinical specimens should be carefully handled and disposed of in accordance with the precautions recommended by the CDC, CLSI, and your facility's recommended procedures.
4. If a urine specimen contains blood, all specimen collection devices must be classified as biohazardous for handling and disposal purposes.
5. Do not squeeze cup.

CAUTION:

1. Practice standard precautions. Use gloves, gowns, eye protection, other personal protective equipment and engineering controls to protect from specimen splatter, leakage, and potential exposure to bloodborne pathogens or other infectious materials.
2. All biologic specimens and devices used to collect or store clinical specimens should be carefully handled and disposed of in accordance with the precautions recommended by the CDC, CLSI, and your facility's recommended procedures.
3. Discard all biologic samples in containers approved for their disposal per your facility's recommended procedures.
4. Transferring a sample to a tube using a syringe and needle is not recommended. Additional manipulation of sharps such as hollow bore needles increases the potential for needlestick injuries.
5. Transferring samples from syringe to an evacuated tube using non-sharps devices should be performed with caution for the reasons described below:
 - Depressing the syringe plunger during transfer can create a positive pressure, forcefully displacing the stopper and sample, causing splatter and potential exposure to infectious material.
 - Evacuated tubes are designed to draw the volume indicated.
 - Using a syringe for specimen transfer may also cause over or under filling of tubes resulting in the incorrect urine to additive ratio. Filling is complete when vacuum no longer continues to draw, though some tubes may partially fill due to plunger resistance when filled from a syringe. Consult your facility's policy and procedure regarding the use of these specimens for diagnostic evaluations.

STORAGE

BD Vacutainer® Plus C&S Boric Acid sodium Borate/Formate Tubes are stable when stored between 4 - 25°C and used prior to expiration date.

METHODS OF COLLECTION

Instructions to transfer the specimen into evacuated tube(s) from BD Vacutainer® Urine Collection Cup:

1. If the urine specimen is collected in a BD Urine Collection Cup, the patient is instructed to give the urine specimen to the healthcare professional immediately after collection. Note: The healthcare professional obtains a cup for the patient and cautions the patient not to remove the cap label to protect against needlestick from the "sharp" contained in the integrated transfer device.
 - a. Place cup upright on clean, flat surface. Container may be tipped at an angle if specimen volume is limited.
 - b. Peel back label on cap to expose the integrated transfer device.
 - c. Place evacuated tube into cavity on cap, stopper down. Advance the tube over puncture point to pierce stopper. BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Tubes should be filled first when collecting multiple tubes.
 - d. Hold tube in position until filled. NOTE: If the sample volume in the BD Vacutainer® Urine Collection Cup is insufficient to fill the urine tube by closed system transfer; follow the instructions below for filling the tubes manually or by transfer straw.
 - e. Remove tube from integrated transfer device.
 - f. For all BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Tubes, mix tubes 8 - 10 times by inversion.
 - g. Repeat steps c - f if another tube is to be filled.
 - h. Replace label over integrated transfer device cavity and reseal. Use caution to avoid contact with needle when replacing label.
2. Label tube(s) for transport to laboratory.
3. Treat the screw cap of the cup as a contaminated sharp and discard in biohazard container approved for sharps disposal as per your facility's recommended procedure.

Instructions to transfer specimen into evacuated tube(s) using BD Vacutainer® Urine Transfer Straw:

1. If the urine specimen is collected into an alternate container, the patient gives the specimen to the healthcare professional.
2. To transfer the specimen into evacuated tube(s):
 - a. Place container on a clean, flat surface.
 - b. Place tip of transfer straw into urine specimen. Container may be tipped at an angle if volume of urine is limited.
 - c. Place evacuated tube into holder, stopper down. Advance the tube over puncture point to pierce stopper.

BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Tubes should be filled first when collecting multiple tubes.

- d. Hold tube in position until filled.
 - e. Remove tube from holder.
 - f. For all BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Tubes, mix tubes 8 - 10 times by inversion.
 - g. If another tube is to be filled, leave transfer straw in container and repeat steps c - f.
3. Lift transfer straw from cup and allow specimen to drain. Discard transfer straw in a biohazard container approved for sharps disposal per your facility's recommended procedure.
4. Label tube(s) for transport to laboratory.

Instructions For Removal of BD Hemogard™ Closure:

1. Grasp the tube with one hand, placing the thumb under the BD Hemogard™ Closure (for added stability, place arm on a solid surface). With the other hand, twist the BD Hemogard™ Closure while simultaneously pushing up with the thumb of the other hand **ONLY UNTIL THE TUBE STOPPER IS LOOSENED**.
2. Move thumb away before lifting closure. **DO NOT** use thumb to push closure off the tube. If the urine specimen contains blood, an exposure hazard exists. To help prevent injury during removal, it is important that the thumb used to push upward on the closure be removed from contact with the tube as soon as the BD Hemogard™ Closure is loosened.
3. Lift closure off tube. In the unlikely event of the plastic shield separating from the rubber stopper, **DO NOT REASSEMBLE CLOSURE**. Carefully remove rubber stopper from tube.

Instructions For Reinsertion of BD Hemogard™ Closure:

1. Place closure over tube.
2. Twist and push down until stopper is fully reseated. Complete reinsertion of stopper is necessary for the closure to remain securely on the tube during handling.

Transport of Urine Specimen:

Properly label and package tubes used to transport specimen to alternate location in accordance with your facility's recommended procedure and applicable local, state and federal requirements.

Materials Provided:

BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Tubes

Equipment Required But Not Supplied for Urine Testing:

Media and supplies for bacterial culture and identification.

LIMITATIONS

1. The quantity of specimen drawn varies with altitude, ambient temperature, barometric pressure, tube age, and filling technique.
2. Urine specimen must be drawn to the minimum fill line in order to maintain the proper additive to urine ratio.
3. It is not recommended to manually fill this tube. Removal of the stopper may compromise the sterility of the tube.
4. The maintenance formula will not inactivate antibiotics.
5. The microbial load in urine from a given patient may be influenced by the time of collection and fluid intake. Symptomatic patients may have counts below 10⁵ microorganisms/mL if specimens are collected late in the day or if diuresis is occurring.³

ANALYTIC EQUIVALENCY

Evaluations of BD Vacutainer® Plus C&S Boric Acid Sodium Borate/Formate Tubes have been performed for selected bacterial strains. The BD Technical Services Department is available to provide information regarding these studies.

Whenever changing any manufacturer's collection tube type, size or storage condition for a particular laboratory assay, the laboratory personnel should review the tube manufacturer's data and their own data to establish/verify the reference range for a specific instrument/reagent system. Based on such information, the laboratory can then decide if a change is appropriate.

REFERENCES

1. Kass, EH. Asymptomatic Infections of the Urinary Tract. *Trans Assoc Amer Phys.* 1956;69:56-64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. *J Lab Clin Med.* 1958;52:463-470.
3. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. *Arch Intern Med.* 1957;100:700-714.
4. Barry AL, et al. *Laboratory diagnosis of urinary tract infections.* Cumitech 2, Washington, DC: American Society for Microbiology, 1975.
5. O'Grady F, Catell WR. Kinetics of urinary tract infections. *Br J Urol.* 1966;38:149-151.
6. Hendman R, et al. Effect of delay on culture of urine. *J Clin Microbiol.* 1976;4:102-103.
7. Jefferson N, et al. Transportation delay and the microbiological quality of clinical specimen. *Am J Clin Pathol.* 1957;64:689-693.

Thank you for using BD Vacutainer® products made only by BD. The BD Vacutainer® brand assures you of high quality in laboratory products. Remember to protect the quality of your specimens, always specify the BD Vacutainer® name.

TECHNICAL SERVICES

U.S. customers please call BD Global Technical Services at 1.800.631.0174

Canada customers please call 1.866.979.9408

For reference material: www.bd.com/vacutainer/referencematerial/

For inquiries: www.bd.com/vacutainer/contact

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1. Kass, EH. Asymptomatic Infections of the Urinary Tract. Trans Assoc Amer Phys. 1956;69:56-64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. J Lab Clin Med. 1958;52:463-470.
3. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. Arch Intern Med. 1957;100:700-714.
4. Barry AL, et al. Laboratory diagnosis of urinary tract infections. Cumitech 2, Washington, DC: American Society for Microbiology, 1975.
5. O'Grady F, Catell WR. Kinetics of urinary tract infections. Br J Urol. 1966;38:149-151.
6. Hendman R, et al. Effect of delay on culture of urine. J Clin Microbiol. 1976;4:102-103.
7. Jefferson N, et al. Transportation delay and the microbiological quality of clinical specimen. Am J Clin Pathol. 1957;64:689-693.

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LITERATUR

1. Kass, EH. Asymptomatic Infections of the Urinary Tract. Trans Assoc Amer Phys. 1956;69:56-64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. J Lab Clin Med. 1958;52:463-470.
3. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. Arch Intern Med. 1957;100:700-714.
4. Barry AL, et al. Laboratory diagnosis of urinary tract infections. Cumitech 2, Washington, DC: American Society for Microbiology, 1975.

5. O'Grady F, Catell WR. Kinetics of urinary tract infections. Br J Urol. 1966;38:149-151.
6. Hendman R, et al. Effect of delay on culture of urine. J Clin Microbiol. 1976;4:102-103.
7. Jefferson N, et al. Transportation delay and the microbiological quality of clinical specimen. Am J Clin Pathol. 1957;64:689-693.

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REFERENCIAS

1. Kass, EH. Asymptomatic Infections of the Urinary Tract. *Trans Assoc Amer Phys.* 1956;69:56-64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. *J Lab Clin Med.* 1958;52:463-470.
3. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. *Arch Intern Med.* 1957;100:700-714.
4. Barry AL, et al. Laboratory diagnosis of urinary tract infections. *Cumitech 2*, Washington, DC: American Society for Microbiology, 1975.
5. O'Grady F, Catell WR. Kinetics of urinary tract infections. *Br J Urol.* 1966;38:149-151.
6. Hendman R, et al. Effect of delay on culture of urine. *J Clin Microbiol.* 1976;4:102-103.
7. Jefferson N, et al. Transportation delay and the microbiological quality of clinical specimen. *Am J Clin Pathol.* 1957;64:689-693.

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VIITED

1. Kass, EH. Asymptomatic Infections of the Urinary Tract. Trans Assoc Amer Phys. 1956;69:56-64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. J Lab Clin Med. 1958;52:463-470.
3. Kass, EH. Bacteriuria and the diagnosis of infections of the urinary tract. Arch Intern Med. 1957;100:700-714.
4. Barry AL, et al. Laboratory diagnosis of urinary tract infections. Cumitech 2, Washington, DC: American Society for Microbiology, 1975.
5. O'Grady F, Catell WR. Kinetics of urinary tract infections. Br J Urol. 1966;38:149-151.
6. Hendman R, et al. Effect of delay on culture of urine. J Clin Microbiol. 1976;4:102-103.
7. Jefferson N, et al. Transportation delay and the microbiological quality of clinical specimen. Am J Clin Pathol. 1957;64:689-693.

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RÉFÉRENCES

1. Kass, EH. Asymptomatic Infections of the Urinary Tract. *Trans Assoc Amer Phys.* 1956;69:56-64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. *J Lab Clin Med.* 1958;52:463-470.
3. Kass, EH. Bacteriuria and the diagnosis of infections of the urinary tract. *Arch Intern Med.* 1957;100:700-714.
4. Barry AL, et al. *Laboratory diagnosis of urinary tract infections.* Cumitech 2, Washington, DC: American Society for Microbiology, 1975.
5. O'Grady F, Catell WR. Kinetics of urinary tract infections. *Br J Urol.* 1966;38:149-151.
6. Hendman R, et al. Effect of delay on culture of urine. *J Clin Microbiol.* 1976;4:102-103.
7. Hendman R, et al. Transportation delay and the microbiological quality of clinical specimen. *Am J Clin Pathol.* 1957;64:689-693.

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HEIMILDIR

1. Kass, EH. Asymptomatic Infections of the Urinary Tract. Trans Assoc Amer Phys. 1956;69:56-64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. J Lab Clin Med. 1958;52:463-470.
3. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. Arch Intern Med. 1957;100:700-714.
4. Barry AL o.fl. Laboratory diagnosis of urinary tract infections. Cumitech 2, Washington, DC: American Society for Microbiology, 1975.
5. O'Grady F, Catell WR. Kinetics of urinary tract infections. Br J Urol. 1966;38:149-151.
6. Hendman R o.fl. Effect of delay on culture of urine. J Clin Microbiol. 1976;4:102-103.
7. Jefferson N o.fl. Transportation delay and the microbiological quality of clinical specimen. Am J Clin Pathol. 1957;64:689-693.

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1. Kass, EH. Asymptomatic Infections of the Urinary Tract. Trans Assoc Amer Phys. 1956;69:56–64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. J Lab Clin Med. 1958;52:463–470.
3. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. Arch Intern Med. 1957;100:700–714.
4. Barry AL, et al. Laboratory diagnosis of urinary tract infections. Cumitech 2, Washington, DC: American Society for Microbiology, 1975.
5. O'Grady F, Catell WR. Kinetics of urinary tract infections. Br J Urol. 1966;38:149–151.
6. Hendman R, et al. Effect of delay on culture of urine. J Clin Microbiol. 1976;4:102–103.
7. Jefferson N, et al. Transportation delay and the microbiological quality of clinical specimen. Am J Clin Pathol. 1957;64:689–693.

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1. Kass, EH. Asymptomatic Infections of the Urinary Tract. *Trans Assoc Amer Phys.* 1956;69:56–64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. *J Lab Clin Med.* 1958;52:463–470.
3. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. *Arch Intern Med.* 1957;100:700–714.
4. Barry AL, et al. Laboratory diagnosis of urinary tract infections. *Cumitech 2, Washington, DC: American Society for Microbiology, 1975.*
5. O'Grady F, Catell WR. Kinetics of urinary tract infections. *Br J Urol.* 1966;38:149–151.
6. Hendman R, et al. Effect of delay on culture of urine. *J Clin Microbiol.* 1976;4:102–103.
7. Jefferson N, et al. Transportation delay and the microbiological quality of clinical specimen. *Am J Clin Pathol.* 1957;64:689–693.

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REFERANSLAR

1. Kass, EH. Asymptomatic Infections of the Urinary Tract. Trans Assoc Amer Phys. 1956;69:56-64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. J Lab Clin Med. 1958;52:463-470.
3. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. Arch Intern Med. 1957;100:700-714.
4. Barry AL ve ark. Laboratory diagnosis of urinary tract infections. Cumitech 2, Washington, DC: American Society for Microbiology, 1975.
5. O'Grady F, Catell WR. Kinetics of urinary tract infections. Br J Urol. 1966;38:149-151.
6. Hendman R ve ark. Effect of delay on culture of urine. J Clin Microbiol. 1976;4:102-103.
7. Jefferson N ve ark. Transportation delay and the microbiological quality of clinical specimen. Am J Clin Pathol. 1957;64:689-693.

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1. Kass, EH. Asymptomatic Infections of the Urinary Tract. *Trans Assoc Amer Phys.* 1956;69:56–64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. *J Lab Clin Med.* 1958;52:463–470.
3. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. *Arch Intern Med.* 1957;100:700–714.
4. Barry AL, et al. Laboratory diagnosis of urinary tract infections. *Cumitech 2*, Washington, DC: American Society for Microbiology, 1975.
5. O'Grady F, Catell WR. Kinetics of urinary tract infections. *Br J Urol.* 1966;38:149–151.
6. Hendman R, et al. Effect of delay on culture of urine. *J Clin Microbiol.* 1976;4:102–103.
7. Jefferson N, et al. Transportation delay and the microbiological quality of clinical specimen. *Am J Clin Pathol.* 1957;64:689–693.

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1. Kass, EH. Asymptomatic Infections of the Urinary Tract. *Trans Assoc Amer Phys.* 1956;69:56-64.
2. Merritt AD, Sanford, JD. Sterile voided urine culture. *J Lab Clin Med.* 1958;52:463-470.
3. Kass EH. Bacteriuria and the diagnosis of infections of the urinary tract. *Arch Intern Med.* 1957;100:700-714.
4. Barry AL, et al. Laboratory diagnosis of urinary tract infections. *Cumitech 2*, Washington, DC: American Society for Microbiology, 1975.

5. O' Grady F, Catell WR. Kinetics of urinary tract infections. Br J Urol. 1966;38:149-151.
6. Hendman R, et al. Effect of delay on culture of urine. J Clin Microbiol. 1976;4:102-103.
7. Jefferson N, et al. Transportation delay and the microbiological quality of clinical specimen. Am J Clin Pathol. 1957;64:689-693.

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