



# KAISER PERMANENTE®

<b>DOCUMENT NUMBER:</b> SCPMG-PPP-0106
<b>DOCUMENT TITLE:</b> ABO Discrepancy
<b>DOCUMENT NOTES:</b>

Added instructions for entering results in Cerner and cold agg info, read highlighted areas

<b>LOCATION:</b> SCPMG-dft	<b>VERSION:</b> 02
<b>DOC TYPE:</b> SCPMG PPP	<b>STATUS:</b> Draft

<b>EFFECTIVE DATE:</b> 02 Nov 2020	<b>NEXT REVIEW DATE:</b>
<b>RELEASE DATE:</b>	<b>EXPIRATION DATE:</b>

<b>AUTHOR:</b> G938509	<b>PREVIOUS NUMBER:</b>
<b>OWNER:</b> G938509	<b>CHANGE NUMBER:</b> SCPMG-CR-0587

## ABO Discrepancy

---

**Purpose** This procedure describes steps and testing to perform when encountering an ABO discrepancy.

- Steps to take when a historical ABO discrepancy (prior result does not agree with current testing) in patient testing is encountered.
  - The step by step procedure on how to resolve ABO discrepancies in red cell (forward) typing and plasma/serum (reverse) typing in patient testing.
  - How to enter results in Cerner when performing additional testing to resolve discrepancy.
  - Transfusion recommendation for patients with ABO discrepancies in testing.
  - Steps to take when donor ABO and/or Rh confirmation testing results are discrepant from the product label.
- 

**Policy**

- **Tube method** is the routine method for ABO discrepancy testing.
- When a forward/reverse type discrepancy is encountered, the discrepant results must be recorded (performed) in Cerner.
  - After resolution the final interpretation of the ABO group can be verified.
- **Before the patient's historical blood type is changed, there must be two matching blood typing results from two different specimen collection events and/or consultation with manager/designee.**
  - Refer to procedure "*Hematopoietic Cell Transplantation and Adoptive Cell Transfer: Patient Management*" for additional considerations.
- **In the event of massive transfusion, the historical blood type should not be changed.**
- **In the event of an unresolved ABO discrepancy and transfusion support is required; type O Red Blood Cells and/or type AB plasma will be prepared and/or dispensed until discrepancy resolved. An emergency waiver is required due to incomplete testing.**
- Cases may be sent out to the Immunohematology Reference Lab (IRL) for resolution.

---

## ABO Discrepancy, Continued

---

### Definitions

**Patient Historical Discrepancy-** An ABO discrepancy exists when prior ABO results do not agree with current testing ABO result.

**Patient Testing Discrepancy -** An ABO testing discrepancy exists when:

- The expected reactions in the forward typing (red cells) do not agree with the expected reactions in the reverse typing (serum/plasma).
- Unexpected weak or mixed-field reaction that cannot be explained with patient condition or history.

**Product RBC Testing Discrepancy-** A discrepancy exists when the unit ABO and/or Rh confirmation test result does not agree with the unit ABORh face label.

### Limitations and procedure notes

- Significant variations in red blood cell suspension (<3% or >5%) may result in false- positive or false negative reactions.
- Reverse grouping tests should never be incubated at 37°C per current manufacturer's instructions for use.
- Anomalous results may be caused by:
  - Incompletely clotted serum
  - Fibrin or particulate matter in serum or plasma or red cells that stick to the sides of the microtube
  - Contaminated specimen
  - Contaminated reagent
  - Wrong/void reagent
  - Specimen error

### Equipment and reagents

- 12 X 75 mm test tubes
- Uniform drop transfer pipettes
- Isotonic saline
- ABO typing reagents
- Anti-A1 lectin
- A2 cells
- Anti-A,B
- Cell washer
- Centrifuge
- Agglutination viewer

## ABO Discrepancy, Continued

- Specimen**
- EDTA anticoagulated plasma (preferred)
  - Serum, completely clotted

**Procedure-  
 Resolving  
 Historical  
 Discrepancy**

Follow the steps below to resolve a patient historical discrepancy in ABO/Rh typing.

NOTE: Do not verify result in Cerner until investigation completed.

Step	Action				
1.	Perform a clerical check - Verify patient information on tube, paper, electronic records, and repeat testing.				
2.	<p>Review patient history (Cerner, Health Connect, contact nursing) for evidence of recent transfusion, stem cell/bone marrow transplant.</p> <ul style="list-style-type: none"> <li>• If it is discovered that there was a prior admission/treatment at a non-SCAL KP facility it may be advisable to contact the facility for patient testing and transfusion history.</li> </ul> <p>Consult with manager/designee if further instruction is needed.</p>				
3.	Repeat initial testing with washed patient cells using manual tube method as false positive forward typing reactions may be easily resolved.				
4.	<table border="1"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>Discrepancy resolved (current results agree with prior results)</td> <td>Follow separate procedure "Result Entry: Resulting Patient Tests in Cerner Millennium" for result entry.</td> </tr> </tbody> </table>	If...	Then...	Discrepancy resolved (current results agree with prior results)	Follow separate procedure "Result Entry: Resulting Patient Tests in Cerner Millennium" for result entry.
If...	Then...				
Discrepancy resolved (current results agree with prior results)	Follow separate procedure "Result Entry: Resulting Patient Tests in Cerner Millennium" for result entry.				

## ABO Discrepancy, Continued

**Resolving  
 Historical  
 Discrepancy  
 con't**

Step	Action						
5.	<p data-bbox="646 352 1242 388">Historical ABO Discrepancy is not resolved...</p> <table border="1" data-bbox="646 388 1385 1421"> <thead> <tr> <th data-bbox="646 388 857 424">If...</th> <th data-bbox="857 388 1385 424">Then...</th> </tr> </thead> <tbody> <tr> <td data-bbox="646 424 857 571">Previous sample is still available</td> <td data-bbox="857 424 1385 571">Repeat ABO/Rh testing on that sample.</td> </tr> <tr> <td data-bbox="646 571 857 1421">Previous sample is not available</td> <td data-bbox="857 571 1385 1421"> <p data-bbox="868 577 1372 724">Review the results on previous sample (if available in system) to ensure forward, reverse, and interpretation agreed.</p> <p data-bbox="868 756 1356 861">Request provider to place an ABORh order to obtain a second sample to retest.</p> <ul data-bbox="933 871 1356 1018" style="list-style-type: none"> <li>For outpatient specimen, inform the provider's office to ask the patient to return for new specimen collection.</li> </ul> <p data-bbox="868 1050 925 1081"><b>OR</b></p> <p data-bbox="868 1123 1339 1270">Obtain an additional sample (example: Hematology specimen) from a different date or time of draw for testing.</p> <ul data-bbox="933 1270 1356 1417" style="list-style-type: none"> <li>Order Double Check to retest.</li> <li>If possible, second sample should be drawn by a different phlebotomist.</li> </ul> </td> </tr> </tbody> </table>	If...	Then...	Previous sample is still available	Repeat ABO/Rh testing on that sample.	Previous sample is not available	<p data-bbox="868 577 1372 724">Review the results on previous sample (if available in system) to ensure forward, reverse, and interpretation agreed.</p> <p data-bbox="868 756 1356 861">Request provider to place an ABORh order to obtain a second sample to retest.</p> <ul data-bbox="933 871 1356 1018" style="list-style-type: none"> <li>For outpatient specimen, inform the provider's office to ask the patient to return for new specimen collection.</li> </ul> <p data-bbox="868 1050 925 1081"><b>OR</b></p> <p data-bbox="868 1123 1339 1270">Obtain an additional sample (example: Hematology specimen) from a different date or time of draw for testing.</p> <ul data-bbox="933 1270 1356 1417" style="list-style-type: none"> <li>Order Double Check to retest.</li> <li>If possible, second sample should be drawn by a different phlebotomist.</li> </ul>
If...	Then...						
Previous sample is still available	Repeat ABO/Rh testing on that sample.						
Previous sample is not available	<p data-bbox="868 577 1372 724">Review the results on previous sample (if available in system) to ensure forward, reverse, and interpretation agreed.</p> <p data-bbox="868 756 1356 861">Request provider to place an ABORh order to obtain a second sample to retest.</p> <ul data-bbox="933 871 1356 1018" style="list-style-type: none"> <li>For outpatient specimen, inform the provider's office to ask the patient to return for new specimen collection.</li> </ul> <p data-bbox="868 1050 925 1081"><b>OR</b></p> <p data-bbox="868 1123 1339 1270">Obtain an additional sample (example: Hematology specimen) from a different date or time of draw for testing.</p> <ul data-bbox="933 1270 1356 1417" style="list-style-type: none"> <li>Order Double Check to retest.</li> <li>If possible, second sample should be drawn by a different phlebotomist.</li> </ul>						

## ABO Discrepancy, Continued

**Resolving  
 Historical  
 Discrepancy  
 con't**

Step	Action						
6.	<p>Retest new sample. Notify manager/designee if unable to obtain a new sample.</p> <table border="1" data-bbox="646 422 1377 1157"> <thead> <tr> <th data-bbox="646 422 857 464">If...</th> <th data-bbox="857 422 1377 464">Then...</th> </tr> </thead> <tbody> <tr> <td data-bbox="646 464 857 831">Repeat specimen sample result matches historical result but does not match 1<sup>st</sup> current sample.</td> <td data-bbox="857 464 1377 831"> <ul style="list-style-type: none"> <li>• Select the correct interpretation in Cerner.</li> <li>• Print unverified results of initial sample and cancel accession (Misabeled Specimen, Test Not Performed)</li> <li>• Submit QIM for possible mis-identified sample-attach printout of unverified results</li> </ul> </td> </tr> <tr> <td data-bbox="646 831 857 1157">New sample result does NOT match the historical result but agrees with the 1<sup>st</sup> (current) sample.</td> <td data-bbox="857 831 1377 1157"> <ul style="list-style-type: none"> <li>• Verify results in Cerner. Refer to ABO discrepancy flowchart in procedure "Result Entry: Resulting Patient Tests in Cerner Millennium"</li> <li>• Submit QIM for further investigation.</li> </ul> </td> </tr> </tbody> </table>	If...	Then...	Repeat specimen sample result matches historical result but does not match 1 <sup>st</sup> current sample.	<ul style="list-style-type: none"> <li>• Select the correct interpretation in Cerner.</li> <li>• Print unverified results of initial sample and cancel accession (Misabeled Specimen, Test Not Performed)</li> <li>• Submit QIM for possible mis-identified sample-attach printout of unverified results</li> </ul>	New sample result does NOT match the historical result but agrees with the 1 <sup>st</sup> (current) sample.	<ul style="list-style-type: none"> <li>• Verify results in Cerner. Refer to ABO discrepancy flowchart in procedure "Result Entry: Resulting Patient Tests in Cerner Millennium"</li> <li>• Submit QIM for further investigation.</li> </ul>
If...	Then...						
Repeat specimen sample result matches historical result but does not match 1 <sup>st</sup> current sample.	<ul style="list-style-type: none"> <li>• Select the correct interpretation in Cerner.</li> <li>• Print unverified results of initial sample and cancel accession (Misabeled Specimen, Test Not Performed)</li> <li>• Submit QIM for possible mis-identified sample-attach printout of unverified results</li> </ul>						
New sample result does NOT match the historical result but agrees with the 1 <sup>st</sup> (current) sample.	<ul style="list-style-type: none"> <li>• Verify results in Cerner. Refer to ABO discrepancy flowchart in procedure "Result Entry: Resulting Patient Tests in Cerner Millennium"</li> <li>• Submit QIM for further investigation.</li> </ul>						
7.	If Rh discrepancy not resolved, follow separate procedure "Rh Testing".						

## ABO Discrepancy, Continued

**Procedure-  
 Resolving  
 Forward  
 typing  
 Discrepancy**

**Problem with forward/red cell typing: Unexpected, extra, mixed field, missing, weaker than expected reactions**

(Expected reaction is less than 2+ reaction grading)

Step	Action										
1.	Repeat initial testing with washed patient cells using manual tube method as false positive forward typing reactions may be easily resolved.										
2.	<p>Review patient history (Cerner, Health Connect, contact nursing) for evidence of massive transfusion, stem cell/bone marrow transplant.</p> <ul style="list-style-type: none"> <li>If it is discovered that there was a prior admission/treatment at a non-SCAL KP facility it may be advisable to contact the facility for patient testing and transfusion history.</li> </ul> <p>Refer to Attachment B “Possible Causes of ABO Typing Discrepancies” and/or consult with manager/designee if further instruction is needed.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 50%;">If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td>Patient received an ABO-mismatched stem cell or bone marrow transplant.</td> <td>                     Incubate forward typing tubes at RT for 30 minutes to enhance reactivity.                     <ul style="list-style-type: none"> <li>Mixed field reactions may be observed.</li> <li>Enter a result note in Cerner “Forward type RT 30”</li> </ul> </td> </tr> <tr> <td>Cold agglutinin suspected</td> <td>                     Forward typing-Wash cells with warmed saline 4-5 times.                     <ul style="list-style-type: none"> <li>Retype using warmed washed red cells.</li> <li>Enter a result note in Cerner “Forward type warm washed cells”</li> </ul> </td> </tr> <tr> <td>Subgroup of type A suspected</td> <td>Perform A1 lectin testing using <b>Anti-A1 lectin</b> procedure below.</td> </tr> <tr> <td>Other subgroup suspected</td> <td>                     Refer to Attachment A: “Serologic Reactions Observed in A and B subgroups”                     <ul style="list-style-type: none"> <li>Consider send out to IRL</li> </ul> </td> </tr> </tbody> </table> <p>Result notes entered above may also be entered into PPI.</p>	If...	Then...	Patient received an ABO-mismatched stem cell or bone marrow transplant.	Incubate forward typing tubes at RT for 30 minutes to enhance reactivity. <ul style="list-style-type: none"> <li>Mixed field reactions may be observed.</li> <li>Enter a result note in Cerner “Forward type RT 30”</li> </ul>	Cold agglutinin suspected	Forward typing-Wash cells with warmed saline 4-5 times. <ul style="list-style-type: none"> <li>Retype using warmed washed red cells.</li> <li>Enter a result note in Cerner “Forward type warm washed cells”</li> </ul>	Subgroup of type A suspected	Perform A1 lectin testing using <b>Anti-A1 lectin</b> procedure below.	Other subgroup suspected	Refer to Attachment A: “Serologic Reactions Observed in A and B subgroups” <ul style="list-style-type: none"> <li>Consider send out to IRL</li> </ul>
If...	Then...										
Patient received an ABO-mismatched stem cell or bone marrow transplant.	Incubate forward typing tubes at RT for 30 minutes to enhance reactivity. <ul style="list-style-type: none"> <li>Mixed field reactions may be observed.</li> <li>Enter a result note in Cerner “Forward type RT 30”</li> </ul>										
Cold agglutinin suspected	Forward typing-Wash cells with warmed saline 4-5 times. <ul style="list-style-type: none"> <li>Retype using warmed washed red cells.</li> <li>Enter a result note in Cerner “Forward type warm washed cells”</li> </ul>										
Subgroup of type A suspected	Perform A1 lectin testing using <b>Anti-A1 lectin</b> procedure below.										
Other subgroup suspected	Refer to Attachment A: “Serologic Reactions Observed in A and B subgroups” <ul style="list-style-type: none"> <li>Consider send out to IRL</li> </ul>										

## ABO Discrepancy, Continued

**Procedure-  
 Resolving  
 Reverse  
 typing  
 Discrepancy**

**Problem with reverse cell typing: Unexpected, extra (due to cold antibodies or other issues), mixed field, missing, weaker than expected reactions.**

Step	Action										
1.	<p>Review patient history. Refer to Attachment B “Possible Causes of ABO Typing Discrepancies” and/or consult with manager/designee if further instruction is needed.</p> <table border="1" data-bbox="649 646 1380 1780"> <thead> <tr> <th data-bbox="649 646 857 682">If...</th> <th data-bbox="857 646 1380 682">Then...</th> </tr> </thead> <tbody> <tr> <td data-bbox="649 682 857 865">Patient is elderly, an infant, or immuno-compromised</td> <td data-bbox="857 682 1380 865">Follow <b>reverse type incubation</b> procedure below.</td> </tr> <tr> <td data-bbox="649 865 857 1159">Patient has disease state such as multiple myeloma, Waldentroms macro-globulinemia,</td> <td data-bbox="857 865 1380 1159">Read A, B, O, or A2 test tubes microscopically for rouleaux appearance (Stack of coins). If rouleaux is observed, perform <b>saline replacement</b> procedure below.</td> </tr> <tr> <td data-bbox="649 1159 857 1600">Patient suspected of being subgroup of A or AB, A2 or A2B with anti-A1 (Forward types as A, reverse type is positive with A1 cells)</td> <td data-bbox="857 1159 1380 1600">Test patient rbc's for A2 (optional). See <b>A2 subgroup</b> testing procedure above.  Refer to Attachment A: “Serologic Reactions Observed in A and B subgroups”</td> </tr> <tr> <td data-bbox="649 1600 857 1780">Patient has unexpected reactivity</td> <td data-bbox="857 1600 1380 1780">Probable cold antibody causing discrepancy. See procedure “Antibody Identification- How to interpret test results” See <b>cold antibody</b> section below</td> </tr> </tbody> </table>	If...	Then...	Patient is elderly, an infant, or immuno-compromised	Follow <b>reverse type incubation</b> procedure below.	Patient has disease state such as multiple myeloma, Waldentroms macro-globulinemia,	Read A, B, O, or A2 test tubes microscopically for rouleaux appearance (Stack of coins). If rouleaux is observed, perform <b>saline replacement</b> procedure below.	Patient suspected of being subgroup of A or AB, A2 or A2B with anti-A1 (Forward types as A, reverse type is positive with A1 cells)	Test patient rbc's for A2 (optional). See <b>A2 subgroup</b> testing procedure above.  Refer to Attachment A: “Serologic Reactions Observed in A and B subgroups”	Patient has unexpected reactivity	Probable cold antibody causing discrepancy. See procedure “Antibody Identification- How to interpret test results” See <b>cold antibody</b> section below
If...	Then...										
Patient is elderly, an infant, or immuno-compromised	Follow <b>reverse type incubation</b> procedure below.										
Patient has disease state such as multiple myeloma, Waldentroms macro-globulinemia,	Read A, B, O, or A2 test tubes microscopically for rouleaux appearance (Stack of coins). If rouleaux is observed, perform <b>saline replacement</b> procedure below.										
Patient suspected of being subgroup of A or AB, A2 or A2B with anti-A1 (Forward types as A, reverse type is positive with A1 cells)	Test patient rbc's for A2 (optional). See <b>A2 subgroup</b> testing procedure above.  Refer to Attachment A: “Serologic Reactions Observed in A and B subgroups”										
Patient has unexpected reactivity	Probable cold antibody causing discrepancy. See procedure “Antibody Identification- How to interpret test results” See <b>cold antibody</b> section below										
2.	Proceed to appropriate testing procedure.										



## ABO Discrepancy, Continued

**Procedure- Reverse type incubation** Follow the steps below to perform reverse type incubation to enhance weak reactions.

Step	Action						
1.	For weak or missing reactivity in reverse typing:  In Cerner (DOE) add the ABO Discrepancy test to the accession for the ABO/Rh test to record results for RT and/or 4°C incubations.						
2.	Label 2 test tubes for A1 cells, B cell.						
3.	Add two drops of patient plasma to each test tube.						
4.	Add one drop of appropriate reagent red cells to appropriate tube.						
5.	Incubate the A1 and B tube at room temperature for 15-30 minutes.						
6.	Mix and centrifuge.						
7.	Re-suspend and read macroscopically. <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 50%;">If...</th> <th style="width: 50%;">Then...</th> </tr> </thead> <tbody> <tr> <td>The reaction is <b>now positive</b> with A and/or B cells.</td> <td>Report result in the ABO discrepancy test with result note of &lt;reverse typing incubated at RT&gt;. Forward and reverse typing agree. Select correct interpretation.</td> </tr> <tr> <td>The reaction is <b>still negative</b> with A and/or B cells.</td> <td>Incubate the A1 and/or B tube at 4°C for 15-30 minutes.                             <ul style="list-style-type: none"> <li>• Mix, centrifuge re-suspend, and read macroscopically.</li> </ul>                             Report result in the ABO discrepancy test with result note of &lt;reverse typing incubated at 4C&gt;. Forward and reverse typing agree. Select correct interpretation.                         </td> </tr> </tbody> </table>	If...	Then...	The reaction is <b>now positive</b> with A and/or B cells.	Report result in the ABO discrepancy test with result note of <reverse typing incubated at RT>. Forward and reverse typing agree. Select correct interpretation.	The reaction is <b>still negative</b> with A and/or B cells.	Incubate the A1 and/or B tube at 4°C for 15-30 minutes. <ul style="list-style-type: none"> <li>• Mix, centrifuge re-suspend, and read macroscopically.</li> </ul> Report result in the ABO discrepancy test with result note of <reverse typing incubated at 4C>. Forward and reverse typing agree. Select correct interpretation.
If...	Then...						
The reaction is <b>now positive</b> with A and/or B cells.	Report result in the ABO discrepancy test with result note of <reverse typing incubated at RT>. Forward and reverse typing agree. Select correct interpretation.						
The reaction is <b>still negative</b> with A and/or B cells.	Incubate the A1 and/or B tube at 4°C for 15-30 minutes. <ul style="list-style-type: none"> <li>• Mix, centrifuge re-suspend, and read macroscopically.</li> </ul> Report result in the ABO discrepancy test with result note of <reverse typing incubated at 4C>. Forward and reverse typing agree. Select correct interpretation.						

## ABO Discrepancy, Continued

<p><b>7.</b> (continued)</p>	<p>Unable to enhance reactions</p>	<p>Forward and reverse typing do <b>NOT</b> agree.</p> <table border="1" data-bbox="846 352 1339 831"> <thead> <tr> <th data-bbox="846 352 1092 390">If</th> <th data-bbox="1092 352 1339 390">Then</th> </tr> </thead> <tbody> <tr> <td data-bbox="846 390 1092 684">Blood products are required to be dispensed</td> <td data-bbox="1092 390 1339 684">Result as "Indeterminate"  Add result comment with forward typing result.</td> </tr> <tr> <td data-bbox="846 684 1092 831">No Blood products are required to be dispensed</td> <td data-bbox="1092 684 1339 831">Leave test as pending</td> </tr> </tbody> </table> <p>Refer to appropriate SOPs for product blood type selection and crossmatch testing.</p> <p><b>Inform manager/designee of "Indeterminate" typing or pending interpretation.</b></p>	If	Then	Blood products are required to be dispensed	Result as "Indeterminate"  Add result comment with forward typing result.	No Blood products are required to be dispensed	Leave test as pending																
If	Then																							
Blood products are required to be dispensed	Result as "Indeterminate"  Add result comment with forward typing result.																							
No Blood products are required to be dispensed	Leave test as pending																							
<p><b>8.</b></p>	<p>Enter patient result(s) for testing completed (with result notes above) and interpretation in Cerner for ABO Discrepancy test (procedure).</p> <table border="1" data-bbox="613 1392 1360 1493"> <thead> <tr> <th>Procedure</th> <th>ID</th> <th>ABO/Rh</th> <th>Comment</th> <th>Anti-A1 Lectin</th> <th>A2 Cells</th> <th>A1 Cells 4C-RT</th> <th>B Cells 4C-RT</th> <th>D Cells</th> <th>Rouleaux?</th> <th>ABO Discrp Intpr</th> </tr> </thead> <tbody> <tr> <td>ABO Discrepancy</td> <td>RIV, TEST PATIENT</td> <td>O POS</td> <td>BP</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>		Procedure	ID	ABO/Rh	Comment	Anti-A1 Lectin	A2 Cells	A1 Cells 4C-RT	B Cells 4C-RT	D Cells	Rouleaux?	ABO Discrp Intpr	ABO Discrepancy	RIV, TEST PATIENT	O POS	BP							
Procedure	ID	ABO/Rh	Comment	Anti-A1 Lectin	A2 Cells	A1 Cells 4C-RT	B Cells 4C-RT	D Cells	Rouleaux?	ABO Discrp Intpr														
ABO Discrepancy	RIV, TEST PATIENT	O POS	BP																					
<p><b>9.</b></p>	<p>Refer to <b>Cerner Entry</b> and <b>Final Interpretation</b> sections below.</p>																							


## ABO Discrepancy, Continued

**Procedure- Anti-A1 Lectin testing** Follow the steps below to perform A1 lectin testing using Anti-A1 lectin reagent and patient cells. **(Perform when unexpected reactivity with A1 cells) Perform and record QC on infrequent QC form.**

Step	Action																				
1.	In Cerner (DOE) add the ABO Discrepancy test to the accession for the ABO/Rh test.																				
2.	Label 3 test tubes for positive control, negative control, and patient testing using anti-A1 lectin reagent. <ul style="list-style-type: none"> <li>• Positive control: Reagent A1 red blood cells.</li> <li>• Negative control: Reagent B red blood cells or 3% selectogen/panel cells (Type O cells).</li> </ul> <table border="1" style="width: 100%; margin-top: 10px;"> <tr> <td style="width: 33%;">Tube 1 Positive Control</td> <td style="width: 33%;">Reagent A1 Blood Cells</td> <td style="width: 33%;">Anti A1 lectin</td> </tr> <tr> <td>Tube 2 Negative Control</td> <td>Reagent B or O Blood Cells</td> <td>Anti A1 lectin</td> </tr> <tr> <td>Tube 3 Patient cells (approximately 3-5%)</td> <td>Patient Cells</td> <td>Anti A1 lectin</td> </tr> </table>	Tube 1 Positive Control	Reagent A1 Blood Cells	Anti A1 lectin	Tube 2 Negative Control	Reagent B or O Blood Cells	Anti A1 lectin	Tube 3 Patient cells (approximately 3-5%)	Patient Cells	Anti A1 lectin											
Tube 1 Positive Control	Reagent A1 Blood Cells	Anti A1 lectin																			
Tube 2 Negative Control	Reagent B or O Blood Cells	Anti A1 lectin																			
Tube 3 Patient cells (approximately 3-5%)	Patient Cells	Anti A1 lectin																			
3.	Add one drop of reagent anti-A1 Lectin to each test tube.																				
4.	Add one drop of reagent red cell to QC tubes.																				
5.	Add one drop of patient 3-5% washed red cell suspension to patient test tube.																				
6.	Mix and centrifuge. Examine for agglutination macroscopically <ul style="list-style-type: none"> <li>• Manually record the required quality control documentation.</li> </ul> Record patient results in Cerner <b>Anti-A1 Lectin</b> result field found in the ABO Discrepancy test (procedure)																				
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 10%;">Procedure</td> <td style="width: 10%;">ID</td> <td style="width: 10%;">ABO/Rh</td> <td style="width: 10%;">Comment</td> <td style="width: 10%;">Anti-A1 Lectin</td> <td style="width: 10%;">A2 Cells</td> <td style="width: 10%;">A1 Cells 4C-RT</td> <td style="width: 10%;">B Cells 4C-RT</td> <td style="width: 10%;">O Cells</td> <td style="width: 10%;">Route</td> </tr> <tr> <td>ABO Discrepancy</td> <td>RIV, TEST PATIENT</td> <td>O POS</td> <td>BP</td> <td style="border: 2px solid red;"></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </table> <p style="text-align: center; margin-top: 5px;">QC group: Not specified</p>		Procedure	ID	ABO/Rh	Comment	Anti-A1 Lectin	A2 Cells	A1 Cells 4C-RT	B Cells 4C-RT	O Cells	Route	ABO Discrepancy	RIV, TEST PATIENT	O POS	BP						
Procedure	ID	ABO/Rh	Comment	Anti-A1 Lectin	A2 Cells	A1 Cells 4C-RT	B Cells 4C-RT	O Cells	Route												
ABO Discrepancy	RIV, TEST PATIENT	O POS	BP																		

## ABO Discrepancy, Continued

**Anti-A1**  
**Lectin testing**  
**con't**

Step	Action						
7.	<p>Interpretation</p> <table border="1" style="width: 100%;"> <thead> <tr> <th style="width: 50%;">If...</th> <th style="width: 50%;">Then...</th> </tr> </thead> <tbody> <tr> <td>The patient red cells tested positive with anti-A1 lectin</td> <td>The patient is group A1 or A1B. Continue discrepancy investigation. Refer to procedure “<i>Antibody Identification - How to interpret test results</i>”</td> </tr> <tr> <td>The patient red cells tested negative with anti-A1 lectin</td> <td>                     Patient is subgroup of A, most likely A2 or A2B.                     <ul style="list-style-type: none"> <li>Assume unexpected reactivity is due to Anti-A1.</li> <li>Test patient plasma against A2 cells (optional). See <b>A2 cells testing</b> below.</li> </ul> </td> </tr> </tbody> </table> <p>Enter ABO Discrepancy Interpretation in Cerner</p> 	If...	Then...	The patient red cells tested positive with anti-A1 lectin	The patient is group A1 or A1B. Continue discrepancy investigation. Refer to procedure “ <i>Antibody Identification - How to interpret test results</i> ”	The patient red cells tested negative with anti-A1 lectin	Patient is subgroup of A, most likely A2 or A2B. <ul style="list-style-type: none"> <li>Assume unexpected reactivity is due to Anti-A1.</li> <li>Test patient plasma against A2 cells (optional). See <b>A2 cells testing</b> below.</li> </ul>
If...	Then...						
The patient red cells tested positive with anti-A1 lectin	The patient is group A1 or A1B. Continue discrepancy investigation. Refer to procedure “ <i>Antibody Identification - How to interpret test results</i> ”						
The patient red cells tested negative with anti-A1 lectin	Patient is subgroup of A, most likely A2 or A2B. <ul style="list-style-type: none"> <li>Assume unexpected reactivity is due to Anti-A1.</li> <li>Test patient plasma against A2 cells (optional). See <b>A2 cells testing</b> below.</li> </ul>						
8.	Refer to <b>Cerner Entry</b> and <b>Final Interpretation</b> sections below.						
9.	<p>Additional Cerner documentation when Anti-A1 detected. (Patient cells are negative with anti-A1 lectin).</p> <p>In Cerner, launch DOE, Task, Accession Add On, Add Antibody ID. Result Anti-A1.</p> <table border="1" style="width: 100%;"> <thead> <tr> <th style="width: 50%;">If patient is</th> <th style="width: 50%;">Then</th> </tr> </thead> <tbody> <tr> <td><b>A subgroup</b> (Forward typing for Anti-B is negative)</td> <td>Add Transfusion Requirement of “Give O RBCs”</td> </tr> <tr> <td><b>A subgroup B</b> (Forward typing for Anti-B is positive)</td> <td>Consult with manager/designee as patient may be transfused with either <b>Type O</b> or Type B RBCs and enter appropriate Transfusion Requirement/patient comment.</td> </tr> </tbody> </table>	If patient is	Then	<b>A subgroup</b> (Forward typing for Anti-B is negative)	Add Transfusion Requirement of “Give O RBCs”	<b>A subgroup B</b> (Forward typing for Anti-B is positive)	Consult with manager/designee as patient may be transfused with either <b>Type O</b> or Type B RBCs and enter appropriate Transfusion Requirement/patient comment.
If patient is	Then						
<b>A subgroup</b> (Forward typing for Anti-B is negative)	Add Transfusion Requirement of “Give O RBCs”						
<b>A subgroup B</b> (Forward typing for Anti-B is positive)	Consult with manager/designee as patient may be transfused with either <b>Type O</b> or Type B RBCs and enter appropriate Transfusion Requirement/patient comment.						

Kaiser Permanente  
Medical Care Program  
California Division South

SCPMG Laboratory Systems  
RL Transfusion Service  
Procedure

## ABO Discrepancy, Continued

	•	
--	---	--

## ABO Discrepancy, Continued

**Procedure-**  
**A2 cells**  
testing

Follow the steps below to perform A2 cell testing using reagent A2 cells and patient plasma, this is additional confirmation of presence of Anti-A1.


**Perform and record QC on infrequent QC form.**

NOTE: **Optional** test as A2 cells may not be available at all medical centers.

Step	Action																									
1.	In Cerner (DOE) add the ABO Discrepancy test to the accession for the ABO/Rh test to record results for A2 cells if not already done when testing with Anti-A1 lectin																									
2.	Label 2 tubes for QC or A2 cell reagent. <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 30%;">Tube 1 Positive Control</td> <td style="width: 30%;">Anti-A reagent</td> <td style="width: 40%;">A2 cells</td> </tr> <tr> <td>Tube 2 Negative Control</td> <td>Anti-B reagent</td> <td>A2 cells</td> </tr> </table>	Tube 1 Positive Control	Anti-A reagent	A2 cells	Tube 2 Negative Control	Anti-B reagent	A2 cells																			
Tube 1 Positive Control	Anti-A reagent	A2 cells																								
Tube 2 Negative Control	Anti-B reagent	A2 cells																								
3.	Label 4 test tubes for A1 cells, A2 cells, O cells, and, auto control (patient cells). (Testing patient plasma)																									
4.	Add two drops of patient plasma/reagent control to each test tube.																									
5.	Add one drop of appropriate cells to each test tube.																									
6.	Mix and centrifuge.																									
7.	Re-suspend and read macroscopically. Record QC results. For patient testing refer to table below. <table border="1" style="width: 100%; border-collapse: collapse; margin-top: 10px;"> <thead> <tr> <th colspan="4">If reaction pattern is...</th> <th>Then...</th> </tr> <tr> <th>A1 cells</th> <th>A2 cells</th> <th>O cells</th> <th>Auto Control</th> <th>Conclusion</th> </tr> </thead> <tbody> <tr> <td>Pos</td> <td>Neg</td> <td>Neg</td> <td>Neg</td> <td>Anti-A1 confirmed</td> </tr> <tr> <td>Neg</td> <td>Neg</td> <td>Neg</td> <td>Neg</td> <td>Unable to confirm Anti-A1</td> </tr> <tr> <td colspan="4">Any other combination of reactions with a positive A2 cell, O cell, and/or positive Auto control</td> <td>Probable cold antibody causing discrepancy. See procedure "<i>Antibody Identification- How to interpret test results</i>"</td> </tr> </tbody> </table>	If reaction pattern is...				Then...	A1 cells	A2 cells	O cells	Auto Control	Conclusion	Pos	Neg	Neg	Neg	Anti-A1 confirmed	Neg	Neg	Neg	Neg	Unable to confirm Anti-A1	Any other combination of reactions with a positive A2 cell, O cell, and/or positive Auto control				Probable cold antibody causing discrepancy. See procedure " <i>Antibody Identification- How to interpret test results</i> "
If reaction pattern is...				Then...																						
A1 cells	A2 cells	O cells	Auto Control	Conclusion																						
Pos	Neg	Neg	Neg	Anti-A1 confirmed																						
Neg	Neg	Neg	Neg	Unable to confirm Anti-A1																						
Any other combination of reactions with a positive A2 cell, O cell, and/or positive Auto control				Probable cold antibody causing discrepancy. See procedure " <i>Antibody Identification- How to interpret test results</i> "																						


## ABO Discrepancy, Continued

**Procedure-  
A2 cells  
testing  
 continued**

Step	Action
8.	Enter patient A2 cells result and interpretation in Cerner for ABO Discrepancy test (procedure).  
9.	Refer to <b>Cerner Entry</b> and <b>Final Interpretation</b> sections below.

## ABO Discrepancy, Continued

**Procedure-  
 Saline  
 Replacement**

Step	Action						
1.	For <b>extra reactivity in reverse typing:</b>  In Cerner (DOE) add the ABO Discrepancy test to the accession for the ABO/Rh test to record results for rouleaux						
2.	Re-centrifuge the A1 and B tube.						
3.	Remove supernatant leaving the red cell button undisturbed.						
4.	Add 2 drops of saline.						
5.	<b>Do not centrifuge. Re-suspend and read microscopically.</b>						
	<table border="1"> <thead> <tr> <th>If...</th> <th>Then...</th> </tr> </thead> <tbody> <tr> <td><b>Positive</b> (The cells no longer appear like a stack of coins, but appear as agglutination)</td> <td>The test result is positive for hemagglutination and further testing may be indicated.</td> </tr> <tr> <td><b>Negative</b> (The agglutination is dispersed- the saline has successfully dispersed the <b>Rouleaux</b>)</td> <td>The test result is negative.  Enter a comment in test results that the sample had rouleaux and was negative with saline replacement.</td> </tr> </tbody> </table>	If...	Then...	<b>Positive</b> (The cells no longer appear like a stack of coins, but appear as agglutination)	The test result is positive for hemagglutination and further testing may be indicated.	<b>Negative</b> (The agglutination is dispersed- the saline has successfully dispersed the <b>Rouleaux</b> )	The test result is negative.  Enter a comment in test results that the sample had rouleaux and was negative with saline replacement.
If...	Then...						
<b>Positive</b> (The cells no longer appear like a stack of coins, but appear as agglutination)	The test result is positive for hemagglutination and further testing may be indicated.						
<b>Negative</b> (The agglutination is dispersed- the saline has successfully dispersed the <b>Rouleaux</b> )	The test result is negative.  Enter a comment in test results that the sample had rouleaux and was negative with saline replacement.						
6.	Enter patient result for testing completed (with result notes above) and interpretation in Cerner for ABO Discrepancy test (procedure).  						
7.	Refer to <b>Cerner Entry</b> and <b>Final Interpretation</b> sections below.						



## ABO Discrepancy, Continued

### Procedure- Cold Antibodies

Follow the steps below to interpret ABO type in the presence of cold auto or allo-antibodies. Unexpected reactivity seen in reverse typing.

The cold antibody must be identified on the sample. If ABO was ordered, you should add the screen/ABID to confirm presence of cold antibodies

Step	Action
1.	Once the presence of cold antibodies are confirmed (refer to “ <i>Antibody Identification - How to interpret test results</i> ” procedure)  NOTE: Testing with cold reagents and/or plasma may result in unexpected reactions in reverse typing.
2.	Record reverse typing reactions at room temperature in the Cerner ABO/Rh test
3.	If only an ABO/Rh was ordered and an ABO discrepancy is encountered, an antibody screen and antibody ID may be added in DOE to resolve the ABO discrepancy.
4.	Result may be interpreted on forward type results only once a cold auto/allo antibody is confirmed. <ul style="list-style-type: none"> <li>Override exception for interpretation with reason “Not applicable”</li> </ul>

### Cerner Entry

**Unit typing ABO/Rh discrepancy:** An ABO/Rh discrepancy exists when the unit ABO confirmation test result does not agree with the unit ABO/Rh face label.

Step	Action
1.	For ABO/Rh test-retain Vision results sent to Cerner
2.	Evaluate additional testing done and recorded with the ABO Discrepancy test.
3.	Enter final ABO/Rh interpretation in Cerner for test <ul style="list-style-type: none"> <li>Select manually as expected result pattern match is not found</li> <li>Override exception for interpretation with reason “Not applicable”</li> </ul>
4.	Result to be flagged on exception report and reviewed.

**ABO Discrepancy**, Continued

**Final Interpretation** When testing does not resolve forward and reverse ABO discrepancies follow the steps below to determine the final ABO/Rh interpretation.

<b>If</b>	<b>Then</b>
Patient received an ABO-mismatched <b>stem cell or bone marrow transplant.</b>	Refer to <i>SOP Hematopoietic Cell Transplantation and Adoptive Cell Transfer: Patient Management</i>
Patient is <b>elderly, an infant, or otherwise immuno-compromised</b>	Interpret ABO/Rh based on forward typing only. <ul style="list-style-type: none"> <li>Add result comment: Unable to enhance reverse typing</li> </ul>
Patient has <b>confirmed cold antibody</b> (auto and/or allo)	Interpret ABO/Rh based on forward typing only. <ul style="list-style-type: none"> <li>Add result comment: "Unexpected reactivity seen with (A1 cells and/or B cells) due to cold antibody."</li> </ul>
Patient has disease state such as <b>multiple myeloma, Waldentroms macro-globulinemia</b>	Interpret ABO/Rh based on forward typing only. <ul style="list-style-type: none"> <li>Add result comment: "Unexpected reactivity seen with (A1 cells and/or B cells) due to rouleaux."</li> </ul>
Patient has been <b>transfused with RBCs outside of blood group</b> and is showing mixed field reactivity in forward typing	Interpret ABO/Rh based on forward typing with mixed field reactions. <ul style="list-style-type: none"> <li>Add result comment: "Mixed field forward typing due to transfusion of (applicable blood product/type)"</li> </ul>
Patient has been <b>transfused with plasma and/or platelets outside of blood group</b> and is showing decreased reactivity in reverse typing	Interpret ABO/Rh based on forward typing with missing expected reverse typing. <ul style="list-style-type: none"> <li>Add result comment: "Unable to enhance reaction with (A1 cells and/or B cells) due to transfusion of (applicable blood product/type)"</li> </ul>

**ABO Discrepancy**, Continued

**Product RBC Testing Discrepancy** Unit typing ABO/Rh discrepancy: An ABO/Rh discrepancy exists when the unit ABO and/or Rh confirmation test result does not agree with the unit ABO/Rh face label.

Step	Action
1.	Perform a clerical check - Verify donor unit information on tube, electronic records, and the blood bag.
2.	Obtain another unit segment.
3.	Repeat testing using washed donor cells.
4.	If the results remain discrepant, immediately quarantine the unit physically and in the Cerner system.
5.	Complete a QIM.
6.	Notify the blood supplier.

**Uncontrolled Documents** Fung, Mark K. Ed. Technical Manual, 19<sup>th</sup> Ed. AABB ,2017  
Ortho Clinical Diagnostics, current manufacture insert (Instructions for Use)

**Controlled Documents** Attachment A: Serologic Reactions Observed in A and B subgroups  
Attachment B: Possible Causes of ABO Typing Discrepancies  
Hematopoietic Cell Transplantation and Adoptive Cell Transfer: Patient Management  
Antibody Identification - How to interpret test results  
Result Entry: Resulting Patient Tests in Cerner Millennium  
Rh Testing

<b>MasterControl History of Change:</b>		
Change type: new, major, minor etc.	Version #	Description of Change
New	1	New SOP
Minor	2	Added instructions for Cerner entry of results. Added sections for cold antibodies and final interpretation. Updated controlled and uncontrolled documents. Updated Attachment A and B (Tables 10-3 and 10-4 from Technical Manual)

**Authors:** Transfusion Service Managers  
Regional Blood Bank Compliance Officer

**ABO Discrepancy**, Continued  
ATTACHMENT ASerologic Reactions Observed in **Selected A and B subgroups**

Red Cell Phenotype	Anti-A	Anti-B	A1 cells	B cells
A <sub>1</sub>	4+	0	0	4+
A <sub>2</sub>	4+	0	0/2+	4+
The following subgroups must be confirmed by Reference Lab				
A <sub>3</sub>	3+/mf	0	0/2+	4+
A <sub>x</sub>	0/wk	0	0/2+	4+
A <sub>el</sub>	0	0	0/2+	4+
A <sub>m</sub>	0/wk	0	0	4+
B	0	4+	4+	0
B <sub>3</sub>	0	3+/mf	4+	0
B <sub>weak</sub>	0	wk/2+	4+	0
B <sub>el</sub>	0	0	4+	0
B <sub>m</sub>	0	0/wk	4+	0

A<sub>1</sub> is most common subgroup of A (80%), A<sub>2</sub> is 2<sup>nd</sup> most common (20%).

**Reporting: Cerner does not distinguish group A<sub>1</sub> or subgroup A<sub>2</sub>. These are reported as group A.**

Acquired B is usually not demonstrable in mono-clonal Anti-B reagents (check manufacturer insert).

Make sure to add "give O RBC" transfusion requirement if patient is A<sub>2</sub>. **Do NOT give A RBC!**

**ABO Discrepancy**, Continued**ATTACHMENT B****Possible Causes of ABO Typing Discrepancies**

<b>Category</b>	<b>Causes</b>
Weak/missing red cell reactivity	<ul style="list-style-type: none"> <li>• ABO subgroup</li> <li>• Leukemia/malignancy</li> <li>• Transfusion (including Intrauterine fetal)</li> <li>• Pregnancy</li> <li>• Transplantation</li> <li>• Excessive soluble blood group substance</li> </ul>
Extra red cell reactivity	<ul style="list-style-type: none"> <li>• Autoagglutinins/excess protein coating rbc's</li> <li>• Unwashed rbc's; plasma proteins</li> <li>• Unwashed rbc's; antibody in patient's serum to reagent constituent</li> <li>• Transplantation</li> <li>• Acquired B antigen or other polyagglutinable conditions</li> <li>• cis AB or B(A) phenomenon</li> <li>• Out of group transfusion</li> </ul>
Mixed-field red cell reactivity	<ul style="list-style-type: none"> <li>• ABO subgroup</li> <li>• Recent transfusion</li> <li>• Transplantation</li> <li>• Fetomaternal hemorrhage</li> <li>• Twin or dispermic (tetragametic) chimerism</li> </ul>
Weak/missing serum reactivity	<ul style="list-style-type: none"> <li>• Age related (&lt;4-6 months old, elderly)</li> <li>• ABO subgroup</li> <li>• Hypogammaglobulinemia</li> <li>• Transplantation</li> <li>• Excessive anti-A or anti-B (prozone)</li> <li>• Hemodilution</li> </ul>
Extra serum reactivity	<ul style="list-style-type: none"> <li>• Cold antibody-allo/auto</li> <li>• Serum antibody to reagent constituent</li> <li>• Excess serum protein</li> <li>• Transfusion of plasma components</li> <li>• Transplantation</li> <li>• Infusion of intravenous immune globulin</li> </ul>