

JOB AID – PROCESSING STONE RISK ANALYSIS ORDER KAISER 24 HOUR URINE COLLECTION CONTAINER

Step	Action
Note	IF KAISER 24 Hour Urine Collection Container is received (see pictures)
	Then: Lab Assistants: MOB Areas (Except Pasadena)
1	MOB lab assistant/s must label specimen received with the KPHC generic label. DO NOT print/use the CERNER accession label for the Stone Risk Analysis order
2	MOB lab assistant/s must log in specimen received in their 24-hour urine specimen log <ul style="list-style-type: none"> • Make sure the start and end date and start and end time of collection are provided
3	MOB lab assistant/s must completely fill out a Quest Manual Requisition form for Test Code 90523
4	MOB lab assistant must create a “Misc Test” order using KRMS <ul style="list-style-type: none"> • Follow the current process for accessioning via CO function <ul style="list-style-type: none"> ○ Enter required patient and demographic information (i.e. patient’s MRN, DOB, ordering location, ordering provider) ○ Include Sort Group, Priority, etc. ○ Select Alpha Search (AS) and select #42 (for tests starting with letter M) ○ Select #06 – Miscellaneous Test ○ Select PR – Process ○ Press the Enter key twice ○ Completely fill out all 3 questions: <ul style="list-style-type: none"> ❖ Test Name: STONE RISK ANALYSIS 24 HOUR U ❖ Test Code: 90523 ❖ Performing Lab: Quest Diagnostics ○ Press the Enter key twice and the CERNER accession label for the Misc Test will print

JOB AID – PROCESSING STONE RISK ANALYSIS ORDER KAISER 24 HOUR URINE COLLECTION CONTAINER

	<ul style="list-style-type: none"> ○ Label the 24-hour urine specimen collection container with the CERNER accession label ○ Log in the specimen into CERNER
5	MOB lab assistant must write the specimen information in the Specimen Tote Manifest
6	MOB lab assistant/s must send entire specimen to the 4867 main laboratory with the completed Quest Manual Requisition
	Then: Lab Assistants: 4867 Main Laboratory
1	4867 Main Laboratory lab assistant/s must bring the entire specimen to the Routine Urinalysis Area together with the Quest manual requisition
	Then: CLS/MLT: 4867 Main Laboratory – Routine Urinalysis Area
1	4867 Main Laboratory lab CLS/MLT must re-print the CERNER accession labels for the Misc Test
2	4867 Main Laboratory lab CLS/MLT must log into the 24 urine specimens' Log
3	4867 Main Laboratory CLS/MLT must measure and record on the log and on the accession label/s the following: <ul style="list-style-type: none"> ▪ Total Volume ▪ pH
4	4867 Main Laboratory CLS/MLT must mix the urine very well by vigorously shaking the container for at least 1 minute
5	4867 Main Laboratory CLS/MLT must make 2 sets of aliquots using sterile urine cups (Each set must have 1 acidified urine and 1 unacidified urine aliquots) <ul style="list-style-type: none"> ▪ 2 acidified ▪ 2 unacidified

JOB AID – PROCESSING STONE RISK ANALYSIS ORDER KAISER 24 HOUR URINE COLLECTION CONTAINER

6	4867 Main Laboratory CLS/MLT must retain 1 set of the aliquots (1 acidified and 1 unacidified) and place these in the walk-in refrigerator
7	<p>4867 Main Laboratory CLS/MLT must bring the aliquots (1 acidified and 1 unacidified) and the Quest manual requisition to the lab assistant processing Outside Lab Specimen Send-outs</p> <ul style="list-style-type: none"> ▪ Do a warm hand-off to the lab assistant assigned in the processing area for outside reference laboratory specimens for send-outs ▪ Instruct lab assistant to ship the specimen refrigerated
	<p>Then: Lab Assistants: 4867 Outside Ref Lab Specimen Send-outs Area</p>
1	4867 lab assistant/s must pack and send the specimen aliquots together with the Quest manual requisition

JOB AID – PROCESSING STONE RISK ANALYSIS ORDER KAISER 24 HOUR URINE COLLECTION CONTAINER



UR99

Type of Container/Media:

Primary Specimen Collection

Description:

24 hour Urine Bottle for non-preservative collections from which an aliquot is submitted

[Urine Chemistry Specimen Guidelines for the Regional Reference Laboratories](#)

