Step	Action
Note	IF KAISER 24 Hour Urine Collection Container is received (see pictures)
1	Then: Lab Assistants: MOB Area (Pasadena)
	Pasadena MOB lab assistant must label specimen received with the KPHC
	generic label. DO NOT print/use the CERNER accession label for the Stone
	Risk Analysis order
2	Pasadena MOB lab assistant must write/document the specimen
	received in their 24-hour urine specimen log
	 Make sure the start and end date and start and end time of
	collection are provided
3	Pasadena MOB lab assistant/s must completely fill out a Quest Manual
	requisition for Test Code 90523
	Pasadena MOB lab assistant must create a "Misc Test" order using KRMS
	 Follow the current process for accessioning via CO function
	 Enter required patient and demographic information (i.e.
	patient's MRN, DOB, ordering location, ordering provider)
	 Include Sort Group, Priority, etc.
	 Select Alpha Search (AS) and select #42 (for tests starting with
4	letter M)
	 Select #06 – Miscellaneous Test
	Select PR – Process
	 Press the Enter key twice
	 Completely fill out all 3 questions:
	❖ Test Name: STONE RISK ANALYSIS 24 HOUR U
	❖ Test Code: 90523
	Performing Lab: Quest Diagnostics
	 Press the Enter key twice and the CERNER accession label for
	the Misc Test will print

	 Label the 24-hour urine specimen collection container with the CERNER accession label Log in the specimen into CERNER
5	Pasadena MOB lab assistant must write the specimen information in the Specimen Tote Manifest
6	Pasadena MOB lab assistant/s must do a warm hand-off the Pasadena
	MOB CLS with the completed Quest Manual Requisition
	Then: CLS: Pasadena MOB Laboratory – Urinalysis Area
1	Pasadena MOB Laboratory lab CLS must re-print the CERNER accession label for the Misc Test
3	Pasadena MOB Laboratory CLS must measure and record on the log and on the accession label/s the following: Total Volume pH
4	Pasadena MOB Laboratory CLS must mix the urine very well by vigorously shaking the container for at least 1 minute
5	Pasadena MOB Laboratory CLS must make 2 sets of aliquots using sterile urine cups (Each set must have 1 acidified urine and 1 unacidified urine aliquots) 2 acidified 2 unacidified
6	Pasadena MOB Laboratory CLS must retain 1 set of the aliquots (1 acidified and 1 unacidified) and keep these in the refrigerator
7	Pasadena MOB Laboratory CLS must give the aliquots (1 acidified and 1 unacidified) with the completed Quest Manual Requisition to the Pasadena MOB lab assistant to send to the 4867 Main Laboratory

	 Do a warm hand-off to the lab assistant assigned in the processing area for outside reference laboratory specimens for send-outs Instruct the lab assistant to send the following to the 4867 Main Laboratory 1 set of aliquots (1 acidified and 1 unacidified) the completed Quest Manual requisition Instruct lab assistant to ship the specimen refrigerated
	Then: Lab Assistants: Outside Ref Lab Specimen Send-outs Area (Pasadena)
1	Pasadena MOB lab assistant/s must pack and send the aliquots and the completed Quest Manual requisition to the 4867 Main Laboratory
	Then: 4867 Lab Assistants: Outside Ref Lab Specimen Send-outs Area
1	4867 lab assistant/s must write/document the specimen information in the Specimen Tote Manifest





UR99

Type of Container/Media:

Primary Specimen Collection

Description:

24 hour Urine Bottle for nonpreservative collections from which an aliquot is submitted

Urine Chemistry Specimen Guidelines for the Regional Reference Laboratories