

**KAISER MEDICAL CARE PROGRAM
ORANGE COUNTY AREA
POLICIES AND PROCEDURES**

TITLE:	COAGULATION	INDEX NO:	03-035-01
SECTION:	PROCEDURES	ORIGIN DATE:	8/05
SUBJECT:	HEPARIN ANTI-Xa ASSAY (UFH)	REVIEW DATE:	9/14
		REVISION DATE:	9/05, 1/06, 2/07, 11/07, 3/08, 3/13, 9/14, 9/15

HEPARIN ANTI-Xa ASSAY (UFH)

PURPOSE This procedure provides instructions for performing the Heparin Anti-Xa Assay. This assay determines the plasma activity levels of unfractionated (UFH).

**DILUENT/
EQUIPMENT**

- STR-R coagulation analyzer
- STA Owren-Koller Buffer
- STA mini Reducer
- Pipettes
- Reagent grade water (NERL)

SPECIMEN

- Blue top tube-citrated plasma.
- Centrifuge within one hour of collection.
- Test must be performed within 2 hours of collection; otherwise FREEZE plasma until analysis can be performed.
- Samples that are less than 80% full, clotted or grossly hemolyzed should be rejected.

**REAGENTS:
CALIBRATORS**

- STA Multi Hep Calibrators for **UFH**:
 - **Reagent 1:** STA®- Multi Hep Calibrator **0**
 - **Reagent 2:** STA®- Multi Hep Calibrator **4**
 - **Reagent 3:** STA®- Multi Hep Calibrator **7**
 - **Reagent 4:** STA®- Multi Hep Calibrator **10**
 - **Reagent 5:** STA®- Multi Hep Calibrator **18**

**REAGENTS:
CONTROLS**

- STA® Quality HNF/UFH Controls for **UFH**:
 - **Reagent 1:** STA® Quality HNF/UFH Control **2**
 - **Reagent 2:** STA® Quality HNF/UFH Control **7**

Continued on next page

SUBJECT: HEPARIN ANTI-Xa ASSAY (UFH)	INDEX NO: 03-035-02
---------------------------------------------	----------------------------

**CONTROL &
CALIBRATOR
RECONSTITUTION**

Reconstitution of reagents and controls is very critical to getting your controls within range. Wait times are very important.

Step	Action
1.	Always use the pipette provided for reconstitution.
2.	Use freshly poured reagent grade water. Do not pipette directly out of the bottle.
3.	Before opening, tap bottle on hard surface to make all the particles fall to the bottom of the vial.
4.	When opening the vials, pop the vacuum seal by slightly lifting the rubber seals. Tap on the rubber stoppers to dislodge the lyophilized reagent on the caps so it goes back into the bottles and not on the counter. Remove the rubber stoppers.
5.	Add one (1) mL reagent grade water. Swirl to mix. Replace rubber stoppers.
6.	Let them sit for five minutes at room temperature, remove and discard the rubber stoppers.
7.	Replace the white caps on the bottles. Swirl again to mix.
8.	Allow controls/calibrators to stand at least another 25 minutes at room temperature.
9.	Load on instrument and let sit 10-15 minutes prior to calibrating or running controls.

**REAGENTS:
STABILITY**

- STA[®] Multi Hep Calibrators for **UFH**:
 - Reconstituted stability on the STA-R is **4 hours**.
- STA[®] Quality HNF/UFH Controls for **UFH**:
 - Reconstituted stability on the STA-R is **4 hours**.
 - **7 days** at 2 – 8 °C in the original **capped** vials.
- STA[®] LIQUID ANTI – Xa
 - **Reagent 1** (Substrate) & **Reagent 2** (F.Xa)
 - **7 days** on STA-R (with Reducer)
 - **3 months** at 2 – 8 °C in the original **capped** vials (without Reducer)

Continued on next page

SUBJECT: HEPARIN ANTI-Xa ASSAY (UFH)	INDEX NO: 03-035-03
---------------------------------------------	----------------------------

REAGENTS: STA[®] LIQUID ANTI - Xa
TESTING

A bar-code insert is provided in each box of reagents. This contains the lot number, kit code number, and reagent code numbers.

- **Reagent 1:** Substrate (ready to use)
 - Allow the reagent to stand at room temperature (18-25 °C) for 30 minutes before use.
 - Then install a mini Reducer in the vial and replace the perforated plastic cap on top.
 - Request the product drawer to open by clicking the Products menu icon, select load products, bar code the reagent and place the reagent in product drawer R1.
 - Equilibrate on instrument 10-15 minutes before using.
- **Reagent 2:** F.Xa (ready to use)
 - Allow the reagent to stand at room temperature (18-25 °C) for 30 minutes before use.
 - Then install a mini Reducer in the vial and replace the perforated plastic cap on top.
 - Request the product drawer to open by clicking the Products menu icon, select “open”, barcode the reagent and place the reagent in product drawer, **R2**.
 - Equilibrate on instrument 10-15 minutes before using.

Continued on next page

SUBJECT: HEPARIN ANTI-Xa ASSAY (UFH)	INDEX NO: 03-035-04
---------------------------------------------	----------------------------

CALIBRATION See P&P 03-060-01 for instructions

QUALITY CONTROL See QC section of Star Startup and Specimen Processing procedure 03-020-04.

PATIENT TESTING See Load Reagents section of Star Startup and Specimen Processing procedure 03-020-03 and running patients 03-020-05.

THERAPEUTIC RANGE See Coagulation P&P 01-100-01

CRITICAL VALUE See Laboratory P&P 05-090-01
Repeat all critical values to verify result.

Continued on next page

SUBJECT: HEPARIN ANTI-Xa ASSAY (UFH)	INDEX NO: 03-035-05
---------------------------------------------	----------------------------

LINEARITY (UFH) The acceptable linear range for this test is 0.10 to 1.10 IU/mL. A result of less or greater than this range must be acted upon. If the result is > 1.10, the STA-R will send a result of 1.10 to merge in CERNER. **All results of 1.10 must be reviewed on the instrument as follows:**

IF UFH RESULT IS 1.10 THEN CLS MUST TAKE ACTION:

Step	Action														
1.	Double click the patient ID. Patient test screen appears.														
2.	Select the UFH test.														
3.	Review the RAW DATA (IU/mL) in lower right corner of window.														
	<table border="1"> <thead> <tr> <th>IF</th> <th>THEN</th> </tr> </thead> <tbody> <tr> <td>Result is 1.10</td> <td>Repeat & if result duplicates, report result.</td> </tr> <tr> <td>Result is >1.10</td> <td>Dilute patient plasma 1:2 with normal pooled plasma¹, then retest.</td> </tr> <tr> <td></td> <td> <table border="1"> <thead> <tr> <th>If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Result is <=1.10</td> <td>Multiply x 2, report result if it is 1.10 or less.</td> </tr> <tr> <td>Result is >1.10</td> <td>If result is >1.10, report ">1.10"</td> </tr> </tbody> </table> </td> </tr> </tbody> </table> <p>IMPORTANT: DO NOT "Convert Results" to "Freetext", just type in < or > and the value, it will default to what is the linearity in Cerner and trigger the critical notification pop up if applicable.</p>	IF	THEN	Result is 1.10	Repeat & if result duplicates, report result.	Result is >1.10	Dilute patient plasma 1:2 with normal pooled plasma ¹ , then retest.		<table border="1"> <thead> <tr> <th>If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Result is <=1.10</td> <td>Multiply x 2, report result if it is 1.10 or less.</td> </tr> <tr> <td>Result is >1.10</td> <td>If result is >1.10, report ">1.10"</td> </tr> </tbody> </table>	If	Then	Result is <=1.10	Multiply x 2, report result if it is 1.10 or less.	Result is >1.10	If result is >1.10, report ">1.10"
IF	THEN														
Result is 1.10	Repeat & if result duplicates, report result.														
Result is >1.10	Dilute patient plasma 1:2 with normal pooled plasma ¹ , then retest.														
	<table border="1"> <thead> <tr> <th>If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Result is <=1.10</td> <td>Multiply x 2, report result if it is 1.10 or less.</td> </tr> <tr> <td>Result is >1.10</td> <td>If result is >1.10, report ">1.10"</td> </tr> </tbody> </table>	If	Then	Result is <=1.10	Multiply x 2, report result if it is 1.10 or less.	Result is >1.10	If result is >1.10, report ">1.10"								
If	Then														
Result is <=1.10	Multiply x 2, report result if it is 1.10 or less.														
Result is >1.10	If result is >1.10, report ">1.10"														
4.	Click on the Bell icon (lower left corner) to view all error codes pertaining to this result.														

¹ Use NPP (normal pooled plasma), located in -20°C freezer.

SUBJECT: HEPARIN ANTI-Xa ASSAY (UFH)	INDEX NO: 03-035-06
---------------------------------------------	----------------------------

IF UFH RESULT IS LESS THAN 0.10, THEN:

Step	Action								
1.	Double click the patient ID on the Star monitor.								
2.	Select the UFH test.								
3.	Review the RAW DATA (lower right corner of window). Verify that a test was completed and not a “QNS” or error. Then report as “<0.10”.								
	<table border="1"> <thead> <tr> <th>IF</th> <th>THEN</th> </tr> </thead> <tbody> <tr> <td>Result is 0.10</td> <td>Report patient result</td> </tr> <tr> <td>Result is <0.10</td> <td>Report the patient result as “<0.10”</td> </tr> <tr> <td>Result is 0.00</td> <td>Make certain the instrument did not get “QNS” or test error before reporting result as “<0.10”</td> </tr> </tbody> </table>	IF	THEN	Result is 0.10	Report patient result	Result is <0.10	Report the patient result as “<0.10”	Result is 0.00	Make certain the instrument did not get “QNS” or test error before reporting result as “<0.10”
IF	THEN								
Result is 0.10	Report patient result								
Result is <0.10	Report the patient result as “<0.10”								
Result is 0.00	Make certain the instrument did not get “QNS” or test error before reporting result as “<0.10”								
4.	Click on the Bell icon (lower left corner) to view all error codes pertaining to this result.								

LIMITATIONS Any release of platelet factor 4 (PF4), which is a potent heparin inhibitor, will lead to an under-estimation of the heparin level in the plasma being tested. Careful and adequate centrifugation is essential

REFERENCE Diagnostica Stago product insert: April 2012
STA-R Operator Manual
