

APTT Mixing Study

Purpose

- This document describes the procedure for performing an activated partial thromboplastin time (APTT) mixing study to screen patient plasmas for inhibitors of clotting in the intrinsic and common pathways.
 - The rationale for performing a mixing study is to differentiate between a factor deficiency and an inhibitor.
-

Policy

In order to perform the APTT mixing study, the baseline APTT result should exceed the upper limit of the laboratory's defined reference range by 5 seconds or more.

Principle

- The APTT mixing study is performed to detect inhibitors of clotting in the intrinsic and common pathways, and to determine whether the prolongation of the APTT is due to deficiency in factor levels or due to a circulating inhibitor, often referred to as a circulating anticoagulant.
 - An APTT is performed on a 1:1 mixture of one part normal pooled plasma (NPP) and one part patient plasma. If the APTT does not correct to within the normal reference range on the immediate mix, the presence of an inhibitor is indicated.
 - If there is correction to within the normal range on the immediate mix, a second APTT must be performed on a timed incubation of the patient and NPP mixture. If the APTT remains corrected to within 3 seconds of the upper limit of the reference range following incubation, a factor deficiency is indicated. If the correction disappears following incubation, the presence of an inhibitor is indicated. This may occur because certain inhibitors, such as Factor VIII inhibitors and about 15% of lupus anticoagulant inhibitors are time-or temperature-dependent.
-

Scope

The intended users of this document include Clinical Laboratory Scientists (CLS) and Laboratory Technical Supervisors handling APTT mixing study samples, issues, or concerns.

Specimen Sources

Plasma from citrated whole blood (blue top) drawn by venipuncture

Continued on next page

APTT Mixing Study, Continued

Specimen Collection and Transport

Citrated whole blood (blue top) should be collected, handled, transported and processed in accordance with CLSI Document H21-A5 *Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-5th Edition*.

- Centrifuge within one hour of collection.
 - Specimens maintained as plasma-based whole blood are stable up to 4 hours.
 - Spun citrated plasma is stable for four hours.
 - If testing cannot be performed within 4 hours of collection, prepare platelet-poor plasma by double centrifugation, then freeze.
 - Refrigeration and transportation of whole blood specimens on ice is not recommended because cold temperatures may lead to a gradual loss of von Willebrand Factor and factor VIII activity.
-

Preanalytical Variables

- The plasma should be evaluated to exclude micro clots or fibrin strands by passing a small wooden stick through the sample or by gently inverting the sample tube for possible clotting. The presence of micro clots or fibrin threads could indicate a difficult venipuncture and pre-activation of some of the factors.
 - A high hematocrit or short draw can result in falsely prolonged APTT due to excess citrate anticoagulant.
 - Grossly hemolyzed specimens should be rejected, if possible. APTT values may increase or decrease because cell lysis products include tissue factors that may activate coagulation.
-

Technical Considerations

- NPP should be made from a pool of donors with normal factor levels, and must be fresh frozen and cell free.
 - Commercial NPP from Precision Biologic must be used in mixing studies. See ordering information in the “Specialty Products Needed” section.
 - Plasma should be stored frozen in 1 or 0.5 ml aliquots. At -70° C, stored NPP will be stable for 6 months. NOTE: Plasma aliquots can alternatively be maintained in a -20°C freezer for up to 6 weeks. If the -20°C freezer has automatic defrost cycles, aliquots must be placed inside a small Styrofoam container inside the freezer.
 - NPP is considered competent for use in the mixing study if the APTT is within 2 seconds of the baseline APTT when the lot was first shipped to the lab (usually 28-31 seconds).
-

Continued on next page

APTT Mixing Study, Continued

**Specialty
Products
Needed**

Description	Vendor	Product Number
Cryocheck Pooled Normal Plasma	Precision Biologic	CCN-10, available in 0.5 mL or 1 mL aliquots

Equipment

- Diagnostica Stago Coagulation Analyzer
 - Pipettes
-

**Materials and
Supplies**

- Pipette tips
 - Micro vials
 - Micro vial adapters
-

Safety

Refer to the safety manual for general safety requirements.

Quality Control

- Refer to Stago Quality Control and Start-up Procedures for specific guidelines for performing quality control for APTT assay.
 - The APTT assay should be performed and documented on the NPP used in mixing study at the beginning of each mixing study run.
 - The APTT result on the NPP used in the mixing study should fall within the established normal reference range of the laboratory, and additionally be within 2 seconds of the baseline APTT performed on the NPP when the lot was first shipped to the laboratory (usually 28-31 seconds).
 - If performing the incubation study, patient plasma and NPP should also be incubated separately for one hour at 37°C without mixing, and then mixed together for the APTT to be performed. This will serve as a control for the timed incubation, which may affect the stability of factors V and VIII.
-

Continued on next page

APTT Mixing Study, Continued

Procedure: Follow the steps below to perform the APTT mixing study.

1:1 Immediate Mixing Study

Step	Action								
1	Perform APTT on patient plasma alone. Record the result on the APTT mixing study worksheet. <table border="1" style="margin-left: 20px;"> <thead> <tr> <th style="text-align: center;">If...</th> <th style="text-align: center;">Then...</th> </tr> </thead> <tbody> <tr> <td>a. patient baseline APTT is normal (within reference range)</td> <td>a mixing study is not indicated and should not be performed. Select [APTT Norm] from comment dropdown to attach canned interpretation message.</td> </tr> <tr> <td>b. patient baseline APTT is minimally prolonged (<5 seconds from the upper limit of reference range)</td> <td>a mixing study is not performed. Select [APTT Min Prolon] from comment dropdown to attach canned interpretation message.</td> </tr> <tr> <td>c. patient baseline APTT is prolonged (≥5 seconds from the upper limit of reference range)</td> <td>proceed to perform a 1:1 immediate mixing study.</td> </tr> </tbody> </table>	If...	Then...	a. patient baseline APTT is normal (within reference range)	a mixing study is not indicated and should not be performed. Select [APTT Norm] from comment dropdown to attach canned interpretation message.	b. patient baseline APTT is minimally prolonged (<5 seconds from the upper limit of reference range)	a mixing study is not performed. Select [APTT Min Prolon] from comment dropdown to attach canned interpretation message.	c. patient baseline APTT is prolonged (≥5 seconds from the upper limit of reference range)	proceed to perform a 1:1 immediate mixing study.
If...	Then...								
a. patient baseline APTT is normal (within reference range)	a mixing study is not indicated and should not be performed. Select [APTT Norm] from comment dropdown to attach canned interpretation message.								
b. patient baseline APTT is minimally prolonged (<5 seconds from the upper limit of reference range)	a mixing study is not performed. Select [APTT Min Prolon] from comment dropdown to attach canned interpretation message.								
c. patient baseline APTT is prolonged (≥5 seconds from the upper limit of reference range)	proceed to perform a 1:1 immediate mixing study.								
2	Gently mix 200 uL of patient plasma and 200 uL of NPP together in a single plastic tube or instrument micro vial.								
3	Immediately after preparation, perform APTT on mixture. <table border="1" style="margin-left: 20px;"> <thead> <tr> <th style="text-align: center;">If after immediate mixing study...</th> <th style="text-align: center;">Then...</th> </tr> </thead> <tbody> <tr> <td>a. the immediate mix APTT corrects to within normal reference range for APTT</td> <td>an additional mixing study with incubation at 37°C should be performed (see Step 4)</td> </tr> <tr> <td>b. the immediate mix APTT does not correct to within the normal reference range for APTT</td> <td>incubation study is not applicable. Results are suggestive of an inhibitor. Select [Does Not Correc] from comment dropdown to attach canned interpretation message.</td> </tr> </tbody> </table>	If after immediate mixing study...	Then...	a. the immediate mix APTT corrects to within normal reference range for APTT	an additional mixing study with incubation at 37°C should be performed (see Step 4)	b. the immediate mix APTT does not correct to within the normal reference range for APTT	incubation study is not applicable. Results are suggestive of an inhibitor. Select [Does Not Correc] from comment dropdown to attach canned interpretation message.		
If after immediate mixing study...	Then...								
a. the immediate mix APTT corrects to within normal reference range for APTT	an additional mixing study with incubation at 37°C should be performed (see Step 4)								
b. the immediate mix APTT does not correct to within the normal reference range for APTT	incubation study is not applicable. Results are suggestive of an inhibitor. Select [Does Not Correc] from comment dropdown to attach canned interpretation message.								

Continued on next page

APTT Mixing Study, Continued

Procedure Follow the steps below to perform the APTT incubation mixing study.

**Mixing Study
 with Incubation**

Step	Action						
4	<p>Test incubation:</p> <ul style="list-style-type: none"> Incubate a 1:1 mixture (e.g. 300µL + 300µL) of patient test plasma and NPP in a single plastic tube for 1 hour at 37 °C. Perform this incubation at the same time as the control incubation step (see Step 5 below). After the 1 hour incubation, run the incubated test APTT on the incubated mixed sample of patient plasma and NPP. 						
5	<p>Control incubation:</p> <ul style="list-style-type: none"> Incubate 300 µL of patient plasma alone and 300 µL NPP alone in separate plastic tubes for 1 hour at 37 °C. Perform this incubation step at the same time as the patient test incubation step (see Step 4 above). After the 1 hour incubation, gently mix the patient plasma and the NPP from their separate tubes into a single plastic tube or instrument microvial. Load plastic tube or micro vial with the mixture onto instrument and run the incubated control APTT. 						
6	<ul style="list-style-type: none"> For the run to be valid, the incubated control APTT will not increase by more than 3 seconds from the immediate mix APTT. If this is the case and the run is valid, use the following to interpret the patient test result. If this is not the case, see Step 7. <table border="1" data-bbox="548 1192 1442 1745"> <thead> <tr> <th data-bbox="548 1192 992 1234">If after incubation ...</th> <th data-bbox="992 1192 1442 1234">Then...</th> </tr> </thead> <tbody> <tr> <td data-bbox="548 1234 992 1451">a. the incubated test APTT remains corrected to within 3 seconds of the upper limit of the reference range (see Interpretation / Results / Alert Values block below)</td> <td data-bbox="992 1234 1442 1451">Mixing studies are suggestive of factor deficiency. Select [Remains Correct] from comment dropdown to attach canned interpretation message.</td> </tr> <tr> <td data-bbox="548 1451 992 1745">b. the incubated test APTT does not remain corrected to within 3 seconds of the upper limit of the reference range (see Interpretation / Results / Alert Values block below)</td> <td data-bbox="992 1451 1442 1745">Mixing studies indicate a time- or temperature-dependent factor inhibitor such as factor VIII inhibitor or some lupus anticoagulants. Select [Does Not Remain] from comment dropdown to attach canned interpretation message.</td> </tr> </tbody> </table>	If after incubation ...	Then...	a. the incubated test APTT remains corrected to within 3 seconds of the upper limit of the reference range (see Interpretation / Results / Alert Values block below)	Mixing studies are suggestive of factor deficiency. Select [Remains Correct] from comment dropdown to attach canned interpretation message.	b. the incubated test APTT does not remain corrected to within 3 seconds of the upper limit of the reference range (see Interpretation / Results / Alert Values block below)	Mixing studies indicate a time- or temperature-dependent factor inhibitor such as factor VIII inhibitor or some lupus anticoagulants. Select [Does Not Remain] from comment dropdown to attach canned interpretation message.
If after incubation ...	Then...						
a. the incubated test APTT remains corrected to within 3 seconds of the upper limit of the reference range (see Interpretation / Results / Alert Values block below)	Mixing studies are suggestive of factor deficiency. Select [Remains Correct] from comment dropdown to attach canned interpretation message.						
b. the incubated test APTT does not remain corrected to within 3 seconds of the upper limit of the reference range (see Interpretation / Results / Alert Values block below)	Mixing studies indicate a time- or temperature-dependent factor inhibitor such as factor VIII inhibitor or some lupus anticoagulants. Select [Does Not Remain] from comment dropdown to attach canned interpretation message.						

Continued on next page

APTT Mixing Study, Continued

Mixing Study with Incubation, continued

Step	Action
7	<ul style="list-style-type: none">• If the incubated control APTT increases by more than 3 seconds from the immediate mix APTT, check the temperature of the heating element to see that it is at 37 °C.• Repeat the mixing study with incubation (Steps 4-6) at the correct temperature.• If the incubated control APTT remains increased by more than 3 seconds from the immediate mix APTT on the repeat study, stop the study and consult a supervisor.• The sample may need to be referred to the Regional Reference Laboratories if additional testing is still required.

Continued on next page

APTT Mixing Study, Continued

**Interpretation /
 Results / Alert
 Values**

Use the following guidelines for interpretation of mixing study results (the results for the incubation study assume that the run was valid, such that the **incubated control APTT** does not differ from the **immediate mix APTT** by more than 3 seconds):

Immediate Mix Result	Incubation Study Result	Interpretation
Complete correction: the immediate mix APTT corrects to within APTT reference range	Correction: the incubated test APTT corrects to within 3 seconds of the upper limit of the normal APTT reference range	Results of these studies are suggestive of factor deficiency. [Remains Correct]
Complete correction: the immediate mix APTT corrects to within APTT reference range	No correction: the incubated test APTT does not correct to within 3 seconds of the upper limit of the normal APTT reference range	Results of these studies indicate a time-or-temperature-dependent factor inhibitor such as factor VIII inhibitor or some lupus anticoagulants. [Does Not Remain]
Partial or no correction: the immediate mix APTT does not correct to within APTT reference range	Not applicable	Results are suggestive of an inhibitor. The presence of anticoagulant inhibitor drugs such as heparin or direct thrombin inhibitors cannot be excluded. [Does Not Correct]

Continued on next page

APTT Mixing Study, Continued

Results

- Result APTT 1:1 Mixing Study manually in LIS using **Accession Result Entry** task module.
- Report clotting times in whole seconds.
- Select appropriate code from comment dropdown to attach canned interpretation message.
- Screenshots of an APTT Mixing Study are presented below:

Procedure	Result	Flags	Status	Reference Range
APTT Mix	48	H	Pending	25 - 37
APTT 1:1 Mix Im	32			25 - 37
APTT 1:1 Mix 1H	34			25 - 37
APTT 1:1 Mix Cn	35			25 - 37
dns APTT Mix Cm	Remains Correct			
APTT Mix Commnt	APTT Norm APTT Mtn Prolon Remains Correct Does Not Remain Does Not Correc		Pending	

- Patient baseline APTT result
- APTT Immediate 1:1 Mix
- APTT 1:1 Mix – 1Hr. Incubation
- APTT 1:1 Mix – 1 Hr. Incubation Control
- Select appropriate message code from dropdown.

Limitations

While this procedure can broadly identify whether a factor deficiency or an inhibitor may be present in the patient sample, it does not identify any one specific factor deficiency or inhibitor by name. Identification of specific factor deficiencies or inhibitors may be performed at the Regional Reference Laboratory if clinically necessary.

Continued on next page

APTT Mixing Study, Continued

Non-Controlled Documents The following non-controlled document supports this policy.

Clinical and Laboratory Standards Institute (CLSI). *Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline-Fifth Edition*. CLSI document H21-A5 (ISBN 1-56238-657-3). Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2008.

Author(s) Ji Yeon Kim, MD, MPH
Bill Brice, MT(ASCP), MBA
Eleanor E. Callasan, MPH, CLS(ASCP)

APTT Mixing Study, Continued

HISTORY PAGE

Type of Change: New Major, Minor	Description of Change(s)	Quality Systems Leader/Date	Operations Director, Area Laboratory Review/Date	CLIA Laboratory Director Review/Date	Date Change Implemented
New					
Major Version 02	<ul style="list-style-type: none"> • Page 4: 1:1 Immediate Mixing Study: Step 1 (rows a and b) and Step 3 (row b): Added instructions for the actions required to attach the appropriate canned interpretation messages. • Page 5: Mixing Study with Incubation: Step 6 (rows a and b): Added instructions for the actions required to attach the appropriate canned interpretation messages. • Page 7: Interpretation / Results / Alert Values: Interpretation column: Added the applicable canned interpretation message that will appear with each result. • Page 8: Results: Added specific instructions for manually entering results and selecting comments in the LIS; replaced the RESULT EXAMPLE table with a screenshot taken from the LIS (includes a description of the contents of each result field and an instruction for using the dropdown). 				

APTT Mixing Study, Continued

HISTORY PAGE

Type of Change: New Major, Minor	Description of Change(s)	Quality Systems Leader/Date	Operations Director, Area Laboratory Review/Date	Laboratory Director Review/Date	Date Change Implemented
Major Version 03	<ul style="list-style-type: none"> • Page 2 Technical Considerations: specified Precision Biologic as the vendor of the of the Normal Pooled Plasma (NPP); deleted instruction to alternatively prepare NPP in-house; defined criteria for determining the competency of the NPP • Page 3 Specialty Products Needed: removed ordering information of Pooled Normal Plasma from George King Biomedical • Page 3 Quality Control: replaced criteria for acceptability of the APTT result of the NPP • Reformatted 				

APTT Mixing Study, Continued

Reviewed and approved by (for Medical Center Area Approval Only):

SIGNATURE	DATE
Name: _____ Operations Director, Area Laboratory	
Name: _____ CLIA Laboratory Director	

Signature Manifest

Document Number: SCPMG-PPP-0051

Revision: 03

Title: APTT Mixing Study - Procedure

All dates and times are in Pacific Standard Time.

Revise APTT procedure

New Document or Change Request

Name/Signature	Title	Date	Meaning/Reason
Vincent Dizon (I713793)			
Mary Anne Umekubo (K076412)	ASST DIR	13 Nov 2014, 04:30:33 PM	Approved

Collaboration

Name/Signature	Title	Date	Meaning/Reason
Mary Anne Umekubo (K076412)			
Vincent Dizon (I713793)	Director of Lab Services, Chem	20 Nov 2014, 04:39:18 PM	Complete

Initial Approval

Name/Signature	Title	Date	Meaning/Reason
Ji Yeon Kim (B727360)	Physician-In-Charge, Chem Svcs	20 Nov 2014, 04:58:39 PM	Approved

Final Approval

Name/Signature	Title	Date	Meaning/Reason
Darryl Palmer-Toy (T188420)	RRL MEDICAL DIRECTOR	22 Dec 2014, 09:38:59 AM	Approved

Set Effective Date

Name/Signature	Title	Date	Meaning/Reason
Vincent Dizon (I713793)	Director of Lab Services, Chem	16 Jan 2015, 10:33:32 AM	Approved

Notify Trainers

Name/Signature	Title	Date	Meaning/Reason
Vincent Dizon (I713793)	Director of Lab Services, Chem	16 Jan 2015, 10:33:32 AM	Email Sent