

Beaumont

Origination 6/21/2022
Last Approved 6/7/2022
Effective 6/21/2022
Last Revised 6/7/2022
Next Review 6/6/2024

Document Contact Kelly Sartor
Area Laboratory-Blood Bank
Applicability Dearborn

Weighing Blood Products - Dearborn Blood Bank

Document Type: Procedure

I. PURPOSE AND OBJECTIVE:

This document will provide Blood Bank staff with the procedure to determine, with precision and accuracy, the weight of blood products.

II. GENERAL INFORMATION:

- A. For the purposes of this document 1 gm is equivalent to 1 mL of blood product.
- B. Each RBC (with the exception of apheresis) units is given a default weight when brought in to the computer system. All other blood components (Fresh Frozen Plasma [FFP], cryoprecipitate, platelets) have their individual weights listed on the face label and are captured in the computer system when being processed.
- C. Occasionally the nurse will request the weight of an RBC unit to be given. The digital scale is used for this purpose.
- D. The standard unit of RBCs including autologous RBCs is collected in 600 ml bags.
- E. Below are the pre-determined empty bag weights for each size of transfer bags available for use in the Blood Bank.
 - 1. 150 ml transfer bag = 15 g
 - 2. 300 ml transfer bag= 25 g
 - 3. 600 ml transfer bag= 30 g

III. SCOPE:

The digital scale is used whenever the estimated volume for a Red blood Cell (RBC) unit is not precise

enough as determined by nursing or when indicated by specific blood bank procedure.

IV. EQUIPMENT:

Ohaus Digital Scale – located in the triage area near the segment sealing device

V. DEFINITIONS/ACRONYMS:

- A. FFP: Thawed Fresh Frozen Plasma
- B. ISBT: International Society of Blood Transfusion

VI. QUALITY CONTROL:

Weigh the standard 200 g weight on the scale and record the weight on the *Daily Temperature and Quality Control Record*. Document the *Daily Temperature and Quality Control Record*, with a “√” to indicate that scale is operating within the expected results.

VII. PROCEDURE:

A. Undivided RBC Components

Each RBC (with exception of apheresis units) is given a default weight when brought into the computer system. When the nursing unit specifically requests the actual weight of the unit be given:

1. Turn on the digital scale or verify the zero setting by using the reset button.
2. Place the RBC component needing to be weighed on a scale, ensuring that the attached segments are not included in the final weight.
3. Subtract 30 grams from the reading to allow for the weight of the empty bag.
4. Record the weight of the product on the *Record of Transfusion* form in grams.

B. Component Aliquots

1. RBC units to be used as parent units for aliquots or divided products must be weighed before any modification is done to them.
 - a. If the unit is already brought into the computer, then the unit must be weighed, and its weight must be modified in the computer. *Example: a unit which needs to be divided for an adult who has severe congestive heart disease and has orders to give one unit over more than 4 hours.*
 - b. If the unit specifically meets the needs for an infant, i.e. quad pack, the weight may be taken before the unit is brought into the computer inventory and modified at that time.
 - c. Obtain the volume of the aliquot(s) by weighing the transfer bag in which the component was aliquoted and subtracting the pre-determine weight of the transfer bag.

- d. The volume for the aliquot(s) will have its weigh recorded by the medical technologist during computer modification. The computer automatically updates the volume of the parent bag. This value should be confirmed with actual weight and modified if necessary. Refer to Transfusion Medicine policy, [Aliquot Preparation](#) for more information.
2. All other blood components (FFP, cryoprecipitate, platelets) have their individual weights listed on the face label and are captured in the computer system.
 3. Obtain the volume of the aliquot by weighing the transfer bag in which the component was aliquoted and subtracting the pre-determine weight of the transfer bag.
 4. Perform the product modification in the computer using the Blood Bank CDM - *Dividing Blood Products*. Record in the computer and record the volume for the prepared aliquot. After entering the volume of the prepared aliquot in the computer the computer automatically updates the volume of the parent bag and allows for the print of new product labels.
 5. Product label with the new volume for the parent bag should be printed after modification or the parent bag's "new" volume for the parent bag should be recorded on the face label if automated ISBT label is not available.
Refer to Transfusion Medicine policy, [Aliquot Preparation](#) and Blood Bank CDM - *Dividing Blood Products* for more information.

VIII. REFERENCES:

1. AABB *Standards for Blood Banks and Transfusion Medicine*, current edition.
2. AABB *Technical Manual*, current edition.
3. College of American Pathologist, *Transfusion Medicine Checklist*, current edition.

Approval Signatures

Step Description	Approver	Date
	Jeremy Powers: Chief, Pathology	6/7/2022
Policy and Forms Steering Committee (if needed)	Kelly Sartor: Supv, Laboratory	6/6/2022
Policy and Forms Steering Committee (if needed)	Gail Juleff: Project Mgr Policy	6/6/2022
	Kimberly Geck: Dir, Lab Operations B	6/4/2022
	Kelly Sartor: Supv, Laboratory	6/3/2022
	Kelly Sartor: Supv, Laboratory	6/3/2022