

Beaumont

Origination 11/9/2021
Last 3/1/2023
Approved
Effective 3/1/2023
Last Revised 3/1/2023
Next Review 2/28/2025

Document Contact Kelly Walewski:
Supv, Laboratory
Area Laboratory-
Chemistry
Applicability All Beaumont
Hospitals

Formulas for Manual Calculations

Document Type: Guideline

I. PURPOSE AND OBJECTIVE:

- A. The purpose of this procedure is to specify formulas needed to manually calculate results.

II. FORMULAS:

- A. Anion Gap = Sodium – (Chloride + Carbon dioxide)
- B. Globulin = Total Protein – Albumin
- C. A/G Ratio = Albumin / (Total Protein – Albumin)
- D. Cholesterol/High Density Lipoprotein Ratio = Chol/HDL
- E. Non HDL Cholesterol = Chol – HDL = mg/dL
- F. LDL Cholesterol = Total Chol – HDL – (Trig/5)
Note: Formula for calculation only if Trig<400 mg/dL
- G. BUN/Creatinine Ratio = BUN/Creatinine
- H. CKD-EPI Estimated GFR (mL/min/1.73 m²) = $142 \times \min(S_{Cr} / \kappa, 1)^\alpha \times \max(S_{Cr} / \kappa, 1)^{-1.200} \times 0.9938^{\text{Age}} \times 1.012$ [if female]
1. S_{Cr} is serum creatinine in mg/dL
 2. κ is 0.7 for females and 0.9 for males
 3. α is -0.241 for females and -0.302 for males
 4. min indicates the minimum of S_{Cr} / κ or 1
 5. max indicates the maximum of S_{Cr} / κ or 1

- a. Example for a female patient: $142 \times \min(S_{Cr}/0.7, 1)^{-0.241} \times \max(S_{Cr}/0.7, 1)^{-1.200} \times 0.9938^{Age} \times 1.012$
- b. Example for a male patient: $142 \times \min(S_{Cr}/0.9, 1)^{-0.302} \times \max(S_{Cr}/0.9, 1)^{-1.200} \times 0.9938^{Age}$
- I. CKD-EPI Cystatin C Equation (2012) Estimated GFR (mL/min/1.73 m²) = $133 \times \min(S_{cys} / 0.8, 1)^{-0.499} \times \max(S_{cys} / 0.8, 1)^{-1.328} \times 0.9936^{Age} \times 0.932$ [if female]
1. S_{cys} is standardized serum cystatin C in mg/L
 2. min indicates the minimum of S_{cys} / 0.8 or 1
 3. max indicates the maximum of S_{cys} / 0.8 or 1
 4. age is years
- J. CKD-EPI Cystatin C Equation Estimated GFR Under 25 (mL/min/1.73 m²) = $K \times (1/S_{cys})$
1. S_{cys} is standardized serum cystatin C in mg/L
 2. Male 1-15 years $K = 87.2 \times 1.011^{(age\ years-15)}$
 3. Male 15-18 years $K = 87.2 \times 0.960^{(age\ years-15)}$
 4. Female 1-12 years $K = 79.9 \times 1.004^{(age\ years-12)}$
 5. Female 12-18 years $K = 79.9 \times 0.974^{(age\ years-12)}$
- K. Estimated GFR (mL/min/1.73 m²) = $133 \times \min(S_{cys} / 0.8, 1)^{-0.499} \times \max(S_{cys} / 0.8, 1)^{-1.328} \times 0.9936^{Age} \times 0.932$ [if female]
- L. % Saturation = $100 \times [\text{serum Iron} / (\text{Transferrin} \times 1.4)]$
- M. Total Iron Binding Capacity = Transferrin \times 1.4
- N. Urine Albumin/Creatinine Ratio = (ualb mg/dL/Urine Creatinine mg/dL) \times 1000
- O. Urine Creatinine Clearance mL/min

$$\frac{\text{Raw Urine Creatinine mg/dL} \times (\text{TV mL} \times 1 \text{ Hr})}{\text{Serum Creatinine mg/dL} \times \text{Hrs collection} \times 60 \text{ min}}$$
- P. Urine (Ur) Creatinine Clearance Pediatric mL/min CC Corrected =

$$\frac{\text{Raw Ur Creat mg/dL} \times (\text{TV mL} \times 1 \text{ Hr}) \times (1.73 \text{ m}^2)}{\text{Serum Creat mg/dL} \times \text{Hrs collection} \times 60 \text{ min} \times \text{surface area m}^2}$$
 NOTE: Surface area m² is determined from child's height and weight.
- Q. Urine Creatinine Timed or 24 Hour mg/Collection (Coll) =

$$\frac{\text{Raw Urine Creatinine mg}}{\text{dL}} \times \text{dilution} \times \frac{(\text{TV mL} \times 10 \text{dL})}{\text{Coll.} \times 1000 \text{mL}}$$
- R. Urine Amylase Timed Ur/Hr =

$$\frac{\text{Raw Urine Amylase U}}{\text{L}} \times \text{dilution} \times \frac{\text{TV mL}}{\text{Hrs of Coll.}} \times \frac{1 \text{L}}{1000 \text{mL}}$$
- S. Urine Calcium Timed or 24 Hour, mg/Collection =

$$\frac{\text{Raw Urine Protein mg/dL}}{\text{Raw Urine Creatinine mg/dL}}$$
 NOTE: reference range 0.0-0.2

AG. Urine Calcium / Creatinine Ratio =
$$\frac{\text{Raw Urine Calcium mg/dL}}{\text{Raw Urine Creatinine mg/dL}}$$
 NOTE: reference range 0.02-0.26

AH. AaPO₂, mmHg = (713mmHg x FiO₂ as decimal) – (PCO₂ / 0.8) - PO₂

AI. PF, no units = PO₂ / FiO₂ (as a decimal value)

AJ. % Free PSA = (Free PSA/PSA) X 100

AK. Albumin gradient; Serum and Ascites fluid, g/dL = Albumin (g/dL) – Albumin fluid (g/dL)
 NOTE: >1.1 = Transudate
 0.2-1.0 = Exudate

Attachments

[Formulas for Manual Calculations Calculation Verification Worksheet Attachment A_1_2023.pdf](#)

Approval Signatures

Step Description	Approver	Date
Medical Directors	Muhammad Arshad: Physician	3/1/2023
Medical Directors	Vaishali Pansare: Chief, Pathology	2/10/2023
Medical Directors	Jeremy Powers: Chief, Pathology	2/8/2023
Medical Directors	Ryan Johnson: OUWB Clinical Faculty	2/2/2023
Medical Directors	Ann Marie Blenc: System Med Dir, Hematopath	2/1/2023
Medical Directors	John Pui: Chief, Pathology	2/1/2023
Policy and Forms Steering Committee Approval (if needed)	Colette Kessler: Mgr, Division Laboratory	2/1/2023

Policy and Forms Steering
Committee Approval (if
needed)

Ilene Hirsch: Project Mgr Policy	2/1/2023
Caitlin Schein: Staff Physician	2/1/2023
Qian Sun: Tech Dir, Clin Chemistry, Path	1/31/2023
Nga Yeung Tang: Tech Dir, Clin Chemistry, Path	1/24/2023
Elzbieta Wystepek: Dir, Lab Operations B	1/24/2023
Kimberly Geck: Dir, Lab Operations B	1/23/2023
Brittnie Berger: Dir, Lab Operations C	1/23/2023
Colette Kessler: Mgr, Division Laboratory	1/23/2023
Colette Kessler: Mgr, Division Laboratory	1/23/2023

Applicability

Dearborn, Farmington Hills, Grosse Pointe, Royal Oak, Taylor, Trenton, Troy, Wayne