

# PROCEDURE

## Corewell Health East - Parallel Testing Using Kit Controls - Dearborn, Taylor, Trenton, Wayne

**This Procedure is Applicable to the following Corewell Health sites:**

Corewell Health Dearborn Hospital, Corewell Health Taylor Hospital, Corewell Health Trenton Hospital, Corewell Health Wayne Hospital

|                                  |                                 |
|----------------------------------|---------------------------------|
| <b>Applicability Limited to:</b> | N/A                             |
| <b>Reference #:</b>              | 34211                           |
| <b>Version #:</b>                | 2                               |
| <b>Effective Date:</b>           | 10/09/2025                      |
| <b>Functional Area:</b>          | Clinical Operations, Laboratory |
| <b>Lab Department Area:</b>      | Lab - Chemistry                 |

### 1. Principle

This document is intended to provide instructions for testing personnel to perform parallel testing for Quality Control (QC) and Reagent Verification for Group A Strep, Mono and C. diff test kits. Clinical laboratory reagents and control materials are exposed to many variables due to conditions during transportation and storage environments in different laboratory settings. The validation of new shipments and/or new lot numbers of test kits with currently in use test kits is performed to ensure that, despite varying environmental conditions, there are no clinically significant differences in the results obtained when different shipments and/or new lot numbers of reagents are used. Control materials are parallel tested to verify that the same results are achieved using reagents from each test kit.

### 2. Responsibility

Personnel who have completed the competency requirements will perform this testing.

### 3. Specimen

- A. Test kit controls, or
- B. Previously positive and previously negative patient samples, or
- C. Previously run positive and negative College of American Pathologist (CAP) survey samples.

### 4. Reagent/Equipment Needed

- A. Acceava® Strep A kit
- B. Sure-View™ Mono Test kit
- C. C. DIFF Quik Check Complete® test kit

### 5. Procedure

- A. See individual test procedures for testing information.
  - 1. [Corewell Health East - Sure-View Mono Test - All Beaumont Hospitals](#)
  - 2. [Corewell Health East - Strep Screen, Group A - All Beaumont Hospitals](#)
  - 3. [Corewell Health East - C. difficile GDH / Toxin - All Beaumont Hospitals](#)

Entities will reference associated Documentation contained within this document as applicable  
Printouts of this document may be out of date and should be considered uncontrolled.

- B. Each new shipment and/or new lot number of test kits are tested in parallel with the “kit in use”.
- C. “In use” controls are tested along with the new lot controls with each new lot or new shipment of reagents. It is also acceptable to run previously positive and negative patient samples or CAP survey samples as parallel controls.
- D. Record results on the appropriate Parallel Testing QC Worksheet.
- E. If acceptable results are obtained, label the kit “ready for use” and initial and date the box.
- F. If unacceptable results are obtained, repeat the test. If the results are still not acceptable, notify the supervisor, lead tech or manager for assistance.
- G. STAT specimens should be sent to another site for testing until acceptable QC can be obtained.
- H. If the test kit contains multiple components, do not mix components with another reagent kit of the same kind unless otherwise advised by the manufacturer.
- I. Additional QC should be performed at a frequency according to the manufacturer's recommendations or based on individual testing procedures.
- J. Parallel testing logs will be reviewed at least monthly by the lead medical technologist, or designee.

## **6. Results/Interpretation**

All qualitative results must agree for the parallel testing to pass.

## **7. Revisions**

Corewell Health reserves the right to alter, amend, modify or eliminate this document at any time without prior written notice.

## **8. Procedure Development and Approval**

### **Document Owner:**

Larah Crawford (Sr Safety Specialist)

### **Writer(s):**

Laura Bellon (Medical Technologist Lead)

### **Reviewer(s):**

Lillian Reid (Medical Technologist Lead), Malak Saad (Medical Technologist Lead), Richard Quick (Medical Technologist Lead), Shani Kastl (Medical Technologist Lead), Tanya Williams (Medical Technologist Lead)

### **Approver:**

Ashley Beesley (Mgr, Laboratory), Brittne Berger (Dir Sr, Lab Operations), Caitlin Schein (Staff Physician), Christopher Ferguson (Dir, Laboratory Services), Helga Groat (Mgr, Laboratory), Jeremy Powers (Chief, Pathology), Kristen DiCicco (Mgr, Laboratory), Leah Korodan (Mgr, Division Laboratory), Muhammad Arshad (Chief, Pathology), Sarah Britton (VP, Laboratory Svcs), Stephanie Mullins (Supv, Laboratory)

## **9. Keywords**

Not Set