

Electrolytes, Hematocrit and Calculated Hemoglobin Using i-STAT Analyzer

Purpose

This procedure provides instructions for performing the i-STAT Chem8⁺ test, an in vitro quantification of electrolytes and hematocrit that uses arterial or venous whole blood. Lytes and H&H testing on the i-STAT 1 Analyzer (“i-STAT”) is classified as a moderately complex, CLIA non-waived test.

Scope

This procedure may be performed by trained personnel approved to perform Point of Care Testing (POCT) as defined in Chapter 510 California Business and Professional Code (CBPC) or per departmental policy, whichever is most restrictive. I-STAT testing is moderately complex non-waived testing performed by RN, CRNA and/or MD.

Policy

- The i-STAT 1 Analyzer will only be customized via the Data Exchange (DE) in Telcor QML by the POCT Coordinator to provide the primary information management capabilities
 - An electronic check is performed automatically every 8 hours by the i-STAT internal simulator.
 - Two levels of External Quality Control (QC) are tested with each new shipment and/or lot number of reagents, and at least every 8 hours of patient testing or at specified intervals found on a laboratory-established Individualized Quality Control Plan (IQCP).
 - The i-STAT is customized to block testing when it fails the internal electronic simulator test. It is also customized to block patient and proficiency testing if scheduled QC requirements are not met.
 - The i-STAT will be downloaded periodically to transmit results.
 - Proficiency testing through CAP is analyzed 2 times a year.
 - Instrument correlation is performed twice a year.
 - Positive patient identification must be made using two unique identifiers prior to patient testing.
 - Thermal probe and CLEW software updated are performed twice per year.
-

Electrolytes, Hematocrit and Calculated Hemoglobin Using i-STAT Analyzer

Specimen sources, collection & handling

Electrolytes / Hematocrit

- Venipuncture:
 - Whole blood collected in balanced heparin, lithium heparin syringes, direct syringes, or
 - Whole blood collected in evacuated tubes containing lithium heparin, as long as the tubes are filled to capacity
- Mix samples by rolling between palms and by gentle inversion for at least 5 seconds before testing.
- Clotted specimens should not be used.

Sample Transfer Device

- A 1cc syringe (such tuberculin) and needle (no smaller than 20 gauge) can be used to withdraw a sample from a blood collection tube.

Test Timing

Analyte	Syringes*	Test Timing	Evacuated Tubes	Test Timing
Sodium Potassium Chloride Glucose BUN/Urea Creatinine Hematocrit	With balanced heparin anticoagulant or lithium heparin anticoagulant (syringe must be filled to labeled capacity) • Remix thoroughly before filling cartridge.	30 minutes	With lithium heparin anticoagulant (tubes must be filled to labeled capacity) • Remix thoroughly before filling cartridge.	30 minutes
	No anticoagulant	3 minutes		

* Do Not Use Heparin lock flush solution syringes

Specimen rejection

- Hemolyzed samples
- Partial-draw blood collection tubes
- Drawing from an IV line.

Electrolytes, Hematocrit and Calculated Hemoglobin Using i-STAT Analyzer

The table below specifies storage requirements for the CHEM8⁺ cartridges.

CHEM8⁺ Cartridges

Description	Storage	
i-STAT Chem8+ Cartridges Vendor: Abbott	Refrigerated storage (2-8°C)	Store until expiration date. Do not use cartridges after the expiration date.
	Room temperature storage	Allow 5 minutes for a single cartridge to equilibrate to room temperature before opening pouches. Allow an entire box of cartridges to equilibrate to room temperature before use, approximately 1 hour. Cartridges may be stored at room temperature for 14 days. Do not expose to temperatures above 30°C. Do not return cartridges to the refrigerator after room temperature equilibration.
	After opening the pouch	Use cartridge immediately after opening pouch. If the pouch has been punctured, the cartridge should not be used.

Equipment



- i-STAT Analyzer 1
 - i-STAT Downloader
 - OMNI Print (portable printer)
-

Electrolytes, Hematocrit and Calculated Hemoglobin Using i-STAT Analyzer

Safety Precautions Handle i-STAT products using the standard safety precautions used when handling any potentially infectious material. Dispose of this product as biohazardous waste according to local, state, and national regulations.

Calibration Internal calibration is automatically performed as part of the test cycle on each cartridge. Should the i-STAT Analyzer become inoperable for any reason, immediately contact the laboratory for a replacement i-STAT.

Electronic Simulator Perform an electronic check on each handheld in use once a day with either the internal or external Electronic Simulator or as needed for regulatory compliance. The internal simulator check is initiated, every 24 hours or according to a customized schedule, when a cartridge is inserted into the cartridge port. If the internal simulator result is PASS, the cartridge test proceeds and the simulator results are stored. If FAIL is displayed for the internal simulator, reinsert the cartridge or use an external simulator.

Electronic Simulator Test	
Step	Action
1	Turn the handheld on by pressing the Power button. 
2	Press Menu to access the Administration Menu. 
3	Press 3 for Quality Tests.
4	Press 4 for Simulator.
5	Scan or enter Operator ID.
6	Enter the Simulator ID (serial number).
7	Insert the simulator into the cartridge port.
8	View results on the handheld's screen.
9	If PASS is displayed, continue to use the handheld.
10	If FAIL is displayed for the external simulator, reinsert the simulator. If FAIL is displayed for a second time, do not use the handheld and contact your Support Services representative.

Electrolytes, Hematocrit and Calculated Hemoglobin Using i-STAT Analyzer

Quality Control

Quality control testing should be done to confirm that the i-STAT Analyzer is working properly and providing reliable results. Only when controls are used routinely and the values are within acceptable ranges can accurate results be assured for patient samples.

Two levels of External QC are tested:

1. At least every 8 hours of patient testing or at specified intervals found on a laboratory established IQCP.
2. With each new shipment and/or lot number of reagents.
 - a. Run QC levels 1 and 3 using the current lot number and on the new lot number of cartridges.
 - b. Record all results on the appropriate log sheet.
 - c. Compare results for each level; results must correlated within Allowable Total Error (20%) before the new lot is put in use.

Description	Storage	
Electrolytes & Hematocrit Quality Control, Level 1 and Level 3	Refrigerated storage (2-8°C)	Store unopened until expiration date.
	Room temperature	Control solutions may also be stored at room temperature for up to 5 days (18 to 30 °C or 64 to 86 °F). Prolonged storage at temperatures greater than 30 °C (86°F) may cause changes in the values of some analytes. Do not use beyond the expiration date on the box and ampule labels.
	For best results, ampules, cartridges, and analyzer should be at the same temperature.	

Electrolytes, Hematocrit and Calculated Hemoglobin Using i-STAT Analyzer

Quality Control Testing Procedure	
Step	Action
1	Program the test: <ol style="list-style-type: none"> Press the ON/OFF key. Press the Menu key. Press 3 to select Quality Tests. Press 1 to select Control. Select fluid type to be run. Select 1-APOC for Abbott Point of Care liquid controls. Enter/scan Operator ID. Scan Control Lot Number. Scan Cartridge Lot Number
2	Immediately before use, shake the ampoule for 5-10 seconds to equilibrate the liquid and gas phase. If necessary, tap the tip of the ampoule to send solution back to the bottom section of ampoule.
3	Protect fingers with gauze, tissue, or glove. Then snap off the tip of the ampoule at the neck.
4	Immediately transfer the solution from the ampoule into a pipet, and then immediately transfer the solution tip into the cartridge, up to the fill line.
5	Immediately seal the cartridge and insert it into the analyzer. The transfer process from ampoule to cartridge must be expedient.
6	The numeric QC results are displayed on the screen with a "Pass" or "Fail" status. An overall determination of "Pass" is made if each analyte in the test panel passes.
Note: If results fail, do not test patients. Check the cartridges and repeat testing. Document on the i-STAT Corrective Action log.	

Electrolytes, Hematocrit and Calculated Hemoglobin Using i-STAT Analyzer




Before you begin

The i-STAT Analyzer 1 must remain on a level, vibration-free surface with the display facing up during testing. A level surface is also required when running the handheld into the downloader/recharger. The analyzer should remain level until a result is obtained.

All cartridges should remain in pouches until time of use. Use cartridges immediately after removing from its protective pouch. Do not contaminate the contact pads on the cartridges. Do not apply pressure to the central area of the label where the calibrant pack could pre-burst.

Procedure

Follow the steps below to perform patient testing on the i-STAT Analyzer 1.

Patient Test Procedure	
Step	Action
1	Turn the handheld on by pressing the Power button. 
2	Press 2 i-STAT Cartridge
3	Follow handheld prompts.
4	<p>Scan the lot number on the cartridge pouch.</p> <ol style="list-style-type: none"> Position barcode 3-9 inches from scanner window on the handheld. Press and hold SCAN to activate the scanner. Align the red laser light so it covers entire barcode. The handheld will beep when it reads the barcode successfully. 
5	<p>Continue normal procedures for preparing the sample, filling, and sealing cartridge.</p> <ol style="list-style-type: none"> Remove cartridge from its pouch. Avoid touching contact pads or exerting pressure over the calibrant pack in the center of the cartridge. Direct the dispensing tip into the sample well and dispense sample until it reaches the fill mark on the cartridge. Close the cover over the sample well until it snaps into place. Do not press over the sample well. 
6	Push the sealed cartridge into the handheld port until it clicks into place. Wait for the test to complete.
7	Review results and record results accordingly.

Electrolytes, Hematocrit and Calculated Hemoglobin Using i-STAT Analyzer

Reporting of Results

- Patient and internal procedural QC results must be recorded on the Patient Logs/Forms or directly in KP Health Connect or per facility requirement.
- Use DocFlow sheets for inpatient results in KP Health Connect.

To transmit results:

1. Dock the analyzer in the Downloader/Recharger.
2. When properly seated, the red light will turn on and the analyzer will automatically transmit all unsent results.
3. Do not move the analyzer while “Communication in progress” is displayed. If transmission is unsuccessful, contact your local Point-of-Care Coordinator (POCC).

To enter patient results in KP Health Connect (KPHC):

1. Log in to KPHC. Select EPIC, then select Hospital Chart.
2. Enter the patient’s medical record number. Select Accept.
3. Select Flowsheet, then select Point of Care Testing.
4. Type “i-STAT” in the search box, then scroll down to find line 544 POCT (Mod Complex according to age group).
5. Type the patient result in. Review results entered. Select Accept.

For unexpected results, when results do not reflect the patient’s condition, repeat the test using a fresh cartridge and sample. Refer to the i-STAT Technical Bulletin Analyzer Coded Messages.

Result Unable to be Displayed	Action
Results outside the system’s reportable range are flagged with a < or >, indicating that the result is below the lower limit or above the upper limit of the reportable range respectively.	Repeat test with new specimen.
Cartridge results that are not reportable based on internal QC rejection criteria are flagged with ***.	Repeat test with new specimen.
A Quality Check message will be reported instead of results if the analyzer detects a problem with the sample, calibrant solution, sensors, or mechanical or electrical functions of the analyzer during the test cycle.	Take the action displayed with the message that identifies the problem.
Call Abbott Technical Support or your local POCC if problem is not able to be resolved. Document all corrective actions taken.	

Reference Ranges

TEST MEASURED	UNITS	REPORTABLE RANGE	REFERENCE RANGES	
			arterial	venous
Na	mmol/L	100-180	138-146	
K	mmol/L	2.0-9.0	3.5-4.9	
Cl	mmol/L	65-140	98-109	
iCA	mmol/L	0.25-2.50	1.12-1.32	
Glu	mg/dL	20-700	70-105	
BUN/Urea	mg/dL	3-140	8-26	
Crea	mg/dL	0.2-20.0	0.6-1.3	
Hematocrit/Hct ***	%	15-75	Male: 43-51 Female: 38-46	

TEST CALCULATED	UNITS	REPORTABLE RANGE	REFERENCE RANGES	
			arterial	venous
AnGap	mmol	(-10) – (+99)	10-20	
Hemoglobin/Hb	g/dL	5.1-25.5	Male: 14-17 Female: 12-15.6	

Preventive Maintenance

Inspect and clean the exterior of the instrument as needed. Clean the display screen and the case using a gauze pad moistened with a mild non-abrasive cleaner, detergent, soap & water, alcohol, 10% bleach solution or PDI Super Sani-Cloth. Wipe the case using another gauze pad moistened with water and dry. Avoid getting excess fluids in the seam between the display screen and the case. Inspect and clean the Downloader as needed.

Limitations The analyte results should be assessed in conjunction with the patient’s medical history, clinical examination, and other findings. If results appear inconsistent with the clinical assessment, the patient sample should be retested using another cartridge.

**Non-
Controlled
Documents** The following non-controlled documents support this procedure.

- Procedure Manual for the iSTAT System
- i-STAT 1 User Guide, Abbott Point of Care
- i-STAT Technical Bulletin Analyzer Coded Messages
- i-STAT Chem8+

The following controlled documents support this procedure.

**Controlled
Documents**

Documents
i-STAT H/H & Lytes Verification Plan
i-STAT ACT Training and Competency Assessment

**Technical
Support** The manufacturer provides a toll-free line for technical support at **1-800-366-8020** or contact your local POCC.