



# KAISER PERMANENTE®

<b>DOCUMENT NUMBER:</b>
<b>DOCUMENT TITLE:</b>
<b>DOCUMENT NOTES:</b>

<b>LOCATION:</b>	<b>VERSION:</b>
<b>DOC TYPE:</b>	<b>STATUS:</b>

<b>EFFECTIVE DATE:</b>	<b>NEXT REVIEW DATE:</b>
<b>RELEASE DATE:</b>	<b>EXPIRATION DATE:</b>

<b>AUTHOR:</b>	<b>PREVIOUS NUMBER:</b>
<b>OWNER:</b>	<b>CHANGE NUMBER:</b>

## ESR Autoverification Rules in Cerner PathNet

---

**Introduction** All results generated from the interfaced Alcor iSED instrument for Erythrocyte Sedimentation Rate (ESR) cross to Cerner Millennium PathNet GenLab (GL), where auto-verification rules are applied. Results within reportable ranges are transmitted to the Cerner Millennium PathNet GL module and verified without CLS intervention.

---

**Scope** The autoverification rules in this document apply to laboratories that have opted and requested for Cerner autoverification to be turned on for the ESR analyzers.

---

**Policy** Autoverification rules for ESR test results from Alcor iSED analyzer:

1. Results from the iSED analyzer will be auto released for any ESR result values that are:
  - within the Alcor iSED analyzer's measurement range, i.e., from 1 to 130, AND
  - without instrument flags present.
2. For any result outside the linearity range (<1 and >130),
  - Cerner will hold for verification of sample volume and result integrity.
  - CLS will rerun samples, preferably on a different instrument if available.
  - CLS will manually release results in Cerner ARE.

---

*Continued on next page*

## ESR Autoverification Rules in Cerner PathNet, Continued

**Rule Testing and Expected Outcome**

- Autoverification table for ESR results

Condition Status	Result Value /Range	Rule	Expected Outcome
Low linearity limit	<1	Low linear	The result is held, requires review and manual verification by CLS
Within iSED s measurement range	1 – 130	Auto release	The result is auto released with appropriate flagging
High linearity limit	>130	High linear	The result is held, requires review and manual verification by CLS

**Rapid Suspension of Autoverification**

In the event of a problem with the test methodology, analytic instrument, or autoverification program, follow the SCPMG-LIS-0083 Procedure on Turning Autoverification On and Off in Cerner PathNet.

**Controlled Documents**

The following controlled documents support this policy.  
 Locally approved versions will have different document numbers.

Document Number	Title
SCPMG-LIS-0083	Procedure Turning Autoverification On and Off in Cerner PathNet

**Author(s)**

- SCPMG Hematology Working Group
- Eleanor E. Callasan
- Julius Salomon

*Continued on next page*

**ESR Autoverification Rules in Cerner PathNet, Continued**

**Attachment A: Annual Validation of ESR Autoverification Form**

---

Location: \_\_\_\_\_

Date: \_\_\_\_\_

Service Resource: \_\_\_\_\_

Performed by: \_\_\_\_\_

Accession Number	Result Value	Condition Status	Expected Outcome	Actual Outcome	Pass or Fail
	<1	Low Linear	The result is held, requires review and manual verification by CLS.		
		Within measurement range	The result is auto released with appropriate flagging.		
		Within measurement range	The result is auto released with appropriate flagging.		
		Within measurement range	The result is auto released with appropriate flagging.		
	>130	High Linear	The result is held, requires review and manual verification by CLS.		

Operations Director Review/Date: \_\_\_\_\_

CLIA Medical Director Review/Date: \_\_\_\_\_

\_\_\_\_\_

---

## Signature Manifest

---

**Document Number:** RIV-PPP-1138

**Revision:** 01

**Title:** ESR Autoverification Rules in Cerner PathNet

**Effective Date:** 14 Nov 2023

---

All dates and times are in Pacific Standard Time.

### ESR Autoverification Rules in Cerner PathNet

---

#### Operations Director Approval

---

Name/Signature	Title	Date	Meaning/Reason
Annaleah Raymond (Q741709)	Laboratory Operations Director	05 Nov 2023, 01:50:02 PM	Approved

#### Medical Director Approval

---

Name/Signature	Title	Date	Meaning/Reason
Mark Taira (P161328)	CLIA Director	08 Nov 2023, 10:11:30 PM	Approved