

## CK – CREATINE KINASE

**Intended Use** CK reagent, when used in conjunction with SYNCHRON LX<sup>®</sup> System(s), UniCel<sup>®</sup> Dx<sub>C</sub> 600/800 System(s), is intended for the quantitative determination of creatine kinase activity in human serum or plasma.

---

**Clinical Significance** Measurements of creatine kinase and its isoenzymes are used in the diagnosis and treatment of myocardial infarction and muscle diseases such as progressive, Duchenne-type muscular dystrophy.

---

**Methodology** CK reagent is used to measure the CK activity by an enzymatic rate method. In the reaction creatine kinase catalyzes the transfer of a phosphate group from the creatine phosphate substrate to adenosine diphosphate (ADP). The subsequent formation of adenosine triphosphate (ATP) is measured through the use of two coupled reactions catalyzed by hexokinase (HK) and glucose-6-phosphate dehydrogenase (G6PDH) which results in the production of reduced β-nicotinamide adenine dinucleotide phosphate (NADPH) from β-nicotinamide adenine dinucleotide phosphate (NADP). The CK assay contains the activator monothioglycerol.

The SYNCHRON<sup>®</sup> System(s) automatically proportions the appropriate sample and reagent volumes into the cuvette. The ratio used is one part sample to 20 parts reagent. The system monitors the change in absorbance at 340 nanometers. This change in absorbance is directly proportional to the activity of CK in the sample and is used by the System to calculate and express CK activity.

---

*Continued on next page*

## CK – CREATINE KINASE, *Continued*

---

### Specimen

---

#### **Type of Specimen**

PST (Lithium Heparin) is the specimen of choice.

---

#### **Specimen Storage and Stability**

1. Tubes of blood are to be kept closed at all times and in a vertical position. It is recommended that the serum or plasma be physically separated from contact with cells within two hours from the time of collection.
  2. Stability of CK activity in sera is not well defined, but is generally poor. Specimens should be assayed as soon after collection as possible since activity loss may occur after specimens have been stored for 4 hours at room temperature, 8 to 12 hours refrigerated or 2 to 3 days when frozen.
- 

#### **Sample Volume**

The optimum volume, when using a 0.5 mL sample cup, is 0.3 mL of sample. For optimum primary sample tube volumes and minimum volumes, refer to the Primary Tube Sample Template for your system.

---

*Continued on next page*


## CK – CREATINE KINASE, *Continued*

---


### Reagents

---

#### Reagent Preparation



For P/N 442635 (2 x 200 tests): Transfer the entire contents of the smallest reagent compartment (C) into the largest reagent compartment (A).



For P/N 476836 (2 x 400 tests): Transfer all the contents of one bottle CK (A-reagent) into the largest reagent compartment (A).

Replace cartridge caps and gently invert the cartridge several times to ensure adequate mixing.

#### Reagent Storage and Stability

CK reagent, when stored unopened at +2°C to +8°C, will remain stable until the expiration date printed on the cartridge label. Once prepared, the reagent cartridge is stable for 30 days at +2°C to +8°C unless the expiration date is exceeded. DO NOT FREEZE.

---

*Continued on next page*

## CK – CREATINE KINASE, *Continued*

---

### **Calibration**

---

### **Calibration Required**

Calibration is not required.

---

### **Quality Control**

See Beckman Policy QC LCHS-5020

---

### **Sample Processing**

See Beckman Policy QC LCHS-5020

---

### **Reference Range**

Male: 15 – 190 IU/L  
Female: 15 – 170 IU/L

---

### **Analytic Range**

<b>Sample Type</b>	<b>Conventional Units</b>
Serum or Plasma	5 - 120000 IU/L
Serum or Plasma (ORDAC)	860–4100 IU/L

---

### **Author**

Benjamin J. Salas, Jr., CLS, MLS (ASCP)CM, CLS(NCA)

---

### **Distributions**

Kaiser Permanente Riverside Medical Center Laboratory

---

*End*

## CK – CREATINE KINASE, *Continued*

Reviewed and approved by:

SIGNATURE	DATE
<b>Denise Topliff, MHA, CLS</b> <b>Director, Clinical Laboratory – Riverside Medical Center</b>	
<b>Mark Taira, M.D.</b> <b>Medical Director, Clinical Laboratory – Riverside Medical Center</b>	

