

BIO-RAD

Autoantibody Atlas



DEVELOPED BY

**The University Of Washington
Department of Laboratory Medicine
Seattle, Washington**

Kathleen W. Hutchinson, MS, MT(ASCP)
Mark H. Wener, MD
Bruce G. Gilliland, MD
Michael L. Astion, MD, PhD

List of Antibodies

[ANA](#)

[Anti Scl-70](#)

[ANCA](#)

[Anti Sm](#)

[Anti Beta-2 Glycoprotein I](#)

[Anti Smooth Muscle](#)

[Anti Cardiolipin](#)

[Anti SS-A](#)

[Anti CCP](#)

[Anti SS-B](#)

[Anti Centromere](#)

[RF](#)

[Anti Chromatin](#)

[Anti dsDNA](#)

[Anti ENA](#)

[Anti GBM](#)

[Anti Histone](#)

[Anti Jo-1](#)

[Anti Mitochondria](#)

[Anti MPO](#)


[Anti PR3](#)

[Anti Parietal](#)

[Anti Ribosomal-P](#)

[Anti RNP](#)

ANA



Sample ANA staining patterns on HEp-2 Substrate

Antibody

ANA
Anti nuclear antibodies

Antigen

Many different antigens that are found in the cell nucleus.

Disease Associations

Dependent upon the antibody specificity and antibody level: systemic autoimmune diseases and autoimmune liver disease. May also be seen in patients with some acute or chronic infections, in association with some medications, and in healthy individuals.

Testing Methods

IFA and ELISA (screening, not specific identification), Multiplexing.


Typical IFA Pattern

Many. See ANA Tutor

Notes

The sensitivity and specificity is very dependent upon the method used for ANA screening. The titer and frequency of anti nuclear antibodies in healthy individuals increases with age.

ANCA



Sample ANCA staining patterns

Antibody

ANCA

Anti neutrophil cytoplasmic antibodies

Antigen

A group of several different antigens that are found in neutrophils. Includes proteinase-3 (PR3) and myeloperoxidase (MPO). Some people also include other neutrophil specific autoantibodies (NSA) such as elastase (EL), and bacterial permeability increase (BPI) protein, Cathepsin G, lysozyme, lactoferrin, and catalase.

Disease Associations

Dependent upon the antibody specificity: ANCA-associated vasculitis, inflammatory bowel disease (Crohn's disease, ulcerative colitis), systemic autoimmune diseases.

Testing Methods

IFA (screening, not specific identification)

Typical IFA Pattern

Dependent on antibody specificity: C-ANCA, P-ANCA, atypical ANCA (see ANCA Tutor)

Notes

Clinical assays are not available for specific identification of many of these antibodies. ANA staining may resemble ANCA staining.

Anti Beta-2 Glycoprotein I



Antibody

IgG anti beta-2 glycoprotein I,
IgA anti beta-2 glycoprotein I,
IgM anti beta-2 glycoprotein I

Alternate Names

Anti β 2GPI, a β 2GPI

Antigen

Glycoprotein that binds to cardiolipin and other anionic phospholipids (cardiolipin cofactor). Also called apolipoprotein H.

Disease Associations

Anti phospholipid syndrome (APS), systemic lupus erythematosus (SLE).

Testing Methods

ELISA


Typical IFA Pattern

Not applicable

Notes

Assays for different isotypes are available (i.e. IgG, IgA, IgM). IgG and IgM antibodies results are included in the 2006 APS Classification Criteria. Thought to be more specific for APS than anti-cardiolipin. APS is characterized by recurrent thrombosis. In pregnant women, antibodies to phospholipids have been associated with repeated, unexplained fetal loss.

Anti Cardiolipin



Phospholipid

Antibody

Anti Cardiolipin

IgG anti cardiolipin, IgA anti cardiolipin, IgM anticardiolipin

Alternate Names

aCL

Antigen

Diphosphatidylglycerol lipid. Originally derived from cow heart tissue and used as the antigen in the Wassermann test for syphilis.

Disease Associations

Anti phospholipid syndrome (APS), systemic lupus erythematosus (SLE), syphilis. IgM assay may be elevated in patients with infections, rheumatoid factor, or cryoglobulins.

Testing Methods

ELISA


Typical IFA Pattern

Not applicable

Notes

Assays are available for different immunoglobulin classes (i.e. IgG, IgA, IgM). IgG and IgM antibodies results are included in the 2006 APS Classification Criteria. One of many aPL (anti phospholipids) that include antibodies to cardiolipin, B2GP1, phosphatidylserine/prothombin, and lupus anticoagulant. May cause a biological false positive syphilis test.

Anti CCP



Antibody

Anti CCP

Anti cyclic citrullinated peptide

Alternate Names

Anti-citrullinated peptide antibody (ACPA), Anti-Sa

Antigen

Antibodies directed against citrullinated proteins and peptides (ACPA) are present in sera from a majority of patients with rheumatoid arthritis. The most common method of testing for ACPA is by detecting antibodies to synthetic cyclic (ring structured) citrullinated peptides (CCP). There are a variety of synthetic peptides that are used for clinical testing. Successive generations of tests are assigned different numbers (e.g. CCP2, CCP3, etc.). Some citrullinated peptides used in testing are derived from the sequences in filaggrin or vimentin.

Disease Associations

Rheumatoid arthritis (RA).

Testing Methods

ELISA, Western blot (anti-Sa)


Typical IFA Pattern

Oral mucosal cells: Anti perinuclear factor (APF), Esophageal epithelial cells: anti keratin (AKA) (Not routinely used for clinical testing)

Notes

Anti CCP has better specificity for RA than rheumatoid factor does. Anti-CCP levels may correlate with disease prognosis. Citrulline is formed by a post translational modification of arginine.

Anti Centromere



Antibody

Anti Centromere

Alternate Names

Anti centromere associated protein (CENP)

Antigen

Proteins present in the centromere (kinetochore) region of the chromosome that includes CENP-A, CENP-B, CENP-C.

Disease Associations

Most common association is with the limited cutaneous form of scleroderma (CREST syndrome). May be present in other systemic autoimmune diseases and primary biliary cirrhosis (PBC). Patients with primary Raynaud's phenomenon and anti-centromere have an increased risk for later development of a systemic autoimmune disease.

Testing Methods

ELISA, Multiplexing, IFA

Typical IFA Pattern


HEp-2: Discrete speckles of the interphase cell nuclei and dividing cell chromatin. May be difficult to identify in the presence of other nuclear or cytoplasmic antibodies.

Notes

Most patient sera with centromere antibodies react with CENP-A, CENP-B and CENP-C.

HEp-2 Centromere Pattern, 100x Objective

Anti Chromatin



Antibody

Anti Chromatin

Alternate Names

Anti Nucleosome

Antigen

Epitopes present on histones that are bound to DNA in the nucleosome structure, i.e. the DNA-histone complex.

Disease Associations

Present in the sera of patients with drug induced lupus, SLE, and several other systemic autoimmune diseases. Some reports indicate that the antibody may play a pathogenic role and be helpful in the diagnosis and establishing prognosis of SLE patients.

Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)

Typical IFA Pattern


HEp-2: Homogeneous ANA, pattern is not unique therefore does not permit specific identification.

Notes

Anti chromatin refers to a family of autoantibodies that react to histones bound to DNA.

Nucleosome structure

Anti dsDNA



Antibody

Anti dsDNA

Anti double stranded deoxyribonucleic acid

Alternate Names

Anti native DNA

Antigen

Double stranded DNA (deoxyribonucleic acid)

Disease Associations

SLE

Testing Methods

ELISA, Multiplexing, IFA on Crithidia, ANA IFA (screening, not specific identification)

Typical IFA Pattern


HEp-2: Homogeneous ANA, pattern is not unique therefore does not permit specific identification. Crithidia: Staining of the kinetoplast.

Notes

Antibody levels often correlate with disease activity, especially high avidity antibodies. Many assays detect both low and high avidity antibodies.

HEp-2 (top), Crithidia (bottom)

Anti ENA



Antibody

Anti ENA

Anti Extractable Nuclear Antigens

Antigen

Group of several different antigens that are easily extracted from the cell nucleus including Sm/RNP, SS-A, SS-B, Scl-70, and Jo-1.

Disease Associations

Systemic autoimmune diseases

Testing Methods

ELISA, Multiplexing, Double Diffusion, CIE (counterimmunoelectrophoresis)

Typical IFA Pattern


HEp-2: most, but not all produce speckled nuclear staining.

Notes

It is common to use a screening method and then proceed to follow-up testing for specific antibodies if the screen is positive. The screening method used will determine which antibodies can be detected.

ELISA plate (top), Double diffusion gel (bottom)

Anti GBM



Antibody

Anti GBM

Anti glomerular basement membrane

Alternate Names

Goodpasture-antibodies

Antigen

Type IV collagen. Most antibodies react with a cryptic epitope located on the C-terminal alpha 3 NC1 domain.

Disease Associations

Goodpasture's syndrome.

Testing Methods

ELISA, Multitplexing, IFA on primate kidney


Typical IFA Pattern

Primate Kidney: Linear staining of the glomerular basement membrane

Notes

Anti-GBM levels may correlate with disease activity. Patients with acute renal failure and/or severe lung hemorrhage are often screened for ANCA and GBM antibodies.

Anti Histone



Antibody

Anti Histone

Antigen

Histone proteins (H1, H2A, H2B, H3, H4, H5)

Disease Associations

Present in the sera of patients with drug induced lupus, spontaneous SLE, and several other systemic autoimmune diseases.

Testing Methods

ELISA, IFA (screening, not specific identification)

Typical IFA Pattern

HEp-2: Homogeneous ANA, pattern is not unique enough for specific identification.


Crithidia: Negative

Notes

Autoantibodies that only react to histones bound to DNA are usually called anti-nucleosome or anti-chromatin (see anti-chromatin section).

HEp-2 Homogeneous Pattern, 40x Objective

Anti Jo-1



Antibody

Anti Jo-1

Alternate Names

Anti histidyl-tRNA synthetase

Antigen

Anti histidyl-tRNA synthetase

Disease Associations

Polymyositis and dermatomyositis

Testing Methods

ELISA, Multiplexing

Typical IFA Pattern


HEp-2: negative ANA, may have weak cytoplasmic staining that is not unique therefore does not permit specific identification.

Notes

Usually considered to be one of the ENA (extractable nuclear antigens) antibodies. Although the antibody is only present in the sera of approximately 20% of patients with myositis, it is the most common of the myositis specific antibodies (MSA).

HEp-2 Negative with weak cytoplasmic staining

Anti Mitochondria



HEP-2 (top), Mouse Stomach Kidney (bottom)

Antibody

Anti - mitochondria

Alternate Names

AMA

Antigen

Nine different subtypes have been identified (M1-M9). Antibodies to M2 are considered to be the most clinically significant. Antibodies to M2 react with the E2 subunit of the pyruvate dehydrogenase complex (PDC-E2).

Disease Associations

Over 90% of patients with primary biliary cirrhosis(PBC) have mitochondrial antibodies.

Testing Methods

ELISA (anti-M2) IFA using mouse stomach kidney substrate (detects multiple subtypes)
Immunoblot (not generally available in clinical labs)

Typical IFA Pattern


MSK: cytoplasmic staining of the renal tubules and stomach parietal cells.

HEp-2: numerous granular speckles in the cytoplasm.

Notes

There are other antibodies that may resemble the mitochondrial staining seen on HEp-2. Different techniques should be used to identify the antibody.

Anti MPO



Ethanol (top) and Formalin fixed neutrophils

Antibody

Anti MPO
Anti myeloperoxidase

Antigen

Myeloperoxidase. A cationic enzyme found in the azurophilic granules of neutrophils.

Disease Associations

ANCA-associated vasculitis. Strongest association is with idiopathic crescentic glomerulonephritis (NCGN) and microscopic polyangiitis (MPA), less common in and Churg-Strauss Syndrome (CSS) and Wegener's granulomatosis (WG).

Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)


Typical IFA Pattern

Ethanol Fixed Neutrophils: P-ANCA (perinuclear) Formalin Fixed Neutrophils: cytoplasmic granules

Notes

One of the anti-neutrophil cytoplasmic antibodies (ANCA). Anti-MPO is the most common antigen associated with the P-ANCA pattern. ANAs may resemble P-ANCA staining (see ANCA Tutor). Anti-MPO levels may correlate with disease activity. Patients with acute renal failure and/or severe lung hemorrhage are often screened for ANCA and GBM antibodies.

Anti PR3



Ethanol (top) and Formalin fixed neutrophils

Antibody

Anti PR3
Anti Proteinase 3

Antigen

Proteinase 3. A cationic enzyme found in the azurophilic granules of neutrophils.

Disease Associations

ANCA-associated vasculitis. Strongest association is with Wegener's granulomatosis (WG), less often seen in microscopic polyangiitis (MPA) and Churg-Strauss Syndrome (CSS).

Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)


Typical IFA Pattern

Ethanol Fixed Neutrophils: C-ANCA (cytoplasmic) Formalin Fixed Neutrophils: cytoplasmic granules

Notes

One of the anti-neutrophil cytoplasmic antibodies (ANCA). Antibody levels often correlate with disease activity. Patients with acute renal failure and/or severe lung hemorrhage are often screened for ANCA and glomerular basement membrane (GBM) antibodies.

Anti Parietal Cell



Antibody

Anti parietal cell

Alternate Names

PCA (parietal cell antibodies), GPC antibodies (gastric parietal cell antibodies)

Antigen

α and β subunits of the H⁺/K⁺ ATPase

Disease Associations

Pernicious anemia

Testing Methods

IFA using mouse stomach kidney substrate

ELISA (anti-GPC)

Typical IFA Pattern


MSK: cytoplasmic staining of the stomach parietal cells. Renal tubules are negative.

Notes

This antibody cannot be detected by IFA if antibodies to mitochondria are present. The antibody frequency in the normal population increases with age, and the antibody is more common in women than in men.

Mouse stomach kidney

Anti Ribosomal-P



Antibody

Anti Ribosomal P

Antigen

Ribosomal proteins P0, P1, and P2

Disease Associations

SLE

Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)

Typical IFA Pattern


HEp-2: cytoplasmic staining that is not unique therefore does not permit specific identification, may have coexisting nucleolar staining. MSK: cytoplasmic staining of the chief cells

Notes

Antigen is present in thymus extract used for ENA (extractable nuclear antigens) antibody testing by double diffusion methods, but is not typically included in ELISA ENA screening kits.

HEp-2 cytoplasmic (top) MSK chief cells (bottom)

Anti RNP



Antibody

Anti RNP
Anti ribonucleoprotein

Alternate Names

Anti U1-RNP, anti U1 snRNP

Antigen

Small ribonucleoproteins in spliceosomes, several proteins including 68-70kd, A, and C

Disease Associations

Several systemic autoimmune diseases including SLE, Sjögren's Syndrome, scleroderma, myositis, occasionally rheumatoid arthritis. Patients with mixed connective tissue disease tend to have very high levels of the antibody, without other autoantibodies.

Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification), immunodiffusion.

Typical IFA Pattern


HEp-2: Speckled ANA, pattern is not unique, therefore does not permit specific identification

Notes

Difficult to prepare the RNP antigens without Sm epitopes also present. Some assays use antigens that have both the Sm and RNP epitopes, these assays are called anti Sm/RNP. One of the ENA (extractable nuclear antigens) antibodies.

HEp-2 Speckled Pattern, 40x Objective

Anti Scl-70



Antibody

Anti Scl-70
Anti Scleroderma 70

Alternate Names

Anti Topo I, Anti topoisomerase I

Antigen

DNA topoisomerase I

Disease Associations

Most common association is with the diffuse cutaneous form of scleroderma. Occasionally present in the serum of patients with limited cutaneous scleroderma.

Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)

Typical IFA Pattern


HEp-2: Homogeneous ANA, often with coexisting nucleolar staining. Pattern is not unique therefore does not permit specific identification.

Notes

One of the ENA (extractable nuclear antigens) antibodies.

HEp-2, Homogenous with nucleolar staining

Anti Sm



Antibody

Anti Sm

Alternate Names

Anti Smith

Antigen

Small ribonucleoproteins in spliceosomes, several proteins including B'/B and D

Disease Associations

SLE

Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification), immunodiffusion

Typical IFA Pattern


HEp-2: Speckled ANA, pattern is not unique therefore does not permit specific identification.

Notes

Almost all samples with anti Sm have coexisting antibodies to RNP. One of the ENA (extractable nuclear antigens) antibodies. Some assays use antigens that have both the Sm and RNP epitopes, these assays are called anti Sm/RNP.

HEp-2 Speckled Pattern, 40x Objective

Anti Smooth Muscle



Antibody

Anti-smooth muscle

Alternate Names

SMA

Antigen

Several potential antigens, including F-actin

Disease Associations

Autoimmune hepatitis (type 1)

Testing Methods

IFA using mouse stomach kidney substrate, ELISA (anti-F-actin)

Typical IFA Pattern


MSK: cytoplasmic staining of the stomach smooth muscle cells and staining of the arterial blood vessel walls.

Notes

Low level staining of the smooth muscle layer is fairly common. Reference range studies should be performed to insure the proper cut-off is used.

Mouse stomach kidney

Anti SS-A



Antibody

Anti SS-A

Anti Sjögren's Syndrome A

Alternate Names

Anti SSA, Anti Ro, Anti SS-A/Ro, Anti- Ro/SS-A

Antigen

Small ribonucleoproteins, two different epitopes: 52kd and 60 kd

Disease Associations

Sjögren's Syndrome, SLE (particularly subacute cutaneous lupus variant), Neonatal Lupus, Congenital Heart Block

Testing Methods

ELISA, Multiplexing

Typical IFA Pattern


Hep-2: Speckled ANA, pattern is not unique therefore does not permit specific identification.

Notes

ANA by IFA may be weak or absent for some anti-SSA antibodies. One of the ENA (extractable nuclear antigens) antibodies.

Hep-2 Speckled Pattern, 40x Objective

Anti SS-B



Antibody

Anti SS-A
Anti Sjögren's Syndrome A

Alternate Names

Anti SSB, Anti La

Antigen

Small ribonucleoproteins

Disease Associations

Sjögren's Syndrome, SLE, Neonatal Lupus, Congenital Heart Block

Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)

Typical IFA Pattern


HEp-2: Speckled ANA, may also have some nucleolar staining, pattern is not unique therefore does not permit specific identification.

Notes

Most samples with anti SS-B have coexisting antibodies to SS-A. One of the ENA (extractable nuclear antigens) antibodies.

HEp-2 Speckled Pattern, 40x Objective

Rh

**Antibody**

RF
Rheumatoid factor

Antigen

Epitopes are in the constant region of the IgG heavy chain.

Disease Associations

Rheumatoid arthritis (RA) and Sjögren's syndrome. Also seen in the sera of patients with other systemic autoimmune disease, mixed cryoglobulinemia, some chronic bacterial and viral infections, and in the sera of healthy individuals (particularly the elderly) in low titer.

Testing Methods

Agglutination, Nephelometry/Turbidometry, ELISA

Typical IFA Pattern

Not applicable

Notes

The name rheumatoid factor was given before the discovery that it was an antibody to IgG. There are several different isotypes of rheumatoid factor, IgM, IgA, and IgG (including the subclasses IgG1, IgG2, IgG3, and IgG4). Most clinical assays in the US detect only IgM RF.

IgM RF binding to normal IgG