
	<b>Coag Specimen Processing</b>  <b>CCL-026</b>	<b>Dept: 324318</b>	<b>Critical Care Labs</b>
		<b>Effective Date:</b>	<b>Feb 2015</b>
		<b>Revised Date:</b>	<b>Feb 2019</b>
		<b>Contact:</b>	<b>Ann Shoffner</b>
<b>Name &amp; Title: Gregory Pomper, MD</b> <b>CLIA Laboratory Director</b>		<b>Date:</b>	<b>2 / 25 / 19</b>
<b>Signature:</b> 			

**1) General Procedure Statement:**

- a. **Purpose:** As of June 2<sup>nd</sup> 2014, the OR Lab no longer performs Coagulation testing. This procedure provides information to the OR Lab staff regarding pre-analytical processing of Coagulation samples for the Core Laboratory.
- b. **Responsible Department/Party/Parties:**
  - i. Procedure owner: Ann Shoffner
  - ii. Procedure: Critical Care Lab (CCL) Staff
  - iii. Supervision: Ann Shoffner
  - iv. Implementation: Ann Shoffner and Critical Care Lab (CCL) Staff

**2) PPE Requirements: Lab Coat, Gloves**

**3) Procedure:**

**A. Coag Sample Requirements:**

1. Citrated plasma should be obtained using a mix of freshly collected blood with one part of 3.2% sodium citrate. Invert the tube gently three or four times immediately after venipuncture to ensure proper mixing of blood to anticoagulant.
2. A syringe or evacuated tubes may be used for collection. If an evacuated tube system is used, the coagulation sample should be the first tube collected.
3. If blood is drawn from a vascular access device (VAD), the line should be flushed with 5.0 mL saline and the first 5.0ml of blood or six dead space volumes of the VAD discarded.
4. Specimens containing clots, collected in the wrong tube or have less than 90% expected fill of the collection tube should be rejected. The plasma volume required for testing is 50uL with a dead volume of 100uL.

**Icteric, Lipemic and/or Turbid samples** – should be forwarded to the Main Coag Lab for verification using airfuge centrifugation, BSC methodology or Fibrometer.

**Hemolyzed samples** – should be noted in the LIS as hemolyzed.

**Hct values >55%** - The citrate concentration of the collection tube must be adjusted in patients who have hematocrit values greater than 55%. Visually inspect the sample. If greater than half of the sample appears to be red cells, suspect an increased Hct. Review the patient's Hct history in the LIS. If a Hct > 55% is suspected, request a special tube from the Core Lab and ask the ordering location to recollect the sample using the special tube.

**B. Coag Sample Preparation, Handling and Processing:**

1. The sample should be transported at room temperature.
2. Coagulation samples should be processed without delay.
3. Check the sample and requisition for proper ID and matching.

4. Check the sample for clot formation using wooden applicator sticks and visual inspection prior to centrifugation.
5. Place the sample in the centrifuge as soon as it arrives in the lab.
6. Centrifuge the blood specimen for a minimum of 3 minutes at 7200 rpm as soon as possible after collection to produce platelet poor plasma (platelet count  $<10 \times 10^3/\mu\text{L}$ ).
7. Patient plasma should be tested within two hours after collection if stored at room temperature, four hours if stored at 2-4°C, within two weeks if stored at -20°C, and within 6 months at -70°C. If immediate testing is to be done, the plasma may remain on the packed cells or separated. To separate plasma, use a plastic transfer pipette; remove the plasma to a plastic tube without disturbing the buffy coat. Do not store the sample on ice. If testing is not complete within 4 hours, the plasma must be frozen. Frozen plasma samples must be rapidly thawed at 37°C while gently mixing and tested immediately after thawing. Mixing is critical before testing.
8. Order/Receive the testing in the LIS.
9. Manually complete a "call to" slip (*CCL-F052 OR Lab Call To Form*). Put a LIS taglet on the call slip and circle the patient's OR Room #. If a CBC is also ordered, complete a separate "call to" slip and give each sample a separate stat bag. This will expedite delivery to the appropriate lab area.
10. Place the call slip in the pocket of a red stat bag.

Place LIS taglet in box below

PLEASE CALL

Operating Suite Phone List

Man Desk	3-2440	PACU	3-4767	Cyto	3-31
Houska/Plt	3-2681	OP PACU	3-2899	Endo	3-77
Ann Burg	3-2944	Ann Room	6-8908		

CPH	3-2700	CPH	3-2716	PREP	CPH	3
CPH	3-2702	CPH	3-2708	CPH	3-27	
CPH	3-2703	CPH	3-2721	CPH	3-27	
CPH	3-2704	CPH	3-2722	CPH	3-27	
CPH	3-2705	CPH	3-2723	CPH	3-27	
CPH	3-2706	CPH	3-2724	CPH	3-27	
CPH	3-2707	CPH	3-2725	CPH	3-27	
CPH	3-2708	CPH	3-2726	CPH	3-27	
CPH	3-2709	CPH	3-2727	CPH	3-27	
CPH	3-2710	CPH	3-2728	CPH	3-27	
CPH	3-2711	CPH	3-2729	CPH	3-27	
CPH	3-2712	CPH	3-2730	CPH	3-27	
CPH	3-2713	CPH	3-2731	CPH	3-27	
CPH	3-2714	CPH	3-2732	CPH	3-27	
CPH	3-2715	CPH	3-2733	CPH	3-27	
CPH	3-2716	CPH	3-2734	CPH	3-27	
CPH	3-2717	CPH	3-2735	CPH	3-27	
CPH	3-2718	CPH	3-2736	CPH	3-27	
CPH	3-2719	CPH	3-2737	CPH	3-27	
CPH	3-2720	CPH	3-2738	CPH	3-27	
CPH	3-2721	CPH	3-2739	CPH	3-27	
CPH	3-2722	CPH	3-2740	CPH	3-27	
CPH	3-2723	CPH	3-2741	CPH	3-27	
CPH	3-2724	CPH	3-2742	CPH	3-27	
CPH	3-2725	CPH	3-2743	CPH	3-27	
CPH	3-2726	CPH	3-2744	CPH	3-27	
CPH	3-2727	CPH	3-2745	CPH	3-27	
CPH	3-2728	CPH	3-2746	CPH	3-27	
CPH	3-2729	CPH	3-2747	CPH	3-27	
CPH	3-2730	CPH	3-2748	CPH	3-27	
CPH	3-2731	CPH	3-2749	CPH	3-27	
CPH	3-2732	CPH	3-2750	CPH	3-27	
CPH	3-2733	CPH	3-2751	CPH	3-27	
CPH	3-2734	CPH	3-2752	CPH	3-27	
CPH	3-2735	CPH	3-2753	CPH	3-27	
CPH	3-2736	CPH	3-2754	CPH	3-27	
CPH	3-2737	CPH	3-2755	CPH	3-27	
CPH	3-2738	CPH	3-2756	CPH	3-27	
CPH	3-2739	CPH	3-2757	CPH	3-27	
CPH	3-2740	CPH	3-2758	CPH	3-27	
CPH	3-2741	CPH	3-2759	CPH	3-27	
CPH	3-2742	CPH	3-2760	CPH	3-27	
CPH	3-2743	CPH	3-2761	CPH	3-27	
CPH	3-2744	CPH	3-2762	CPH	3-27	
CPH	3-2745	CPH	3-2763	CPH	3-27	
CPH	3-2746	CPH	3-2764	CPH	3-27	
CPH	3-2747	CPH	3-2765	CPH	3-27	
CPH	3-2748	CPH	3-2766	CPH	3-27	
CPH	3-2749	CPH	3-2767	CPH	3-27	
CPH	3-2750	CPH	3-2768	CPH	3-27	
CPH	3-2751	CPH	3-2769	CPH	3-27	
CPH	3-2752	CPH	3-2770	CPH	3-27	
CPH	3-2753	CPH	3-2771	CPH	3-27	
CPH	3-2754	CPH	3-2772	CPH	3-27	
CPH	3-2755	CPH	3-2773	CPH	3-27	
CPH	3-2756	CPH	3-2774	CPH	3-27	
CPH	3-2757	CPH	3-2775	CPH	3-27	
CPH	3-2758	CPH	3-2776	CPH	3-27	
CPH	3-2759	CPH	3-2777	CPH	3-27	
CPH	3-2760	CPH	3-2778	CPH	3-27	
CPH	3-2761	CPH	3-2779	CPH	3-27	
CPH	3-2762	CPH	3-2780	CPH	3-27	
CPH	3-2763	CPH	3-2781	CPH	3-27	
CPH	3-2764	CPH	3-2782	CPH	3-27	
CPH	3-2765	CPH	3-2783	CPH	3-27	
CPH	3-2766	CPH	3-2784	CPH	3-27	
CPH	3-2767	CPH	3-2785	CPH	3-27	
CPH	3-2768	CPH	3-2786	CPH	3-27	
CPH	3-2769	CPH	3-2787	CPH	3-27	
CPH	3-2770	CPH	3-2788	CPH	3-27	
CPH	3-2771	CPH	3-2789	CPH	3-27	
CPH	3-2772	CPH	3-2790	CPH	3-27	
CPH	3-2773	CPH	3-2791	CPH	3-27	
CPH	3-2774	CPH	3-2792	CPH	3-27	
CPH	3-2775	CPH	3-2793	CPH	3-27	
CPH	3-2776	CPH	3-2794	CPH	3-27	
CPH	3-2777	CPH	3-2795	CPH	3-27	
CPH	3-2778	CPH	3-2796	CPH	3-27	
CPH	3-2779	CPH	3-2797	CPH	3-27	
CPH	3-2780	CPH	3-2798	CPH	3-27	
CPH	3-2781	CPH	3-2799	CPH	3-27	
CPH	3-2782	CPH	3-2800	CPH	3-27	
CPH	3-2783	CPH	3-2801	CPH	3-27	
CPH	3-2784	CPH	3-2802	CPH	3-27	
CPH	3-2785	CPH	3-2803	CPH	3-27	
CPH	3-2786	CPH	3-2804	CPH	3-27	
CPH	3-2787	CPH	3-2805	CPH	3-27	
CPH	3-2788	CPH	3-2806	CPH	3-27	
CPH	3-2789	CPH	3-2807	CPH	3-27	
CPH	3-2790	CPH	3-2808	CPH	3-27	
CPH	3-2791	CPH	3-2809	CPH	3-27	
CPH	3-2792	CPH	3-2810	CPH	3-27	
CPH	3-2793	CPH	3-2811	CPH	3-27	
CPH	3-2794	CPH	3-2812	CPH	3-27	
CPH	3-2795	CPH	3-2813	CPH	3-27	
CPH	3-2796	CPH	3-2814	CPH	3-27	
CPH	3-2797	CPH	3-2815	CPH	3-27	
CPH	3-2798	CPH	3-2816	CPH	3-27	
CPH	3-2799	CPH	3-2817	CPH	3-27	
CPH	3-2800	CPH	3-2818	CPH	3-27	
CPH	3-2801	CPH	3-2819	CPH	3-27	
CPH	3-2802	CPH	3-2820	CPH	3-27	
CPH	3-2803	CPH	3-2821	CPH	3-27	
CPH	3-2804	CPH	3-2822	CPH	3-27	
CPH	3-2805	CPH	3-2823	CPH	3-27	
CPH	3-2806	CPH	3-2824	CPH	3-27	
CPH	3-2807	CPH	3-2825	CPH	3-27	
CPH	3-2808	CPH	3-2826	CPH	3-27	
CPH	3-2809	CPH	3-2827	CPH	3-27	
CPH	3-2810	CPH	3-2828	CPH	3-27	
CPH	3-2811	CPH	3-2829	CPH	3-27	
CPH	3-2812	CPH	3-2830	CPH	3-27	
CPH	3-2813	CPH	3-2831	CPH	3-27	
CPH	3-2814	CPH	3-2832	CPH	3-27	
CPH	3-2815	CPH	3-2833	CPH	3-27	
CPH	3-2816	CPH	3-2834	CPH	3-27	
CPH	3-2817	CPH	3-2835	CPH	3-27	
CPH	3-2818	CPH	3-2836	CPH	3-27	
CPH	3-2819	CPH	3-2837	CPH	3-27	
CPH	3-2820	CPH	3-2838	CPH	3-27	
CPH	3-2821	CPH	3-2839	CPH	3-27	
CPH	3-2822	CPH	3-2840	CPH	3-27	
CPH	3-2823	CPH	3-2841	CPH	3-27	
CPH	3-2824	CPH	3-2842	CPH	3-27	
CPH	3-2825	CPH	3-2843	CPH	3-27	
CPH	3-2826	CPH	3-2844	CPH	3-27	
CPH	3-2827	CPH	3-2845	CPH	3-27	
CPH	3-2828	CPH	3-2846	CPH	3-27	
CPH	3-2829	CPH	3-2847	CPH	3-27	
CPH	3-2830	CPH	3-2848	CPH	3-27	
CPH	3-2831	CPH	3-2849	CPH	3-27	
CPH	3-2832	CPH	3-2850	CPH	3-27	
CPH	3-2833	CPH	3-2851	CPH	3-27	
CPH	3-2834	CPH	3-2852	CPH	3-27	
CPH	3-2835	CPH	3-2853	CPH	3-27	
CPH	3-2836	CPH	3-2854	CPH	3-27	
CPH	3-2837	CPH	3-2855	CPH	3-27	
CPH	3-2838	CPH	3-2856	CPH	3-27	
CPH	3-2839	CPH	3-2857	CPH	3-27	
CPH	3-2840	CPH	3-2858	CPH	3-27	
CPH	3-2841	CPH	3-2859	CPH	3-27	
CPH	3-2842	CPH	3-2860	CPH	3-27	
CPH	3-2843	CPH	3-2861	CPH	3-27	
CPH	3-2844	CPH	3-2862	CPH	3-27	
CPH	3-2845	CPH	3-2863	CPH	3-27	
CPH	3-2846	CPH	3-2864	CPH	3-27	
CPH	3-2847	CPH	3-2865	CPH	3-27	
CPH	3-2848	CPH	3-2866	CPH	3-27	
CPH	3-2849	CPH	3-2867	CPH	3-27	
CPH	3-2850	CPH	3-2868	CPH	3-27	
CPH	3-2851	CPH	3-2869	CPH	3-27	
CPH	3-2852	CPH	3-2870	CPH	3-27	
CPH	3-2853	CPH	3-2871	CPH	3-27	
CPH	3-2854	CPH	3-2872	CPH	3-27	
CPH	3-2855	CPH	3-2873	CPH	3-27	
CPH	3-2856	CPH	3-2874	CPH	3-27	
CPH	3-2857	CPH	3-2875	CPH	3-27	
CPH	3-2858	CPH	3-2876	CPH	3-27	
CPH	3-2859	CPH	3-2877	CPH	3-27	
CPH	3-2860	CPH	3-2878	CPH	3-27	
CPH	3-2861	CPH	3-2879	CPH	3-27	
CPH	3-2862	CPH	3-2880	CPH	3-27	
CPH	3-2863	CPH	3-2881	CPH	3-27	
CPH	3-2864	CPH	3-2882	CPH	3-27	
CPH	3-2865	CPH	3-2883	CPH	3-27	
CPH	3-2866	CPH	3-2884	CPH	3-27	
CPH	3-2867	CPH	3-2885	CPH	3-27	
CPH	3-2868	CPH	3-2886	CPH	3-27	
CPH	3-2869	CPH	3-2887	CPH	3-27	
CPH	3-2870	CPH	3-2888	CPH	3-27	
CPH	3-2871	CPH	3-2889	CPH	3-27	
CPH	3-2872	CPH	3-2890	CPH	3-27	
CPH	3-2873	CPH	3-2891	CPH	3-27	
CPH	3-2874	CPH	3-2892	CPH	3-27	
CPH	3-2875	CPH	3-2893	CPH	3-27	
CPH	3-2876	CPH	3-2894	CPH	3-27	
CPH	3-2877	CPH	3-2895	CPH	3-27	
CPH	3-2878	CPH	3-2896	CPH	3-27	
CPH	3-2879	CPH	3-2897	CPH	3-27	
CPH	3-2880	CPH	3-2898	CPH	3-27	
CPH	3-2881	CPH	3-2899	CPH	3-27	

4) **Related Procedures:** none

5) **Attachments:** none

6) **Related Forms:**

- CCL-F052 OR Lab Call To Form found under G:\Lab\_Shared\ICU\_ORLab\FORMS and LOGS\REQUISITIONS\CCL-F052 OR Call To Form.xls
- CCL-F075 Platelet Poor Plasma Log found under G:\Lab\_Shared\Houska\_Shoffner\QA Plt Poor Plasma Log.xlsx.

7) **References:**

- Clinical Laboratory Standards Institute. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays: Approved Guideline-Fifth Edition. CLSI Publication H21-A5. Wayne, PA, January, 2008.
- Pittiglio, DH. Clinical Hematology and Fundamentals of Hemostasis. USA, 1987.

8) **Review/Revision/Implementation:**

- Review Cycle: All procedures must be reviewed at least every 2 years.
- Office of Record: Department of Pathology, Critical Care Laboratory

9) **Previous Revision Date(s):** 2/17

10) **Revised/Reviewed Dates and Signatures:**

Reviewed/Revision Date: \_\_\_\_\_

Signature: \_\_\_\_\_

Reviewed/Revision Date: \_\_\_\_\_

Signature: \_\_\_\_\_

Reviewed/Revision Date: \_\_\_\_\_

Signature: \_\_\_\_\_