

	<b>CP 15 Spin Procedures</b>	<b>Dept:</b>	<b>324317 Central Processing</b>
		<b>Effective Date:</b>	<b>10/2000</b>
		<b>Revised Date:</b>	<b>9/16/19</b>
		<b>Contact:</b>	<b>Julie H Simmons</b>
<b>CLIA Medical Director Signature:</b>		<b>Approved Date:</b>	
See paper copy in the Central Processing procedure manual			

1. **General Procedure Statement:** To obtain a standard process for the Spin Rotation
  - A. **Purpose:** To provide guidelines for Central Processing staff to process patient samples in a consistent manner in order to produce a quality specimen for testing
  - B. **Responsible Department/Scope:**
    - i. Procedure owner/Implementer: Central Processing Lab
    - ii. Procedure prepared by: Julie H Simmons
    - iii. Who performs procedure: Central Processing staff
  
2. **Procedure:**
  1. The person in spin should rotate through the accessioning area to pick up pre-sorted samples in racks. The spin person has been provided a cart in order to carry the workload and distribute to the various lab sections.
    - a. Review tests to be performed
    - b. Note any special handling requirements. Some samples may not require centrifugation. ( Refer to the charts for details)
  
  2. Centrifuge Samples
    - a. Routine Samples:
      - Placed in a centrifuge carrier for the large centrifuges
      - Place cover on carrier to minimize aerosols in the event of leakage or breakage of any tubes
      - Spin all fluids ( including urines) that are to be tested in Chemistry
    - b. STAT samples
      - Spin using designated centrifuges
    - c. Spin samples at 5000 rpm for 5 minutes with the lid locked
    - d. Remove samples when centrifuge comes to a complete stop and the lid can be opened.

3. Hematology will check samples for clots
4. Aliquot Samples.
  - a. Initial aliquot verifying ID of primary tube with aliquot label as this is the last opportunity to assure the ID of the sample with the primary tube.
5. Distribute Samples:
  - a. All STAT Samples should be placed in the front of the rack for the Automation line.
  - b. Automation line Samples
    - a. All Chemistry samples are racked on the Chemistry automation line with the exception of the tests listed in the spin.
    - b. Baby samples, short samples and other samples should be placed after spinning in the designated Chemistry tray in Chemistry.
  - c. Extra samples are ordered in LIS and placed on the automation line for tracking and storage except for 2-b above which will be tracked and assigned a storage manually.
  - d. Tests that are not performed at WFBMC should be delivered to Send Out area. Chemistry personnel should be notified on 2<sup>nd</sup>, 3<sup>rd</sup>, weekends and holidays when sample is in send out. Chemistry personnel will direct CP personnel to perform any processing needed to the sample.
6. Urine samples:
  - Chemistry samples**
    - a. the person accessioning a 24 hour urine is responsible for measuring the volume and recording it on requisition with their initials  
*Refer to [Section g \(Determining Total Volume of Urines\)](#) to determine volume.*
    - b. Spin should pour off and label an aliquot off in a 10ml conical tube, labeled and placed in the extra rack in the spin refrigerator.
    - c. Acid, as a preservative for a 24hour urine collection should be added to only an orange acid resistant jug. Any collection that does not require preservative, use appropriate jug form the store room (item# 15-044407)
    - d. 24 hour urine jugs are not provided for WFUP clinics. Jugs are provided for inpatient locations and OPD clinics. Acid is provided to all locations.
    - e. An instruction sheet for collection 24 hour urine is provided with each jug that acid is added.

### **Urinalysis samples**

- a. Pour-up an aliquot in a conical centrifuge tube
- b. Label tube with an LIS Barcode and deliver to urinalysis
- c. Retain all ED urine samples and place in spin refrigerator

### **Microbiology samples**

- a. Tube staff should identify urine for culture only and deliver to microbiology immediately.
- b. If other testes are ordered in addition to micro deliver to STAT container. Urine for cultures must be expedited as urine cultures must be set up within two hours of collection.
- c. Accession should order the appropriate test and place barcode labels in the outer pocket of the bag and send to spin with a print screen or a paper order for the culture.
- d. Spin staff will aliquot the tubes appropriate for testing and forward the remaining urine sample in the original container with requisition to micro.

## **7. Determining Total Volume of Urines**

### **With a Matching container:**

- a. Choose a matching container from the containers available
- b. Place the empty container on the scale. PRESS and HOLD the tare button until the scale reads “00”. Return the empty container.
- c. Place the patient container on the scale and record the Total volume on the document that came with it or record patient information on a blank sheet of paper and record volume if nothing came with urine.  
Note: 1 gram – 1 ml
- d. Place document with information in the routine accession box.

### **With a Unique container:**

- a. PRESS and HOLD the tare button on the scale until it reads “0”.
- b. Place the patient container on the scale and record the total volume on the provided/or created document.
- c. Label a urine collection cup with patient information and fill the cup from the 24 hour container.
- d. Empty the 24 hour container down the sink and follow with water to completely cleanse the sink.
- e. Weigh the empty 24 hour container.
- f. Subtract the weight of the empty container from the total weight of the filled container. The result is the Total volume of urine.
- g. Write the Total Volume on the document AND on the lid of the urine cup.
- h. Place document with information in the routine accession box.

**3. Review/Revision/Implementation:**

All procedures must be reviewed at least every 2 years.

- All new procedures and procedures that have major revisions must be signed by the CLIA Laboratory Medical Director.
- All reviewed procedures and procedures with minor revisions can be signed by the designated section medical director.

**4. Related Procedures: NA**

**5. References: NA**

**6. Attachments: NA**

**7. Revised/Reviewed Dates and Signatures:**

<u>Review/Revision</u> Date: <u>3/2015</u>	Signature: <u>Tami Bradley</u>
<u>Review/Revision</u> Date: <u>8/11/16</u>	Signature: <u>Jennifer Hausmann</u>
<u>Review/Revision</u> Date: <u>12/6/18</u>	Signature: <u>Jennifer Hausmann</u>
<u>Review/Revision</u> Date: _____	Signature: _____
<u>Review/Revision</u> Date: _____	Signature: _____
<u>Review/Revision</u> Date: _____	Signature: _____
<u>Review/Revision</u> Date: _____	Signature: _____
<u>Review/Revision</u> Date: _____	Signature: _____
<u>Review/Revision</u> Date: _____	Signature: _____