



PANOCELL - 16 MASTER LIST



IMMUCOR, INC. Norcross, GA 30071 USA

US LICENSE NO: 886
LOT NO: 41815
EXPIRES: 24 Oct 02

NAME Mike Blue SEX M
NO. 2613 INSTITUTION _____
PREGNANCY _____ TRANSFUSION _____
BLOOD GROUP _____
ANTIBODY IDENTITY _____

TECH _____ DATE _____

Main data table with columns: Special Type, Donor, Rh-Hr (D, C, E, e, f, V, C*), Kell (K, k, K2, K3, J2, J3, J4), Duffy (Fy1, Fy2, J1, J2, J3, J4), Kidd (Ki1, Ki2, Jk1, Jk2, Jk3, Jk4), Lewis (Ls1, Ls2, Lo, Lp), MN (M, N, S, s), Luth. (L1, L2, L3, L4, L5, L6, L7, L8), X4, and Patient's Test Results (IS, S, P, A, M, G, C, C, S, E, L).

Some examples of anti U may react weakly with the S+ cells above.
All cells are positive for L, Gc, V1a, J1a, V1c, V1d, and V1e except where noted. * indicates phase antigens whose presence or absence may have been determined using only a single sample of a specific antibody.
In those instances where a patient's serum is known to contain anti U, it may be desirable to perform antibody screening first with D- red cells. The panel cells whose red numbers are set off by brackets (B) can be used together in a three- or four- cell D negative antibody screening reagent. All bracketed cells must be

Table with columns: DIRECT ANTIGLOBULIN TEST, ELUATE RESULT, POLY, %G, C3, LOT, REF: S11 T

3

Lot No. 8RA151 Exp. Date 14 JAN 2003

PATIENT NAME: Brown, Charlie

PATIENT ID: 10010

DATE: 1/28/03

CONCLUSION:

Ortho-Clinical Diagnostics

a Johnson & Johnson company

Panel Reagent Red Blood Cells
0.8% Resolve® Panel A
Antigram® Antigen Profile

A

Shaded columns indicate those antigens which are destroyed or depressed by enzyme treatment.

Cells	Rh-ir	Rh-ir				KELL				DUFFY			KIDD			LEWIS			MNS			P			Special Antigen Typing	Test Results						
		D	C	e	i	K	k	Kp ^a	Kp ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Xg ^a	Xg ^b	Le ^a	Le ^b	S	s	M	N	P ₁	P ₂	P ₃			Lu ^a	Lu ^b				
1	R1wR1	+	0	0	+	0	+	0	+	+	+	0	+	+	+	0	+	0	0	+	+	+	0	+	+	0	+		+		IS 31 ^a AK	AK
2	R1R1	+	0	0	+	+	+	0	+	+	+	0	+	+	+	0	+	0	0	+	+	+	0	+	+	0	+		+		AK 31	AK
3	R2R2	+	0	+	0	0	+	0	+	0	+	0	+	0	0	0	+	0	0	+	+	+	0	+	+	0	+		+		AK 31	AK
4	Ror	+	0	+	0	0	+	0	+	0	+	0	+	0	0	0	+	0	0	+	+	+	0	+	+	0	+		+		-	-
5	r1	0	+	0	+	0	+	0	+	0	+	0	+	0	0	0	+	0	0	+	+	+	0	+	+	0	+		+		AK 31	AK
6	r1	0	+	0	+	0	+	0	+	0	+	0	+	0	0	0	+	0	0	+	+	+	0	+	+	0	+		+		-	-
7	r1	0	0	+	+	0	+	0	+	0	+	0	+	0	0	0	+	0	0	+	+	+	0	+	+	0	+		+		-	-
8	r1	0	0	+	+	0	+	0	+	0	+	0	+	0	0	0	+	0	0	+	+	+	0	+	+	0	+		+		-	-
9	r1	0	0	+	+	0	+	0	+	0	+	0	+	0	0	0	+	0	0	+	+	+	0	+	+	0	+		+		-	-
10	r1	0	0	+	+	0	+	0	+	0	+	0	+	0	0	0	+	0	0	+	+	+	0	+	+	0	+		+		-	-
11	R1R1	+	0	0	+	0	+	0	+	0	+	0	+	0	0	0	+	0	0	+	+	+	0	+	+	0	+		+		AK 31	AK

* i antigen status may have been determined presumptively based on Rh-ir phenotype

Cells	Rh-ir	Rh-ir				KELL				DUFFY			KIDD			LEWIS			MNS			P			Special Antigen Typing	Test Results						
		D	C	e	i	K	k	Kp ^a	Kp ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Xg ^a	Xg ^b	Le ^a	Le ^b	S	s	M	N	P ₁	P ₂	P ₃			Lu ^a	Lu ^b				
Cell																																

