

16

Ortho-Clinical Diagnostics, Inc.  
a Johnson & Johnson company  
Raritan, New Jersey 08859

PATIENT NAME: SKINN, Hariv  
PATIENT ID: 123456  
DATE: 01/25/03 TECH: DR  
CONCLUSION:

Lot No. RA357 Exp. Date 18 FEB 2003  
Panel **A** Reagent Red Blood Cells  
Resolve Panel A  
Antigram Antigen Profile

Shaded columns indicate those antigens which are destroyed or depressed by enzyme treatment.

Cell#	Rh-ir	Rh-ir					KELL					DUFFY					KIDD					LEWIS					MNS					P					Special Antigen Typing	Test Results
		D	C	E	c	e	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Js <sup>a</sup>	Js <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Xg <sup>a</sup>	Xg <sup>b</sup>	Le <sup>a</sup>	Le <sup>b</sup>	S	s	M	N	P <sub>1</sub>	P <sub>2</sub>	P <sub>3</sub>	P <sub>4</sub>										
1	R1wR1	+	0	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		15 37 A B C C 6cd					
2	R1R1	+	0	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		488 - 2+					
3	R2R2	+	0	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		1+ - 1+					
4	Ror	+	0	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		2+ - 2+					
5	r1	0	+	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		1+ - 2+					
6	r1	0	+	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		1+ - 2+					
7	r1	0	+	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		1+ - 1+					
8	r1	0	+	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		2+ - 2+					
9	r1	0	+	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		1+ - 1+					
10	r1	0	+	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		2+ - 2+					
11	R1R1	+	0	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0		2+ - 2+					
Patient Cells																																						
Mode of Reactivity		37°C/Antiglobulin																																				

\* I antigen status may have been determined presumptively based on Rh-ir phenotype.

Cell#	Rh-ir	Rh-ir					KELL					DUFFY					KIDD					LEWIS					MNS					P					Special Antigen Typing	Test Results
		D	C	E	c	e	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Js <sup>a</sup>	Js <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Xg <sup>a</sup>	Xg <sup>b</sup>	Le <sup>a</sup>	Le <sup>b</sup>	S	s	M	N	P <sub>1</sub>	P <sub>2</sub>	P <sub>3</sub>	P <sub>4</sub>										
Additional Cells																																						
Donor Number																																						
Rh-ir																																						
Cell#																																						

17



Ortho-Clinical Diagnostics  
a Johnson & Johnson company

PATIENT NAME: Sammy Sneed  
PATIENT ID: 1094074

Lot No. BRA142 Exp. Date 07 MAY 2002

Panel Reagent Red Blood Cells  
0.8% Resolve® Panel A  
Antigram® Antigen Profile

**A**

DATE: \_\_\_\_\_ TECH: \_\_\_\_\_  
CONCLUSION: \_\_\_\_\_

Shaded columns indicate those antigens which are destroyed or depressed by enzyme treatment.

Cell #1 of BRA142 is designed with an A.

Cell #	Rh-ir	RH-ir				KELL				DUFFY				KIDDO				LEWIS				MNS				P				Special Antigen Typing	Test Results		
		D	C	E	f	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Xg <sup>a</sup>	Xg <sup>b</sup>	Le <sup>a</sup>	Le <sup>b</sup>	S	s	M	N	P <sub>1</sub>	P <sub>2</sub>	L <sup>a</sup>	L <sup>b</sup>	Cell #	Anticell	Result					
1	R1WR1	+	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	1	-	-	
2	R1R1	+	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	2	-	-	
3	R2R2	+	0	+	0	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	3	-	-	
4	For	+	0	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	4	-	-	
5	r1	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	5	-	-	
6	r1	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	6	-	-	
7	r	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	7	-	-	
8	r	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	8	-	-	
9	r	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	9	-	-	
10	r	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	10	-	-	
11	R1R1	+	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	11	-	-	

\* All antigen status may have been determined presumptively based on Rh-ir phenotype.

Cell #	Rh-ir	RH-ir				KELL				DUFFY				KIDDO				LEWIS				MNS				P				Special Antigen Typing	Test Results		
		D	C	E	f	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Xg <sup>a</sup>	Xg <sup>b</sup>	Le <sup>a</sup>	Le <sup>b</sup>	S	s	M	N	P <sub>1</sub>	P <sub>2</sub>	L <sup>a</sup>	L <sup>b</sup>	Cell #	Anticell	Result					
1	R1WR1	+	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	1	-	-	
2	R1R1	+	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	2	-	-	
3	R2R2	+	0	+	0	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	3	-	-	
4	For	+	0	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	4	-	-	
5	r1	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	5	-	-	
6	r1	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	6	-	-	
7	r	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	7	-	-	
8	r	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	8	-	-	
9	r	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	9	-	-	
10	r	0	0	+	+	0	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	10	-	-	
11	R1R1	+	0	0	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	+	11	-	-	



PATIENT NAME: Stein, Ben  
 PATIENT ID: 141 35 26  
 DATE: 1-29-03 TECH.  
 CONCLUSION: Anti. C-1k

Lot No. 8RA144 Exp. Date 02 JUL 2002



**Panel** Reagent Red Blood Cells  
 0.8% Resolve® Panel A  
 Antigram® Antigen Profile

**A**

Shaded columns indicate those antigens which are destroyed or depressed by enzyme treatment.

Cell#	Rh-ir	Donor Number	Rh-ir										MNS				LEWIS			BUFFY			KIDNEY			KELL			SPECIAL ANTIGEN TYPING			Test Results		
			D	C	E	c	e	f	C <sup>w</sup>	V	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Jsb	Jsa	Fyb	Fya	Jkb	Jka	Xgb	Xga	Lgb	Lga	S	s	M	N	P <sub>1</sub>	P <sub>2</sub>	Lub		Lub	Cell#
1	R1wR1	108737	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	3+	3+
2	R1R1	102827	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	3+	3+
3	R2R2	112808	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
4	Ror	114283	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
5	r'	107129	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	+	
6	r'	107444	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	+	
7	r'	112004	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	+	
8	r'	105893	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	+	
9	r'	82113	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	+	
10	r'	113948	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	+	
11	R1R1	111771	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
Mode of Reactivity			37°C/Antiglobulin										Antiglobulin			Variable			Cold			Var.												

\* If antigen status may have been determined presumptively based on Rh-ir phenotype.

Cell#	Rh-ir	Donor Number	Rh-ir										MNS				LEWIS			BUFFY			KIDNEY			KELL			SPECIAL ANTIGEN TYPING			Test Results		
			D	C	E	c	e	f	C <sup>w</sup>	V	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Jsb	Jsa	Fyb	Fya	Jkb	Jka	Xgb	Xga	Lgb	Lga	S	s	M	N	P <sub>1</sub>	P <sub>2</sub>	Lub		Lub	Cell#

