


Job Aide CP 1.3 Eyewash Maintenance

	DOCUMENT TYPE: <input checked="" type="checkbox"/> Form	ORIGIN DATE IN TITLE 21 2/11/2019
CLIA Lab Director: Dr. Gregory Pomper	LAB DEPARTMENT: Central Processing Lab	CONTACT: Central Processing Manager

Procedure:

Note: The eyewash station shall be flushed weekly to ensure the station is functioning properly.

1. Pull the eyewash bar counter-clockwise over the sink to start the flow of water
2. Verify the water pressure freely removed the caps
3. Allow the water to flow from the eyewash station for at least three (3) minutes
4. Verify good water flow is maintained
5. Turn off eyewash station by pushing clockwise over the counter
6. If the eyewash station fails to operate properly, notify the hospital engineering department by submitting an online work order or calling the Service Response Center at 336-716-1111 and document corrective action on the Central Processing Spin Rotation Maintenance Checklist.