

Title: Liquid Nitrogen Use, Precautions, and Exposure Policy (NCBH)		Document Number: 54926
Document Type: <input checked="" type="checkbox"/> Policy <input type="checkbox"/> Procedure <input type="checkbox"/> Guideline <input type="checkbox"/> Other		Last Review/Revision Date: 02/21/2023
Content Applies to Patient Care: (Select all that apply)	Content Applies to: (Select One)	Effective Date: 02/21/2023
<input checked="" type="checkbox"/> Adults <input type="checkbox"/> Pediatrics (Under 18)	<input checked="" type="checkbox"/> Clinical <input type="checkbox"/> Administrative	
Scope: <input type="checkbox"/> Enterprise <input type="checkbox"/> MW Region <input type="checkbox"/> SE Region <input type="checkbox"/> WI <input type="checkbox"/> IL <input type="checkbox"/> Greater Charlotte Market <input type="checkbox"/> Navicent Market <input type="checkbox"/> Wake Market <input type="checkbox"/> Floyd Market <input checked="" type="checkbox"/> Entity Only (Entity Name): NCBH <input checked="" type="checkbox"/> Department Only (Department Name): SCTCT		

I. PURPOSE

Liquid Nitrogen (LN₂) is used in the Stem Cell Transplant and Cellular Therapy (SCTCT) Lab to cryopreserve stem cells. LN₂ is a cryogenic liquid, or a liquefied gas that has a boiling point of less than -73°C. LN₂ therefore has hazardous properties, and specific instructions for its handling and dealing with exposures.

II. SCOPE

This document applies to NCBH SCTCT staff and management

III. DEFINITIONS/ABBREVIATIONS

LN₂: Liquid Nitrogen
SCTCT: Stem Cell Transplant and Cellular Therapy
PPE: Personal Protective Equipment

IV. POLICY

SECTION I: HANDLING LN₂

- Access to LN₂ sources shall be limited to authorized personnel.
- Employees of the SCTCT lab must be trained in the use of LN₂ equipment. Refer to [Microbulk Liquid Nitrogen Supply](#) and [Changing Liquid Nitrogen Tanks](#).
- Never store cryovials in the liquid phase of nitrogen as they might explode
- Store products in the vapor phase of liquid nitrogen when possible to avoid splashing and possible exposure.
- Always use a lab coat, cryogloves, and a face shield when handling products stored in liquid nitrogen, or when changing liquid nitrogen hoses on dewars.
- Maintain LN₂ supplied equipment in well-ventilated areas.
- An oxygen monitor should be employed in a space containing LN₂ supplied equipment.

Liquid Nitrogen Use, Precautions, and Exposure Policy (NCBH)

- Equipment containing cryogenic fluid should be stored in specific, validated containers which are not above eye level.
- Never touch uninsulated pipes or vessels without PPE.
- Never remove product identification labels from LN² dewars.
- LN₂ containers and equipment should be stored away from corrosives and combustible materials.
- LN₂ dewars are **NEVER** to be transported in an elevator with persons on the elevator.
- In case of breach of LN₂ vacuum-jacketed piping, immediately shut down the bulk system and evacuate area; Notify engineering, management, and medical director immediately.
- There are 2 emergency shutoff buttons for the bulk LN₂ system; One is next to the door in the Cryo Room and one is in the Blood Bank on the wall on the way from the Blood Bank to the Blood Bank Conference Room.
- If you must evacuate the SCTCT Lab prior to shutting down the bulk system during a large LN₂ breach, utilize the Blood Bank shutoff valve to disable the system.
- The LN₂ system can be reenabled by pulling out on the emergency stop button.

SECTION II: PPE AND CLOTHING CONSIDERATIONS TO BE USED WITH NITROGEN

- Do not wear open or canvas shoes while handling LN₂.
- Do not wear pants with cuffs while handling LN₂.
- Always wear cryogloves, a lab coat and a face shield when handling LN₂.
- Cryogloves should be loose fitting so they can be easily removed in case of exposure.
- Shoulder-length cryogloves are available for some tasks.
- Never touch uninsulated pipes or vessels without PPE.

SECTION III. EXPOSURE TO LN₂

- Treatment for LN₂ exposure should be given pending the care of a physician **ONLY** if there is no immediate risk to the caregiver.
- For skin contact, remove any clothing that may restrict circulation to the area.
- Do not rub frozen parts, as tissue damage may result.
- Do not walk on frostbitten feet.
- As soon as possible, place affected area in warm (NOT HOT) water.
- Never use dry heat to treat LN₂ exposure.
- If part of body is frozen, thaw with warm (NOT HOT) water, and cover with a dry sterile dressing pending medical care.
- In case of massive exposure, shower the victim with warm water while removing clothing.
- Victims suffering from lack of oxygen should be moved immediately to fresh air.
- If the victim is not breathing, administer artificial respiration.

Liquid Nitrogen Use, Precautions, and Exposure Policy (NCBH)

- If breathing is difficult, administer oxygen.



V. CROSS REFERENCES

- Microbulk Liquid Nitrogen Supply Procedure
- Changing Liquid Nitrogen Tanks Procedure
- Material Safety Data Sheet, LN2 Form
- Airgas Safe Handling and Use of Specialty Gases Appendix Form (NCBH)
- Airgas Cryogenic Liquid Containers Appendix Form (NCBH)
- Pathology Laboratory Safety Manual
- Precautions for Mechanical and Thermal Hazards

VI. RESOURCES AND REFERENCES

- Safe Handling and Use of Specialty Gases and Cryogenic Liquid Containers, Appendix, Airgas.com

VII. ATTACHMENTS

Not Applicable