

**BIO-RAD**

# Autoantibody Atlas



## DEVELOPED BY

**The University Of Washington  
Department of Laboratory Medicine  
Seattle, Washington**

Kathleen W. Hutchinson, MS, MT(ASCP)  
Mark H. Wener, MD  
Bruce G. Gilliland, MD  
Michael L. Astion, MD, PhD

## List of Antibodies

[ANA](#)

[ANCA](#)

[Anti Beta-2 Glycoprotein I](#)

[Anti Cardiolipin](#)

[Anti CCP](#)

[Anti Centromere](#)

[Anti Chromatin](#)

[Ani DFS70](#)

[Anti dsDNA](#)

[Anti ENA](#)

[Anti GBM](#)

[Anti Golgi](#)

[Anti Histone](#)

[Anti Jo-1](#)

[Anti Mitochondria](#)

[Anti MPO](#)

[Anti NuMA](#)

[Anti PR3](#)

[Anti Parietal](#)

[Anti PCNA](#)

[Anti Ribosomal-P](#)

[Anti RNP](#)

[Anti Scl-70](#)

[Anti Sm](#)

[Anti Smooth Muscle](#)


[Anti SS-A](#)

[Anti SS-B](#)

[RF](#)

[Anti RNA Pol III](#)

## ANA



Sample ANA staining patterns on HEp-2 Substrate.  
Staining patterns clockwise from top-left:  
Centromere, Nucleolar, Homogeneous, Speckled.

### Antibody

ANA  
Anti nuclear antibodies

### Antigen

Many different antigens that are found in the cell nucleus.

### Disease Associations

Dependent upon the antibody specificity and antibody level: systemic autoimmune diseases and autoimmune liver disease. May also be seen in patients with some acute or chronic infections, in association with some medications, and in healthy individuals.

### Testing Methods

IFA and ELISA (screening, not specific identification), Multiplexing.


### Typical IFA Pattern

Many. See ANA Tutor

### Notes

The sensitivity and specificity is very dependent upon the method used for ANA screening. The titer and frequency of anti nuclear antibodies in healthy individuals increases with age.

## ANCA



Sample ANCA staining patterns.; Clockwise from Top-left: C-ANCA, P-ANCA, ANA, Atypical ANCA. Ethanol Fixed ANCA substrate.

### Antibody

ANCA

Anti neutrophil cytoplasmic antibodies

### Antigen

A group of several different antigens that are found in neutrophils. Includes proteinase-3 (PR3) and myeloperoxidase (MPO). Some people also include other neutrophil specific autoantibodies (NSA) such as elastase (EL), and bacterial permeability increase (BPI) protein, Cathepsin G, lysozyme, lactoferrin, and catalase.

### Disease Associations

Dependent upon the antibody specificity: ANCA-associated vasculitis, inflammatory bowel disease (Crohn's disease, ulcerative colitis), systemic autoimmune diseases.

### Testing Methods

IFA (screening, not specific identification)

### Typical IFA Pattern

Dependent on antibody specificity: C-ANCA, P-ANCA, atypical ANCA (see ANCA Tutor)

### Notes

Clinical assays are not available for specific identification of many of these antibodies. ANA staining may resemble ANCA staining.

## Anti Beta-2 Glycoprotein I



### Antibody

IgG anti beta-2 glycoprotein I,  
IgA anti beta-2 glycoprotein I,  
IgM anti beta-2 glycoprotein I

### Alternate Names

Anti  $\beta$ 2GPI, a $\beta$ 2GPI

### Antigen

Glycoprotein that binds to cardiolipin and other anionic phospholipids (cardiolipin cofactor). Also called apolipoprotein H.

### Disease Associations

Anti phospholipid syndrome (APS), systemic lupus erythematosus (SLE).

### Testing Methods

ELISA


### Typical IFA Pattern

Not applicable

### Notes

Assays for different isotypes are available (i.e. IgG, IgA, IgM). IgG and IgM antibodies results are included in the 2006 APS Classification Criteria. Thought to be more specific for APS than anti-cardiolipin. APS is characterized by recurrent thrombosis. In pregnant women, antibodies to phospholipids have been associated with repeated, unexplained fetal loss.

## Anti Cardiolipin



Phospholipid

### Antibody

Anti Cardiolipin

IgG anti cardiolipin, IgA anti cardiolipin, IgM anticardiolipin

### Alternate Names

aCL

### Antigen

Diphosphatidylglycerol lipid. Originally derived from cow heart tissue and used as the antigen in the Wassermann test for syphilis.

### Disease Associations

Anti phospholipid syndrome (APS), systemic lupus erythematosus (SLE), syphilis. IgM assay may be elevated in patients with infections, rheumatoid factor, or cryoglobulins.

### Testing Methods

ELISA


### Typical IFA Pattern

Not applicable

### Notes

Assays are available for different immunoglobulin classes (i.e. IgG, IgA, IgM). IgG and IgM antibodies results are included in the 2006 APS Classification Criteria. One of many aPL (anti phospholipids) that include antibodies to cardiolipin, B2GP1, phosphatidylserine/prothombin, and lupus anticoagulant. May cause a biological false positive syphilis test.

## Anti CCP



### Antibody

Anti CCP

Anti cyclic citrullinated peptide

### Alternate Names

Anti-citrullinated peptide antibody (ACPA), Anti-Sa

### Antigen

Antibodies directed against citrullinated proteins and peptides (ACPA) are present in sera from a majority of patients with rheumatoid arthritis. The most common method of testing for ACPA is by detecting antibodies to synthetic cyclic (ring structured) citrullinated peptides (CCP). There are a variety of synthetic peptides that are used for clinical testing. Successive generations of tests are assigned different numbers (e.g. CCP2, CCP3, etc.). Some citrullinated peptides used in testing are derived from the sequences in filaggrin or vimentin.

### Disease Associations

Rheumatoid arthritis (RA).

### Testing Methods

ELISA, Western blot (anti-Sa)


### Typical IFA Pattern

Oral mucosal cells: Anti perinuclear factor (APF), Esophageal epithelial cells: anti keratin (AKA) (Not routinely used for clinical testing)

### Notes

Anti CCP has better specificity for RA than rheumatoid factor does. Anti-CCP levels may correlate with disease prognosis. Citrulline is formed by a post translational modification of arginine.

## Anti Centromere



### Antibody

Anti Centromere

### Alternate Names

Anti centromere associated protein (CENP)

### Antigen

Proteins present in the centromere (kinetochore) region of the chromosome that includes CENP-A, CENP-B, CENP-C.

### Disease Associations

Most common association is with the limited cutaneous form of scleroderma (CREST syndrome). May be present in other systemic autoimmune diseases and primary biliary cirrhosis (PBC). Patients with primary Raynaud's phenomenon and anti-centromere have an increased risk for later development of a systemic autoimmune disease.

### Testing Methods

ELISA, Multiplexing, IFA

### Typical IFA Pattern

HEp-2: Discrete speckles of the interphase cell nuclei and dividing cell chromatin. May be difficult to identify in the presence of other nuclear or cytoplasmic antibodies.


### Notes

Most patient sera with centromere antibodies react with CENP-A, CENP-B and CENP-C.

HEp-2 Centromere Pattern, 100x Objective



## Anti Chromatin



### Antibody

Anti Chromatin

### Alternate Names

Anti Nucleosome

### Antigen

Epitopes present on histones that are bound to DNA in the nucleosome structure, i.e. the DNA-histone complex.

### Disease Associations

Present in the sera of patients with drug induced lupus, SLE, and several other systemic autoimmune diseases. Some reports indicate that the antibody may play a pathogenic role and be helpful in the diagnosis and establishing prognosis of SLE patients.

### Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)

### Typical IFA Pattern


HEp-2: Homogeneous ANA, pattern is not unique therefore does not permit specific identification.

### Notes

Anti chromatin refers to a family of autoantibodies that react to histones bound to DNA.

Nucleosome structure

## Anti DFS70



### Antibody

Anti DFS70, Anti Dense Fine Speckled 70

### Alternate Names

Anti LEDGF, Anti LEDGFp75 (LEDGF = lens epithelium-derived growth factor),

### Antigen

75 kd protein

### Disease Associations

Variety of non-autoimmune conditions including cancer, asthma, eye diseases, atopic dermatitis, chronic fatigue syndrome, and in apparently healthy individuals. In adults, very low frequency in systemic autoimmune rheumatic diseases (SARD).

### Testing Methods

ELISA, IFA

### Typical IFA Pattern


Dense fine speckles of interphase cell nuclei. Intense staining of the chromatin of mitotic cells.

### Notes

ASome studies indicate that for adults, if the only ANA antibody present is anti-DFS70, the likelihood that the patient has SARD is low.

Hep-2 DFS70 pattern, 40x Objective

## Anti dsDNA



### Antibody

Anti dsDNA  
Anti double stranded deoxyribonucleic acid

### Alternate Names

Anti native DNA

### Antigen

Double stranded DNA (deoxyribonucleic acid)

### Disease Associations

SLE

### Testing Methods

ELISA, Multiplexing, IFA on Crithidia, ANA IFA (screening, not specific identification)

### Typical IFA Pattern


HEp-2: Homogeneous ANA, pattern is not unique therefore does not permit specific identification. Crithidia: Staining of the kinetoplast.

### Notes

Antibody levels often correlate with disease activity, especially high avidity antibodies. Many assays detect both low and high avidity antibodies.

HEp-2 (top), Crithidia (bottom)

## Anti ENA



### Antibody

Anti ENA

Anti Extractable Nuclear Antigens

### Antigen

Group of several different antigens that are easily extracted from the cell nucleus including Sm/RNP, SS-A, SS-B, Scl-70, and Jo-1.

### Disease Associations

Systemic autoimmune diseases

### Testing Methods

ELISA, Multiplexing, Double Diffusion, CIE (counterimmunoelectrophoresis)

### Typical IFA Pattern


HEp-2: most, but not all produce speckled nuclear staining.

### Notes

It is common to use a screening method and then proceed to follow-up testing for specific antibodies if the screen is positive. The screening method used will determine which antibodies can be detected.

ELISA plate (top), Double diffusion gel (bottom)

## Anti GBM



### Antibody

Anti GBM

Anti glomerular basement membrane

### Alternate Names

Goodpasture-antibodies

### Antigen

Type IV collagen. Most antibodies react with a cryptic epitope located on the C-terminal alpha 3 NC1 domain.

### Disease Associations

Goodpasture's syndrome.

### Testing Methods

ELISA, Multitplexing, IFA on primate kidney


### Typical IFA Pattern

Primate Kidney: Linear staining of the glomerular basement membrane

### Notes

Anti-GBM levels may correlate with disease activity. Patients with acute renal failure and/or severe lung hemorrhage are often screened for ANCA and GBM antibodies.

## Anti Golgi



### Antibody

Anti Golgi

### Antigen

Several different autoantigens have been found in the Golgi apparatus including several different molecular weight golgins, giantin/macrogolgin.

### Disease Associations

Variety of autoimmune rheumatic diseases including SLE, Sjogren's syndrome, Rheumatoid Arthritis. Also seen in patients with liver disease and cancer.

### Testing Methods

IFA

### Typical IFA Pattern


Coarse cytoplasmic speckles that tend to cluster on one side of the nucleus.

### Notes

Rare pattern. May be a transient finding for some patients.

HEp-2 Golgi pattern, 40x Objective

## Anti Histone



### Antibody

Anti Histone

### Antigen

Histone proteins (H1, H2A, H2B, H3, H4, H5)

### Disease Associations

Present in the sera of patients with drug induced lupus, spontaneous SLE, and several other systemic autoimmune diseases.

### Testing Methods

ELISA, IFA (screening, not specific identification)

### Typical IFA Pattern

HEp-2: Homogeneous ANA, pattern is not unique enough for specific identification.


Crithidia: Negative

### Notes

Autoantibodies that only react to histones bound to DNA are usually called anti-nucleosome or anti-chromatin (see anti-chromatin section).

HEp-2 Homogeneous Pattern, 40x Objective

## Anti Jo-1



### Antibody

Anti Jo-1

### Alternate Names

Anti histidyl-tRNA synthetase

### Antigen

Anti histidyl-tRNA synthetase

### Disease Associations

Polymyositis and dermatomyositis

### Testing Methods

ELISA, Multiplexing

### Typical IFA Pattern

HEp-2: negative ANA, may have weak cytoplasmic staining that is not unique therefore does not permit specific identification.


### Notes

Usually considered to be one of the ENA (extractable nuclear antigens) antibodies. Although the antibody is only present in the sera of approximately 20% of patients with myositis, it is the most common of the myositis specific antibodies (MSA).

HEp-2 Negative nuclear staining with weak cytoplasmic staining



## Anti Mitochondria



HEp-2 (top), Mouse Stomach Kidney (bottom)

### Antibody

Anti - mitochondria

### Alternate Names

AMA

### Antigen

Nine different subtypes have been identified (M1-M9). Antibodies to M2 are considered to be the most clinically significant. Antibodies to M2 react with the E2 subunit of the pyruvate dehydrogenase complex (PDC-E2).

### Disease Associations

Over 90% of patients with primary biliary cirrhosis(PBC) have mitochondrial antibodies.

### Testing Methods

ELISA (anti-M2) IFA using mouse stomach kidney substrate (detects multiple subtypes)  
Immunoblot (not generally available in clinical labs)

### Typical IFA Pattern


MSK: cytoplasmic staining of the renal tubules and stomach parietal cells.

HEp-2: numerous granular speckles in the cytoplasm.

### Notes

There are other antibodies that may resemble the mitochondrial staining seen on HEp-2. Different techniques should be used to identify the antibody.

## Anti MPO



Ethanol (top) and Formalin fixed neutrophils

### Antibody

Anti MPO  
Anti myeloperoxidase

### Antigen

Myeloperoxidase. A cationic enzyme found in the azurophilic granules of neutrophils.

### Disease Associations

ANCA-associated vasculitis. Strongest association is with idiopathic crescentic glomerulonephritis (NCGN) and microscopic polyangiitis (MPA), less common in and Churg-Strauss Syndrome (CSS) and Wegener's granulomatosis (WG).

### Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)


### Typical IFA Pattern

Ethanol Fixed Neutrophils: P-ANCA (perinuclear) Formalin Fixed Neutrophils: cytoplasmic granules

### Notes

One of the anti-neutrophil cytoplasmic antibodies (ANCA). Anti-MPO is the most common antigen associated with the P-ANCA pattern. ANAs may resemble P-ANCA staining (see ANCA Tutor). Anti-MPO levels may correlate with disease activity. Patients with acute renal failure and/or severe lung hemorrhage are often screened for ANCA and GBM antibodies.

## Anti NuMA



### Antibody

Anti NuMA (nuclear mitotic apparatus))

### Alternate Names

Anti NuMA1

### Antigen

236 kDa nuclear matrix protein that is distributed to the spindle poles during mitosis

### Disease Associations

SLE, Sjogren's syndrome, occasionally associated with active infections and cancer..

### Testing Methods

IFA

### Typical IFA Pattern


Speckled nuclear staining of the interphase (non-dividing) cells and staining of the mitotic spindle fibers during mitosis.

### Notes

Rare pattern. Anti NuMA-2 is a different, but related, mitotic spindle pattern.

HEp-2 NuMA pattern, 40x Objective

## Anti PR3



Ethanol (top) and Formalin fixed neutrophils

### Antibody

Anti PR3  
Anti Proteinase 3

### Antigen

Proteinase 3. A cationic enzyme found in the azurophilic granules of neutrophils.

### Disease Associations

ANCA-associated vasculitis. Strongest association is with Wegener's granulomatosis (WG), less often seen in microscopic polyangiitis (MPA) and Churg-Strauss Syndrome (CSS).

### Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)


### Typical IFA Pattern

Ethanol Fixed Neutrophils: C-ANCA (cytoplasmic) Formalin Fixed Neutrophils: cytoplasmic granules

### Notes

One of the anti-neutrophil cytoplasmic antibodies (ANCA). Antibody levels often correlate with disease activity. Patients with acute renal failure and/or severe lung hemorrhage are often screened for ANCA and glomerular basement membrane (GBM) antibodies.

## Anti Parietal Cell



### Antibody

Anti parietal cell

### Alternate Names

PCA (parietal cell antibodies), GPC antibodies (gastric parietal cell antibodies)

### Antigen

$\alpha$  and  $\beta$  subunits of the H<sup>+</sup>/K<sup>+</sup> ATPase

### Disease Associations

Pernicious anemia

### Testing Methods

IFA using mouse stomach kidney substrate

ELISA (anti-GPC)

### Typical IFA Pattern


MSK: cytoplasmic staining of the stomach parietal cells. Renal tubules are negative.

### Notes

This antibody cannot be detected by IFA if antibodies to mitochondria are present. The antibody frequency in the normal population increases with age, and the antibody is more common in women than in men.

Mouse stomach kidney

## Anti PCNA



Hep-2 PCNA Pattern, 40x and 100x Objective

### Antibody

Anti PCNA (Anti Proliferating Cell Nuclear Antigen)

### Alternate Names

Anti PCNA/Cyclin

### Antigen

34 kDA factor of an auxiliary protein for DNA polymerase  $\delta$ . Cell cycle dependent, expressed during the S-phase of mitosis.

### Disease Associations

Historically believed to be highly specific for SLE, not confirmed by later studies. Has been identified in patients with Sjogren's syndrome autoimmune thyroiditis, hepatitis B, hepatitis C, and some cancers.

### Testing Methods

IFA. Specific assays (e.g. double immunodiffusion, immunoblot, ELISA) are usually lab developed or research use only reagents.


### Typical IFA Pattern

Variable staining. Approximately 2/3 of the interphase (non-dividing) cells are negative. Positive cells have nuclear speckled staining that varies from coarse, discrete speckles to fine speckles.

### Notes

Rare pattern. There are several other rare antibodies that have cell cycle dependent staining that resembles anti-PCNA. Many have been associated with cancer.

## Anti Ribosomal-P



### Antibody

Anti Ribosomal P

### Antigen

Ribosomal proteins P0, P1, and P2

### Disease Associations

SLE

### Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)

### Typical IFA Pattern


HEp-2: cytoplasmic staining that is not unique therefore does not permit specific identification, may have coexisting nucleolar staining. MSK: cytoplasmic staining of the chief cells

### Notes

Antigen is present in thymus extract used for ENA (extractable nuclear antigens) antibody testing by double diffusion methods, but is not typically included in ELISA ENA screening kits.

HEp-2 cytoplasmic (top) MSK chief cells (bottom)

## Anti RNP



### Antibody

Anti RNP  
Anti ribonucleoprotein

### Alternate Names

Anti U1-RNP, anti U1 snRNP

### Antigen

Small ribonucleoproteins in spliceosomes, several proteins including 68-70kd, A, and C

### Disease Associations

Several systemic autoimmune diseases including SLE, Sjögren's Syndrome, scleroderma, myositis, occasionally rheumatoid arthritis. Patients with mixed connective tissue disease tend to have very high levels of the antibody, without other autoantibodies.

### Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification), immunodiffusion.

### Typical IFA Pattern

HEp-2: Speckled ANA, pattern is not unique, therefore does not permit specific identification


### Notes

Difficult to prepare the RNP antigens without Sm epitopes also present. Some assays use antigens that have both the Sm and RNP epitopes, these assays are called anti Sm/RNP. One of the ENA (extractable nuclear antigens) antibodies.

HEp-2 Speckled Pattern, 40x Objective



## Anti Scl-70



### Antibody

Anti Scl-70  
Anti Scleroderma 70

### Alternate Names

Anti Topo I, Anti topoisomerase I

### Antigen

DNA topoisomerase I

### Disease Associations

Most common association is with the diffuse cutaneous form of scleroderma. Occasionally present in the serum of patients with limited cutaneous scleroderma.

### Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)

### Typical IFA Pattern


HEp-2: Homogeneous ANA, often with coexisting nucleolar staining. Pattern is not unique therefore does not permit specific identification.

### Notes

One of the ENA (extractable nuclear antigens) antibodies.

HEp-2, Homogenous with nucleolar staining

## Anti Sm



### Antibody

Anti Sm

### Alternate Names

Anti Smith

### Antigen

Small ribonucleoproteins in spliceosomes, several proteins including B'/B and D

### Disease Associations

SLE

### Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification), immunodiffusion

### Typical IFA Pattern


HEp-2: Speckled ANA, pattern is not unique therefore does not permit specific identification.

### Notes

Almost all samples with anti Sm have coexisting antibodies to RNP. One of the ENA (extractable nuclear antigens) antibodies. Some assays use antigens that have both the Sm and RNP epitopes, these assays are called anti Sm/RNP.

HEp-2 Speckled Pattern, 40x Objective

## Anti Smooth Muscle



### Antibody

Anti-smooth muscle

### Alternate Names

SMA

### Antigen

Several potential antigens, including F-actin

### Disease Associations

Autoimmune hepatitis (type 1)

### Testing Methods

IFA using mouse stomach kidney substrate, ELISA (anti-F-actin)

### Typical IFA Pattern


MSK: cytoplasmic staining of the stomach smooth muscle cells and staining of the arterial blood vessel walls.

### Notes

Low level staining of the smooth muscle layer is fairly common. Reference range studies should be performed to insure the proper cut-off is used.

Mouse stomach kidney

## Anti SS-A



### Antibody

Anti SS-A

Anti Sjögren's Syndrome A

### Alternate Names

Anti SSA, Anti Ro, Anti SS-A/Ro, Anti- Ro/SS-A

### Antigen

Small ribonucleoproteins, two different epitopes: 52kd and 60 kd

### Disease Associations

Sjögren's Syndrome, SLE (particularly subacute cutaneous lupus variant), Neonatal Lupus, Congenital Heart Block

### Testing Methods

ELISA, Multiplexing

### Typical IFA Pattern


HEp-2: Speckled ANA, pattern is not unique therefore does not permit specific identification.

### Notes

ANA by IFA may be weak or absent for some anti-SSA antibodies. One of the ENA (extractable nuclear antigens) antibodies.

HEp-2 Speckled Pattern, 40x Objective

## Anti SS-B



### Antibody

Anti SS-A  
Anti Sjögren's Syndrome A

### Alternate Names

Anti SSB, Anti La

### Antigen

Small ribonucleoproteins

### Disease Associations

Sjögren's Syndrome, SLE, Neonatal Lupus, Congenital Heart Block

### Testing Methods

ELISA, Multiplexing, IFA (screening, not specific identification)

### Typical IFA Pattern


HEp-2: Speckled ANA, may also have some nucleolar staining, pattern is not unique therefore does not permit specific identification.

### Notes

Most samples with anti SS-B have coexisting antibodies to SS-A. One of the ENA (extractable nuclear antigens) antibodies.

HEp-2 Speckled Pattern, 40x Objective

## RF

**Antibody**

RF  
Rheumatoid factor

**Antigen**

Epitopes are in the constant region of the IgG heavy chain.

**Disease Associations**

Rheumatoid arthritis (RA) and Sjögren's syndrome. Also seen in the sera of patients with other systemic autoimmune disease, mixed cryoglobulinemia, some chronic bacterial and viral infections, and in the sera of healthy individuals (particularly the elderly) in low titer.

**Testing Methods**

Agglutination, Nephelometry/Turbidometry, ELISA

**Typical IFA Pattern**


Not applicable

**Notes**

The name rheumatoid factor was given before the discovery that it was an antibody to IgG. There are several different isotypes of rheumatoid factor, IgM, IgA, and IgG (including the subclasses IgG1, IgG2, IgG3, and IgG4). Most clinical assays in the US detect only IgM RF.

IgM RF binding to normal IgG

## Anti RNA Pol III



HEp-2 RNA Pol III Pattern, 40x Objective

### Antibody

Anti RNA Pol III (RNA polymerase III)

### Alternate Names

Anti RNAP3

### Antigen

RNA Polymerase III, some ELISAs may use a recombinant immunodominant fragment of the antigen.

### Disease Associations

Systemic sclerosis (scleroderma)

### Testing Methods

IFA, ELISA

### Typical IFA Pattern

Typically nucleolar staining of interphase (non-dividing) cells, but some positive samples may have speckled nuclear staining. The IFA pattern may not be a sensitive marker for the presence of the antibody.

### Notes

Included as one of the ACR/EULAR criteria for the classification of systemic sclerosis.