

# GAS Storage and Stability of Samples, Controls and Reagents

## PURPOSE

- This procedure provides instructions for storage and stability of sample buffer tubes and reagents.

## SAFETY CONSIDERATIONS

- Standard precautions. Refer to [MB 2.02](#) Biohazard Containment
- Use of engineering controls: Refer to [MB 3.01](#) Engineering Controls to Prevent Nucleic Acid Contamination

## ABBREVIATIONS

- BSC: BioSafety Cabinet
- BSL: BioSafety level
- GAS: Group A Strep
- GASD: Group A Strep Detection
- MM: master mix
- NEGC: negative control
- NFW: nuclease free water
- PCTL: process control
- POSC: positive control
- PP: primer – pair
- RT: room temperature
- SEAC: Simplexa extraction and amplification control
- TE buffer: Tris – EDTA buffer
- Area/Room 1: Clean room
- Area/Room 2: Processing room
- Area/Room 3: Amplification room

## MATERIALS REQUIRED

Equipment	Reagents	Supplies
<b>Room 1</b> <ul style="list-style-type: none"> <li>Refrigerator 2 – 8° C</li> <li>-10 to -30° C freezer</li> <li>Mini-centrifuge</li> <li>Laminar flow hood</li> <li>Eppendorf Repeater pipette</li> </ul>	GAS primer pair ( 50 µl)	Orange barrier wipes
	TA master mix (2 vials, 200 µl ea)	Nitrile gloves (powder-free)
	GAS positive control (POSC) 100 µl	Cryovial storage box
	NFW (NEGC)	Test tube rack
	TE buffer 1X pH 8.0 (100 ml)	Scissors
<b>Room 2</b> <ul style="list-style-type: none"> <li>Refrigerator 2 – 8° C</li> <li>BSC BSL-2</li> <li>-70° C freezer</li> </ul>	SEAC <ul style="list-style-type: none"> <li>Amplification Control DNA</li> <li>Amplification Control primer pair</li> </ul>	Eppendorf pipette tip, 5 ml

**PROCEDURE A:** Follow the activity below for the proper storage of sample buffers

### Storage and Stability of Processed Specimens

Activity	Step	Action	Related Doc																
<b>Specimen Processing</b> <b>Room 2</b>	1	Prepare swabs for testing	<a href="#">MB 1.01</a> Specimen Management																
		<table border="1"> <thead> <tr> <th>Step</th> <th>Action</th> </tr> </thead> <tbody> <tr> <td>a</td> <td>Number patients on GASD worksheet in consecutive order</td> </tr> <tr> <td>b</td> <td>Number primary container and associated label with assigned test number on worksheet</td> </tr> <tr> <td>c</td> <td>Number cap of a 250 µl TE tube according to assigned number on worksheet</td> </tr> <tr> <td>d</td> <td>Properly label TE tube with patient aliquot label matching the number on the cap to the number on the label</td> </tr> <tr> <td>e</td> <td>Verify number on primary container before transfer</td> </tr> <tr> <td>f</td> <td>Using a barrier wipe, break s swab off in the sample TE buffer tube with corresponding number on cap</td> </tr> <tr> <td>g</td> <td>Vortex 2 min, vortex setting 9</td> </tr> </tbody> </table>		Step	Action	a	Number patients on GASD worksheet in consecutive order	b	Number primary container and associated label with assigned test number on worksheet	c	Number cap of a 250 µl TE tube according to assigned number on worksheet	d	Properly label TE tube with patient aliquot label matching the number on the cap to the number on the label	e	Verify number on primary container before transfer	f	Using a barrier wipe, break s swab off in the sample TE buffer tube with corresponding number on cap	g	Vortex 2 min, vortex setting 9
		Step		Action															
		a		Number patients on GASD worksheet in consecutive order															
		b		Number primary container and associated label with assigned test number on worksheet															
		c		Number cap of a 250 µl TE tube according to assigned number on worksheet															
		d		Properly label TE tube with patient aliquot label matching the number on the cap to the number on the label															
		e		Verify number on primary container before transfer															
f	Using a barrier wipe, break s swab off in the sample TE buffer tube with corresponding number on cap																		
g	Vortex 2 min, vortex setting 9																		

Activity	Step	Action	Related Doc	
Sample storage	2	Store sample buffer tubes as follows:		
		Temperature		Stability
		Room temp		4 hr
		2 – 8° C		48 h
		-70° C		90 days

**PROCEDURE B:** Follow the activity below for proper storage of reagents. Refer to Tables 1 – 3.

**Information for Reagent Storage**

Activity	Step	Action	Related Doc
General Information	1	<b><i>Always change gloves prior to handling new reagents and cartridges</i></b>	<a href="#">MB 5.02</a> Standards of Practice
		GAS PCR reagents are shipped frozen on dry ice <ul style="list-style-type: none"> <li>▪ Do not use reagents if thawed upon arrival</li> <li>▪ Do not use reagents if vials have been damaged</li> <li>▪ Contact <b>Focus Customer Service at 1.800.838.4548</b> for shipping issues</li> </ul>	
	2	Store reagents at -10 to -30° C until expiration date located on the vial unless otherwise noted. Refer to Table 1.	
	3	Discard reagents that have not been stored properly or have expired according to lab safety policy	
	4	Remove only the required amount of reagents from storage needed for testing. <ul style="list-style-type: none"> <li>▪ Clean gloves required</li> </ul>	
	5	Protect from excess heat and light; store in dark	
	6	Thaw reagents at room temperature before use	
	7	Once thawed, store reagents at 2 – 8° C up to 30 days <ul style="list-style-type: none"> <li>▪ <b>Do not refreeze</b></li> </ul>	
	8	Do not allow contact with reactive vapors from bleach or Extran or dust as these may affect the performance.	
9	Do not interchange the reagent tube caps		

**Table 1: Simplexa GAS Reagents**

Reagent	Unopened Reagent		Stability	Opened Reagent		Stability
	Temp (° C)	Location		Temp (° C)	Location	
GAS POSC (red)	-70	Room 2	expiry date	2 – 8	Room 2	30 days
TA MM (green)	-10 to -30	Room 1	expiry date	2 – 8	Room 1	30 days
GAS pp, conc. 15 µM (brown)	-10 to -30	Room 1	expiry date	2 – 8	Room 1	30 days
SEAC (blue)	-10 to -30	Room 1	expiry date	2 – 8	Room 1	30 days

**Table 2: Molecular Grade Water (RNase and DNase free)**

Reagent	Unopened/Opened		Aliquot Storage		In Use Aliquots	
	Temp	Location	Temp (° C)	Location	Temp (° C)	Location
Nuclease free water (NFW)	RT	Room 1	2 – 8	Room 1	2 – 30	Room 2

**Table 3: TE Buffer and Aliquot Storage**

Reagent	Unopened/Opened temp		Aliquot Storage		In Use Aliquots, temp (° C)	
	Temp	Location	Temp (° C)	Location	Temp (° C)	Location
TE buffer 1X	RT	Room 1	2 – 8	Room 1	2 – 30	Room 2

**Table 4: Process Control Storage**

Reagent	Temp (° C)	Location	Stability	Temp (° C)	Location	Stability
GAS Process Control (matrix)	≤ 70	Room 2	1 year	2 – 8	Room 2	7 days
GAS Process Control (swab)	NA	NA	NA	2 – 8	Room 2	7 days

## REFERENCES

1. GAS PCR Clinical Verification and Validation Study performed at Children’s Hospitals and Clinics of MN August 2014
2. Simplexa™ Group A Strep Molecular Control Circular PI.MOL8033.IVD, Rev. C, 16-April-2013, Focus Diagnostics, Cypress, CA 90630
3. Group A Primer Pair (50 µl) ASR, Circular PI.MOL9033 Rev. C, 21 May 2012, Focus Diagnostics, Cypress, CA 90630
4. Simplexa™ Extracton & Amplification Control Set, Circular PI.MOL9000, Rev. D, CE, 7 Mar 2013, Focus Diagnostics, Cypress, CA 90630

## Historical Record

Version	Written/Revised by:	Effective Date:	Summary of Revisions
1	P. Ackerman	9.9.14	Initial Version
2	P. Ackerman	07.26.16	Reformatttted for CMS upload