

Bordetella Equipment and Room Decontamination

PURPOSE:

- This procedure provides instructions for decontamination of equipment and molecular rooms used for *Bordetella* testing

DOCUMENTATION/RECORDS

- Daily Maintenance Log

SAFETY CONSIDERATIONS

- Standard precautions. Refer to [MB 2.02](#) Biohazard Containment
- Use of engineering controls: Refer to [MB 3.01](#) Engineering Controls to Prevent Nucleic Acid Contamination

MATERIALS REQUIRED

Equipment	Reagents	Supplies
BioSafety Cabinet	Sani-Cloth Bleach Wipes (10%)	Nitrile gloves (powder-free)
Pipettes	70% alcohol	Tacky mats
Test tube racks	5% Extran	Lint free absorbent cloths
Disc cooling block		Cotton tip swabs
Liaison MDX		
Stainless steel scissors		

ABBREVIATIONS

- PPE: personal protective equipment
- BSC: biosafety cabinet
- UV: ultraviolet light
- Area/Room 1: Clean room
- Area/Room 2: Processing room
- Area/Room 3: Amplification room

PROCEDURE A: Follow the activities in the table below for equipment and room decontamination **Equipment and Room Decontamination**

Activity	Step	Action (refer to Table 1 for recommended schedule)
PPE	1	Gloves and lab coat required
Tube racks Room 2	2	Soak tube racks in 5% Extran for minimum of 5 min after each run set-up and when visibly contaminated <ul style="list-style-type: none">▪ Rinse with water▪ Air dry

Activity	Step	Action (refer to Table 1 for recommended schedule)																								
Cooling blocks	3	Wipe down Nalgene cooling block (rm 1) with Sani-Cloth Bleach Wipe at the end of the day and when visibly contaminated <ul style="list-style-type: none"> Allow bleach to sit for 4 – 5 min Remove bleach with water followed by 70% alcohol 																								
	4	Wipe down aluminum disc cooling block with Sani-Cloth Bleach Wipe after each use. <ul style="list-style-type: none"> Allow bleach to sit for 4 – 5 min Rinse with water Air dry <p>Caution: Do not leave bleach on aluminum equipment more than 5 min – the metal will corrode.</p>																								
Scissors	5	Wipe down scissors with Sani-Cloth Bleach Wipe after each use <ul style="list-style-type: none"> Allow bleach to sit for 4 – 5 min Rinse with water Air dry Store individually in clean Zip lock bag until use 																								
Hoods, Pipettes and Bench tops Room 1, 2, 3	5	Wipe down hoods, pipettes and benchtops after each procedure and when visibly contaminated																								
		<table border="1"> <thead> <tr> <th>Location</th> <th>Step</th> <th>Action</th> </tr> </thead> <tbody> <tr> <td rowspan="3"> Room 1 <ul style="list-style-type: none"> 5% Extran 70% alcohol </td> <td>1</td> <td>Allow 5% Extran to sit for approximately 5 min</td> </tr> <tr> <td>2</td> <td>Remove Extran from surfaces with 70% alcohol</td> </tr> <tr> <td>3</td> <td>Lay pipettes flat in hood</td> </tr> <tr> <td rowspan="3"> Room 2 and 3 <ul style="list-style-type: none"> Sani-Cloth Bleach Water 70% alcohol </td> <td>1</td> <td>Allow bleach to sit 4 – 5 min</td> </tr> <tr> <td>2</td> <td>Remove bleach residue with water followed by 70% alcohol</td> </tr> <tr> <td>3</td> <td>Lay pipettes flat in hood</td> </tr> <tr> <td rowspan="4"> Pipettes (additional information) </td> <td>1</td> <td>Carefully clean the pipettor handle and barrel</td> </tr> <tr> <td>2</td> <td>Electronic pipettes: Avoid electronic buttons and AC charger socket when cleaning</td> </tr> <tr> <td>3</td> <td>Manual pipettes: Do not allow cleaning solutions to enter the plunger and gear mechanism</td> </tr> <tr> <td>4</td> <td>The above solutions are corrosive and over time can contribute to the electronic components of the pipettes and increased resistance when adjusting volumes on manual pipettes</td> </tr> </tbody> </table>	Location	Step	Action	Room 1 <ul style="list-style-type: none"> 5% Extran 70% alcohol 	1	Allow 5% Extran to sit for approximately 5 min	2	Remove Extran from surfaces with 70% alcohol	3	Lay pipettes flat in hood	Room 2 and 3 <ul style="list-style-type: none"> Sani-Cloth Bleach Water 70% alcohol 	1	Allow bleach to sit 4 – 5 min	2	Remove bleach residue with water followed by 70% alcohol	3	Lay pipettes flat in hood	Pipettes (additional information)	1	Carefully clean the pipettor handle and barrel	2	Electronic pipettes: Avoid electronic buttons and AC charger socket when cleaning	3	Manual pipettes: Do not allow cleaning solutions to enter the plunger and gear mechanism
Location	Step	Action																								
Room 1 <ul style="list-style-type: none"> 5% Extran 70% alcohol 	1	Allow 5% Extran to sit for approximately 5 min																								
	2	Remove Extran from surfaces with 70% alcohol																								
	3	Lay pipettes flat in hood																								
Room 2 and 3 <ul style="list-style-type: none"> Sani-Cloth Bleach Water 70% alcohol 	1	Allow bleach to sit 4 – 5 min																								
	2	Remove bleach residue with water followed by 70% alcohol																								
	3	Lay pipettes flat in hood																								
Pipettes (additional information)	1	Carefully clean the pipettor handle and barrel																								
	2	Electronic pipettes: Avoid electronic buttons and AC charger socket when cleaning																								
	3	Manual pipettes: Do not allow cleaning solutions to enter the plunger and gear mechanism																								
	4	The above solutions are corrosive and over time can contribute to the electronic components of the pipettes and increased resistance when adjusting volumes on manual pipettes																								
UV Hoods	7	UV-irradiate hoods at the end of the day <table border="1"> <thead> <tr> <th>Step</th> <th>Action</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Turn off lights</td> </tr> <tr> <td>2</td> <td>Lower sash</td> </tr> <tr> <td>3</td> <td>Turn on UV light for 15 min</td> </tr> </tbody> </table>	Step	Action	1	Turn off lights	2	Lower sash	3	Turn on UV light for 15 min																
Step	Action																									
1	Turn off lights																									
2	Lower sash																									
3	Turn on UV light for 15 min																									
Tacky mats	8	Change tacky mattes in rooms 1, 2 and 3 daily; more frequently if needed																								
Record	9	Fill out daily maintenance log; initial																								

Table 1: Routine Decontamination Schedule (*increase frequency if contamination/spills occur*)

Room	Step	Action	Frequency
Reagent Prep Room 1	1	Clean hood and UV-irradiate hood	Daily
	2	Clean benchtops, pipettes, racks, and cooling blocks	Daily
	3	Replace lab coats	Weekly
	4	Replace tacky matt	Daily; more frequently if needed
	5	Discard waste	As needed
Specimen Processing Room 2	1	Clean hood	After each procedure
	2	Clean benchtops, pipettes and cooling blocks	After each procedure
	3	Clean specimen racks and scissors	After each use
	4	Clean centrifuges and rotors	As needed <ul style="list-style-type: none"> ▪ Possible contamination ▪ Spill clean-up
	5	Replace lab coats	Weekly
	6	Replace tacky matt	Daily; more frequently if needed
	7	Discard waste	Daily
Amplification Room 3	1	Clean benchtops	After each procedure; more frequently if required
	2	Clean Liaison MDX	As needed <ul style="list-style-type: none"> ▪ Possible contamination ▪ When returning instrument for service
	3	Replace lab coats	Weekly
	4	Discard amplification waste	After each run (Zip lock)
	5	Replace tacky matt	Daily; more frequently if needed

REFERENCES

1. Simplexa™ 3M™ Integrated Cycler Studio 5.0 , 3M™ Integrated Cycler Operator Manual Reference 34-8710-8382-9, PI.MOL1101.UD_REV. F for use with user defined assays, Focus Diagnostics 2009-2012, Focus Diagnostics, Inc. Cypress, CA [Simplexa Operator Manual](#)

Historical Record

Version	Written/Revised by:	Effective Date:	Summary of Revisions
1	P. Ackerman	1.23.16	Initial Version
2	P. Ackerman	07.20.16	Reformatted for CMS upload
3	P. Ackerman	03.29.17	Instrument name change from Focus Integrated Cycler to DiaSorin Liaison MDX; fixed hyperlinks for SharePoint upload