

# RVP Sample Handling and Storage

## PURPOSE

- This procedure provides instructions for sample acceptability, handling and storage

## SAFETY CONSIDERATIONS

- Standard precautions. Refer to [MB 2.02](#) Biohazard Containment
- Use of engineering controls: Refer to [MB 3.01](#) Engineering Controls to Prevent Nucleic Acid Contamination

## MATERIALS REQUIRED

Equipment	Media	Supplies
BioSafety Cabinet (BSC)	Viral Transport Media (VTM) with glass beads	Dispo Pipettes
Refrigerator, 2 - 8° C		2 ml cryovials
Vortex mixer		

## SAMPLE

### A. Acceptable specimens

Specimen type	Specimen code	Volume	Transport container
Nasopharyngeal swab	<b>NP</b>	2 wire NP swabs	<ul style="list-style-type: none"> <li>NP CultureSwab™ Rayon wire mini-tip swab in Liquid Stuart's, <b>green</b> top</li> </ul>
Nasopharyngeal aspirate	<b>NASP</b>	1 – 2 mL (0.5 ml minimum)	Sterile, plastic leak proof container
Nasopharyngeal washing	<b>NW</b>		
Bronchoscopy	<b>BR</b>		
Bronchial alveolar lavage	<b>BAL</b>		
Bronch aspirate	<b>BASP</b>		
Bronch wash	<b>BRW</b>		

### B. Transport and Storage: For additional information refer to [Specimen Collection Manual](#)

Temperature	Sample Stability in VTM	Purified Nucleic Acid Stability
Room temperature	1 hr	1 hr
Refrigerated 2 - 8° C	7 days	7 days
Frozen at ≤ - 70° C	1 year, thawed up to 2 times	1 year, thawed up to 2 times

## PROCEDURE A: Follow the activity below for sample handling

Activity	Step	Action	Related Doc
<b>Identification</b> Location: Microbiology	1	Verify that the patient identification on the primary container corresponds to the accompanying order	<a href="#">MB 1.01</a> Specimen Management <a href="#">Organizational policy 630.00</a> Laboratory Specimen Labeling
	2	Receive sample in Sunquest and generate label <ul style="list-style-type: none"> <li>Sunquest location: MC</li> <li>Sunquest code: RVP</li> </ul>	

Activity	Step	Action	Related Doc
Identification	3	Confirm the name on the label is the same as the name on primary container	
	4	Affix LIS accession label to corresponding primary container	
Quality	5	Evaluate the quality of the sample.	<a href="#">MB 1.02</a> Specimen Rejection Criteria
	6	Vortex nasal/bronchial washes to obtain an even suspension	
Process	7	Transfer 1 – 2 ml of nasal/ bronch sample or 2 NP swabs into VTM <ul style="list-style-type: none"> <li>▪ Note: If NP swabs are used, leave swabs in VTM; do not discard</li> </ul>	
	8	Vortex for 30 – 60 s to break up mucus and release virus from cells	
Store	9	Store specimen refrigerator at 2 - 8° C	
	10	Place sample label in molecular designated area	
Aliquots	11	If additional testing is necessary, refer to Procedure B: <b>Sample Aliquots</b>	

**PROCEDURE B:** Follow the activity below for aliquoting samples and preventing cross-contamination

**Sample Aliquots**

Activity	Step	Action	Related Doc														
Identification of secondary container	1	Sample identification of all aliquots must be traceable to the primary specimen	<a href="#">Organizational policy 630.00</a> Laboratory Specimen Labeling														
	2	Confirm the name and accession number on the aliquot label is the same as on the primary container															
	3	Affix LIS aliquot label with corresponding accession number on secondary container															
Avoiding cross-contamination	4	Handle specimens to avoid cross contamination of primary sample and aliquots as follows:	Refer to assay specific procedures for additional information														
		<table border="1"> <thead> <tr> <th>Step</th> <th>Action</th> </tr> </thead> <tbody> <tr> <td>a</td> <td>Deliver primary specimens unopened to the molecular laboratory when possible or perform the following steps</td> </tr> <tr> <td>b</td> <td>Properly label secondary container with patient aliquot label</td> </tr> <tr> <td>c</td> <td>Verify name on primary and secondary container before transfer</td> </tr> <tr> <td>d</td> <td>Use sterile pipettes and technique when transferring samples</td> </tr> <tr> <td>e</td> <td>Aliquot one specimen at a time with only one tube open at a time</td> </tr> <tr> <td>f</td> <td>Never return the aliquot to the original container</td> </tr> </tbody> </table>		Step	Action	a	Deliver primary specimens unopened to the molecular laboratory when possible or perform the following steps	b	Properly label secondary container with patient aliquot label	c	Verify name on primary and secondary container before transfer	d	Use sterile pipettes and technique when transferring samples	e	Aliquot one specimen at a time with only one tube open at a time	f	Never return the aliquot to the original container
		Step		Action													
		a		Deliver primary specimens unopened to the molecular laboratory when possible or perform the following steps													
		b		Properly label secondary container with patient aliquot label													
		c		Verify name on primary and secondary container before transfer													
		d		Use sterile pipettes and technique when transferring samples													
e	Aliquot one specimen at a time with only one tube open at a time																
f	Never return the aliquot to the original container																

**REFERENCES**

1. eSensor® Respiratory viral Panel, PI1032 REV:D, December 2013, Clinical Micro Sensors, Inc. dba GenMark Diagnostics, Inc., 5964 La Place Court, Carlsbad, CA 92008, 1-800-373-6767, ww.genmarkdx.com

**Historical Record**

Version	Written/Revised by:	Effective Date:	Summary of Revisions
1	P. Ackerman	05.06.2015	Initial Version
2	P. Ackerman	08.27.2016	Reformatted for CMS upload; changed logo