

# RVP New Lot and/or New Shipment Quality Control

## PURPOSE

- This procedure provides instructions for verifying RVP reagent performance

## ABBREVIATIONS

- CLSI: Clinical Laboratory and Standards Institute
- EXC: extraction control
- MM: master mix
- NEGC: negative control
- NFW: nuclease free water
- RT-PCR: reverse transcription polymerase chain reaction
- PCTL: process control
- POSC: positive control
- QC: quality control
- RT: room temperature
- RVP: Respiratory Viral Panel
- VTM: viral transport media
- Area/Room 1: Clean room
- Area/Room 2: Processing room
- Area/Room 3: Amplification room

## MATERIALS

Equipment	Reagents	Supplies
<b>Room 1</b> <ul style="list-style-type: none"> <li>▪ Adjustable pipettes</li> <li>▪ Cold block</li> <li>▪ Freezer, -20° C</li> <li>▪ Laminar air-flow hood</li> <li>▪ Refrigerator 2 – 8° C</li> <li>▪ Vortex mixer</li> </ul> <b>Room 2</b> <ul style="list-style-type: none"> <li>▪ Adjustable pipettes</li> <li>▪ BioHit 8 channel pipette</li> <li>▪ Bio-Safety Cabinet (BSC)</li> <li>▪ Cold Block</li> <li>▪ Freezer, -70° C</li> <li>▪ Magnetic rack</li> <li>▪ Mini-centrifuge</li> <li>▪ NucliSens easyMag</li> <li>▪ Refrigerator 2 – 8° C</li> <li>▪ Tube racks, 1.5 – 2 ml</li> <li>▪ Vortex mixer</li> </ul> <b>Room 3</b> <ul style="list-style-type: none"> <li>▪ Adjustable pipettes</li> <li>▪ Cold Block</li> <li>▪ Freezer, -20° C</li> <li>▪ GenMark eSensor XT-8 instrument</li> <li>▪ Mini-centrifuge</li> <li>▪ PCR thermocycler</li> <li>▪ PCR workstation</li> <li>▪ Vortex mixer</li> </ul>	eSensor RVP kit: Product No. MT005102	Sterile filtered 10 µl pipette tips
	easyMAG Lysis buffer, 2 ml	Sterile filtered 30 µl pipette tips
	easyMAG Buffer 1	Sterile filtered 100 µl pipette tips
	easyMAG Buffer 2	Sterile filtered 200 µl pipette tips
	easyMAG Buffer 3	Sterile filtered 1000 µl pipette tips
	MagSil	Micro tubes 1.5 ml, RNase/DNase free
	Molecular grade water, nuclease free	Nitrile gloves (powder-free)
	Viral transport media (VTM)	PCR 8 tube strips with caps
	Viral isolates: H1, H3, RSV, Flu B	easyMag disposable vessel strips and tips
	Patient hMPV sample	BioHit pipette tips
	Sani-Cloth Bleach Wipes (10%)	BioHazard wipes
	70% alcohol	Gripper rack
	Household bleach	Sharps disposal container
	MMQCI RVP Control Panel	

**PROCEDURE A:** Follow the instructions for use of the MMQCI RVP control panel  
**RVP control panel**

Activity	Step	Action	Related Doc
<b>MMQCI RVP Control Panel</b>  1X use	1	The RVP control panel consists of 2 vials M244 and M245, single use only	
	2	Allow the vials to warm to RT	
	3	Vortex each vial for 5 s prior to use	<a href="#">MB 11.08.F2</a> Reagent QC  <a href="#">MB 11.08.F3</a> MMQCI Inventory Record
	4	Spin for 5 s to pull down the matrix	
	5	Extract both vials on the easyMag, 200 µl supernatant adding 10 µl IC to each; final elution 60 µl	
	6	Vortex the eluate for 5 sec; allow to sit in magnetic rack for 10 min	
	7	Analyze extracted controls; proceed to Procedure B for testing	

**PROCEDURE B:** Follow the activities for testing reagent reactivity in the table below  
**RVP New reagent lot and/or new shipment verification**

Activity	Step	Action	Related Doc
Testing requirements	1	Test new reagent lots and/or shipments against old reagent lots before or concurrently with being placed in service.	<a href="#">MB 5.02</a> Standards of Practice
MMQCI Control Panel	2	Each analyte detected in the RVP assay must be challenged with the MMQCI RVP panel to verify that the reagents are working properly	<a href="#">MB 11.08.F2</a> Reagent QC  <a href="#">MB 11.08.F1</a> RVP Inventory Record
	3	The MMQCI panel must have been tested on the previous RVP lot prior to testing the new lot/shipment	
	4	Include a POSC/EXC and NEGC with known performance in each test run	
Results	5	Record results on RVP QC worksheet; staple QC worksheet to the eSensor reports	
Records	6	Check off RVP inventory form, date and initial	
	7	Archive result forms in <i>New Lot Inventory and QC</i> manual	
	8	The MMQCI results must be equivalent to the previous RVP lot/shipment to confirm acceptability of the new RVP lot/shipment	
	9	Verify that all reagents and materials meet expiration date and QC parameters as per CLSI document MM3-A2.	

**PROCEDURE C:** Follow the activities for troubleshooting verification failures in the table below  
**Performance Failures**

Activity	Step	Action	Related doc														
Troubleshooting Failures	1	Verify that the reagent performance is acceptable before implementation of a new lot and/or shipment	<a href="#">RVP Retest</a> Recommendations by Report Type  <a href="#">RVP Technical Support</a> and Troubleshooting														
		<table border="1"> <thead> <tr> <th>If</th> <th>Then</th> </tr> </thead> <tbody> <tr> <td>Any Control fails</td> <td> <ul style="list-style-type: none"> <li>Document observation/corrective action on QC log</li> <li>Do not implement new lot/shipment</li> <li>Repeat all testing, extraction, PCR and XT-8 analysis</li> </ul> </td> </tr> <tr> <td>IC fails</td> <td> <ul style="list-style-type: none"> <li>Poor amplification of sample</li> <li>Poor recovery from extracted sample</li> <li>System error</li> <li>If one or more targets are positive, retest not necessary</li> <li>No targets positive, repeat extraction, PCR and XT-8 analysis</li> </ul> </td> </tr> <tr> <td>POSC/EXC fails</td> <td> <ul style="list-style-type: none"> <li>Amplification failure</li> <li>Poor recovery from extracted sample</li> <li>Possible reagent or system failure: Review MM preparation and assay set-up</li> <li>Repeat run extraction, PCR and XT-8 analysis</li> </ul> </td> </tr> <tr> <td>NEGC fails Contamination</td> <td> <ul style="list-style-type: none"> <li>Possible carryover or reagent contamination: Review pipetting technique, glove contamination, possible aerosol creation, and MM preparation</li> <li>Repeat run extraction, PCR and XT-8 analysis</li> </ul> </td> </tr> <tr> <td>NEGC fails System error</td> <td> <ul style="list-style-type: none"> <li>Report message reads "error" for any target</li> <li>Rerun cartridge; call for password, GenMark technical service at <b>1.800.373.6767</b></li> <li>Retest NEGC only; repeat PCR and XT-8 analysis</li> </ul> </td> </tr> <tr> <td>Problem unresolved</td> <td> <ul style="list-style-type: none"> <li>Call GenMark technical service at <b>1.800.373.6767</b></li> <li>Notify section technical director or designee</li> </ul> </td> </tr> </tbody> </table>		If	Then	Any Control fails	<ul style="list-style-type: none"> <li>Document observation/corrective action on QC log</li> <li>Do not implement new lot/shipment</li> <li>Repeat all testing, extraction, PCR and XT-8 analysis</li> </ul>	IC fails	<ul style="list-style-type: none"> <li>Poor amplification of sample</li> <li>Poor recovery from extracted sample</li> <li>System error</li> <li>If one or more targets are positive, retest not necessary</li> <li>No targets positive, repeat extraction, PCR and XT-8 analysis</li> </ul>	POSC/EXC fails	<ul style="list-style-type: none"> <li>Amplification failure</li> <li>Poor recovery from extracted sample</li> <li>Possible reagent or system failure: Review MM preparation and assay set-up</li> <li>Repeat run extraction, PCR and XT-8 analysis</li> </ul>	NEGC fails Contamination	<ul style="list-style-type: none"> <li>Possible carryover or reagent contamination: Review pipetting technique, glove contamination, possible aerosol creation, and MM preparation</li> <li>Repeat run extraction, PCR and XT-8 analysis</li> </ul>	NEGC fails System error	<ul style="list-style-type: none"> <li>Report message reads "error" for any target</li> <li>Rerun cartridge; call for password, GenMark technical service at <b>1.800.373.6767</b></li> <li>Retest NEGC only; repeat PCR and XT-8 analysis</li> </ul>	Problem unresolved	<ul style="list-style-type: none"> <li>Call GenMark technical service at <b>1.800.373.6767</b></li> <li>Notify section technical director or designee</li> </ul>
	If	Then															
	Any Control fails	<ul style="list-style-type: none"> <li>Document observation/corrective action on QC log</li> <li>Do not implement new lot/shipment</li> <li>Repeat all testing, extraction, PCR and XT-8 analysis</li> </ul>															
	IC fails	<ul style="list-style-type: none"> <li>Poor amplification of sample</li> <li>Poor recovery from extracted sample</li> <li>System error</li> <li>If one or more targets are positive, retest not necessary</li> <li>No targets positive, repeat extraction, PCR and XT-8 analysis</li> </ul>															
	POSC/EXC fails	<ul style="list-style-type: none"> <li>Amplification failure</li> <li>Poor recovery from extracted sample</li> <li>Possible reagent or system failure: Review MM preparation and assay set-up</li> <li>Repeat run extraction, PCR and XT-8 analysis</li> </ul>															
	NEGC fails Contamination	<ul style="list-style-type: none"> <li>Possible carryover or reagent contamination: Review pipetting technique, glove contamination, possible aerosol creation, and MM preparation</li> <li>Repeat run extraction, PCR and XT-8 analysis</li> </ul>															
	NEGC fails System error	<ul style="list-style-type: none"> <li>Report message reads "error" for any target</li> <li>Rerun cartridge; call for password, GenMark technical service at <b>1.800.373.6767</b></li> <li>Retest NEGC only; repeat PCR and XT-8 analysis</li> </ul>															
Problem unresolved	<ul style="list-style-type: none"> <li>Call GenMark technical service at <b>1.800.373.6767</b></li> <li>Notify section technical director or designee</li> </ul>																

**REFERENCES**

1. eSensor® Respiratory viral Panel, PI1032 REV:D, December 2013, Clinical Micro Sensors, Inc. dba GenMark Diagnostics, Inc., 5964 La Place Court, Carlsbad, CA 92008, 1-800-373-6767, [www.genmarkdx.com](http://www.genmarkdx.com)
2. eSensor XT-8 RVP Control Panel package insert; circular M243 102914.001, Maine Molecular Quality Controls, Inc. [www.mmqci.com](http://www.mmqci.com)
3. CLSI *Molecular Diagnostic Methods for Infectious Diseases*; Approved Guideline – Second Edition, CLSI document MM3-A2, Wayne, PA, Clinical and Laboratory Standards Institute; 2006
4. CLSI *Establishing Molecular Testing in Clinical Laboratory Environments*; Approved Guideline, MM19-A, Vol. 31. No. 21, Wayne, PA, Clinical and Laboratory Standards Institute; 2011
5. CAP All Common Checklist COM.30450 New reagent Lot Confirmation of Acceptability revised 7/29/2013

**Historical Record**

Version	Written/Revised by:	Effective Date:	Summary of Revisions
1	P. Ackerman	07.25.15	Initial Version
2	P. Ackerman	08.28.2016	Reformatted for CMS upload; changed logo