

Situation: We recently had a patient that had a discrepancy between the WBC-BF (White Blood Cell Count) and the TC-BF# (Total Nucleated Cell Count) on a CSF that was analyzed on the Sysmex. In most cases these numbers will always match or at least be very close together. In our laboratory we will only report a WBC count from the analyzer when these numbers match.

Background: In the attached example, (#1) I have highlighted the discrepancy between the WBC count and the TNC count. In the WDF scattergram there is a large population of cells that is being excluded from the WBC count.

In example #2 there is obviously an abnormal scattergram in the WDF channel, the cell counts have been flagged and the WBC Abn Scattergram message is present.

Assessment: In example #1 there is no "cue" for the tech that would alert them to this discrepancy. I have enabled the WBC Abn Scattergram message to appear whenever there is a discrepancy of 5% or more between "normal" WBC's and other nucleated cells that would be counted in this "High Fluorescent (HF) Region" I have also included the rule "use manual method for counting".

Example #3 shows the same patient analyzed again, note the highlighted codes.

Recommendation: Any body fluid counts that appear with asterisks, WBC Abn Scattergram flag, or the message "use manual method for counting" should be counted manually using the hemocytometer method.


Sample No.: W27575
Patient ID: ████████
Name:
Sample Comment:

Rack:
Ward:

Position: 06/22/2017 11:14:31
Doctor:
Birth:
Sex:
Instrument ID: XN-10^24041

EXAMPLE # 1

WBC	[10 ³ /uL]	
RBC	[10 ⁶ /uL]	
HGB	[g/dL]	
HCT	[%]	
MCV	[fL]	
MCH	[pg]	
MCHC	[g/dL]	
RDW-SD	[fL]	
RDW-CV	[%]	
PLT	[10 ³ /uL]	
MPV	[fL]	
NEUT	[10 ³ /uL]	[%]
LYMPH	[10 ³ /uL]	[%]
MONO	[10 ³ /uL]	[%]
EO	[10 ³ /uL]	[%]
BASO	[10 ³ /uL]	[%]
IG	[10 ³ /uL]	[%]
NRBC	[10 ³ /uL]	[%]
RET	[%]	[10 ⁶ /uL]
IRF	[%]	
RET-He	[pg]	
IPF	[%]	
WBC-BF	0.124	
RBC-BF	0.000	
MN	0.117	94.3 [%]
PMN	0.007	5.7 [%]
TC-BF#	0.798	



NEUT	
BAND	
LYMPH	
MONO	
EOS	
BASO	

META	
MYELO	
PRO	
BLAST	
ATYP	
NRBC	

MACRO	
MICRO	
HYPO	
POLY	
OVAL	
FRAG	

PLT	
OTHER	
TECH:	

WBC IP Message

RBC IP Message

PLT IP Message

Sample No.: HH7273
Patient ID: 2351232
Name: ~~ELSON SAMUEL JOSEPH~~
Sample Comment:

Rack:
Ward: PRIV

Position: 06/08/2017 15:11:58
Doctor: MCNAMARA, JOHN J MD
Birth: 06/22/2000 Sex: Male
Instrument ID: XN-10^24042

Positive

#1: Make smear and perform WBC est/MAND


Morph.

WBC [10³/uL]
RBC [10⁶/uL]
HGB [g/dL]
HCT [%]
MCV [fL]
MCH [pg]
MCHC [g/dL]
RDW-SD [fL]
RDW-CV [%]
PLT [10³/uL]
MPV [fL]

NEUT [10³/uL] [%]
LYMPH [10³/uL] [%]
MONO [10³/uL] [%]
EO [10³/uL] [%]
BASO [10³/uL] [%]
IG [10³/uL] [%]
NRBC [10³/uL] [%]
RET [%] [10⁶/uL]
IRF [%]
RET-He [pg]
IPF [%]

WBC-BF 16.927 * [10³/uL]
RBC-BF 0.026 [10⁶/uL]
MN 11.877 * [10³/uL] 70.2 * [%]
PMN 5.050 * [10³/uL] 29.8 * [%]
TC-BF# 18.192 * [10³/uL]

EXAMPLE #2



NEUT	
BAND	
LYMPH	
MONO	
EOS	
BASO	

META	
MYELO	
PRO	
BLAST	
ATYP	
NRBC	

MACRO	
MICRO	
HYPO	
POLY	
OVAL	
FRAG	

PLT	
OTHER	
TECH:	

WBC IP Message

RBC IP Message

PLT IP Message

WBC Abn Scattergram

Body Fluid Analysis IP Message

Abnormal, WBC Abn Scattergram


The WBC Abn Scattergram IP message is generated during body fluid analysis whenever clustering in the WDF scattergram is abnormal. This flag can also be generated when the HF-BF#* or HF-BF%* result exceeds the HF-BF user defined limits in the WBC Abnormal Flag (Body Fluid Analysis) setting.

*These parameters are not reportable, and are used only in the algorithm for the WBC Abn Scattergram flag during body fluid analysis.

Dashes may appear in place of data that was not calculated.

NOTE: Body fluid analysis is not available on all XN-series analyzer configurations.

XN-Series Results



WDF Scatter
(Close-Up)



Normal WDF Scatter
(Reference)



Abnormal, WBC Abn Scattergram (continued)

Suggested Action Steps:

1. Dashes (— —) in place of numeric data:
 - Verify WBC-BF, TNC-BF and differential results according to your laboratory's policy. Possible actions may include:
 - repeating the sample
 - diluting the sample with CELLPACK DCL diluent to minimize interference.
 - performing a manual cell count
 - performing a manual differential
2. Asterisk (*) next to results:
 - Verify WBC-BF, TNC-BF and differential results according to your laboratory's policy. Possible actions may include:
 - performing a manual cell count
 - performing a manual differential
 - If the WBC ABN Scattergram IP Message is present due to HF-BF results exceeding the HF-BF user defined limits— scan the slide for the presence of mesothelial and/or abnormal cells. Report results according to your laboratory protocol.
 - If no abnormalities are found when reviewing the smear and the WBC and TNC estimate matches the analyzer reported WBC-BF and TNC-BF, the results with asterisks (*) may be reported.

NOTE: Decisions to report with a comment, perform a dilution or perform an alternate method should be based on your local laboratory protocol.

Sample No.: W27575
Patient ID: XXXXXXXXXX
Name:
Sample Comment:

Ward: Rack:

Position: 06/22/2017 12:09:50
Doctor: BF
Birth: Sex:
Instrument ID: XN-10^24041


Positive

#1: Make smear and perform WBC est/MAND
#89: Use manual method for counting

Morph.

WBC	[10 ³ /uL]
RBC	[10 ⁶ /uL]
HGB	[g/dL]
HCT	[%]
MCV	[fL]
MCH	[pg]
MCHC	[g/dL]
RDW-SD	[fL]
RDW-CV	[%]
PLT	[10 ³ /uL]
MPV	[fL]
NEUT	[10 ³ /uL]
LYMPH	[10 ³ /uL]
MONO	[10 ³ /uL]
EO	[10 ³ /uL]
BASO	[10 ³ /uL]
IG	[10 ³ /uL]
NRBC	[10 ³ /uL]
RET	[%]
IRF	[%]
RET-He	[pg]
IPF	[%]
WBC-BF	0.125 [10 ³ /uL]
RBC-BF	0.000 [10 ⁶ /uL]
MN	0.108 [10 ³ /uL]
PMN	0.017 [10 ³ /uL]
TC-BF#	0.765 [10 ³ /uL]

EXAMPLE #3



[%]
[10⁶/uL]

86.4 [%]
13.6 [%]

NEUT	
BAND	
LYMPH	
MONO	
EOS	
BASO	

META	
MYELO	
PRO	
BLAST	
ATYP	
NRBC	

MACRO	
MICRO	
HYPO	
POLY	
OVAL	
FRAG	

PLT	
OTHER	
TECH:	

WBC IP Message


RBC IP Message

PLT IP Message


WBC Abn Scattergram

XN-3000 Quick Guide


WNR




WDF




RET



PLT-F



WDF—Body Fluid Mode



(NORMAL SCATTERGRAM)

