
FORMULAS FOR MANUAL CALCULATIONS (LIS Downtime)

RC.CH.LOP.LIS.RG.001.r03

Anion Gap = Sodium – (Cl + CO₂)

Globulin = Total Protein – Albumin

A/G Ratio = Albumin / (Total Protein – Albumin)

Chol/HDL Ratio = Chol/HDL

Non HDL Cholesterol = Chol – HDL = mg/dL

LDL Cholesterol = Total Chol – HDL – (Trig/5)

Note: Formula for calculation only if Trig < 400 mg/dL

BUN/Creatinine Ratio = BUN/Creatinine

CKD-EPI Estimated GFR (mL/min/1.73 m²) =

$141 \times \min(S_{Cr}/\kappa, 1)^{\alpha} \times \max(S_{Cr}/\kappa, 1)^{-1.209} \times 0.993^{Age} \times 1.018 \text{ [if female]} \times 1.159 \text{ [if black]}$

S_{Cr} is serum creatinine in mg/dL,

κ is 0.7 for females and 0.9 for males,

α is -0.329 for females and -0.411 for males,

min indicates the minimum of S_{Cr}/κ or 1, and

max indicates the maximum of S_{Cr}/κ or 1.

Example for a female patient:

$141 \times \min(S_{Cr}/0.7, 1)^{-0.329} \times \max(S_{Cr}/0.7, 1)^{-1.209} \times 0.993^{Age} \times 1.018 \times 1.159 \text{ [if black]}$

Example for a male patient:

$141 \times \min(S_{Cr}/0.9, 1)^{-0.411} \times \max(S_{Cr}/0.9, 1)^{-1.209} \times 0.993^{Age} \times 1.159 \text{ [if black]}$

(Use www.kidney.org to calculate)

% Saturation = $100 \times [\text{serum Iron} / (\text{Transferrin} \times 1.4)]$

TIBC = Transferrin $\times 1.4$

CK-MB Index = $(\text{CKMB} / \text{CK}) \times 100$

Formulas for Manual Calculations (LIS Downtime)

Urine MicroAlbumin/Creatinine Ratio = (ualb mg/dL/Urine Creatinine mg/dL) x 1000

Urine Creatinine Clearance mL/min

$$\frac{\text{Raw Urine Creatinine mg/dL}}{\text{Serum Creat mg/dL}} \times \left(\frac{\text{TV mL}}{\text{Hrs collection}} \times \frac{1 \text{ Hr}}{60 \text{ min}} \right)$$

Urine Creatinine Clearance Pediatric mL/min CC Corrected =

$$\frac{\text{Raw Ur Creat mg/dL}}{\text{Serum Creat mg/dL}} \times \left(\frac{\text{TV mL}}{\text{Hrs collection}} \times \frac{1 \text{ Hr}}{60 \text{ min}} \right) \times \left(\frac{1.73 \text{ m}^2}{\text{surface area m}^2} \right)$$

NOTE: Surface area m² is determined from child's height and weight.

Urine Creatinine Timed or 24 Hour mg/Coll =

$$\text{Raw Urine Creatinine } \frac{\text{mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10 \text{ dL}}{1000 \text{ mL}} \right)$$

Urine Amylase Timed U/Hr =

$$\text{Raw Urine Amylase } \frac{\text{U}}{\text{L}} \times \text{dilution} \times \frac{\text{TV mL}}{\text{Hrs of Coll.}} \times \frac{1 \text{ L}}{1000 \text{ mL}}$$

Urine Calcium Timed or 24 Hour, mg/Collection =

$$\text{Raw Urine Calcium } \frac{\text{mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10 \text{ dL}}{1000 \text{ mL}} \right)$$

Urine Glucose Timed and 24 Hour, g/Collection =

$$\text{Raw Urine Glucose } \frac{\text{mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10 \text{ dL}}{1000 \text{ mL}} \right) \times \frac{\text{g}}{1000 \text{ mg}}$$

Urine Magnesium Timed or 24 Hour, mg/Collection =

$$\text{Raw Magnesium } \frac{\text{mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10 \text{ dL}}{1000 \text{ mL}} \right)$$

Urine Phosphorous Times or 24 Hour, g/Collection =

$$(\text{Raw Urine Phos } \frac{\text{mg}}{\text{dL}} \times \frac{1 \text{ g}}{1000 \text{ mg}}) \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10 \text{ dL}}{1000 \text{ mL}} \right)$$

Printed copies of this document are not considered up-to-date. Please verify current version date with online document.

Formulas for Manual Calculations (LIS Downtime)

Urine Protein Timed or 24 Hour, mg/Collection =

$$\text{Raw Urine Protein } \frac{\text{mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right)$$

Urine Urea Nitrogen Timed or 24 Hour, mg/Collection =

$$\text{Raw Urine Urea Nitrogen } \frac{\text{mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right)$$

Urine Uric Acid Timed or 24 Hour, mg/Collection =

$$\text{Raw Urine Uric Acid } \frac{\text{mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right)$$

Urine Sodium Timed or 24 Hour, mmol/Collection =

$$\text{Raw Urine Sodium } \frac{\text{mmol}}{\text{L}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{1\text{L}}{1000\text{mL}} \right)$$

Urine Potassium Timed or 24 Hour, mmol/Collection =

$$\text{Raw Urine Potassium } \frac{\text{mmol}}{\text{L}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{1\text{L}}{1000\text{mL}} \right)$$

Urine Chloride Timed or 24 Hour, mmol/Collection =

$$\text{Raw Urine Chloride } \frac{\text{mmol}}{\text{L}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{1\text{L}}{1000\text{mL}} \right)$$

% Urine Amylase/Creatinine Clearance Ratio =

$$\frac{\text{Raw Urine Amylase U/L}}{\text{Serum Amylase U/L}} \times \frac{\text{Serum Creatinine mg/dL}}{\text{Raw Urine Creatinine mg/dL}} \times 100$$

NOTE: reference range 1.3-4.3%

Urine Protein/ Creatinine Ratio =

$$\frac{\text{Raw Urine Protein mg/dL}}{\text{Raw Urine Creatinine mg/dL}}$$

NOTE: reference range 0.0-0.2

Printed copies of this document are not considered up-to-date. Please verify current version date with online document.

Formulas for Manual Calculations (LIS Downtime)

Urine Calcium/ Creatinine Ratio =

$$\frac{\text{Raw Urine Calcium mg/dL}}{\text{Raw Urine Creatinine mg/dL}}$$

NOTE: reference range 0.02-0.26

AAPO₂, mmHg =

$$((713\text{mmHg} \times \text{FiO}_2 \text{ as decimal}) - (\text{PCO}_2 / 0.8)) - \text{PO}_2$$

PF, no units =

$$\text{PO}_2 / \text{FiO}_2 \text{ (as a decimal value)}$$

% Free PSA =

$$(\text{Free PSA/PSA}) \times 100$$

Albumin gradient; Serum and Ascites fluid, g/dL =

$$\text{Albumin (g/dL)} - \text{Albumin fluid (g/dL)}$$

NOTE: >1.1 = Transudate

0.2-1.0 = Exudate

Authorized Reviewers: Section Medical or Technical Director

Printed copies of this document are not considered up-to-date. Please verify current version date with online document.

Formulas for Manual Calculations (LIS Downtime)

Document Location of Master:

Control Master electronic file stored on the Clinical Pathology server under
S:/ Automated Chemistry/ Document Control Library/ LOP/LIS/Masters
Master printed document stored in the General Policy and Procedure Manual, Core Lab

Number of Controlled Copies posted for educational purposes: 0

Number of circulating Controlled Copies: 1

Location of circulating Controlled Copies:Automated Chemistry General Policy & Procedure Manual, Stat Lab

Uncontrolled Copies posted in Stat Lab (red) downtime binders

Document History

Signature	Date			
Prepared by: M. Landskroener MT(ASCP)	12/05/2008			
Approved by:				
Reviewed by: (Signature)	Date	Revision	Modification	Related Documents Reviewed/Updated
Raymond E. Karcher, PhD	12/05/2008			
Raymond E. Karcher, PhD	12/04/2009			
Vivek Kumar, PhD	12/07/2010			
Vivek Kumar, PhD	02/14/2012	r01	Added the following equations: Non HDL Chol, %CPSA, % Free PSA, Albumin gradient: serum and ascites fluid. Changed microalbumin units to mg/dL,	
Kenneth Simkowski	11/19/2014			
Steven Truscott, PhD	09/28/2015	r02	Updated GFR equation to CKD-EPI (2009 version)	
Kenneth Simkowski, PhD	11/03/2017			
Elizabeth Sykes, MD	02/02/2018			
Peter Millward, MD	09/17/2018		New medical director	
Peter Millward, MD	11/19/2018			
Updated by: Robin Carey-Ballough MT(ASCP)	10/31/19	r03	Updated Calculations to Abbott methods. Removed calculations not in use	
Approved: Qian Sun PhD	11/1/2019			

Printed copies of this document are not considered up-to-date. Please verify current version date with online document.