

Albumin	The quantitative determination of albumin can aid in the diagnosis and management of numerous diseases including those involving the liver and kidneys. It may also be used to assess nutritional status although prealbumin is a better indicator of malnutrition.
A1-Antitrypsin	More than 70 genetic variants of alpha-1-antitrypsin (AAT) have been described. Not all of these are associated with decreased AAT levels or with clinical disease. An individual homozygous for PiZ has about 15% normal AAT and Pi null has no AAT. Such individuals are at significantly increased risk for development of pulmonary emphysema at an earlier age than individuals with a normal AAT phenotype; this process is accelerated by smoking. Development of liver disease may occur in infants (hepatitis and cirrhosis) and in older individuals (chronic hepatitis and cirrhosis). In individuals with a decreased level of AAT, Alpha 1 Antitrypsin Phenotyping is recommended.
ACE	The angiotensin converting enzyme assay is used to aid in the diagnosis of active sarcoidosis. It may be useful for confirmation of Gaucher's disease.
Alk Phos	The quantitative determination of alkaline phosphatase activity aids in the diagnosis and management of liver and bone diseases.
ALT	The quantitative determination of ALT (alanine aminotransferase) aids in the diagnosis and management of liver disease.
Amylase	Amylase measurements aid in the diagnosis and management of pancreatitis (inflammation of the pancreas).
ASO	Streptolysin O is one of several immunogenic exoenzymes produced by Group A, Beta-hemolytic streptococci. An elevated anti-streptolysin O (ASO) titer is usually indicative of a recent infection with a group A streptococci and is a routine part of the diagnosis and management of acute rheumatic fever and acute glomerulonephritis. In the absence of complications or reinfection, antibody levels usually fall to preinfection levels 6-12 months following infection.
AST	The quantitative determination of AST (aspartate aminotransferase) aids in the diagnosis and management of certain types of liver disease. AST is no longer recommended for diagnosis of myocardial infarction.
B2-Microglobulin	This assay is used to evaluate renal disease and to assess prognosis and monitor lymphoproliferative disorders, such as multiple myeloma and chronic lymphocytic leukemia.
BUN	The BUN assay is frequently requested in conjunction with the serum creatinine test for the differential diagnosis of prerenal (cardiac decompensation, water depletion, increases protein catabolism), renal (glomerulonephritis, chronic nephritis, polycystic kidney, nephrosclerosis, tubular necrosis), and postrenal (obstructions of the urinary tract) hyperuremia.
C3 Complement	C3 quantitation is used to detect individuals with congenital C3 deficiencies, or individuals who have depleted their complement levels due to an immunological process. C3 has the highest serum concentration of any complement component. Complement C3 is used as a screening test because it is consumed by activation of either the classical or alternative pathway. Individual assays for C3 and C4 are most useful in monitoring patients with immunologic diseases. Functional assays (e.g., the CH50 test) measure the activity of the entire complement cascade and are more likely to detect inherited deficiencies.
C4 Complement	C4 measurements should be performed whenever a complement activating disease is suspected or whenever hyposynthesis due to inherited deficiency is a possibility. C4 is the second most abundant complement protein in serum. C4 is only used in the classical pathway. Conditions affecting only the alternate pathway will not affect C4 levels. Individual assays for C3 and C4 are most useful in monitoring patients with immunologic diseases. Functional assays (e.g., the CH50 test) measure the activity of the entire complement cascade and are more likely to detect inherited deficiencies.
Calcium	The quantitative determination of calcium aids in the diagnosis and management of a variety of diseases including those involving the parathyroid glands, bone and kidneys.
Ceruloplasmin	Wilson's Disease is an autosomal recessive trait resulting in a copper metabolism disorder. It affects males and females equally. The onset of the disease is commonly seen in late childhood and early adult life. Affected individuals usually have ceruloplasmin levels less than 20 mg/dL. In these patients, free copper accumulates in selected areas of the body and may result in cirrhosis of the liver and central nervous system dysfunctions. These symptoms can improve with treatment. In untreated patients, the disease progresses and is usually fatal. Menke's Disease (also known as "kinky hair" disease) is a sex-linked disease that produces hypoceruloplasminemia. The disease affects only males and is characterized by steely hair, defective cross-linking of collagen and elastin, and neurologic findings. Menke's Disease is usually fatal within 3 years.
Cholesterol	Total cholesterol is used to assess the risk of ASCVD. It is recommended that HDL Cholesterol, non-LDL Cholesterol, LDL Cholesterol and Triglycerides also be obtained in initial screening. Several organizations have issued guidelines for management of dyslipidemias, all aiming to standardize and optimize patient care. The recent ACC/AHA guidelines aim to reduce/prevent heart disease, peripheral vascular disease and stroke by taking into account life style and lipid levels (1). Based on this information an estimate of ASCVD risk can be calculated and a decision on whether or not to treat (e.g. with statins) and modify life style can be made. The ACC/AHA guidelines do not recommend specific cholesterol set-points, but aim for a particular percent decrease in LDL cholesterol. Our lab will continue to use the ATP guideline cut-points in lipid reporting (2). The National Lipid Association also has recommendations that are similar to the ATP III guidelines (3). Total cholesterol may be decreased after acute myocardial infarction (AMI). Assessment of lipid status should therefore be determined within 24 hours of chest pain or 12 weeks following the AMI.
CK	The quantitative determination of CK (creatinine kinase) and its isoenzymes aid in the diagnosis and management of myocardial, skeletal, and muscle diseases.
Cl	The quantitation of chloride aids in the diagnosis and treatment of electrolyte and metabolic disorders such as acidosis or alkalosis, dehydration, renal failure and hormone imbalance.
CO2	The CO ₂ assay aids in the evaluation of acid-base balance.
Creatinine	The quantitative determination of creatinine aids in the diagnosis and treatment of renal diseases, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes. It is most sensitive in detecting renal impairment when used as part of a Creatinine Clearance test.

CRP	C-reactive protein (CRP) is an acute phase reactant and is used as a sensitive and quantitative measure of the body's acute phase response. CRP is not diagnostic for any specific disease. Increased CRP levels are consistently found in patients with acute bacterial and viral infections, rheumatoid arthritis, acute myocardial infarction, and widespread malignant disease. CRP levels respond to inflammation within 8 hours of onset and peak levels are reached within 24-48 hours. Levels may rise to 2000 times normal levels. CRP levels associated with viral infection, rheumatoid arthritis, and neoplasia are usually 10-40 mg/L. CRP levels of 40 to greater than 300 mg/L are usually associated with acute bacterial infections. Monitoring serum CRP levels aids in the detection and evaluation of post-operative complications associated with inflammation and/ or tissue necrosis. CRP levels evaluated 48-72 hours postoperatively may be 250-350 mg/L. These levels return to normal within one week.
CRP High Sensitivity	In an apparently healthy adult or in the absence of a known inflammatory process, hs-CRP assesses an individual's risk of developing a coronary event, stroke or peripheral vascular disease.
Dbili	The quantitative determination of direct bilirubin is used in the evaluation of liver and biliary disease.
GGT	The quantitative determination of GGT (gamma-glutamyl transferase) aids in the diagnosis and monitoring of hepatobiliary disease. GGT is currently the most sensitive enzymatic indicator of liver disease. Normal GGT values are rarely found in the presence of clinically significant hepatic disease.
Glucose	Glucose measurements are used in the diagnosis and management of disorders of carbohydrate metabolism; these include diabetes mellitus, neonatal hypoglycemia, idiopathic hypoglycemia, and pancreatic islet cell tumors.
Haptoglobin	Haptoglobin binds to hemoglobin released into the circulation by intravascular hemolysis. Haptoglobin is an acute phase reactant. Serial assays are used to detect and monitor hemolytic states. Haptoglobin is decreased in diseases associated with intravascular hemolysis. In severe hemolysis, haptoglobin may be totally depleted, requiring up to 1 week to return to normal. In chronic hemolytic states such as hemoglobinopathies and mechanical heart valves, there may be a steady decline in haptoglobin levels. In these conditions, serial assays provide a better index of ongoing hemolysis than a single haptoglobin value. Increased serum haptoglobin levels are present in infection, neoplasia, and other inflammatory conditions characterized by tissue injury and repair
HDL Chol	The quantitative determination of HDL (high-density lipoprotein) cholesterol aids in the diagnosis and treatment of lipid disorders. It should be used in the risk assessment of coronary artery disease (CAD). Low levels of HDL cholesterol are associated with an increased risk of CAD.
IgA	Selective IgA deficiency is a primary immunodeficiency disorder characterized by reduced production of IgA with recurrent respiratory and gastrointestinal infections. Selective IgA deficiency can result from congenital intrauterine infection with rubella virus, <i>Toxoplasma gondii</i> , or cytomegalovirus. A transient IgA deficiency may result following the treatment with penicillamine of Wilson's disease. Most patients with selective IgA deficiencies are asymptomatic. Symptomatic patients usually present with recurrent ear infections, sinusitis, pneumonia, diarrhea, asthma, autoimmune diseases and/or allergies. Administration of blood products containing IgA can cause some IgA deficient patients to develop antibodies against IgA. If an anti-IgA antibody develops, a massive allergic reaction can result during blood or plasma transfusions.
IgE	Hyperimmunoglobulinemia E syndrome is an autosomal dominant disorder of unknown cause. Patients present with recurrent staphylococcal abscesses involving the skin, lungs, joints, and soft tissues. Job's syndrome is closely related to this disorder. Patient's typically have normal serum immunoglobulin concentrations except for IgE, which is typically markedly elevated to levels in excess of 2000 IU/ml and as high as 20,000 to 50,000 IU/ml.
IgG	IgG levels can be used to evaluate humoral immunity and aids in the diagnosis of conditions associated with IgG excess or depression. IgG can cross the placenta. IgG antibodies are the most important and persistent antibodies of the secondary immune response.
IgM	IgM levels can be used to evaluate humoral immunity and to assist in the diagnosis of conditions associated with IgM excess or depression. IgM is the first antibody to appear in a primary antibody response. IgM does not cross the placenta. Increased IgM levels in the newborn are associated with intrauterine infections.
Iron	This assay is used in the evaluation of iron deficiency, hemochromatosis and to verify acute iron poisoning.
K	The quantitation of potassium is used to monitor electrolyte balance. Pseudohyperkalemia: If an increase in platelets or leukocytes is suspected as a cause of hyperkalemia, a heparin tube should be obtained for plasma potassium.
LD	Lactate dehydrogenase is present in multiple cells and tissue types and therefore its utility in patient diagnosis is questionable. It is useful in the assessment of in-vivo hemolysis and hematologic disorders (benign and malignant) conditions in which LD is often increased. LD is also increased following myocardial infarction, pulmonary and renal cortical infarction, liver disease, skeletal muscle disease, and many other conditions. Use of more specific enzymes and protein markers is preferable to LD in the diagnosis of myocardial infarction (use Troponin or CK-MB), skeletal muscle (use CK) and liver disease.
LDL Chol (Direct)	This direct quantitative determination of LDL cholesterol aids in the diagnosis and management of coronary atherosclerosis. Several organizations have issued guidelines for management of dyslipidemias, all aiming to standardize and optimize patient care. The recent ACC/AHA guidelines aim to reduce/prevent heart disease, peripheral vascular disease and stroke by taking into account life style and lipid levels. Based on this information an estimate of ASCVD risk can be calculated and a decision on whether or not to treat (e.g. with statins) and modify life style can be made. The ACC/AHA guidelines do not recommend specific cholesterol set-points, but aim for a particular percent decrease in LDL cholesterol. Our lab will continue to use the ATP guideline cut-points in lipid reporting. The National Lipid Association also has recommendations that are similar to the ATP III guidelines.
Lipase	This assay aids in the diagnosis of patients suspected of having acute pancreatitis.

Lipoprotein A1	Elevated levels of Lp(a) indicate a major risk for the development of atherosclerosis and coronary heart disease, independent of LDL-cholesterol levels. The wide differences in lipoprotein(a) levels seen among individuals are largely due to hereditary factors and cannot be controlled by dietary or lifestyle changes. Nevertheless, the identification of individuals at risk can be useful in alerting them to the need to eliminate or control other high risk factors.
Magnesium	The quantitation of magnesium aids in the investigation of unexplained hypocalcemia, in the management of patients following cardiac surgery or those with cardiac arrhythmias, and in the management of patients being treated for pre-eclampsia or eclampsia. An association between severe hypomagnesemia and aminoglycoside therapy has been described.
Na	The quantitation of sodium is used to monitor electrolyte balance.
Phosphorus	Measurement of inorganic phosphorus aids in the diagnosis and management of various disorders, including those involving the parathyroid gland, kidney, bone, and vitamin D metabolism.
PLACA	Lipoprotein-Associated Phospholipase A2 (Lp-PLA2) has been shown to be an independent inflammatory marker of cardiovascular risk and events. It is produced by macrophages in response to the presence of oxidized lipids and circulates primarily in association with low-density lipoprotein particles (LDL). Whereas hsCRP detects inflammation that is either part of atherosclerosis or some other systemic or localized process, Lp-PLA2 is much more specific for vascular inflammation and appears to be a marker of unstable atherosclerotic plaques. In the West of Scotland Coronary Prevention Study (WOSCOPS) (2), there was a two-fold risk of CHD in individuals in the highest quintile compared to the lowest quintile and in the Atherosclerosis Risks in Communities Study (ARIC) (3), there was almost a two-fold risk of ischemic stroke in individuals with an increased Lp-PLA2 level.
Plasma Hgb	This assay is useful for determining whether hemolysis is occurring <i>in vivo</i> . Plasma hemoglobin is increased with intravascular hemolysis, ABO incompatible transfusion, falciparum malaria, burns and march hemoglobinuria. Slight increases may occur with extravascular hemolysis, delayed transfusion reactions, and sickle cell anemia.
Pre-Albumin	Serum prealbumin levels are used as an index of subclinical or marginal protein-calorie malnutrition, as an indicator of response to total parenteral nutrition (TPN), a marker of nutritional status in premature infants, and as an index of liver function in hepatobiliary disease.
Rheumatoid Factor	Rheumatoid factor assay is one of the most frequently requested tests in the clinical investigation of patients with joint symptoms. Rheumatoid factor (RF) is usually an IgM autoantibody that reacts with the Fc portion of IgG. In the presence of pathogen-specific IgG antibodies, IgM RF can produce false-positive results in IgM assays. RF is not a screening test. The test performs poorly when applied to the general population.
Tbili	Total Bilirubin, Serum. The quantitative determination of bilirubin aids in the evaluation of liver disease, in the detection of hemolytic anemia, and in the evaluation of jaundice. Bilirubin, Cord Blood: Cord blood bilirubin may be a useful indicator of developing jaundice in newborns and a useful predictor of significant hyperbilirubinemia in the neonate.
Total Protein	Total protein measurements aid in the assessment of nutritional status (see Prealbumin) and may be useful in the diagnosis and management of a variety of diseases involving the liver, kidney, and bone marrow.
Transferrin	Transferrin functions as the principal plasma protein responsible for the transport of iron. Transferrin binds and transports iron and serves as a negative acute phase reactant (levels decrease during inflammatory processes). Transferrin levels in serum aid in the diagnosis of iron deficiency anemia, malnutrition, acute inflammation, infection, assessment of renal function, and red blood cell disorders. Serum transferrin concentration increases in iron deficiency anemia, pregnancy, and patients taking estrogens. Decreased transferrin levels are associated with chronic infections, malignancy, iron overload, hemolytic anemia, hemochromatosis, sickle cell disease, atransferrinemia, renal disease, and hepatic failure. Atransferrinemia is a rare congenital disorder. Patients with this disorder have very low levels of plasma transferrin. They also have iron overload and severe anemia that results from their inability to mobilize the body's iron stores.
Triglycerides	Triglyceride measurements aid in the diagnosis and management of patients with primary and secondary lipid disorders (e.g., diabetes mellitus, renal disease, liver obstruction, hypothyroidism).
Uric Acid	Uric acid measurements aid in the diagnosis and management of numerous renal and metabolic disorders, including renal failure, gout, leukemia, psoriasis, starvation or other wasting conditions, and of patients receiving cytotoxic drugs. It is recommended that patients being treated for gout maintain a uric acid level of less than 6 mg/dL.
Ammonia	The diagnostic utility of ammonia measurements is limited. The test is used mainly in the diagnosis of urea cycle defects and in the detection of Reye's syndrome. Levels of ammonia do not correlate well with CNS changes in end-stage liver disease.
B-OH Butyrate	This assay aids in the diagnosis and monitoring of therapy for diabetic ketoacidosis. This assay may also aid in the diagnosis of any patient presenting to the emergency room with documented hypoglycemia, acidosis, alcohol ingestion, or an unexplained increase in the anion gap and in the investigation of inborn errors of metabolism.
Ethanol	This assay is used to document prior consumption or administration of ethanol. Ethanol is the single most important abused substance in the U.S. Ethanol is found in beer, wine, and other liquors. Ethanol depresses cerebral functions similar to general anesthetics. Symptoms of ethanol abuse may include impaired thought, clouded judgment, and changed behavior. Blood ethanol levels correlate directly with the degree of intoxication.
Lactate	This assay aids in the evaluation of metabolic acidosis, regional or diffuse tissue hypoperfusion, hypoxia, shock, congestive heart failure, dehydration, complicated postoperative state, ketoacidosis or nonketotic acidosis in diabetes mellitus, patients with infections, inflammatory states, certain myopathies, acute leukemia and other neoplasia, enzyme defects, glycogen storage disease (type1), thiamine deficiency and hepatic failure.

Uric Acid Rasburicase	In patients being treated with Rasburicase (Elitek), "Uric Acid Rasburicase" should be ordered for 3 days following drug administration. This recommendation is based on the drug's half-life and the expected duration of activity. Rasburicase catalyzes the oxidation of uric acid into an inactive and soluble metabolite. It is used for prophylaxis and treatment of chemotherapy-induced or spontaneous acute tumor lysis syndrome. Rasburicase can cause spuriously low plasma uric acid levels if the specimen is transported or stored at room temperature (20-26°C or 68-78.8°F) or if processing is delayed. Uric Acid Rasburicase differs from Uric Acid in the collection and handling of the specimen. Transporting the specimen on ice and running it STAT helps to minimize (but not necessarily completely eliminate) the artifactual effects of rasburicase. Uric acid levels should be ordered in this way only when rasburicase has been administered to a patient.
Peritoneal Fluid	
Albumin	A serum ascites albumin gradient greater than or equal to 1.1 g/dL is consistent with portal hypertension from causes such as cirrhosis, congestive heart failure, or portal vein thrombosis. A low gradient (less than 1.1 g/dL) occurs in conditions such as peritoneal carcinomatosis, peritoneal tuberculosis, pancreatitis, serositis, and nephrotic syndrome. In the past a value of greater than or equal to 1.1 g/dL was interpreted as a transudate, and if less than 1.1 g/dL the fluid was interpreted as an exudate.
Amylase	Useful initially, in the classification of an effusion as an exudate or a transudate.
Tbili	Elevated body fluid bilirubin is suggestive of an exudative fluid. This testing should be performed in conjunction with other testing including serum bilirubin analysis, body fluid/serum protein ratio, body fluids/serum lactate dehydrogenase ratio, and serum lactate dehydrogenase. Determination of body fluid bilirubin concentration can aid in the distinction between a transudative and an exudative fluid. Bilirubin values tend to be higher in exudates than in transudates, although there is some overlap between groups. However, a ratio of body fluids to serum bilirubin has been reported to identify exudative body fluids with sensitivity, specifically, positive predictive accuracy, and absolute accuracy equivalent to that obtained using Light's criteria for an exudative pleural fluid (pleural/serum protein ratio greater than 0.5, pleural/serum lactate dehydrogenase ratio greater than 0.6, and serum lactate dehydrogenase greater than 200 U/L).
Cholesterol	Peritoneal fluid cholesterol determination can distinguish cirrhotic versus malignant ascites.
Creatinine	Measurement of creatinine is useful to differentiate between peritoneal fluid and urine. Elevated peritoneal fluid creatinine, in association with normal serum creatinine, suggests urinary bladder leakage or rupture.
Glucose	Glucose measurement in body fluid may be useful with other laboratory tests to evaluate effusions. Decreased concentrations are associated with bacterial infections, inflammation such as rheumatoid arthritis, and occasionally malignancy.
LD	Measuring LD in fluid aspirated from a pleural effusion (or pericardial effusion) can help in the distinction between exudates (actively secreted fluid, e.g. due to inflammation) and transudates (passively secreted fluid, due to a high hydrostatic pressure or a low oncotic pressure). The most reliable method for differentiating between transudates and exudates is the simultaneous analysis of fluid and serum for lactic dehydrogenase (LD) and total protein level.
Total Protein	Total Protein Interpretation: Measurement of total protein in body fluids other than blood, urine, or cerebrospinal fluid is usually done to differentiate an inflammatory fluid collection (exudate) from one that is non-inflammatory (transudate). In pericardial, peritoneal, pleural, and synovial fluids, 3 g/dL is usually taken as the cut-off value for differentiating transudates from exudates. Some authors use a lower cut-off of 2.5 g/dL. Some references suggest using a ratio of fluid to serum protein to differentiate transudate from exudate. Protein is just one of several markers that can be used for differentiating transudates from exudates. Low total protein is seen in patients with cirrhosis of the liver when ascites develops late in the disease. Patients with a low value, below 1.5 g/dL, are at greater risk of developing spontaneous bacterial peritonitis. Knowing that the concentration is low has some prognostic value, although it should not be a reason for beginning prophylactic antibiotic therapy.
Triglycerides	Peritoneal fluid triglyceride determination can distinguish cirrhotic versus malignant ascites
Amniotic Fluid	
Glucose	The following criteria is/are considered to be predictive of a positive amniotic fluid culture: 1) Positive gram stain for bacteria or 2) WBC count >30 cells/mm ³ and 3) Low amniotic fluid glucose less than 15 mg/dL. Laboratory studies on amniotic fluid should be performed on non-bloody fluid obtained by amniocentesis only.
Pleural Fluid	
Albumin	Serum albumin/pleural fluid albumin gradient of less than or equal to 1.2 g/dL is consistent with exudate.
Amylase	Useful initially, in the classification of an effusion as an exudate or a transudate.
Tbili	Pleural fluid bilirubin/serum bilirubin ratio of greater than or equal to 0.60 is consistent with exudate.
Cholesterol	Fluid cholesterol/serum cholesterol ratio greater than or equal to 0.3 or fluid cholesterol greater than 45 mg/dL is consistent with exudate.
Glucose	Pleural fluid glucose of less than 60 mg/dL or pleural fluid glucose/serum glucose ratio of less than 0.5 is abnormal. Abnormal pleural fluid glucose is seen in rheumatoid pleuritis and parapneumonic exudates. It can also be seen in malignancy, tuberculosis, SLE, and esophageal rupture.
LD	Pleural fluid LD/serum LD ratio of greater than or equal to 0.60 or pleural fluid LD greater than or equal to 2/3 rd s the upper limit of normal serum LD level is consistent with exudate. Pleural fluid LD/serum LD ratio of less than 0.60 or pleural fluid LD less than or equal to 2/3 rd s the upper limit of normal serum LD level is consistent with transudate.

Total Protein	Pleural fluid TP/serum TP ratio of greater than 0.5 or pleural fluid total protein level greater than 3.0 g/dL is consistent with exudate. Pleural fluid TP/serum TP ratio of less than or equal to 0.5 or pleural fluid total protein level less than or equal to 3.0 g/dL is consistent with transudate. Using total protein alone misclassifies exudates and transudates by approximately 30%. Sensitivity and specificity increases to 98% and 80%, respectively, when using both total protein and LD criteria for classifying exudates or transudates.
Triglycerides	Pleural fluid triglyceride levels greater than or equal to 110 mg/dL is indicative of chylous effusion. Pleural fluid triglyceride level less than 50 mg/dL is indicative of non-chylous effusion. Levels between 50-109 mg/dL are equivocal.
Dialysate Fluid	
Creatinine	Results of Creatinine, Glucose, and Urea Nitrogen are used by clinical staff in calculations to assess the adequacy of peritoneal dialysis.
Glucose	Results of Creatinine, Glucose, and Urea Nitrogen are used by clinical staff in calculations to assess the adequacy of peritoneal dialysis.
Urea Nitrogen	Results of Creatinine, Glucose, and Urea Nitrogen are used by clinical staff in calculations to assess the adequacy of peritoneal dialysis.
Pancreatic Fluid	
Amylase	Testing is used to determine whether a pancreatic cyst is likely to be benign or malignant. However these results cannot be used in isolation and should be used in conjunction with clinical information, imaging studies, and cytology.
Breast Milk	
Sodium	The composition of breast milk varies with the stage of lactation. It is usually not necessary to analyze breast milk electrolytes, however if a neonate is experiencing an electrolyte disorder, it may be appropriate to measure sodium or other electrolyte concentrations.
K	
Chloride	
Calcium	
Magnesium	
Stool	
K	Fecal osmolality is useful in cases of chronic diarrhea. It may be helpful to the physician to know whether the diarrhea is: 1. A secretory type caused by either an organism or an abnormality of water or electrolyte transport across the cell wall of the gut. 2. An osmotic type caused by malabsorption of non-electrolyte substances, most commonly carbohydrates or certain laxatives (e.g., magnesium). Fecal osmolality should be similar to serum osmolality. If the fecal osmolality is significantly lower than the serum or plasma osmolality (< 220 mOsmol/kg), factitious diarrhea (i.e., addition of water or liquid to stool by patient) should be suspected. If the fecal sample was not refrigerated immediately after collection, and if necessary frozen, the measured osmolality may be inappropriately elevated (> 330 mOsmol/kg). This change is due to bacterial metabolism which results in production of osmotically active substances. The Osmotic Gap is equal to the measured osmolality (mOsmol/kg) minus the calculated osmolality (in mOsmol/kg, equal to 2 times the fecal sodium plus fecal potassium). An Osmotic Gap > 125 mOsmol/kg with a fecal sodium < 60 mmol/L suggests an osmotic cause of the diarrhea. An Osmotic Gap < or = 50 mOsmol/kg with a fecal sodium > 90 mmol/L suggests a secretory cause of the diarrhea. The test results should be integrated into the clinical context for interpretation. Fecal chloride concentration is markedly elevated > 60 mmol/L in infants and > 100 mmol/L in adults associated with congenital and secondary chloridorrhea. Fecal chloride may be elevated (> 35 mmol/L) in phenolphthalein (or phenolphthalein plus magnesium hydroxide) induced diarrhea. Fecal chloride may be low (< 20 mmol/L) in sodium sulfate induced diarrhea.
Chloride	
Sodium	

Urine	
Amylase	The urinary amylase assay aids in the diagnosis of pancreatitis, pancreatic pseudocyst and macroamylasemia.
Calcium	24 hour urine calcium reflects intake, rates of intestinal calcium absorption, bone resorption and renal loss. Those processes relate to parathyroid hormone and vitamin D levels. Urinary calcium measurements are most useful for evaluation of patients with renal stones or a possible diagnosis of familial hypocalciuric hypercalcemia.
Chloride	Urinary chloride measurements aid in the differentiation of causes of metabolic alkalosis and help classify them as chloride responsive or unresponsive.
Creatinine	This assay aids in the diagnosis and management of renal diseases (when done as part of a Creatinine Clearance test), in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes.
Glucose	This assay aids in the evaluation of glucosuria and renal tubular defects. It is rarely needed in the management of diabetes mellitus.
Magnesium	Magnesium measurements aid in the diagnosis and management of hypomagnesemia (abnormally low plasma levels of magnesium) and hypermagnesemia (abnormally high plasma levels of magnesium). Urine magnesium is measured as part of the Stone Former Panel and is useful in assessing the likelihood of stone formation.
Micro Albumin	Microalbuminuria is an amount of albumin in the urine above normal (10 mg/L) but below that detected by dipsticks for urinary protein (greater than 30 mg/dL). Microalbuminuria has an important predictive value in determining diabetic patients at risk of developing nephropathy. Microalbuminuria may also be caused by poor metabolic regulation, physical exercise, newly diagnosed diabetes, hypertension, and non-diabetic renal or systemic disease.
Phosphorus	Inorganic phosphorus measurements aid in the diagnosis and management of various disorders, including parathyroid gland and kidney diseases, vitamin D imbalance and kidney stones.

Potassium	Urinary potassium measurements are used to monitor electrolyte balance in the diagnosis and treatment of disease conditions characterized by low or high blood potassium levels.
Sodium	Urinary sodium measurements aid in the evaluation of patients with acute oliguria (low urine output), hyponatremia, and volume depletion.
Total Protein	Urinary protein measurements aid in the diagnosis and management of primary diseases involving the kidney and diseases which may secondarily involve the kidney, such as collagen-vascular disease, multiple myeloma, amyloidosis, metal poisoning, diabetes mellitus, pre-eclampsia, or eclampsia.
Urea (UUN)	Urine urea nitrogen is performed as part of the Stone Former Workup to evaluate a patient for the likelihood of renal calculi formation.
Uric Acid	Uric acid measurements aid in the diagnosis and treatment of numerous renal and metabolic disorders, including renal failure, gout, leukemia, psoriasis, starvation, or other wasting conditions, and of patients receiving cytotoxic drugs. It is also part of the Stone Former Workup and is useful in assessing the likelihood of stone formation.
CSF	
CSF Glucose	Evaluation of meningitis, neoplastic involvement of meninges, and other neurological disorders.
CSF Lactate	Evaluation of meningitis, neoplastic involvement of meninges, and other neurological disorders.
CSF LD	Evaluation of meningitis, neoplastic involvement of meninges, and other neurological disorders.
CSF Protein	Evaluation of meningitis, neoplastic involvement of meninges, and other neurological disorders.
Therapeutic Drugs	
Acetaminophen	Acetaminophen is an analgesic, antipyretic drug lacking in significant anti-inflammatory activity. This assay is used to monitor the therapeutic drug level and evaluate the toxicity of acetaminophen. Serum concentration and half-life are the only way to assess the degree of intoxication in the early stages since other liver function studies (bilirubin and liver function enzymes) will not show clinically significant increases until after tissue damage has occurred, at which point therapy is ineffective.
Amikacin	Amikacin is a semisynthetic aminoglycoside antibiotic with a broad spectrum of activity against gram-negative bacteria. This assay is used to monitor the therapeutic drug level and evaluate the toxicity of amikacin.
Carbamazepine	Carbamazepine is an anticonvulsant drug used in the treatment of generalized and partial seizures. This assay aids in monitoring carbamazepine levels to ensure appropriate therapy.
Digoxin	This assay is used to monitor the therapeutic drug level and evaluate the toxicity of digoxin. Digoxin is a cardiac glycoside that is commonly prescribed to treat congestive heart failure.
Gentamicin	Gentamicin is an aminoglycoside antibiotic which exhibits high potency and a broad spectrum bacterial action primarily against gram-negative organisms. Gentamicin is associated with renal and ototoxicity upon extended use. Therapeutic monitoring is advantageous particularly in patients with diminished renal function.
Lithium	This assay is used to monitor the therapeutic drug level and evaluate the toxicity of lithium.
Phenobarbital	Phenobarbital is an anti-convulsant and sedative-hypnotic drug. It is used for the treatment of epilepsy, particularly for controlling focal motor or sensory seizures and grand mal seizures. It is frequently co-administered with phenytoin for the control of complex seizure disorders and with valproic acid for complex partial seizures. Monitoring the serum concentrations of phenobarbital has been shown to improve patient therapy by aiding the physician in adjusting their dosage. Phenobarbital has a narrow therapeutic index and wide inter-individual variability in the rate of metabolism and clearance. This assay is used to monitor the therapeutic drug level and evaluate the toxicity of phenobarbital.
Phenytoin	Phenytoin is an anticonvulsant drug. It is occasionally used as an antiarrhythmic. In the treatment of epilepsy, phenytoin is indicated for grand mal epilepsy, cortical focal seizures and temporal lobe epilepsy. Phenytoin has a narrow therapeutic index and a wide interindividual variability in the rate of metabolism and clearance necessitating the determination of blood levels for patients receiving therapy. This assay is used to monitor the therapeutic drug level and evaluate the toxicity of phenytoin.
Salicylate	The salicylate assay aids in the diagnosis and treatment of salicylate overdose and in monitoring salicylate levels to insure appropriate therapy. Salicylates are a group of compounds used as analgesics, antipyretics and anti-inflammatory agents. Acetylsalicylic acid (aspirin) is the most commonly used salicylate. Salicylates are readily available over-the-counter and most salicylate therapy is the result of patient self-medication. For this reason, salicylates are often seen in overdose cases. Salicylate poisoning is seldom fatal, but causes side effects ranging from nausea, vomiting and tinnitus to fever, lethargy and coma. Prompt recognition and management are necessary to avoid serious consequences.
Theophylline	Theophylline is a naturally occurring compound with bronchodilator effects. It is used in the treatment of bronchospasm associated with bronchial asthma, chronic bronchitis and pulmonary emphysema. This assay is used to monitor the therapeutic drug level and evaluate the toxicity of theophylline. Due to the drug's narrow therapeutic range and the wide personal variation in the rate of metabolism and clearance, essentially every patient taking theophylline should have their serum levels monitored.
Tobramycin	Tobramycin is an aminoglycoside antibiotic. Tobramycin has a narrow therapeutic index which makes its use hazardous, especially in patients with impaired renal function. Accurate monitoring of the serum level in such patients is mandatory. This assay aids in the diagnosis and treatment of tobramycin overdose and in monitoring levels of tobramycin to ensure appropriate therapy.
Valproic Acid	Valproic acid is a broad-spectrum anticonvulsant drug used alone to treat simple and complex absence seizures and in combination with other anticonvulsant drugs for control of generalized seizures that include absence seizures. This assay is used to monitor the therapeutic drug level and evaluate the toxicity of valproic acid.

Vancomycin	Vancomycin is a glycopeptide antibiotic which is bactericidal against many gram-positive and some gram-negative cocci. It is used to treat severe staphylococcal infections in patients who cannot receive or who have failed to respond to the penicillins and cephalosporins. Vancomycin should be used with care due to its potential nephrotoxic and ototoxic effects. This assay aids in the diagnosis and treatment of vancomycin overdose and in monitoring levels of vancomycin to ensure appropriate therapy.
Methotrexate	This assay is used to monitor the therapeutic drug level and evaluate the toxicity of methotrexate. Methotrexate is an antineoplastic agent. It is used in the treatment of malignancies with rapid cell proliferation such as acute lymphoblastic leukemia, choriocarcinoma, trophoblastic tumors in women, and carcinomas of the breast, tongue, pharynx, and testis. It is also indicated in the treatment of severe psoriasis and adult rheumatoid arthritis.

DAU

Amphetamine/ Methamphetamine	Amphetamine and methamphetamine are central nervous system stimulants. They are usually administered orally or by intravenous injection. These drugs are prescribed for treatment of obesity, narcolepsy, hypotension, attention deficit disorder, and hyperactivity disorder. Because of their stimulant effects, the drugs are commonly sold illicitly and abused. This assay is used to aid in the diagnosis and treatment of amphetamine and/or methamphetamine use or abuse. This assay provides a preliminary analytical result. Positive results are confirmed by chromatographic methods.
Barbiturate	This assay is used to aid in the diagnosis and treatment of barbiturate use or abuse. This assay provides a preliminary analytical result. Positive results are confirmed by chromatographic methods. Barbiturates are central nervous system depressants. They are used therapeutically as sedatives, hypnotics, and anticonvulsants. The legal availability of barbiturates has declined but they remain frequently abused sedatives or hypnotic drugs. The most commonly abused barbiturates are the short-acting compounds such as secobarbital, pentobarbital and amobarbital. Tolerance to these drugs can develop from chronic use and death may occur from overdose.
Benzodiazepine	This assay is used to aid in the diagnosis and treatment of benzodiazepine use or abuse. This assay provides a preliminary analytical result. Positive results are confirmed by chromatographic methods. Benzodiazepines are central nervous system depressants composed of over 20 compounds. They are used clinically as sedatives, hypnotics, anxiolytics, muscle relaxants, antiepileptics, and in the treatment of alcohol withdrawal. Benzodiazepines are extensively metabolized with half-lives ranging from 1 to 100 hours. Urine levels vary due to each patient's metabolic and excretion rates.
Cannabinoid	This assay is used to aid in the diagnosis and treatment of cannabinoid use or abuse. This assay provides a preliminary analytical result. Positive results are confirmed by chromatographic methods. The primary psychoactive component of marijuana is delta-9-tetrahydrocannabinol (THC). The THC concentration determines the potency of the marijuana. THC primarily affects the cardiovascular and central nervous systems. There is no reliable method for predicting the degree of impairment from cannabinoid concentrations measured in urine at this time.
Cocaine	Cocaine is a frequently abused drug. The drug is administered by nasal insufflation, intravenous injection or in the free base form as smoke inhalation. The urinary elimination of cocaine and its metabolite begins within 20 minutes of its intranasal administration. This assay is used to aid in the diagnosis and treatment of cocaine use or abuse. This assay provides a preliminary analytical result. Positive results are confirmed by chromatographic methods.
Methadone	Methadone is a synthetic narcotic analgesic that is similar to morphine. Methadone is commonly prescribed as the drug of choice in the maintenance treatment of heroin addicts. Patients on methadone therapy are routinely screened for methadone as a measure of compliance. This assay is used to aid in the diagnosis and treatment of methadone use or abuse. This assay provides a preliminary analytical result. Positive results are confirmed by chromatographic methods.
Opiate	Opiates act on several sites of the central nervous system. Their use results in analgesia, drowsiness, mood changes and some mental clouding. This assay is used to aid in the diagnosis and treatment of opiate use or abuse. This assay provides a preliminary analytical result. Positive results are confirmed by chromatographic methods.
Phencyclidine	Phencyclidine (PCP) is a drug of abuse. PCP acts as a stimulant, depressant, hallucinogenic and analgesic. It can be self-administered by smoking, nasal insufflation, intravenous injection or by oral ingestion. This assay is used to aid in the diagnosis and treatment of phencyclidine use or abuse. It provides a preliminary analytical result. Positive results are confirmed by chromatographic methods.

Assay	Method	Sample Type	Reagent On-board Stability	Reagent Prep	Calibrator	Calibration Frequency	Reagent Storage
General Chemistry							
Albumin BCG	Bromocresol Green	Serum or Li Heparin	42 days	Liquid Ready to Use	Multiconstituent	41 days	RT
Alkaline Phosphatase	Para-nitrophenol Phosphate Enzymatic	Serum or Li Heparin	8 days	Liquid Ready to Use	Internal Water Blank Cal Factor	8 days	2°-8°C
ALT	NADH	Serum or Li Heparin	30 days	Liquid Ready to Use	Internal Water Blank Cal Factor	27 days	2°-8°C
Ammonia	Glutamate Dehydrogenase	EDTA Plasma	15 days	Liquid Ready to Use	Calibrator supplied with reagent	24 hours	2°-8°C
Amylase	CNPG3 Substrate	Serum, Li Heparin, Urine	19 days	Liquid Ready to Use	Internal Water Blank Cal Factor	19 days	2°-8°C
AST	NADH	Serum or Li Heparin	30 days	Liquid Ready to Use	Internal Water Blank Cal Factor	30 days	2°-8°C
B2-Microglobulin	Immunoturbidimetric	Serum or EDTA Plasma	30 days	Ready to Use Invert to mix	B2 Microglobulin Calibrator	31 days	2°-8°C
B-OH Butyrate	Stanbio	Serum or Li Heparin	30 days	Ready to Use	BHBT Calibrator	14 days	2°-8°C
Bilirubin, Direct	Diazo Reaction	Serum or Li Heparin	28 days	Liquid Ready to Use	Bilirubin Calibrator	14 days	2°-8°C
Bilirubin, Total	Diazonium Salt	Serum or Li Heparin	21 days	Liquid Ready to Use	Bilirubin Calibrator	14 days	2°-8°C
Calcium	Arsenazo III Dye	Serum, Li Heparin, Urine	30 days	Liquid Ready to Use	Multiconstituent	30 days	RT
Chloride	indirect potentiometry	Serum, Li Heparin, Urine	ICT Dil 14 days	Liquid Ready to Use	ICT Serum and Urine Calibrator	24 hours	RT
Potassium	indirect potentiometry	Serum, Li Heparin, Urine	ICT Dil 14 days	Liquid Ready to Use	ICT Serum and Urine Calibrator	24 hours	RT
Sodium	indirect potentiometry	Serum, Li Heparin, Urine	ICT Dil 14 days	Liquid Ready to Use	ICT Serum and Urine Calibrator	24 hours	RT
Cholesterol	Enzymatic	Serum or Li Heparin	30 days	Liquid Ready to Use	Multiconstituent Calibrator	30 days	2°-8°C
CK	NAC (N-acetyl-L-cysteine)	Serum or Li Heparin	30 days	Liquid Ready to Use	Internal Water Blank Cal Factor	30 days	2°-8°C
CO2	PEP Carboxylase	Serum or Li Heparin	14 days	Liquid Ready to Use	Carbon Dioxide Calibrator	14 days	2°-8°C
Creatinine	Kinetic Alkaline Picrate	Serum, Li Heparin, Urine	5 days	Liquid Ready to Use	Multiconstituent Calibrator	5 days	RT
Creatinine Ez	Enzymatic	Serum, Li Heparin	30 days	Liquid Ready to Use	Clin Chem Cal	7 days	2°-8°C
GGT	L-Gamma-glutamyl-3-carboxy-4-nitroanilide Substrate	Serum or Li Heparin	27 days	Liquid Ready to Use	Internal Water Blank Cal Factor	27 days	2°-8°C
Glucose	Hexokinase/G-6-PDH	Serum, Li Heparin, Urine, CSF	30 days	Liquid Ready to Use	Multiconstituent Calibrator	30 days	2°-8°C
HDL, Ultra	Accelerator Selective Detergent	Serum or Li Heparin	28 days	Liquid Ready to Use	Lipid Multiconstituent Calibrator	28 days	2°-8°C
Iron	Ferene	Serum or Li Heparin	60 days	Liquid Ready to Use	Multiconstituent Calibrator	14 days	2°-8°C
Lactic Acid	Lactic Acid to Pyruvate	sodium fluoride plasma	30 days	Liquid Ready to Use	Multiconstituent Calibrator	30 days	2°-8°C
LD	Lactate to Pyruvate	Serum or Li Heparin	30 days	Liquid Ready to Use	Internal Water Blank Cal Factor	30 days	2°-8°C
LDL Direct	Measured, Liquid Selective Detergent	Serum or Li Heparin	28 days	Liquid Ready to Use	Lipid Multiconstituent Calibrator	28 days	2°-8°C

Lipase	Quinone Dye	Serum or Li Heparin	11 days	1. Pour contents of R1 into R1A container and mix by gentle inversion until dissolved. 2. Pour contents of R1A back into R1 cartridge and mix again by gentle inversion. R2 reagent is liquid ready to use	Lipase Calibrator	11 days	2°-8°C
Magnesium	Enzymatic	Serum, Li Heparin, Urine	30 days	Liquid Ready to Use	Multiconstituent Calibrator	30 days	2°-8°C
Phosphorus	Phosphomolybdate	Serum, Li Heparin, Urine	65 days	Liquid Ready to Use	Multiconstituent Calibrator	41 days	RT
Total Protein	Biuret	Serum or Li Heparin	23 days	Liquid Ready to Use	Multiconstituent Calibrator	23 days	RT
Triglyceride	Glycerol Phosphate Oxidase	Serum or Li Heparin	42 days	Liquid Ready to Use	Multiconstituent Calibrator	41 days	2°-8°C
Urea Nitrogen	Urease	Serum, Li Heparin, Urine	25 days	Liquid Ready to Use	Multiconstituent Calibrator	7 days/ required with each new cartridge	2°-8°C
Uric Acid	Uricase	Serum, Li Heparin, Urine	60 days	Liquid Ready to Use	Multiconstituent Calibrator	60 days	2°-8°C
Urine/CSF Protein	Benzethonium chloride	Urine, CSF	41 days	Liquid Ready to Use	Urine/CSF Protein Calibrator	41 days	RT

Special Proteins

A1-Antitrypsin	Immunoturbidimetric	Serum or Li Heparin	35 days	Ready to Use	Proteins	30 days	2°-8°C
ASO	Immunoturbidimetric	Serum	35 days	Ready to Use Invert to Mix	ASO Calibrator	35 days	2°-8°C
Beta 2- Microglobulin	Immunoturbidimetric	Serum	30 days	Ready to Use Invert to Mix	Proteins	31 days	2°-8°C
Ceruloplasmin	Immunoturbidimetric	Serum or Li Heparin	60 days	Liquid Ready to Use	Multigent Plasmaproteins Cal	10 days	2°-8°C
Complement 3	Immunoturbidimetric	Serum or Li Heparin	57 days	Liquid Ready to Use	Specific Proteins Multiconstituent Calibrator	57 days	2°-8°C
Complement 4	Immunoturbidimetric	Serum or Li Heparin	57 days	Liquid Ready to Use	Specific Proteins Multiconstituent Calibrator	57 days	2°-8°C
CRP Vario (Wide Range)	Immunoturbidimetric	Serum or Li Heparin	60 days	Liquid Ready to Use	Multigent CRP Calibrator WR	15 days	2°-8°C
CRP Vario (Ultrasensitive)	Immunoturbidimetric	Serum or Li Heparin	60 days	Liquid Ready to Use	Multigent CRP Calibrator HS	15 days	2°-8°C
Haptoglobin	Immunoturbidimetric	Serum or Li Heparin	57 days	Liquid Ready to Use	Specific Proteins Multiconstituent Calibrator	57 days	2°-8°C
Immunoglobulin A	Immunoturbidimetric	Serum or Li Heparin	28 days	Liquid Ready to Use	Specific Proteins Multiconstituent Calibrator	25 days	2°-8°C
Immunoglobulin E	Immunoturbidimetric	Serum or Li Heparin	20 days	Ready to Use Invert to Mix	Quantia IgE Calibrator	30 days	2°-8°C
Immunoglobulin G	Immunoturbidimetric	Serum or Li Heparin	23 days	Liquid Ready to Use	Specific Proteins Multiconstituent Calibrator	23 days	2°-8°C
Immunoglobulin M	Immunoturbidimetric	Serum or Li Heparin	57 days	Liquid Ready to Use	Specific Proteins Multiconstituent Calibrator	57 days	2°-8°C
Lipoprotein A1	Immunoturbidimetric	Serum, Li Heparin	35 days	Liquid Ready to Use	Quantia Lp(a) Calibrator	30 days	2°-8°C

Microalbumin	Immunoturbidimetric	Urine	28 days	Liquid Ready to Use	Microalbumin Calibrators	28 days	2°-8°C
Prealbumin	Immunoturbidimetric	Serum	57 days	Liquid Ready to Use	Prealbumin Calibrator	57 days	2°-8°C
Rheumatoid Factor (RF)	Immunoturbidimetric	Serum	30 days	Ready to Use Invert to Mix	Rheumatoid Factor Calibrator	60 days	2°-8°C
Transferrin	Immunoturbidimetric	Serum or Li Heparin	57 days	Liquid Ready to Use	Specific Proteins Multiconstituent Calibrator	57 days	2°-8°C

Therapeutic Drugs

Amikacin	PETINIA	Non SST Serum	54 days	Ready to Use Invert to Mix	TDM Multiconstituent Calibrator	54 days	2°-8°C
Carbamazepine	PETINIA	Non SST Serum	45 days	Ready to Use Invert to Mix	TDM Multiconstituent Calibrator	7 days	2°-8°C
Digoxin	PETINIA	Non SST Serum	60 days	Ready to Use Invert to Mix	TDM Multiconstituent Calibrator	60 days	2°-8°C
Gentamicin	PETINIA	Non SST Serum	55 days	Ready to Use Invert to Mix	TDM Multiconstituent Calibrator	28 days	2°-8°C
Lithium	Colorimetric	Non SST Serum	18 days	Liquid Ready to Use	Clin Chem Cal	5 days	2°-8°C
Phenobarbital	PETINIA	Non SST Serum	40 days	Ready to Use Invert to Mix	TDM Multiconstituent Calibrator	14 days	2°-8°C
Phenytoin	EIA	Non SST Serum	40 days	Liquid Ready to Use	TDM Multiconstituent Calibrator	7 days	2°-8°C
Theophylline	EIA	Non SST Serum	40 days	Liquid Ready to Use	TDM Multiconstituent Calibrator	7 days	2°-8°C
Tobramycin	PETINIA	Non SST Serum	32 days	Ready to Use Invert to Mix	Tobramycin Calibrator	7 days	2°-8°C
Valproic Acid	PETINIA	Non SST Serum	54 days	Ready to Use Invert to Mix	TDM Multiconstituent Calibrator	27 days	2°-8°C
Vancomycin	PETINIA	Non SST Serum	45 days	Ready to Use Invert to Mix	TDM Multiconstituent Calibrator	45 days	2°-8°C

Drugs of Abuse/Toxicology

Acetaminophen	Enzymatic/Colorimetric	Non SST Serum/Urine	8 days	Liquid Ready to Use R1 reagent one bottle, R2 reagent 2 bottles poured into 20 mL wedges	Acetaminophen Calibrator	24 hours	2°-8°C
Amphet/Meth	EIA	Urine	56 days	Liquid Ready to Use	DOA Neg Cal and DOA Cal #1	13 days	2°-8°C
Barbiturates	EIA	Urine	56 days	Liquid Ready to Use	DOA Neg Cal and DOA Cal #2	13 days	2°-8°C
Benzodiazepines (Urine)	EIA	Urine	56 days	Liquid Ready to Use	DOA Neg Cal and DOA Cal #1, 2, 3, 4	7 days	2°-8°C
Cannabinoids	EIA	Urine	30 days	Liquid Ready to Use	DOA Neg Cal and THC 50	13 days	2°-8°C
Cocaine	EIA	Urine	56 days	Liquid Ready to Use	DOA Neg Cal and DOA Cal #2	13 days	2°-8°C
Methadone	EIA	Urine	56 days	Liquid Ready to Use	DOA Neg Cal and DOA Cal #2	10 days	2°-8°C
Opiates	EIA	Urine	56 days	Liquid Ready to Use	DOA Neg Cal and Opiate 300	13 days	2°-8°C
Phencyclidine	EIA	Urine	56 days	Liquid Ready to Use	DOA Neg Cal and DOA Cal #3	14 days	2°-8°C
Salicylate	Enzymatic/Colorimetric	Non SST Serum/Urine	43 days	Liquid Ready to Use	Salicylate Calibrator	43 days	2°-8°C

Bulk Solutions

	<u>Purpose</u>	<u>Onboard Stability</u>	<u>Preparation</u>	<u>Expiration</u>
ICT Reference	ICT Analysis	Expiration date on bottle	Ready to use	
Alkaline Wash	Wash cuvettes	Expiration date on bottle	Ready to use	
Acid Wash	Wash cuvettes	Expiration date on bottle	Ready to use	
0.5% Acid Wash	Wash sample probe	1 days	5 mL Acid Wash 995 mL reagent grade water	Stability listed on Acid Wash Bulk
10% Detergent B	Wash Reagent Probes	14 days	50 mL Detergent B 450 mL reagent grade water	14 days
Water Bath Additive	Daily Maintenance	N/A	Ready to use- do not top off bottle	
Detergent A	Wash Sample Probe	Expiration date on bottle	Ready to use	
Saline	Sample Dilution	Expiration date on bottle	Ready to use	
ICT Cleaning Fluid	Daily Maintenance	N/A	12 mL ICT Cleaning fluid Added to the ICT Lyophilized bottle	14 days at 2° to 8°C

Calibrator	Abbott Name	Preparation	Open Stability (Days)	Storage Temp after opening	Number of Levels	Requires Set Point Change with New LOT	Assays Calibrated
ICT Serum Cal	ICT Low/High	Ready to Use	7	2-8° C	2	N	CI-C K-C Na-C
ICT Urine Cal	ICTU Low/High	Ready to Use	7	2-8° C	2	N	CI-CU K-CU Na-CU
Multiconstituent Cal	MCC	Ready to Use	7	2-8° C	2	Y	AlbG CaC CaCU Chol CreaC CreaCU Glu GluCU Fe-PI LactA MAG MAGU Phos PhosU TP Trig UA Urea UreaU
BHBT Cal	BHOB	Ready to Use	Until Exp Date	2-8° C	1	N	BHBT
Bilirubin Cal	Bil	Ready to Use	7	2-8° C	2	Y	BiID BilIT
CO2 Cal	CO2C	Ready to Use	30	2-8° C	2	Y	CO2C
Lipid Multiconstituent Cal	LIPIDMCC	ADD 1mL H2O	7	2-8° C	1	Y	UHDl DLDL
Urine/ CSF Protein Cal	Upro	Ready to Use	Until Exp Date	2-8° C	5	N	Upro
Ammonia Cal (included with Reagent)	Amm	Ready to Use	120	2-8° C	1	N	Amm
Lipase Cal	Lipase	ADD 3mL H2O	28	2-8° C	1	Y	Lip
Specific Protein Cal	SP	Ready to Use	30	2-8° C	5	Y	C3 C4 Hapt IgA IgG IgM TRF
Microalbumin Cal	uAlbCal	Ready to Use	6 months	2-8° C	5	N	uAlb
Prealbumin Cal	PreAlb	Ready to Use	30	2-8° C	5	N	Palb
CRP Cal Set *	32CRP	Ready to Use	90	2-8° C	5	N	CRP32 * See below for details on CRP Calibration
CRP Cal Set *	16CRP	Ready to Use	90	2-8° C	5	N	CRP16 * See below for details on CRP Calibration
Quantia Protein Cal	PROT	Ready to Use	Until Exp Date	2-8° C	5	Y	A1AT
Plasma Protein Cal	PPCS3	Ready to Use	30	2-8° C	1	Y	Cerul
Quantia IgE	IgE	Ready to Use	Until Exp Date	2-8° C	5	N	IgE
Rheumatoid Factor Cal	RF	Ready to Use	90	2-8° C	5	N	RF
Quantia ASO Cal	ASO	ADD 1mL H2O	15	2-8° C	1	N	ASO
Ethanol Cal	EtohM	Ready to Use	Until Exp Date	2-8° C	2	N	ETOH
B2 Microglobulin Cal	B2M	Ready to Use	Until Exp Date	2-8° C	1	Y	B2M
Quantia Lp(a) Cal	Lpa	ADD 1 mL H2O	15	2-8° C	5	Y	Lpa
Internal Water							AlkP ALT Amy AmyU AST CK GGT LD
TDM Multiconstituent Cal	TDM MCC	Ready to Use	60	2-8° C	6	Y	AMIK CARB DIG GENT PHENO PHENY THEO VALP VANCO
Acetaminophen Cal	ACET	Ready to Use	Until Exp Date	2-8° C	1	N	ACETA
Salicylate Cal		Ready to Use	Until Exp Date	2-8° C	1	N	SALIC
Tobramycin Cal	TOBRA	Ready to Use	Until Exp Date	2-8° C	6	N	TOBRA
Clin Chem Cal	CCC-S	ADD 3mL H2O***	2 days #	2-8° C	1	Y	LITHI CRENZ #Can be frozen for 14 days *** Wait 30 Minutes
DOA Neg	DOAQL1	Ready to Use	Until Exp Date	2-8° C	1	N	Amphet Barb Benzo Cannab Cocaine Methadone Opiate PCP
DOA MC 1	See below	Ready to Use	Until Exp Date	2-8° C	1	N	Amphet Benzo
DOA MC 2	See below	Ready to Use	Until Exp Date	2-8° C	1	N	Barb Benzo Cocaine Methadone
DOA MC 3	See below	Ready to Use	Until Exp Date	2-8° C	1	N	PCP Benzo
DOA MC 4	See below	Ready to Use	Until Exp Date	2-8° C	1	N	Benzo
THC 50	See below	Ready to Use	Until Exp Date	2-8° C	1	N	Cannab
Opiate 300	See below	Ready to Use	Until Exp Date	2-8° C	1	N	Opiate

CRP Details for Calibration

	CAL 1	CAL 2	CAL 3	CAL 4	CAL 5	CAL 6
CRP 16	2.5 Bright Yellow	5.0 White	10 Light Yellow	20 Green	80 Pink	160 Red
CRP 32	5 White	10 Light Yellow	20 Green	40 Blue	160 Red	320 Brown

Drugs of Abuse Details

Analyte	Calibrators		Abbott Name Negative	Abbott Name for Cut off Cal
Amphetamines	DOA Neg	DOA Cal #1	DOAQL1	DOAQL 2
Barbituates	DOA Neg	DOA Cal #2	DOAQ 1	DOAQ 2
Benzodiazepine	DOA Neg	DOA Cal #1,2,3,4	DOASQ1	DOASQ2,3,4,5
Cannabinoids	DOA Neg	THC 50	THCQ1	THCQ2
Cocaine	DOA Neg	DOA Cal #2	DOAQ1	DOAQ2
Opiates	DOA Neg	Opiate 300	OPQ1	OPQ2
Methadone	DOA Neg	DOA Cal # 2	DOAQ1	DOAQ2
PCP	DOA Neg	DOA Cal #3	DOAQ1	DOAQ3

	Age	Reference Range		Units	
ALBUMIN	0 - 6 DAYS	2.5	3.4	g/dL	
	7 DAYS - 5 YEARS	3.9	5.0		
	6 - 18 YEARS	4.0	5.3		
	19 - 60 YEARS	3.5	5.1		
	> 60 YEARS	3.5	4.9		
ACE		12	60	U/L	
ALT	Male	0 - < 1 YEAR	5	33	U/L
		1 - 12 YEARS	9	25	
		13 - 18 YEARS	9	24	
		19-Adult	9	47	
	Female	0 - < 1 YEAR	5	33	
		1 - 12 YEARS	9	25	
		13 - 18 YEARS	8	22	
		19-Adult	8	37	
ALP	Male	0 - 14 DAYS	90	273	U/L
		15 DAYS - < 1 YEAR	134	518	
		1 - 9 YEARS	156	369	
		10 - 12 YEARS	141	460	
		13 - 14 YEARS	127	517	
		15 - 16 YEARS	89	365	
		17 - 18 YEARS	59	164	
		> 19 YEARS	33	120	
	Female	0 - 14 DAYS	90	273	
		15 DAYS - < 1 YEAR	134	518	
		1 - 9 YEARS	156	369	
		10 - 12 YEARS	141	460	
		13 - 14 YEARS	62	280	
		15 - 16 YEARS	54	128	
		17 - 18 YEARS	48	95	
		19 - 60 YEARS	33	120	
		> 60 YEARS	37	135	
		ALPHA-1-ANTITRYPSIN		90	
AMMONIA		15	50	umol/L	
AMYLASE			< 100	U/L	
ASO	0-4 YEARS	<100		IU/mL	
	5-17 YEARS	<250			
	ADULT	<200			
AST (w/o P5P)		< 35		U/L	
BETA-2_MICROGLOBULIN		<2.5		mg/L	
BETA-OH BUTYRATE		0.02	0.27	mmol/L	
BILIRUBIN, DIRECT	0 - 27 DAYS	0.0	1.0	mg/dL	
	28 DAYS - ADULT	0.0	0.4		

BILIRUBIN, TOTAL	< 1 DAY	0.3	7.9	mg/dL
	1 DAY	0.3	12.0	
	2 - 13 DAYS	0.3	15.0	
	14 - 20 DAYS	0.3	8.0	
	21 - 27 DAYS	0.3	6.0	
	28 DAYS - ADULT	0.3	1.2	
C3		82	193	mg/dL
C4		10	43	mg/dL
CALCIUM	0 - < 1 YEAR	8.5	11.0	mg/dL
	1 - 18 YEARS	9.2	10.7	
	19 YEARS - ADULT	8.5	10.5	
CERULOPLASMIN		17	40	mg/dL
CHOLESTEROL	DESIRABLE		< 170	mg/dL
	BORDERLINE	170	199	
	HIGH		>= 200	
	DESIRABLE		< 200	
	BORDERLINE	200	239	
	HIGH		>= 240	
CK, TOTAL				U/L
<i>Male</i>		40	230	
<i>Female</i>		30	150	
CHLORIDE		98	111	mmol/L
CARBON DIOXIDE	< 1 YEAR	20	28	mmol/L
	1 - 13 YRS	17	28	
	> 13 YRS	20	29	
Creatinine	<i>Male</i> 0 - 14 DAYS	0.42	1.05	mg/dL
	15 DAYS - < 1 YEAR	0.32	0.53	
	1 - 3 YEARS	0.38	0.54	
	4 - 6 YEARS	0.44	0.64	
	7 - 11 YEARS	0.52	0.69	
	12 - 14 YEARS	0.57	0.8	
	15 - 16 YEARS	0.66	1.04	
	17 - 18 YEARS	0.69	1.1	
	19 - ADULT	0.60	1.30	
	<i>Female</i> 0 - 14 DAYS	0.42	1.05	
	15 DAYS - < 1 YEAR	0.32	0.53	
	1 - 3 YEARS	0.38	0.54	
	4 - 6 YEARS	0.44	0.64	
	7 - 11 YEARS	0.52	0.69	
	12 - 14 YEARS	0.57	0.8	
	15 - 16 YEARS	0.59	0.86	
	17 - 18 YEARS	0.6	0.88	
	19 - ADULT	0.50	1.10	

CRP		<8.0		mg/L
ETOH		0	10	mg/dL
GGT <i>MALE & FEMALE</i>	0 - 14 DAYS	23	219	U/L
	15 DAYS - < 1 YEAR	8	127	
	1 - 10 YEARS	6	16	
	11 - 18 YEARS	7	21	
	<i>Male</i> ADULT MALE	13	60	
	<i>Female</i> ADULT FEMALE	8	39	
GLUCOSE	FASTING	60	99	mg/dL
	RANDOM	60	139	
GLUCOSE, CSF		50	80	mg/dL
HS CRP	LOW RISK		< 1.0	mg/L
	AVERAGE RISK	1.0	3.0	
	HIGHER RISK	3.1	10.0	
	ASSOC. WITH INFECTION OR INFLAMMATION		> 10.0	
HAPTOGLOBIN		40	250	mg/dL
HDLC	0 - 19 YEARS	35	90	mg/dL
	> 19 YEARS	40	90	
IGA	< 1 YEAR	0	30	mg/dL
	1 - 2 YEARS	0	90	
	3 - 5 YEARS	30	150	
	6 - 13 YEARS	50	220	
	14 - 18 YEARS	50	209	
	Adult	70	365	
IGE			< 100	IU/mL
IGG	0 - 14 DAYS	320	1400	mg/dL
	15 DAYS - 11 MO	110	700	
	1 - 3 YEARS	320	1150	
	4 - 9 YEARS	540	1360	
	10 - 18 YEARS	660	1530	
	ADULT	550	1650	
IGM	0 - 14 DAYS	10	40	mg/dL
	15 DAYS - 12 WEEKS	10	70	
	13 WEEKS - 51 WKS	20	90	
	1 - 18 YEARS	40	150	
	ADULT	40	293	
IRON	<i>Male</i>	65	175	ug/dL
	<i>Female</i>	50	170	
POTASSIUM	0 - 1 DAY	5.0	7.5	mmol/L
	2 DAYS - 3 MONTHS	4.0	6.0	
	4 MONTHS - ADULT	3.5	5.2	
LACTATE		0.5	2.2	mmol/L

LD	Male	0 - 30 DAYS	125	735	U/L		
		31 DAYS - 12 MONTHS	170	450			
		13 MONTHS - 3 YEARS	155	345			
		4 - 6 YEARS	155	345			
		7 - 9 YEARS	145	300			
		10 - 12 YEARS	120	325			
		13 - 15 YEARS	120	290			
		16 - 18 YEARS	105	235			
		19 - ADULT	100	240			
	Female	0 - 30 DAYS	145	765			
		31 DAYS - 12 MONTHS	190	420			
		13 MONTHS - 3 YEARS	165	395			
		4 - 6 YEARS	135	345			
		7 - 9 YEARS	140	280			
		10 - 12 YEARS	120	260			
		13 - 15 YEARS	100	275			
		16 - 18 YEARS	105	230			
		19 - ADULT	100	240			
		LDL, DIRECT		50		129	mg/dL
		LP(a)		< = 30		mg/dL	
LIPASE		7	60	U/L			
MAGNESIUM		1.7	2.5	mg/dL			
SODIUM	0 - 1 DAY	126	166	mmol/L			
	2 - 29 DAYS	134	144				
	30 DAYS - 1YR	139	146				
	2 - 12 YRS	138	145				
	13 YRS - ADULT	135	145				
PHOSPHORUS	0 - 14 DAYS	5.6	10.5	mg/dL			
	15 DAYS - 1 YEAR	4.8	8.4				
	1 - 4 YEARS	4.3	6.8				
	5 - 12 YEARS	4.1	5.9				
	13 - 15 YEARS	3.2	6.2				
	16 - 18 YEARS	2.9	5.0				
	19 - ADULT	2.3	4.4				
PLACA		225		nmol/m/ in/mL			
PREALBUMIN		18	44	mg/dL			
PROTEIN, TOTAL	0 - 14 DAYS	5.3	8.3	g/dL			
	15 DAYS - < 1 YEAR	4.4	7.1				
	1 - 5 YEARS	6.1	7.5				
	6 - 8 YEARS	6.4	7.7				
	9 - 18 YEARS	6.5	8.1				
	19 - 60 YEARS	6.4	8.3				
	> 60 YEARS	6.2	8.1				
RHEUMATOID FACTOR		< 15		IU/mL			

TIBC		250	425	
	Male			
	Female			
TRANSFERRIN		163	382	mg/dL
% TRANSFERRIN SAT		15	50	%
TRIGLYCERIDES			< 150	mg/dL
UREA (BUN)	0- 14 DAYS	3	23	mg/dL
	15 DAYS - < 1 YEAR	3	17	
	1 - 10 YEARS	9	22	
	11 - 18 YEARS	7	21	
	19 - ADULT	7	25	
URIC ACID				mg/dL
	Male	0 - 13 YRS	1.5	7.6
		14 YRS - ADULT	3.5	7.2
	Female	0- ADULT	2.6	6.0

24 hour URINE

AMYLASE, URINE	2	17	U/Hour
CALCIUM, URINE	100	300	mg/Collection
CALCIUM/CREATININE RATIO	0.02	0.26	
CHLORIDE, URINE	110	250	mmol/collection
CREATININE, URINE	800	2500	mg/collection
GLUCOSE, URINE		< 0.9	g//24 hours
POTASSIUM, URINE	25	125	mmol/collection
MAGNESIUM, URINE	50	200	mg/collection
MICROALBUMIN		<30	mg/24 hours
MICROALBUMIN Excretion Rate		</=20	mcg/min
MICROALBUMIN/Creatinine Ratio		<30	mg/g
SODIUM, URINE	40	220	mmol/collection
PHOSPHORUS, URINE	0.3	1.3	g//24 hours
PROTEIN, URINE		< 150	mg/collection
PROTEIN/CREATININE RATIO	0	0.2	
UREA, URINE	12	20	g//24 hours
URIC ACID, URINE	250	750	mg/Collection

CSF

GLUCOSE (ADULT)	50	80	mg/dL
0 - 9 YEARS			
LACTATE	0.5	2.0	mmol/L
LDH	5	30	U/L
PROTEIN	15	45	g/dL

TDM

Acetaminophen	0	30	mcg/mL
Amikacin Peak			
	Peak	20	30
	Trough	4	8
Carbamazepine	6	12	mcg/mL

Digoxin		0.8	2	ng/mL
Gentamicin				
	Peak	5	10	mcg/mL
	Trough	>2		mcg/mL
Lithium		0.6	1.2	mmol/L
Phenobarbital		20	40	mcg/mL
Phenytoin				
	0-3 years	6	15	mcg/mL
	>3 years	10	20	mcg/mL
Salicylate		0	25	mg/dL
Theophylline		10	20	mcg/mL
Tobramycin				
	Peak	5	10	mcg/mL
	Trough		<2	mcg/mL
Valproic Acid				
	<15 years	50	150	mcg/mL
	>15 years	50	100	mcg/mL
Vancomycin				
	Peak	30	45	mcg/mL
	Trough	<=20		mcg/mL

		LOW	HIGH	Onboard Dilution	EXTENDED RANGE	Diluent/ Dilution	Max Reportable
ACE	U/L	1	120				>120
ALBUMIN	g/dL	0.4	10.5				>10.5
ALBUMIN, Fluid	g/dL	0.4	10.5				>10.5
ALP	U/L	5	4555				>4555
ALT	U/L	6	4113				>4113
Ammonia	umol/L	8	997	X1.85	1844	Saline/X2	>3688
AMYLASE	U/L	3.0	6554				>6554
AMYLASE, Fluid	U/L	3.0	6554				>6554
AST	U/L	3	4202				>4202
BHBT	mmol/L	0.01	4.5	X10	45		>45
BILIRUBIN, DIRECT	mg/dL	0.1	15.0				>15.0
BILIRUBIN, TOTAL	mg/dL	0.1	25.0	X10	250		>250
BILIRUBIN, TOTAL, Fluid	mg/dL	0.1	25.0	X10	250		>250
CALCIUM	mg/dL	2.0	24.0				>24
CHOLESTEROL	mg/dL	7	705	X4	2820		>2820
CHOLESTEROL, Fluid	mg/dL	7	705	X4	2820		>2820
CK	U/L	7	4267	X10	42670	Saline	42670**
CHLORIDE	mmol/L	50	150				>150
CO2	mmol/L	5	50	X2	100		>100
CREATININE	mg/dL	0.20	37.00				>37
CREATININE, Fluid	mg/dL	0.20	37.00				>37
CREATININE ENZ	mg/dL	0.20	37.00				>37
ETOH	mg/dL	10	600				>600
GGT	U/L	4	9256				>9256
GLUCOSE	mg/dL	5	800	X5	4000		>4000
GLUCOSE, Fluid	mg/dL	5	800	X5	4000		>4000
HDLC	mg/dL	5	180				>180
IRON	ug/dL	5	1000	X6.55	6550		>6550
POTASSIUM	mmol/L	1.0	10.0				>10.0
Lactate	mmol/L	0.2	13.3	X4	53		>53
Lactate, CSF	mmol/L	0.2	13.3	X4	53		>53
LD (Plasma/Fluid OHS)	U/L	10	1995	X5	3325	Saline	>4500
LD	U/L	10	4500				>4500
LD, Fluid	U/L	10	4500				>4500
LDLC - DIRECT	mg/dL	1	800				>800
LIPASE	U/L	4	1200				>1200
Lp(a)	mg/dL	1.3	90	X4			>360
MAGNESIUM	mg/dL	0.6	9.5				>9.5
SODIUM	mmol/L	100	200				>200
PHOSPHORUS	mg/dL	0.7	25.3				>25.3
PLACA	nmol/m in/mL	10	382				>382
Plasma Hemoglobin	mg/dL	0	800				>800
PROTEIN, TOTAL	g/dL	0.8	18.4				>18.4
PROTEIN, TOTAL, Fluid	g/dL	0.8	18.4				>18.4
TRIGLYCERIDES	mg/dL	7	1420	X4	5680	Saline X10	>7500
TRIGLYCERIDES, Fluid	mg/dL	7	1420	X4	5680	Saline X10	>7500
UREA, BLOOD	mg/dL	2	125	X5	625		>625
UREA, Fluid	mg/dL	2	125	X5	625		>625
URIC ACID	mg/dL	1.0	33.1				>33.1

**Samples above the CK stated maximum dilution are diluted in duplicate with 2 different dilution factors before reporting. (ex. X20 and x40) Dilutions should agree within 10%.

		LOW	HIGH	Onboard Dilution	EXTENDED RANGE	Diluent/ Dilution	Max Reportable
AMYLASE, URINE	U/L	3	6554				>6554
CALCIUM, URINE	mg/dL	2.0	24.0			Saline X5	>120
CALCIUM, Breast Milk	mg/dL	2.0	24.0			Saline X5	>120
CHLORIDE, URINE	mmol/L	20	300				>300
CHLORIDE, Fecal /Breast Milk	mmol/L	20	300				>300
CREATININE, URINE	mg/dL	5.0	740				>740
GLUCOSE, CSF	mg/dL	1	800				>800
GLUCOSE, URINE	mg/dL	1	800				>800
LD, CSF	U/L	10	2000		4500		>4500
MAGNESIUM, URINE	mg/dL	1.8	26.7	X9	80.1		>80.1
MAGNESIUM, Breast Milk	mg/dL	1.8	26.7	X9	80.1		>80.1
MICROALBUMIN	mg/dL	0.5	50	X4	200	Saline X4	>800
SODIUM, URINE	mmol/L	20	400				>400
SODIUM, Fecal/Breast Milk	mmol/L	20	400				>400
PHOS, URINE	mg/dL	5	186				>186
POTASSIUM, URINE	mmol/L	1.0	300				>300
POTASSIUM, Fecal/Breast Milk	mmol/L	1.0	300				>300
UREA, URINE	mg/dL	40	1980				>1980
URIC ACID, URINE	mg/dL	5.0	250				>250
Protein, CSF	mg/dL	7	200	X10	2000		>2000
PROTEIN, URINE	mg/dL	7	200	X10	2000		>2000
A1AT	mg/dL	25.0	300.0	X5	1500		>1500
ASO	IU/mL	50	850	X13.88	4247		>4247
B2 Microglobulin*	mg/L	0.25	16	X6	96		>96
PREALBUMIN	mg/dL	3	50	X4	200		>200
TRANSFERRIN	mg/dL	19	525	X2	1050		>1050
RF	IU/mL	15	200	X10	2000		>2000
C3	mg/dL	11	340	X2	680		>680
C4	mg/dL	3	60	X2	120		>120
CERULOPLASMIN	mg/dL	2	65				>65
HS CRP	mg/L	0.1	160	X10	1600		>1600
CRP	mg/L	0.2	320	X10	3200		>3200
HAPTOGLOBIN	mg/dL	8	255	X4	1020		>1020
IGA *	mg/dL	5	3650	X10	7300		>7300
IGE	IU/mL	25	1000	X10	10000		>10000
IGG*	mg/dL	109	4150	X4	16000		>16000
IGM*	mg/dL	5	1600	X10	3200	Saline X20	>7000

* Verify dilution condition is undiluted before reporting IGA and IGM as "Less Than"

Verify dilution condition is 3:1 before reporting IGG as "Less Than"

		LOW	HIGH	Onboard Dilution	EXTENDED RANGE	Diluent/ Dilution	Max Reportable
Acetaminophen	ug/mL	3	377		377	Saline	Endpoint
Amikacin	ug/mL	2	50	X4	200	Saline	Endpoint
Carbamazepine	ug/mL	1.9	20	X2	40	Saline	Endpoint
Digoxin	ng/mL	0.2	5.0	X4 or X8	40	Saline	Endpoint
Gentamicin	ug/mL	0.5	10	X2	20	Saline	Endpoint
Lithium	mmol/L	0.10	3.51	X40	7.02	Saline	Endpoint
Phenobarbital	ug/mL	2.0	80	X2	160	Saline	Endpoint
Phenytoin	ug/mL	1.8	40	X4	160	Saline	Endpoint
Salicylate	mg/dL	5.0	100	X5	500	Saline	Endpoint
Theophylline	ug/mL	2.0	40	X4	160	Saline	Endpoint
Tobramycin	ug/mL	0.2	10			Tobra Cal 1	Endpoint
Valproic Acid	ug/mL	13	150	X4 or X8	1200	Saline	Endpoint
Vancomycin	ug/mL	1.1	100	X4	400	Saline	Endpoint

Due to the need to access the potential for toxicity and to assist in monitoring a patient after a toxic ingestion it is necessary to manually dilute until a concentration within the AMR is obtained. Laboratory policy for TDM is to dilute to a concentration unless otherwise indicated.

		Cut off	Dilution
Amphet/Methamphet	ng/mL	500	DO NOT DILUTE
Barbituate	ng/mL	200	DO NOT DILUTE
Benzodiazapine	ng/mL	300	DO NOT DILUTE
Cannabannoid	ng/mL	50	DO NOT DILUTE
Cocaine	ng/mL	300	DO NOT DILUTE
Methadone	ng/mL	300	DO NOT DILUTE
Opiates	ng/mL	300	DO NOT DILUTE
PCP	ng/mL	50	DO NOT DILUTE

Dilution Guide			
Dilution	DILUENT	SAMPLE	Program at
X2	100uL	100uL	2
X4	300uL	100uL	4
X5	400uL	100uL	5
X10	900uL	100uL	10
X16	750uL	50uL	16
X20	950uL	50uL	20

HEMOLYSIS		
Cancel and request redraw at 500		
ABBOTT	Comment	Hemolysis Value
Ammonia	CANCEL	50
CPK	INCREASED	50
L.D.H.	INCREASED	50
Potassium	INCREASED	50
Protein, Total	INCREASED	50
Iron	DECREASED/INCREASED	50
Lactic Acid	INCREASED	100-199
ALT	INCREASED	100
AST	INCREASED	100
Magnesium	INCREASED	100
Phosphorous	INCREASED	100
Lactic Acid	CANCEL	200
Acetaminophen	CANCEL	200
Albumin	INCREASED	200
Alkaline Phos.	DECREASED	250

LIPEMIA		
Sample will be airfuged at 200		
ABBOTT	Comment	Lipemia Value
Magnesium	DECREASED	50
Ammonia	RESULT MAY BE COMPROMISED	100-199
Calcium	DECREASED	125
Urea	DECREASED	125
Alkaline Phos.	DECREASED	250
Iron	DECREASED	250
Protein, Total	INCREASED	250
Lipase	DECREASED	400
ALT	DECREASED	500
AST	DECREASED	500
G.G.T.	DECREASED	750

ICTERUS		
ABBOTT	Comment	Icterus Value
Iron	RESULT MAY BE COMPROMISED	2.5
Protein, Total	DECREASED	10
Ammonia	Program Auto Dilution 1:1.85 Manually dilute X2 and program with dilution if needed to resolve	20 or absorbance error
Phosphorous	INCREASED	25
Creatinine	DECREASED	30

Assay	Dearborn	Royal Oak	Grosse Pointe	Troy	Canton	Taylor	Trenton	Wayne	Farmington Hills
Acetaminophen	X	X	X	X	X	X	X	X	X
Amphet/Methamphetamine	X	X	X	X	X	X	X	X	X
Barbiturates	X	X	X	X	X	X	X	X	X
Benzodiazepines	X	X	X	X	X	X	X	X	X
B-OH Butyrate	X	X	X	X	X	X	X	X	X
Cannabinoids	X	X	X	X	X	X	X	X	X
Cocaine	X	X	X	X	X	X	X	X	X
Ethanol	X	X	X	X	X	X	X	X	X
Methadone	X	X	X	X	X	X	X	X	X
Opiates	X	X	X	X	X	X	X	X	X
Phencyclidine (PCP)	X	X	X	X	X	X	X	X	X
Salicylate	X	X	X	X	X	X	X	X	X
Albumin BCG	X	X	X	X	X	X	X	X	X
Alkaline Phosphatase	X	X	X	X	X	X	X	X	X
ALT	X	X	X	X	X	X	X	X	X
Ammonia	X	X	X	X		X	X	X	X
Amylase	X	X	X	X	X	X	X	X	X
AST	X	X	X	X	X	X	X	X	X
Calcium	X	X	X	X	X	X	X	X	X
Cholesterol	X	X	X	X		X	X	X	X
CK	X	X	X	X	X	X	X	X	X
CO2	X	X	X	X	X	X	X	X	X
Creatinine	X	X	X	X	X	X	X	X	X
Direct Bilirubin	X	X	X	X	X	X	X	X	X
Direct LDL	X	X							
GGT	X	X		X					X
Glucose	X	X	X	X	X	X	X	X	X
HDL Ultra	X	X	X	X		X	X	X	X
ICT Module (Cl, NA, K)	X	X	X	X	X	X	X	X	X
Iron	X	X	X	X			X	X	X
Lactic Acid	X	X	X	X	X	X	X	X	X
LDH	X	X		X	X	X	X	X	X
Lipase	X	X	X	X	X	X	X	X	X
Magnesium	X	X	X	X	X	X	X	X	X
Phosphorus	X	X	X	X	X	X	X	X	X
Total Bilirubin	X	X	X	X	X	X	X	X	X
Total Protein	X	X	X	X	X	X	X	X	X
Triglyceride	X	X	X	X		X	X	X	X
Urea Nitrogen	X	X	X	X	X	X	X	X	X
Uric Acid	X	X	X	X	X	X	X	X	X
Urine/CSF Protein	X	X	X	X	X	X	X	X	X
Alpha-1-Antitrypsin	X	X							
ASO	X	X							
Beta2 Microglobulin	X	X							
Ceruloplasmin	X	X							
Complement 3	X	X							
Complement 4	X	X							
CRP	X	X	X	X			X		X
hsCRP	X	X							X
Haptoglobin	X	X							
Immunoglobulin A	X	X							
Immunoglobulin E	X	X							
Immunoglobulin G	X	X							
Immunoglobulin M	X	X							
Lp(a)		X							
Microalbumin	X	X	X	X					X
Prealbumin	X	X	X	X			X		X

RF	X	X							X
Transferrin	X	X	X	X			X	X	X
Amikacin		X							
Carbamazepine	X	X				X		X	
Digoxin	X	X	X	X	X	X	X	X	X
Gentamicin	X	X	X	X				X	
Lithium	X	X				X			
Phenobarbital	X	X							
Phenytoin	X	X	X	X		X		X	X
Theophylline	X	X							
Tobramycin	X	X	X	X					
Valproic Acid	X	X	X	X		X	X	X	X
Vancomycin	X	X	X	X		X	X	X	X

	Centrifuged SST Tubes			Red Top Tubes without Gel barrier			Plasma					
	RT	2° - 8° C	-20° C pourover	RT	2° - 8° C	-20° C pourover	RT	2° - 8° C	-20° C pourover			
Albumin	4 hours	7 days	7 days	4 hours			48 hours	7 days	7 days			
Alkaline Phosphatase	4 hours	7 days	7 days	4 hours			72 hours	7 days	7 days			
ALT	4 hours	7 days	7 days	4 hours			72 hours	7 days	7 days			
Ammonia (EDTA)								2 hours	14 days			
Amylase	4 hours	7 days	7 days	4 hours			72 hours	7 days	7 days			
AST	4 hours	7 days	7 days	4 hours			72 hours	7 days	7 days			
BHBT	4 hours	7 days	7 days	4 hours			4 hours	7 days	7 days			
Calcium	4 hours	7 days	7 days	4 hours			72 hours	7 days	7 days			
Chloride	4 hours	7 days	7 days	4 hours			72 hours	7 days	7 days			
Cholesterol	4 hours	7 days	7 days	4 hours								
CK	4 hours	7 days	7 days	4 hours			72 hours	7 days	7 days			
CO2 (closed tube)	4 hours	72 hours	7 days	4 hours			4 hours	72 hours	7 days			
Creatinine	4 hours	7 days	7 days	4 hours			72 hours	7 days	7 days			
Direct Bilirubin (light protected)	4 hours	7 days	7 days	4 hours			4 hours	7 days	7 days			
Direct LDL	4 hours	7 days	7 days	4 hours								
Ethanol	24 hours	72 hours	7 days	24 hrs	72 hours	7 days	24 hours	72 hours	7 days			
GGT	4 hours	7 days	7 days	4 hours								
Glucose (SST or FIOx)	4 hours	7 days	7 days	4 hours						4 hours (PST) 24 hours(FIOx)	72 hours (PST) 7 days(FIOx)	7 days
HDL Ultra	4 hours	7 days	7 days	4 hours								
Iron	4 hours	7 days	7 days	4 hours								
Lactic Acid (FIOx)							72 hours	7 days	7 days			
LD	4 hours	48 hours	7 days	4 hours								
Lipase	4 hours	7 days	7 days	4 hours								
Magnesium	4 hours	7 days	7 days	4 hours								
Phosphorus	4 hours	7 days	7 days	4 hours								
PLACA	24 hours	14 days	18 months	4 hours	30 hours							
Plasma Hemoglobin							2-4 hours	7 days				
Potassium	4 hours	7 days	7 days	4 hours								
Potassium, PI (Li Hep)							72 hours	5 days	7 days			
							2 hours	7 days				
Sodium	4 hours	7 days	7 days	4 hours								
							72 hours	7 days	7 days			

Total Bilirubin	4 hours	7 days	7 days	4 hours		4 hours	7 days	7 days
Total Protein	4 hours	7 days	7 days	4 hours		72 hours	7 days	7 days
Triglyceride	4 hours	7 days	7 days	4 hours				
Urea Nitrogen	4 hours	7 days	7 days	4 hours		4 hours	7 days	7 days
Uric Acid	4 hours	7 days	7 days	4 hours		72 hours	7 days	7 days
Alpha-1-Antitrypsin	12 hours	7 days	3 months	12 hours				
ASO	8 hours	7 days	3 months	2 hours				
Beta2 Microglobulin	2 hours	7 days	14 days	2 hours				
Ceruloplasmin	12 hours	7 days	3 months	12 hours				
Complement 3	12 hours	7 days	3 months	12 hours				
Complement 4	12 hours	7 days	3 months	12 hours				
CRP	24 hours	7 days	3 months	12 hours				
hsCRP	24 hours	7 days	3 months	12 hours				
Haptoglobin	4 hours	7 days	7 days	4 hours				
Immunoglobulin A	12 hours	7 days	3 months	12 hours				
Immunoglobulin E	4 hours	7 days	7 days	4 hours				
Immunoglobulin G	12 hours	7 days	3 months	12 hours				
Immunoglobulin M	12 hours	7 days	3 months	12 hours				
Lp(a)	4 hours	3 days	2 months	4 hours				
Prealbumin	4 hours	7 days	7 days	4 hours				
RF	4 hours	7 days	3 months	4 hours				
Transferrin	12 hours	7 days	3 months	12 hours		12 hours	3 days	6 months

	Red Top Tubes without Gel barrier			Without Gel Barrier Pourovers		
	RT	2° - 8° C	-20° C	RT	2° - 8° C	-20° C
Acetaminophen	2 hours	24 hours		2 hours	7 days	3 months
Amikacin	2 hours	24 hours		4 hours	7 days	1 month
Carbamazepine	2 hours	24 hours		2 hours	7 days	3 months
Digoxin	2 hours	24 hours		3 hours	8 days	4 months
Gentamicin	2 hours	24 hours		2 hours	7 days	3 months
Lithium	2 hours	8 hours		24 hours	7 days	3 months
Phenobarbital	2 hours	24 hours		2 hours	7 days	3 months
Phenytoin	2 hours	24 hours		2 hours	7 days	3 months

Salicylate	2 hours	24 hours		2 hours	7 days	3 months
Theophylline	2 hours	24 hours		2 hours	7 days	3 months
Tobramycin	2 hours	24 hours		2 hours	7 days	3 months
Valproic Acid	2 hours	24 hours		2 hours	7 days	3 months
Vancomycin	2 hours	24 hours		2 hours	7 days	3 months

	Urine		
	RT	2° - 8° C	-20° C pourover
Amphet/Methamphet	2 hours	24 hours	3 months
Barbiturates	2 hours	24 hours	3 months
Benzodiazepines	2 hours	24 hours	3 months
Cannabinoids	2 hours	24 hours	3 months
Cocaine	2 hours	24 hours	3 months
Methadone	2 hours	24 hours	3 months
Opiates	2 hours	24 hours	3 months
Phencyclidine (PCP)	2 hours	24 hours	3 months
Amylase, URINE	4 hours	4 days	Allowed
Calcium, URINE	4 hours	7 days	Allowed
Chloride, URINE	4 hours	7 days	Allowed
Creatinine, URINE	4 hours	4 days	Allowed
Glucose, URINE	4 hours	4 days	Allowed
Magnesium, URINE	4 hours	7 days	Allowed
Microalbumin	4 hours	7 days	Allowed
Sodium, URINE	4 hours	7 days	Allowed
Phosphorus, URINE	4 hours	7 days	Allowed
Potassium, URINE	4 hours	7 days	Allowed
Urea, URINE	4 hours	4 days	Allowed
Uric Acid, URINE	4 hours	7 days	Allowed
Protein, URINE	4 hours	7 days	Allowed

	Centrifuged SST Tubes			Red Top Tubes without Gel barrier			Plasma					
	RT	2° - 8° C	-20° C pourover	RT	2° - 8° C	-20° C pourover	RT	2° - 8° C	-20° C pourover			
Albumin	4 hours	7 days	7 days	4 hours			48 hours	7 days	7 days			
Alkaline Phosphatase	4 hours	7 days	7 days	4 hours			2 days	4 days	7 days			
ALT	4 hours	7 days	7 days	4 hours			24 hours	7 days	7 days			
Ammonia (EDTA)								2 hours	14 days			
Amylase	4 hours	7 days	7 days	4 hours			7 days	7 days	7 days			
AST	4 hours	7 days	7 days	4 hours			2 days	7 days	7 days			
BHBT	4 hours	7 days	7 days	4 hours						7 days		
Calcium	4 hours	7 days	7 days	4 hours			2 days	7 days	7 days			
Chloride	4 hours	7 days	7 days	4 hours			2 days	7 days	7 days			
Cholesterol	4 hours	7 days	7 days	4 hours			2 days	7 days				
CK	4 hours	7 days	7 days	4 hours			24 hours	7 days	7 days			
CO2 (closed tube)	4 hours	72 hours	7 days	4 hours			24 hours	3 days	7 days			
Creatinine	4 hours	7 days	7 days	4 hours			2 days	7 days	7 days			
Direct Bilirubin (light protected)	4 hours	7 days	7 days	4 hours			24 hours	7 days	7 days			
Direct LDL	4 hours	7 days	7 days	4 hours								
Ethanol (FIOx)	24 hours	72 hours	7 days	24 hrs	72 hours	7 days		3 days	7 days			
GGT	4 hours	7 days	7 days	4 hours			2 days	7 days				
Glucose (SST or FIOx)	4 hours	7 days	7 days	4 hours			4 hours (PST) 24 hours(FIOx)	72 hours (PST) 7 days(FIOx)	7 days			
HDL Ultra	4 hours	7 days	7 days	4 hours			24 hours	7 days				
Iron	4 hours	7 days	7 days	4 hours			24 hours	7 days	7 days			
Lactic Acid (FIOx)							8 hours	7 days	30 days			
LD	4 hours	48 hours	7 days	4 hours			24 hours	7 days				
Lipase	4 hours	7 days	7 days	4 hours					7 days			
Magnesium	4 hours	7 days	7 days	4 hours			7 days	7 days	7 days			
Phosphorus	4 hours	7 days	7 days	4 hours			24 hours	7 days	7 days			
PLACA	24 hours	14 days	18 months	4 hours	30 hours							
Plasma Hemoglobin							2-4 hours	7 days				
Potassium	4 hours	7 days	7 days	4 hours						24 hours	7 days	7 days
Potassium, PI (Li Hep)							2 hours	7 days				
Sodium	4 hours	7 days	7 days	4 hours						2 days	7 days	7 days

Total Bilirubin	4 hours	7 days	7 days	4 hours	Pearson, Henton, Taylor, Wayne, Canton	24 hours	7 days	7 days
Total Protein	4 hours	7 days	7 days	4 hours		24 hours	7 days	7 days
Triglyceride	4 hours	7 days	7 days	4 hours		2 days	7 days	
Urea Nitrogen	4 hours	7 days	7 days	4 hours		24 hours	7 days	7 days
Uric Acid	4 hours	7 days	7 days	4 hours		2 days	7 days	7 days
Alpha-1-Antitrypsin	12 hours	7 days	3 months	12 hours				
ASO	8 hours	7 days	3 months	2 hours				
Beta2 Microglobulin	2 hours	7 days	14 days	2 hours				
Ceruloplasmin	12 hours	7 days	3 months	12 hours				
Complement 3	12 hours	7 days	3 months	12 hours				
Complement 4	12 hours	7 days	3 months	12 hours				
CRP	24 hours	7 days	3 months	12 hours				
hsCRP	24 hours	7 days	3 months	12 hours				
Haptoglobin	4 hours	7 days	7 days	4 hours				
Immunoglobulin A	12 hours	7 days	3 months	12 hours				
Immunoglobulin E	4 hours	7 days	7 days	4 hours				
Immunoglobulin G	12 hours	7 days	3 months	12 hours				
Immunoglobulin M	12 hours	7 days	3 months	12 hours				
Lp(a)	4 hours	3 days	2 months	4 hours				
Prealbumin	4 hours	7 days	7 days	4 hours				
RF	4 hours	7 days	3 months	4 hours				
Transferrin	12 hours	7 days	3 months	12 hours				

	Red Top Tubes without Gel barrier			Without Gel Barrier Pourovers		
	RT	2° - 8° C	-20° C	RT	2° - 8° C	-20° C
Acetaminophen	2 hours	24 hours		2 hours	7 days	3 months
Amikacin	2 hours	24 hours		4 hours	7 days	1 month
Carbamazepine	2 hours	24 hours		2 hours	7 days	3 months
Digoxin	2 hours	24 hours		3 hours	8 days	4 months
Gentamicin	2 hours	24 hours		2 hours	7 days	3 months
Lithium	2 hours	8 hours		24 hours	7 days	3 months
Phenobarbital	2 hours	24 hours		2 hours	7 days	3 months
Phenytoin	2 hours	24 hours		2 hours	7 days	3 months

Dearborn, Trenton, Taylor, Wayne, Canton

Salicylate	2 hours	24 hours		2 hours	7 days	3 months
Theophylline	2 hours	24 hours		2 hours	7 days	3 months
Tobramycin	2 hours	24 hours		2 hours	7 days	3 months
Valproic Acid	2 hours	24 hours		2 hours	7 days	3 months
Vancomycin	2 hours	24 hours		2 hours	7 days	3 months

	Urine		
	RT	2° - 8° C	-20° C pourover
Amphet/Methamphet	2 hours	24 hours	3 months
Barbiturates	2 hours	24 hours	3 months
Benzodiazepines	2 hours	24 hours	3 months
Cannabinoids	2 hours	24 hours	3 months
Cocaine	2 hours	24 hours	3 months
Methadone	2 hours	24 hours	3 months
Opiates	2 hours	24 hours	3 months
Phencyclidine (PCP)	2 hours	24 hours	3 months
Amylase, URINE	4 hours	4 days	Allowed
Calcium, URINE	4 hours	7 days	Allowed
Chloride, URINE	4 hours	7 days	Allowed
Creatinine, URINE	4 hours	4 days	Allowed
Glucose, URINE	4 hours	4 days	Allowed
Magnesium, URINE	4 hours	7 days	Allowed
Microalbumin	4 hours	7 days	Allowed
Sodium, URINE	4 hours	7 days	Allowed
Phosphorus, URINE	4 hours	7 days	Allowed
Potassium, URINE	4 hours	7 days	Allowed
Urea, URINE	4 hours	4 days	Allowed
Uric Acid, URINE	4 hours	7 days	Allowed
Protein, URINE	4 hours	7 days	Allowed

Approved List of fluids for Abbott Architect Chemistry Analyzers								
	CSF	Peritoneal	Pleural	Dialysate	Pancreatic	Fecal	Amniotic	Breast Milk
Albumin		X	X					
Amylase		X	X		X			
Total Bilirubin		X	X					
Calcium								X
Chloride						X		X
Cholesterol		X	X					
Creatinine		X		X				
Glucose		X	X	X			X	
Glucose, CSF	X							
Lactic Acid	X							
LD	X	X	X					
Magnesium								X
Potassium						X		X
Protein, CSF	X							
Protein, Total		X	X					
Sodium						X		X
Triglycerides		X	X					
Urea Nitrogen				X				

CSF has been approved by the FDA for Glucose and Protein testing on the analyzer and will be tested using the analyzer designated method. (Glucose, CSF and Protein, CSF). The Peritoneal, Pleural, Dialysate, Pancreatic and Amniotic fluid have been laboratory validated using the Serum Method. Fecal and Breast Milk have been validated using the Urine Method. Refer to **Attachment C** for the reportable ranges.

Reference Ranges

Albumin	
Peritoneal Fluid	Because of the wide range of albumin levels seen in peritoneal fluid, results are best evaluated using the serum-ascites albumin gradient. A serum to fluid gradient >1.1 g/dL is seen in transudates. A serum to fluid gradient <= 1.1 g/dL is seen in exudates.
Pleural Fluid	A serum to fluid gradient <= 1.2 g/dL is consistent with exudate.

Amylase	
Peritoneal Fluid	Values greater than or equal to 3X a simultaneously analyzed serum value are considered abnormal
Pleural Fluid	<104 U/L or a fluid to serum ratio less than 1.5-2.0

Pancreatic Fluid	Pancreatic pseudocysts generally contain significantly elevated amylase (e.g.) >250 IU/L, whereas mucinous and serous cysts and adenocarcinoma usually have lower concentrations. Results should be used in conjunction with clinical information, imaging studies, cytology and other pancreatic cyst markers.
------------------	---

Total Bilirubin	
Peritoneal Fluid	Ascitic fluid bilirubin < 6mg/dL and a ascitic fluid to serum bilirubin <1.0 mg/dL
Pleural Fluid	Fluid to serum ratio >= 0.6 is consistent with exudate

Cholesterol	
Peritoneal Fluid	<48 mg/dL
Pleural Fluid	>45 mg/dL = exudate

Creatinine	
Peritoneal Fluid	0.5-2.0 mg/dL

Glucose	
Amniotic Fluid	>15mg/dL
CSF	50-80 mg/dL
Peritoneal Fluid	approximates that found in serum
Pleural Fluid	>60 mg/dL

Lactic Acid	
CSF	0.5-2.0 mmol/L

LD	
CSF	5- 30 U/L
Peritoneal Fluid	<63 U/L
Pleural Fluid	>2/3

Protein	
CSF	15-45 mg/dL
Peritoneal Fluid	Transudate <3.0 g/dL, Exudate >=3.0 g/dL
Pleural Fluid	Transudate <3.0 g/dL, Exudate >=3.0 g/dL

Triglycerides	
Peritoneal	<65 mg/dL
Pleural Fluid	>= 110 mg/dl indicative for chylous effusion, <50 mg/dL indicative of non-chylous effusion, 50-109 mg/dL are equivocal

Fecal Electrolytes	
Sodium, Fecal	Not Established
Potassium, Fecal	Not established

Chloride, Fecal	Fecal chloride may be elevated (>35mmol/L) in phenolphthalein (or phenolphthalein plus magnesium hydroxide) induced diarrhea. Fecal chloride may be low(<20 mmol/L) in sodium sulfate induced diarrhea. Fecal chloride concentration is markedly elevated >60 mmol/L in infants and > 100 mmol/L in adults associated with congenital and secondary
-----------------	---

Breast Milk	
Sodium	4-11 mmol/L
Potassium	8.7-17.5 mmol/L
Chloride	6-16 mmol/L
Calcium	20.8-48.0 mg/dL
Magnesium	2.2-5.1 mg/dL