

---

## RETICULOCYTE COUNTS MILLER DISK METHOD

RC.HM.PR.032.r04

---

### Principle

Reticulocytes are immature red cells which contain remnants of RNA. When blood is mixed with a supravital stain, the precipitated remnants of RNA are seen as granules or filaments within the cell. Generally, the most mature retics are those with the least granulation.

The Miller disk method has a greater degree of precision than the standard method because of ease of counting a restricted field.

**NOTE:** *This method is to be employed when automated reticulocyte counts need verification.*

### Specimen Collection and Handling

---

|                                      |  |  |
|--------------------------------------|--|--|
| Type:                                | Whole blood collected in a 4 mL vacutainer. This is the preferred sample.<br>OR<br>Capillary blood collected in a microtainer. |  |
| Anticoagulant:                       | K <sub>2</sub> EDTA  |  |
| Amount:                              | Whole blood  | - Minimum sample size is 2.0 mL<br>- Optimum sample size is 4.0 mL   |
|                                      | Capillary blood  | - Minimum sample size is 300 mcL<br>- Optimum sample size is 500 mcL |
| Special Handling:                    | Specimen must be well mixed for minimum of two minutes before mixing with stain.   |  |
| Timing:                              | Specimen is stable for 8h at room temperature;<br>24h at 4°C.  |  |
| Criteria for Unacceptable Specimens: | Specimens containing clots or inappropriate volumes are unacceptable and must be redrawn.                                      |  |

---

## RETICULOCYTE COUNTS – MILLER DISK METHOD

### Supplies

#### Reagent

1. Thermo Scientific Richard-Allan Scientific Reticulocyte Stain

#### Solution:

#### Ingredients:

| Component                       | Weight %   |
|---------------------------------|------------|
| Methylene blue trihydrate       | <1%        |
| Sodium chloride                 | <1%        |
| Water                           | 95 - 100 % |
| Oxalate, potassium, monohydrate | <1%        |

**Storage:** Store at 15-30°C.

**Stability:** Make certain that product has been capped immediately after each use and it will remain stable for the stated expiration date. Do not use product past expiration date printed on the reagent label. Record date opened and expiration date on container.

**Filter before use.**

#### Equipment

1. Heat Block
2. Ocular containing a Miller disk (Scientific Instruments)

### Quality Control

Manual retics are usually only done to verify the automated retic values. If the manual retic smear estimate/count is  $\pm 25\%$  of the automated retic count, report the automated count and attached comment "verified". If  $>25\%$  difference, report the manual retic smear count. Remember to enter (.) for IRF if reporting manual retic.

### Procedure

1. Combine 3 drops of well mixed EDTA anticoagulated blood with 2 drops of filtered reticulocyte stain in a 10 x 75 mm glass test tube and mix well. If patient is obviously anemic, add more blood to stain.
2. Cover with parafilm and warm to 37°C in a dry heating block.
3. Allow to stand for 15 minutes at 37°C (see Note #1).
4. **Mix** suspension **well** and make a thin, evenly spread smear (see Note #2). Allow to air dry for 15 minutes.
5. Counting: Use the Miller Ocular (contains Miller disk). (See Figure A.)
  - a. The Miller disk consists of a large square with a small square in one corner - the ratio of the large square area to small square area is 9:1.




Figure A.

- b. Use 100X oil objective.
- c. Find appropriate counting area in thin portion of smear where cells are **evenly** distributed and not overlapping each other nor widely separated.

*Printed copies of this document are not considered up-to-date. Please verify current version date with online document.*

## RETICULOCYTE COUNTS – MILLER DISK METHOD

- d. Count all **RBCs (mature and reticulated)** in the **small square** and at the same time count the **reticulocytes in the total large square**. Count consecutive fields until the RBCs counted equal 111. (See Notes 3 and 6).

**NOTE:** Only cells that touch the left and top lines of the Miller squares are to be counted. Do NOT count cells touching the right and bottom lines.

---

### Calculations and Interpretations

1. % reticulocytes = # of retics counted

$$\text{Example A: } \frac{23 \text{ retics}}{10} = 2.3\%$$

2. Absolute number of reticulocytes = % of retics x RBC count

$$\begin{aligned} \text{Example B: } & 2.3\% \times 4.00 \times 10^{12}/\text{L} = \\ & .023 \times 4.00 \times 10^{12}/\text{L} = \\ & .092 \times 10^{12}/\text{L} = \\ & 92 \times 10^9/\text{L} = \underline{\underline{92 \text{ bil/L}}} \end{aligned}$$

### Reporting Results

Refer to Manual Resulting in SoftLAB procedure.

### Expected Values

#### **NORMAL RANGE:**

For current normal values, refer to Hematology Normal Values procedure.

---

### Limitations

The standard error in the retic count varies, depending on the reticulocyte level and number of cells counted. At the 1% retic level, the error is approximately  $\pm 60\%$  and decreases to  $\pm 30\%$  at the 4% retic level.

Possible sources of error: Precipitated stain, siderocytes, Heinz bodies, Howell-Jolly bodies, superimposed platelets, inappropriate amounts of stain or blood.

---

### Notes

1. The suspension may also be incubated for 15 minutes at room temperature.
2. It is extremely important that the blood and stain be mixed **well** prior to making smears. Reticulocytes have a lower specific gravity than mature red blood cells and, therefore, settle on top of the red blood cells in the mixture.
3. The College of American Pathologists' (CAP) definition of a reticulocyte is that it must have at least two or more clumps of granules visible without fine focus.

*Printed copies of this document are not considered up-to-date. Please verify current version date with online document.*

## RETICULOCYTE COUNTS – MILLER DISK METHOD

4. To preserve retic smear from fading in storage, counterstain with Wright's stain.
  5. Alternate counting method: The number of reticulocytes in 1000 RBC's counted in consecutive fields. Calculate percentage as described above.
  6. When counting 111 RBCs on the Miller ocular, we are actually counting the number of retics per 1000 RBCs.
  7. When retic percentage is "0", report as "0".
  8. Enter (.) for IRF if reporting a manual retic that does not agree with instrument IRF count.
- 

### References

1. Brecher G. New methylene blue as a reticulocyte stain. Am J Clin Path 1949; 19: 895.
  2. Miale J. Laboratory medicine - hematology. 6th Ed. St Louis: CV Mosby. 1982: 865.
  3. NCCLS Document H44-P Reticulocyte counting by flow cytometry; proposed guideline. November, 1993: Vol. 13, No. 18.
- 

### Authorized Reviewers

Medical Director, Hematology

---

*Printed copies of this document are not considered up-to-date. Please verify current version date with online document.*

## RETICULOCYTE COUNTS – MILLER DISK METHOD

### Document Control

**Location of Master:** Hematology Procedure Manual

**Master electronic file stored on the Clinical Pathology server:**

S:\HEMACOAG\Document Control\Hematology\Procedure\Master Documents\reticman.doc

**Number of Controlled Copies posted for educational purposes: 0**

**Number of circulating Controlled Copies: 0**

**Location of circulating Controlled Copies: NA**

### Document History

| Signature                              | Date       | Revision # |   | Related Documents Reviewed/ Updated |
|--|------------|------------|---|-------------------------------------|
| Prepared by: Nancy Ramirez, MT(ASCP)SH | 10/1985    |            |   |                                     |
| Approved by: Joan C. Mattson, MD       |            |            |   |                                     |
|  |            |            |   |                                     |
| Reviewed by: (Signature)               | Date       | Revision # | Modification  | Related Documents Reviewed/ Updated |
| Joan C Mattson, MD                     | 12/29/1987 |            | OK  |                                     |
| Joan C Mattson, MD                     | 02/20/1989 |            | No change   |                                     |
| Joan C Mattson, MD                     | 03/12/1990 |            | Retyped in NCCLS format   |                                     |
| Joan C Mattson, MD                     | 12/24/1991 |            | No change   |                                     |
| Joan C Mattson, MD                     | 12/05/1992 |            | OK  |                                     |
| Joan C Mattson, MD                     | 12/27/1993 |            | OK  |                                     |
| Joan C Mattson, MD                     | 12/12/1994 |            | OK  |                                     |
| Joan C Mattson, MD                     | 12/22/1995 |            | Updated procedure step 6 to clarify counting protocol as per NCCLS.   |                                     |
| Joan C Mattson, MD                     | 02/07/1997 |            | Pg. 1 updated min vol.  |                                     |
| Noelle Procopio, MT(ASCP)SH            | 01/05/1998 |            | No change   |                                     |
| Noelle Procopio, MT(ASCP)SH            | 01/04/1999 |            | No change   |                                     |
| Joan C Mattson, MD                     | 01/24/2000 |            | No change   |                                     |
| Joan C Mattson, MD                     | 01/28/2000 |            | Updated R3000 reference to R3500 references and when to do manual retic; pg. 1, updated STKS references to Sysmex (Reporting Results) Pg. 3 |                                     |
| Noelle Procopio, MT(ASCP)SH            | 12/04/2001 |            | No change   |                                     |
| Noelle Procopio, MT(ASCP)SH            | 12/30/2002 |            | No change   |                                     |
| Joan C Mattson, MD                     | 12/29/2003 |            | Updated K <sub>2</sub> EDTA, pg. 1; added reagent stability, pg. 1  |                                     |
| Noelle Procopio, MT(ASCP)SH            | 12/15/2004 |            | No change   |                                     |
| Joan C Mattson, MD                     | 02/16/2005 | 00         | Removed references to Sysmex R3000; standardized procedure format.  |                                     |
| Noelle Procopio, MT(ASCP)SH            | 10/24/2006 |            | No change   |                                     |
| Ann Marie Blenc, MD                    | 04/19/2007 |            | No change; new director   |                                     |
| Ann Marie Blenc, MD                    | 09/05/2007 | 01         | Updated specimen tube   |                                     |

**Printed copies of this document are not considered up-to-date. Please verify current version date with online document.**

## RETICULOCYTE COUNTS – MILLER DISK METHOD

|                     |            |    |  |    |
|---------------------|------------|----|--|----|
|                     |            |    | volume, pg. 1  |    |
| Ann Marie Blenc, MD | 03/11/2008 |    | No change  |    |
| Ann Marie Blenc, MD | 02/23/2009 |    | No change  |    |
| Ann Marie Blenc, MD | 02/16/2010 | 02 | Referred expected values to Normal Range procedure.  |    |
| Ann Marie Blenc, MD | 07/27/2011 | 03 | Removed references to MIsys approach to resulting; clarified counting RBCs in small square (from “retics” to “reticulated”); added comment re (.) for IRF if not reporting instrument retic count; changed reporting 0% retics from “none seen” to reporting as to “0”.  | OK |
| Ann Marie Blenc, MD | 11/07/2013 |    | Logo update only.  | OK |
| Ann Marie Blenc, MD | 04/16/2015 |    | No change  | OK |
| Ann Marie Blenc, MD | 03/20/2017 |    | Logo update only.  | OK |
| Elizabeth Sykes, MD | 02/02/2018 |    |  |    |
| Peter Millward, MD  | 01/30/2019 |    | New Medical Director   |    |
| Ann Marie Blenc, MD | 05/26/2020 | 04 | Updated reagent from Basic Blue 24 (New Methylene Blue) to Thermo Scientific Richard-Allan Scientific Reticulocyte Stain. Updated reagent stability from 3 months to the expiration date printed on the reagent label. Removed initial procedure step of warming stain prior to adding blood. In procedure step 1 changed number of drops of blood from 2 to 3 and changed “10 x 75 mm test tube” to “10 x 75 mm glass test tube.” In procedure step 3 changed “Allow to stand for 10-15 minutes at 37 °C” to “Allow to stand for 15 minutes at 37 °C.” In procedure step 4 changed “Allow to air dry” to “Allow to air dry for 15 minutes.” Removed Note 1: “The time allowed for staining of reticulocytes is not critical. However, it should never be less than 5 minutes.” Updated Note 1 to state “The suspension may also be incubated for 15 minutes at room temperature.” | OK |
|                     |            |    |  |    |
|                     |            |    |  |    |
|                     |            |    |  |    |

*Printed copies of this document are not considered up-to-date. Please verify current version date with online document.*