



Current Status: Active

PolicyStat ID: 9127737

# Beaumont

Origination: 5/1/2021  
Effective: 5/1/2021  
Last Approved: 5/1/2021  
Last Revised: 5/1/2021  
Next Review: 5/1/2023  
Document Contact: *Tamara Sabih: Medical  
Technologist Lead*  
Area: *Laboratory-Hematology*  
Key Words:  
Applicability: *Royal Oak*

## Coagulation Autoverification-RO

Document Type; Policy

### I. PURPOSE AND OBJECTIVE:

The purpose of this policy is to identify how to detect any failure of the software system to release or not release sample results for auto-verification based on special criteria and reported to Information Technology (IT) for correction / modification of the software. This is a zero-tolerance system for errors in the Laboratory Information System (LIS) software application.

### II. POLICY STATEMENT:

The Coagulation laboratory initiated auto-verification on 10/6/01. It is the department policy that all normal results and selected abnormal results from the IL ACL-Top coagulation analyzers will be auto-verified in the laboratory's LIS.

### III. ANNUAL REVIEW:

Validation of auto-verification will be performed annually and whenever there is a change to the system that could affect the autoverification logic.

### IV. ACRONYMS:

- A. Activated Partial Thromboplastin Time (aPTT)
- B. Anticoagulant Management Services (AMS)
- C. Fibrinogen Equivalent Units (FEU)
- D. Fibrinogen (FIB)
- E. International Normalized Ratio (INR)
- F. Seconds (sec)
- G. Thrombin Time (TT)

### V. PROCEDURE:

- A. The following criteria will be used for autoverification:

#### Autoverification Criteria

### IL ACL TOP Coagulation Analyzers

B. When all of the following criteria evaluate as "True", then the LIS system will auto-verify INR, aPTT, FIB, TT, or D-dimer.

Criteria Name	Range
INR	0.9 – 3.0
INR (AMS)	0.9 – 3.0
aPTT	25 - 95 sec
FIB	101 - 449 mg/dL
TT	16 - 25 sec
D-dimer	250 – 10,000 ng/mL FEU

### Delta Check Rules

C. In addition, an assay not passing the following delta check rules assigned into the LIS will not be auto-verified:

Criteria Name:	Range:
All INR's	Difference Between Previous INR and Current INR $\geq 2.0$

D. Validation of auto-verification done by using previously tested patient samples, which are selected in order to challenge each of the previously stated criteria as listed below.

Assay	Criteria
INR/INR AMS	0.9-3.0 sec, passed delta check, autoverified
INR/INR AMS	0.9-3.0 sec failed delta check, did not autoverify
INR/INR AMS	> 3.0 sec, did not autoverify
aPTT	< 25 sec, did not autoverify
aPTT	25-95 sec, autoverified
aPTT	> 95 sec, did not autoverify
FIB	0-100 mg/dL, did not autoverify
FIB	101-449 mg/dL, autoverified
FIB	450-1000 mg/dL, did not autoverify
TT	<16 sec, did not autoverify
TT	16-25 sec, autoverified
TT	> 25 sec, did not autoverify
D-Dimer	<250 ng/mL FEU, did not autoverify
D-Dimer	250 – 10,000 ng/mL FEU, autoverified
D-Dimer	>10,000 ng/mL FEU, did not autoverify

## Attachments

[Attachment A-Coagulation Autoverification Validation Worksheet.pdf](#)

## Approval Signatures

Step Description	Approver	Date
CP Chief Medical Director	Peter Millward: Chief, Pathology Service Line	5/1/2021
Coagulation Medical Director Designee	Marc Smith: System Med Dir, Coagulation	4/30/2021
Policy and Forms Steering Committee Approval (if needed)	Tamara Sabih: Medical Technologist Lead	4/22/2021
Policy and Forms Steering Committee Approval (if needed)	Gail Juleff: Project Mgr Policy	4/22/2021
System Manager	Rebecca Bacarella: Mgr Laboratory	4/22/2021
	Tamara Sabih: Medical Technologist Lead	4/7/2021

## Applicability

Royal Oak