



Current Status: Active

PolicyStat ID: 9586407

Beaumont

Origination: 5/13/2021
Effective: 5/13/2021
Last Approved: 5/13/2021
Last Revised: 5/13/2021
Next Review: 5/13/2023
Document Contact: *Tamara Sabih: Medical
Technologist Lead*
Area: *Laboratory-Hematology*
Key Words:
Applicability: *Royal Oak*

Coagulation Equipment Quality Control- RO

Document Type: Procedure

I. PURPOSE AND OBJECTIVE:

The purpose of this procedure is to guide the medical technologist with their daily and quarterly quality controls.

II. ACRONYMS:

- A. Clinical and Laboratory Standards Institute (CLSI)
- B. Platelet Poor Plasma (PPP)
- C. Quality Control (QC)

III. PROCEDURE:

- A. Daily Quality Control:
 - 1. Record Temperatures (Temperature Documentation Log)
 - a. Freezers
 - b. Centrifuges
 - c. Waterbath
 - d. Refrigerators
 - 2. Temperatures must fall within posted tolerance limits on the daily temperature documentation log. If temperatures are outside specified limits for:
 - a. Freezers: Temperature may be adjusted manually. Document corrective action on the Coag QC Corrective Action form. If unable to adjust, notify hematology management or call Facility Management (248-551-6300) to place a service call. Move items if necessary.
 - b. Centrifuges: Temperatures may be adjusted manually. Document corrective action on the Coag QC Corrective Action form. If unable to adjust, notify hematology management or Facility Management (248-551-6300) to place a service call and place the centrifuge out of service.
 - c. Waterbath: Temperature may be adjusted manually. Document corrective action on the Coag QC Corrective Action form. If unable to adjust, notify hematology management or call Facility

Management (248-551-6300) to place a service call. Use alternate 37°C water bath (special testing) to thaw specimens.

- d. Refrigerators: Notify hematology management or call Facility Management (248-551-6300) to place a service call. Document corrective action on the Coag QC Corrective Action form. Move items if necessary.

B. Quarterly Quality Control:

1. **Centrifuge Platelet Checks:** The integrity of coagulation specimens must be checked every three months. Platelet-poor-plasma (PPP) with a platelet count of <10 bill/L should be obtained after centrifuging routine coagulation specimens in all coagulation centrifuges.
 - a. Choose five fresh (non BRL) specimens that have not been centrifuged. Number them 1-5 without defacing the bar code. Invert them gently 2 or 3 times.
 - b. Spin them in the selected centrifuge at correct speed and time, making note of the Asset tag # of the centrifuge.
 - c. Remove from the centrifuge immediately upon stopping and perform platelet check within ten minutes.
 - d. Position sample probe of hematology analyzer vertically into plasma directly in the middle of the tube (avoid the sides of the tube) and approximately ¼" from the top surface of the plasma. Refer to Iris "Guideline for Platelet Poor Plasma Testing" (Attachment A).
 - e. Record the platelet counts on the Quarterly Platelet Count Check log sheet.
 - f. All specimens must achieve a platelet count of <10 bill/L.
 - g. Specimens can then be processed on the coagulation instruments for requested testing.
 - h. If two or more specimens have a platelet count above 10 bill/L go directly to steps a, b and c. If any one specimen has a platelet count above 10 bill/L, another specimen must be used. Mix the new fresh specimen gently, centrifuge it (using the centrifuge which gave the unacceptable platelet count). Run it on a hematology analyzer for a platelet count. If, on the new specimen, a platelet count above 10 bill/L is achieved:
 - i. Fill out the "Equipment Corrective Action Log" on the back of the Platelet Check log sheet.
 - ii. Notify hematology management or call Facility Management (248-551-6300) to place a service call for the centrifuge.
 - iii. Do not use the centrifuge until the problem is resolved.
2. **Centrifuges:** Facility Management (248-551-6300) performs and documents the following centrifuge checks quarterly:
 - a. Timer check.
 - b. Speed check using a tachometer.

IV. REFERENCES:

- A. Guideline for Platelet Poor Plasma Testing: Clinical Laboratory Standards Institute. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays: Approved Guideline-Fourth Edition, CLSI document #H21-A4. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003

Attachments

No Attachments

Approval Signatures

Step Description	Approver	Date
CP Chief Medical Director	Peter Millward: Chief, Pathology Service Line	5/13/2021
Coagulation Medical Director Designee	Marc Smith: System Med Dir, Coagulation	5/13/2021
Policy and Forms Steering Committee Approval (if needed)	Tamara Sabih: Medical Technologist Lead	5/11/2021
Policy and Forms Steering Committee Approval (if needed)	Gail Juleff: Project Mgr Policy	5/11/2021
System Manager	Rebecca Bacarella: Mgr Laboratory	5/11/2021
	Tamara Sabih: Medical Technologist Lead	5/10/2021

Applicability

Royal Oak