

Beaumont

Origination: 11/5/2021
Effective: 11/5/2021
Last Approved: 11/5/2021
Last Revised: 11/5/2021
Next Review: 11/5/2023
Document Contact: *Colette Kessler: Mgr*
Laboratory
Area: *Laboratory-Chemistry*
Key Words:
Applicability: *Royal Oak*

T.S. Refractometer - Royal Oak

Document Type: Procedure

I. PURPOSE AND OBJECTIVE:

The T.S. Meter is a Goldberg Refractometer, which has been designed specifically for use in the medical and paramedical fields. Although the actual measurement is refractive index, the scales of the instrument have been calibrated in terms of specific gravity for urine, and protein concentrations of plasma or serum (gms/dL). The value is read on the appropriate scale as seen through the eyepiece where the sharp boundary between the dark and the light fields crosses the scale.

II. SPECIMEN COLLECTION AND HANDLING:

Determinations require only a drop of fluid sample.

III. QUALITY CONTROL (QC):

- A. A STANDARD SOLUTION OF NACL WILL BE CHECKED MONTHLY AND RECORDED ON QC SHEETS and in UNITY. Dissolve 2.8 grams of dessicated NaCl in a 100 mL volumetric flask containing approximately 50 mL deionized water. Quantity sufficient (QS) to 100 mL with deionized water. Stable indefinitely. Refrigerate.
- B. Expected specific gravity (s.g.) for this QC is 1.014-1.016.

IV. PROCEDURE:

- A. Rotate the instrument to a horizontal position, and place the sample liquid on the exposed portion on top or bottom of the measuring prism, so that the liquid will be drawn into the space between the prisms by capillary action.
- B. To take a reading, keep the cover plate in contact with the prism and point the instrument toward the illuminating light source. The instrument should be tilted with respect to the light source until best results are obtained.
- C. Bring the scale seen in the eyepiece into best focus by rotation of the eyepiece.
- D. Make the reading on the appropriate scale at the point where the dividing line between bright and dark fields crosses the scale.
- E. Use a soft cloth or soft tissue moistened with distilled water for wiping the prism. Dry the prism with a soft

cloth or tissue.

V. CALCULATIONS AND INTERPRETATIONS:

A. If the reading is greater than 1.035, then the specimen should be diluted 1:2 with distilled H₂O. Reread the diluted sample and multiply the last two digits of the specific gravity by 2.

1. Example:

- a. Undiluted specimen = >1.035
- b. Add 0.5 ml distilled H₂O to 0.5 ml specimen.
- c. Mix diluted specimen = $1.020 \times 2 = 1.040$.

VI. REFERENCE RANGE:

1.001-1.030

VII. REPORTABLE RANGE:

Specific Gravity 1.000-1.035

VIII. REFERENCES:

1. American Optical Manual – T.S. Refractometer

Attachments

No Attachments

Approval Signatures

Step Description	Approver	Date
Medical Director	Ann Marie Blenc: System Med Dir, Hematopath	11/5/2021
Policy and Forms Steering Committee Approval (if needed)	Gail Juleff: Project Mgr Policy	11/5/2021
Policy and Forms Steering Committee Approval (if needed)	Colette Kessler: Mgr Laboratory	10/5/2021
Lab Chemistry Best Practice Committee	Qian Sun: Tech Dir, Clin Chemistry, Path	10/1/2021
Lab Chemistry Best Practice Committee	Elizabeth Sykes: System Med Dir, Chemistry	9/24/2021
	Colette Kessler: Mgr Laboratory	9/24/2021

Applicability

Royal Oak