

Beaumont

Origination: 9/29/2021
Effective: 2/3/2022
Last Approved: 2/3/2022
Last Revised: 2/3/2022
Next Review: 2/3/2024
Document Contact: *Colette Kessler: Mgr
Laboratory*
Area: *Laboratory-Chemistry*
Key Words:
Applicability: *All Beaumont Hospitals*

Iricell Shutdown and Startup Procedure

Document Type: Procedure

I. PURPOSE AND OBJECTIVE:

The Iricell 2000 may require a system restart or shutdown to correct errors encountered while analyzing specimens. This procedure describes the appropriate way to restart or shutdown. It has the steps for the Velocity and the iQ200

II. PROCEDURE:

- A. **Velocity and iQ 200 Shutdown:** When it is necessary to restart the instrument, shut down the instrument in the following order:
1. Select Instrument on the top right of the main screen
 2. Select **Go Off Line**. The system will prompt: *Continue going off line?*
 3. Select **Yes**
 4. Select **Maintenance** located at the bottom of the Instrument Screen
 5. Select **Shutdown**. The system will prompt: *Do you want to shut down the instrument?*
 6. Select **Yes**. Windows will close
 7. Turn OFF the iChem VELOCITY and/or iQ200 by pressing the ON/OFF button
 8. **NOTE:** Restart the Velocity and iQ200 Modules as soon as possible. Do not allow the instrument to cool down or the measurements will be affected.
 9. To restart the computer, open the flap on the front of the main computer
 10. Press the Power button.
 11. The software will automatically load after about 1 minute
 12. Log onto the system.
 13. Select **Go online**.
- B. **Velocity and iQ200 Restart:** This button is used to shut down and to restart only the user interface software, without any effect on the systems. This function must be used after software settings have been modified, communication settings, or error message requiring a restart of the instrument
1. Select Instrument on the top right of the main screen

2. Select **Go Off Line**.
3. Select **Maintenance** located at the bottom of the Instrument Screen
4. From the Maintenance screen, select **Restart**. The system will prompt: *Do you want to restart the instrument?*
5. Select **Yes**
6. Log onto the system.
7. Select **Go online**.

Attachments

No Attachments

Approval Signatures

Step Description	Approver	Date
CLIA Medical Directors	Jeremy Powers: Chief, Pathology	2/3/2022
CLIA Medical Directors	Vaishali Pansare: Chief, Pathology	2/1/2022
CLIA Medical Directors	Muhammad Arshad: Chief, Pathology	2/1/2022
CLIA Medical Directors	John Pui: Chief, Pathology	2/1/2022
Policy and Forms Steering Committee Approval (if needed)	Colette Kessler: Mgr Laboratory [RC]	2/1/2022
Policy and Forms Steering Committee Approval (if needed)	Gail Juleff: Project Mgr Policy	2/1/2022
	Ann Marie Blenc: System Med Dir, Hematopath	2/1/2022
	Elizabeth Sykes: System Med Dir, Chemistry	1/31/2022
	Qian Sun: Tech Dir, Clin Chemistry, Path	1/31/2022
	Kimberly Geck: Dir, Lab Operations B	1/31/2022
	Amy Conners: Dir, Lab Operations A	1/31/2022
	Colette Kessler: Mgr Laboratory [RC]	1/31/2022
	Colette Kessler: Mgr Laboratory [RC]	1/31/2022

Applicability

Dearborn, Farmington Hills, Grosse Pointe, Royal Oak, Taylor, Trenton, Troy, Wayne