

Beaumont

Origination 4/1/2021
 Last Approved 1/26/2023
 Effective 1/26/2023
 Last Revised 1/26/2023
 Next Review 1/25/2025

Document Contact Kelly Sartor:
 Supv, Laboratory
 Area Laboratory-Blood
 Bank
 Applicability All Beaumont
 Hospitals

ORTHO VISION Analyzer Manual Card Review

Document Type: Policy

I. PURPOSE AND OBJECTIVE:

The ORTHO VISION™ has been configured to send gel cards to the service rack to help the technologists identify those gel cards that need to be reviewed before results are sent across the interface to the laboratory information system (LIS). The purpose of this document is to provide guidelines for extracting, identifying, and managing gel cards that are sent to the Manual Review Rack by the ORTHO VISION™.

II. POLICY STATEMENT:

Any gel cards sent to the Manual Review Rack by the ORTHO VISION™ must be reviewed before results are sent across the interface to the Blood Bank computer system. The ORTHO VISION™ sends gel cards to the Manual Review Rack for several reasons, including but not limited to:

- A. ABO/Rh cards with a 1+ or 2+ graded reaction in the forward ABO typing, or a 1+ graded reaction in the reverse typing.
- B. ABO/Rh cards with a 1+, 2+, or 3+ graded reaction in the Rh(D) typing.
- C. Antibody screen (ABSC) cards with any positive graded reaction.
- D. Any gel card with a well discrepancy which are not interpreted by the ORTHO VISION™.
- E. Any gel card with a positive crossmatch of any strength.
- F. All cards will be sent to the manual review rack if a required maintenance is due but has not been performed.
- G. Gel cards that the ORTHO VISION™ cannot read; the *Card Reader Error* warning message appears.
- H. All panel cards.

III. DEFINITIONS / ACRONYMS:

- A. **LIS:** Laboratory information system.

B. **MT Lead:** Medical Technologist Lead.

C. **Manual Review Rack:** Also referred to by the ORTHO VISION™ as the service rack.

IV. INSTRUCTIONS:

A. **Technologist Must Review / Modify Graded Reactions of all Gel Cards Sent to the Manual Review Rack.**

1. All gel cards that are sent to the Manual Review rack shall be reviewed by a technologist; the gel card itself shall be visually reviewed. If the technologist agrees with the graded reaction they may accept the result and send to the LIS. If the technologist grades a reaction differently than the ORTHO VISION™, then they shall modify the graded reaction in the ORTHO VISION™ accordingly before the results are sent across the interface to Soft. Reactions are graded as described in Transfusion Medicine policy, [Reading, Grading, and Recording Test Reactions](#).
2. Note: Reaction grading can only be modified **once** on the ORTHO VISION™, it is unable to be modified again. If addition modification is required, you must modify the reaction grading in the Blood Bank computer system prior to verifying the results.

B. **Valid Graded Reactions and Interpretations.**

1. Based on Blood Bank procedures, graded ABO and Rh(D) reactions must be valid in order to interpret the results. The graded reactions must be of specific strengths to be considered valid by the Blood Bank procedures. Note that the ORTHO VISION™, Blood Bank computer system, and the Blood Bank procedures interpret results differently. All tests must be interpreted from the actual gel card based on the Blood Bank procedures, not necessarily how the Blood Bank computer system or the ORTHO VISION™ interprets them.

C. **Gel Cards with Positive Crossmatch Results.**

1. The ORTHO VISION™ is configured to send gel cards with any positive crossmatch reaction strength to the Manual Review Rack. As with all gel cards that are sent to the service rack, the graded reactions shall be reviewed by a medical technologist. If the technologist grades a reaction differently than the ORTHO VISION™ then the technologist shall modify the graded reaction in the ORTHO VISION™ accordingly before the results are sent across the interface to the Blood Bank computer system. Refer to Transfusion Medicine policy, [Investigation of Incompatible Crossmatches](#).

D. **ABO/Rh Gel Cards that the ORTHO VISION™ Interprets as NRD (No Results Determined).**

ABO Example						
Anti-A	Anti-B	Anti-D	Control	A Cells	B Cells	Interpretation
0	0	4+	0	1+	0	? Pos
Rh(D) Example						
Anti-A	Anti-B	Anti-D	Control	A Cells	B Cells	Interpretation
0	0	4+	1+	4+	4+	O ?

1. After the technologist extracts the gel card from the Manual Review Rack, the technologist will review the gel cards and, if applicable, modify the results.
2. The technologist will attempt to resolve the ABO/Rh discrepancy as described in

Transfusion Medicine policy, [Resolution of ABO and Rh\(D\) Discrepancies](#).

3. If the technologist is unable to resolve the discrepancy, then a supervisor, MT Lead or a technologist with appropriate access who can interpret the results as GND (Group Not Determined) or RND (Rh Not Determined).
4. If the technologist is able to resolve the discrepancy, then:
 - a. Print a copy of the results by Show Order Report.
 - b. Reject the results on the ORTHO VISION™.
 - c. The printed copy of results from the ORTHO VISION™ will be used to document the results in the Blood Bank computer system.

E. ABO/Rh Gel Cards with Graded Reactions that are Interpreted by the ORTHO VISION™, but that should not be Interpreted Based on Blood Bank Procedures.

ABO Example						
Anti-A	Anti-B	Anti-D	Control	A Cells	B Cells	Interpretation
1+	0	4+	0	0	1+	A Pos
Rh(D) Example (see note below)*						
Anti-A	Anti-B	Anti-D	Control	A Cells	B Cells	Interpretation
0	0	2+	0	4+	4+	O Pos

1. Although considered invalid by the Blood Bank procedures, these results are interpreted by the ORTHO VISION™. The technologist should not send these results across the interface to the Blood Bank computer system.
The technologist will:
 - a. Print a copy of the results by Show Order Report.
 - b. Reject the results on the ORTHO VISION™.
 - c. Investigate the discrepancy as described in Transfusion Medicine policy, [Resolution of ABO and Rh\(D\) Discrepancies](#).
 - d. The Results by Show Order Report will be used to document the results in the Blood Bank computer system.
2. If the technologist is unable to resolve the discrepancy, or if the technologist does not have the required computer access to interpret the results in the Blood Bank computer system, a supervisor, MT Lead or a technologist who has this access will interpret the results as GND or RND.
3. Note: The Rh(D) Example, above:
 - a. Based on Transfusion Medicine policy, [Resolution of Rh\(D\) Discrepancies](#), the 2+ graded reaction for the Rh(D) is not a valid graded reaction. The ORTHO VISION™ will interpret this result as Rh pos, but the Blood Bank computer system will interpret this result as weak D pos. This Rh discrepancy must be investigated as described in Transfusion Medicine policy, [Resolution of Rh\(D\) Discrepancies](#). Although the most likely cause for this 2+ reaction is that the patient is weak D or partial D positive, there are many other potential causes. Examples of other causes include recent Rh(D) dissimilar RBC transfusions, a bone marrow transplant, a mistyped sample, and a wrong blood in tube event

(WBIT).

F. Gel Cards with Mixed-Field Reactions.

1. Mixed-field reactions are referred to as MF reactions by the ORTHO VISION™. MF reactions are considered a well discrepancy; however they can still be interfaced into the Blood Bank computer system. MF reactions may occur in patients who have been recently transfused; refer to Transfusion Medicine policy, [Resolution of ABO and Rh\(D\) Discrepancies for Recently Transfused Patients](#).
 - a. Accept the MF results on the ORTHO VISION™ and click “Send to LIS.”
 - b. Print a copy of the results by “Show Order Report” under the Results tab.
 - c. Give the results to a supervisor, MT Lead or technologist who has access to investigate the discrepancy.
 - d. The results are interpreted as described in Transfusion Medicine policy, [Resolution of ABO and Rh\(D\) Discrepancies for Recently Transfused Patients](#).

G. Gel Cards with Well Discrepancies (Codes).

1. The ORTHO VISION™ is configured to send gel cards with well discrepancies to the Manual Review Rack. Results with well discrepancies are not interpreted by ORTHO VISION™ and do not interface to the Blood Bank computer system.
 - a. As with all gel cards that are sent to the Manual Review Rack, the graded reactions shall be reviewed by a medical technologist. If the technologist grades a reaction differently than the ORTHO VISION™ then the technologist shall modify the graded reaction in the ORTHO VISION™ accordingly.
 - b. If the technologist is unable to resolve the discrepancy, or if the technologist does not have the required computer access to interpret the results in Soft, the technologist will print the Results by Show Order Report from the ORTHO VISION™ and submit it to a supervisor, MT Lead or technologist who has this access.
 - c. See attachment *ORTHO VISION™ Analyzer System Flags and Codes*, for a listing of codes along with the conditions and suggested actions.

V. PROCEDURE:


A. To Change Results:

1. Retrieve the card from the Manual Review Rack. Note that when opening the Dual Purpose Drawer, all cards that are in the Manual Review portion must be removed at that time.
 - a. Resources.
 - b. Manual Load Review.
 - c. Load-Unload.
2. Edit the Reaction Grade.
 - a. Touch Results.
 - b. Select the test to be edited.

c. Show Details.

d. If there is more than one test ordered on that sample, Select Test box on left side of screen.

Identify the card needed for review by matching unique card ID number to card ID number on the ORTHO VISION™ screen.



Unique card ID number is the last 5 digits of the card

e. Touch Edit Grades.

f. Scan card barcode.

g. Select reaction grade to be edited.

h. Change the reaction strength. Note that all reactions are graded by inspecting the gel card, not using the ORTHO VISION™ screen.

i. Touch Next.

j. Add comment (if necessary).

k. Touch Next.

l. Input password.

m. Touch Confirm Password.


3. Accept the results.

a. The results will not cross the interface to the LIS until they have been accepted.

4. Send the results to the LIS.
 - a. If a Type & Screen is ordered on the specimen, both card results (the type and the antibody screen card) will need to be sent to the LIS.
5. Verify the results in the Blood Bank computer system.

B. To Review Results That Don't Need Editing:

1. Retrieve the card from Manual Review Rack. Note that when opening the Dual Purpose Drawer, all cards that are in the Manual Review portion must be removed at that time.
 - a. Touch Resources.
 - b. Touch Manual Load Review.
 - c. Touch Load-Unload.
2. View the Reaction / Well.
 - a. Touch Results.
 - b. Select the test to be reviewed.
 - c. Touch Show Details.
 - d. If there is more than one test ordered on that sample, Select Test box on left side of screen.
 - e. Identify the card needed for review by matching unique card ID number to card ID number on the ORTHO VISION™ screen.



Unique card ID number is the last 5 digits of the card

- f. If there is agreement with the results in the well, touch accept.
 - g. If there is non-agreement with the wells in question, edit the wells. See the *To Change Results* procedure section above.
3. Accept the results.
 - a. The results will not cross the interface to the LIS until they have been accepted.
4. Send the results to the LIS.
 - a. If a Type & Screen is ordered on the specimen, both card results (the type and the antibody screen card) will need to be sent to the LIS.
5. Verify the results in the Blood Bank computer system.

VI. REFERENCES:

- A. AABB *Standards for Blood Banks and Transfusion Services*, current edition.
- B. Cohn, C.S., Delaney, M, Johnson, S.T., Katz, L.M. (2020) *Technical Manual*.(19th ed.). AABB.
- C. Ortho Clinical Diagnostics, Rochester, NY, Ortho Vision General Operator Training Manual, Publication J56102.
- D. ORTHO VISION® Analyzer ID-MTS Gel Cards Reference Guide J40050.
- E. ORTHO VISION® Analyzer ID-MTS Gel Cards Self-Service Customer Procedures Guide J40055ENNA.
- F. ORTHO VISION® Analyzer Electronic Library, Software version 5.13.2.0
- G. Ortho Clinical Diagnostics, Rochester, NY, Electronic Publication number J56102.
- H. ID-Micro Typing System® Implementation Guide 6902200.

Attachments

[ORTHO VISION Analyzer System Flags and Codes \(02/09/2021\)](#)

Approval Signatures

Step Description	Approver	Date
Policy and Forms Steering Committee (if needed)	Kristina Davis: Staff Physician	1/26/2023
	Jeremy Powers: Chief, Pathology	12/30/2022
	Vaishali Pansare: Chief, Pathology	12/29/2022
	Ryan Johnson: OUWB Clinical Faculty	12/28/2022
	Muhammad Arshad: Physician	12/27/2022
	Ann Marie Blenc: System Med Dir, Hematopath	12/21/2022
	John Pui: Chief, Pathology	12/21/2022
	Kelly Sartor: Supv, Laboratory	12/21/2022
Policy and Forms Steering Committee (if needed)	Gail Juleff: Project Mgr Policy	12/21/2022

Rebecca Thompson: Medical Technologist Lead	12/21/2022
Katherine Persinger: Mgr Laboratory	12/12/2022
Teresa Lovins: Supv, Laboratory	12/12/2022
Hilary Morey: Medical Technologist Lead	12/9/2022
Kelly Sartor: Supv, Laboratory	12/8/2022
Michele Ferla: Medical Technologist Lead	12/8/2022
Kristen Lafond: Mgr Laboratory	12/7/2022
Michael Rasmussen: Supv, Laboratory	12/7/2022
Ashley Beesley: Mgr Laboratory	12/7/2022
Karrie Torgerson: Supv, Laboratory	12/7/2022
Brooke Klapatch: Medical Technologist Lead	12/7/2022
Kelly Sartor: Supv, Laboratory	12/7/2022

