

Beaumont

Origination 11/9/2021
 Last Approved 3/1/2023
 Effective 3/1/2023
 Last Revised 3/1/2023
 Next Review 2/28/2025

Document Contact Colette Kessler:
 Mgr, Division Laboratory
 Area Laboratory-Chemistry
 Applicability All Beaumont Hospitals

Formulas for Manual Calculations

Document Type: Guideline

I. PURPOSE AND OBJECTIVE:

- A. The purpose of this procedure is to specify formulas needed to manually calculate results.

II. FORMULAS:

- A. Anion Gap = Sodium – (Chloride + Carbon dioxide)
 B. Globulin = Total Protein – Albumin
 C. A/G Ratio = Albumin / (Total Protein – Albumin)
 D. Cholesterol/High Density Lipoprotein Ratio = Chol/HDL
 E. Non HDL Cholesterol = Chol – HDL = mg/dL
 F. LDL Cholesterol = Total Chol – HDL – (Trig/5)
 Note: Formula for calculation only if Trig<400 mg/dL
 G. BUN/Creatinine Ratio = BUN/Creatinine
 H. CKD-EPI Estimated GFR (mL/min/1.73 m²) = $142 \times \min(S_{Cr}/\kappa, 1)^{\alpha} \times \max(S_{Cr}/\kappa, 1)^{-1.200} \times 0.9938^{Age} \times 1.012$ [if female]
1. S_{Cr} is serum creatinine in mg/dL
 2. κ is 0.7 for females and 0.9 for males
 3. α is -0.241 for females and -0.302 for males
 4. min indicates the minimum of S_{Cr}/κ or 1
 5. max indicates the maximum of S_{Cr}/κ or 1

- a. Example for a female patient: $142 \times \min(S_{Cr}/0.7, 1)^{-0.241} \times \max(S_{Cr}/0.7, 1)^{-1.200} \times 0.9938^{Age} \times 1.012$
- b. Example for a male patient: $142 \times \min(S_{Cr}/0.9, 1)^{-0.302} \times \max(S_{Cr}/0.9, 1)^{-1.200} \times 0.9938^{Age}$
- I. CKD-EPI Cystatin C Equation (2012) Estimated GFR (mL/min/1.73 m²) = $133 \times \min(S_{cys}/0.8, 1)^{-0.499} \times \max(S_{cys}/0.8, 1)^{-1.328} \times 0.9936^{Age} \times 0.932$ [if female]
1. S_{cys} is standardized serum cystatin C in mg/L
 2. min indicates the minimum of $S_{cys}/0.8$ or 1
 3. max indicates the maximum of $S_{cys}/0.8$ or 1
 4. age is years
- J. CKD-EPI Cystatin C Equation Estimated GFR Under 25 (mL/min/1.73 m²) = $K \times (1/S_{cys})$
1. S_{cys} is standardized serum cystatin C in mg/L
 2. Male 1-15 years $K = 87.2 \times 1.011^{(age\ years-15)}$
 3. Male 15-18 years $K = 87.2 \times 0.960^{(age\ years-15)}$
 4. Female 1-12 years $K = 79.9 \times 1.004^{(age\ years-12)}$
 5. Female 12-18 years $K = 79.9 \times 0.974^{(age\ years-12)}$
- K. Estimated GFR (mL/min/1.73 m²) = $133 \times \min(S_{cys}/0.8, 1)^{-0.499} \times \max(S_{cys}/0.8, 1)^{-1.328} \times 0.9936^{Age} \times 0.932$ [if female]
- L. % Saturation = $100 \times [\text{serum Iron} / (\text{Transferrin} \times 1.4)]$
- M. Total Iron Binding Capacity = Transferrin $\times$ 1.4
- N. Urine Albumin/Creatinine Ratio = (ualb mg/dL/Urine Creatinine mg/dL) $\times$ 1000
- O. Urine Creatinine Clearance mL/min

$$\frac{\text{Raw Urine Creatinine mg/dL} \times (\text{TV mL} \times 1 \text{ Hr})}{\text{Serum Creatinine mg/dL} \times \text{Hrs collection} \times 60 \text{ min}}$$
- P. Urine (Ur) Creatinine Clearance Pediatric mL/min CC Corrected =

$$\frac{\text{Raw Ur Creat mg/dL} \times (\text{TV mL} \times 1 \text{ Hr}) \times (1.73 \text{ m}^2)}{\text{Serum Creat mg/dL} \times \text{Hrs collection} \times 60 \text{ min} \times \text{surface area m}^2}$$
 NOTE: Surface area m² is determined from child's height and weight.
- Q. Urine Creatinine Timed or 24 Hour mg/Collection (Coll) =

$$\frac{\text{Raw Urine Creatinine mg}}{\text{dL}} \times \text{dilution} \times \frac{(\text{TV mL} \times 10 \text{dL})}{\text{Coll.} \times 1000 \text{mL}}$$
- R. Urine Amylase Timed Ur/Hr =

$$\frac{\text{Raw Urine Amylase U}}{\text{L}} \times \text{dilution} \times \frac{\text{TV mL} \times 1 \text{L}}{\text{Hrs of Coll.} \times 1000 \text{mL}}$$
- S. Urine Calcium Timed or 24 Hour, mg/Collection =

$$\frac{\text{Raw Urine Calcium mg}}{\text{dL}} \times \text{dilution} \times \frac{(\text{TV mL} \times 10 \text{dL})}{\text{Coll.} \times 1000 \text{mL}}$$
- T. Urine Glucose Timed and 24 Hour, g/Collection =

- Raw Urine Glucose $\frac{\text{mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right) \times \frac{\text{g}}{1000\text{mg}}$
- U. Urine Magnesium Timed or 24 Hour, mg/Collection =
 $\frac{\text{Raw Magnesium mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right)$
- V. Urine Phosphorous Times or 24 Hour, g/Collection =
 $\left(\frac{\text{Raw Urine Phos mg}}{\text{dL}} \times \frac{1\text{g}}{1000\text{mg}} \right) \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right)$
- W. Urine Protein Timed or 24 Hour, mg/Collection =
 $\frac{\text{Raw Urine Protein mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right)$
- X. Urine Urea Nitrogen Timed or 24 Hour, mg/Collection =
 $\frac{\text{Raw Urine Urea Nitrogen mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right)$
- Y. Urine Uric Acid Timed or 24 Hour, mg/Collection =
 $\frac{\text{Raw Urine Uric Acid mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right)$
- Z. Urine Sodium Timed or 24 Hour, mmol/Collection =
 $\frac{\text{Raw Urine Sodium mmol}}{\text{L}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{1\text{L}}{1000\text{mL}} \right)$
- AA. Urine Potassium Timed or 24 Hour, mmol/Collection =
 $\frac{\text{Raw Urine Potassium mmol}}{\text{L}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{1\text{L}}{1000\text{mL}} \right)$
- AB. Urine Chloride Timed or 24 Hour, mmol/Collection =
 $\frac{\text{Raw Urine Chloride mmol}}{\text{L}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{1\text{L}}{1000\text{mL}} \right)$
- AC. Urine Albumin excretion, mcg/min =
 $\frac{\text{Raw Urine Albumin mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right) \times 1000 \text{ mcg} \times \left(\frac{\text{Collection X 1 hr}}{\text{X 60 min}} \right)$
- AD. Urine Albumin timed or 24 hour, mg/Collection=
 $\frac{\text{Raw Urine Albumin mg}}{\text{dL}} \times \text{dilution} \times \left(\frac{\text{TV mL}}{\text{Coll.}} \times \frac{10\text{dL}}{1000\text{mL}} \right)$
- AE. % Urine Amylase / Creatinine Clearance Ratio =
 $\frac{\text{Raw Urine Amylase U/L}}{\text{Serum Amylase U/L}} \times \frac{\text{Serum Creatinine mg/dL}}{\text{Raw Urine Creatinine mg/dL}} \times 100$
 NOTE: reference range 1.3-4.3%
- AF. Urine Protein / Creatinine Ratio =
 $\frac{\text{Raw Urine Protein mg/dL}}{\text{Raw Urine Creatinine mg/dL}}$
 NOTE: reference range 0.0-0.2

AG. Urine Calcium / Creatinine Ratio =

$$\frac{\text{Raw Urine Calcium mg/dL}}{\text{Raw Urine Creatinine mg/dL}}$$

NOTE: reference range 0.02-0.26

AH. AAPO₂, mmHg = ((713mmHg x FiO₂ as decimal) – (PCO₂ / 0.8)) - PO₂

AI. PF, no units = PO₂ / FiO₂ (as a decimal value)

AJ. % Free PSA = (Free PSA/PSA) X 100

AK. Albumin gradient; Serum and Ascites fluid, g/dL = Albumin (g/dL) – Albumin fluid (g/dL)

NOTE: >1.1 = Transudate

0.2-1.0 = Exudate

Attachments

[Formulas for Manual Calculations CalculationVerification Worksheet Attachment A_1_2023.pdf](#)

Approval Signatures

Step Description	Approver	Date
Medical Directors	Muhammad Arshad: Physician	3/1/2023
Medical Directors	Vaishali Pansare: Chief, Pathology	2/10/2023
Medical Directors	Jeremy Powers: Chief, Pathology	2/8/2023
Medical Directors	Ryan Johnson: OUWB Clinical Faculty	2/2/2023
Medical Directors	Ann Marie Blenc: System Med Dir, Hematopath	2/1/2023
Medical Directors	John Pui: Chief, Pathology	2/1/2023
Policy and Forms Steering Committee Approval (if needed)	Colette Kessler: Mgr, Division Laboratory	2/1/2023
Policy and Forms Steering Committee Approval (if needed)	Ilene Hirsch: Project Mgr Policy	2/1/2023
	Caitlin Schein: Staff Physician	2/1/2023

Qian Sun: Tech Dir, Clin Chemistry, Path	1/31/2023
Nga Yeung Tang: Tech Dir, Clin Chemistry, Path	1/24/2023
Elzbieta Wysteppek: Dir, Lab Operations B	1/24/2023
Kimberly Geck: Dir, Lab Operations B	1/23/2023
Brittnie Berger: Dir, Lab Operations C	1/23/2023
Colette Kessler: Mgr, Division Laboratory	1/23/2023
Colette Kessler: Mgr, Division Laboratory	1/23/2023

COPY