

Beaumont

Origination 2/8/2023
Last Approved 4/17/2023
Effective 4/17/2023
Last Revised 4/17/2023
Next Review 4/16/2025

Document Contact Colette Kessler:
Mgr, Division Laboratory
Area Laboratory-Chemistry
Applicability Royal Oak

Bilirubin In Urine (Multistix 10 SG) - Royal Oak

Document Type: Procedure

I. PURPOSE AND OBJECTIVE:

- A. Bilirubin is derived from the degradation of hemoglobin and other heme containing proteins. It is conjugated in the liver by the enzyme glucuronyl transferase to a water-soluble form. This conjugated bilirubin can be excreted via the kidneys, and thus can appear in the urine. Bilirubinuria occurs in intra-and extra hepatic obstructive jaundice, acute/chronic hepatitis, cirrhosis, and Dubin-Johnson and Rotor syndromes. The urine is dark (amber) and may have a yellow foam. Bilirubinuria is associated with elevated serum bilirubin (conjugated).
- B. The Multistix test is based on the coupling of bilirubin with diazotized dichloroaniline in a strongly acid medium. The intensity of the tan color produced increases with increasing concentration of bilirubin.
- C. This document describes the steps for this procedure to assist technologists.

II. SPECIMEN COLLECTION AND HANDLING:

Fresh, well-mixed, uncentrifuged urine. It is recommended that testing be done within one hour after voiding. Otherwise immediately refrigerate the specimen and return to room temperature before testing.

III. REAGENTS:

- A. Siemens Multistix 10 SG (#2161):
- B. 0.4% weight for weight (w/w) 2,4-dichloroaniline diazonium salt
- C. 37.3% w/w buffer
- D. 62.3% w/w nonreactive ingredients

IV. QUALITY CONTROL (QC):

- A. Both Normal and Abnormal Kova-Trols are run and results are recorded:
 - 1. at the beginning of each shift
 - 2. whenever a new lot number of reagent strips is opened
 - 3. whenever a new shipment of reagent strips is received
 - 4. whenever troubleshooting warrants it

V. PROCEDURE:

- A. Briefly dip the test area of the strip in fresh, **well-mixed** uncentrifuged urine.
- B. While removing the strip, run the edge against the rim of the urine container to remove excess urine. Hold the strip in a horizontal position to prevent mixing of chemicals from adjacent reagent areas and/or contaminating the hands with urine.
 - 1. If reading visually, compare the **BILIRUBIN** reagent area to the corresponding Color Chart on the bottle label at **30 seconds**. Hold strip close to color blocks and match carefully.
 - 2. If reading instrumentally, follow directions given in the Clinitek Advantus procedure.
- C. NOTE: The Clinitek Advantus appears to over-read the Bilirubin test pad giving false positives especially in the presence of brown colored urine.

VI. REFERENCE RANGE:

- A. Negative
- B. Normal adult urine contains about 0.02 mg bilirubin/dL. This is not detectable in urine by the usual methods. Even trace amounts of detectable bilirubin are sufficiently abnormal to require further investigation.

VII. SENSITIVITY:

The test pad on the Multistix reagent strip is sensitive from 0.4 – 0.8 mg/dL for bilirubin.

VIII. REPORTABLE RANGE:

Multistix 10SG (backup) has a color comparison chart with **FOUR** color blocks of various tan shades. These represent bilirubin as negative or present in increasing amounts. Results are reported as follows:

Negative	Negative
1+	Positive
2+	Positive
3+	Positive

IX. LIMITATIONS/INTERFERING SUBSTANCES:

- A. Drug metabolites that color the urine red/orange (e.g. Indican) either give false positive results or mask the presence of bilirubin.
- B. Metabolites of Lodine (etodolac) may cause positive or atypical results.
- C. Ascorbic acid concentrations > 15 mg/dL may cause false negatives.
- D. Oxidation of bilirubin in specimens that have stood too long, especially when exposed to light, will yield false negative findings.

X. REFERENCES:

- 1. Multistix 10 SG. Miles, Inc. Diagnostic Division, Elkhart, IN 46515, rev. 04/99
- 2. Henry, J.B., Clinical Diagnosis and Management by Laboratory Methods, 20th edition, Philadelphia, W.B. Saunders Co., 2001 p. 376-378.
- 3. Hundley, J.M. and Fleming, J.K., Urine Analysis American Society of Clinical Pathologists Workshop, Dearborn MI, 1991.

Approval Signatures

Step Description	Approver	Date
Medical Director	Ann Marie Blenc: System Med Dir, Hematopath	4/17/2023
Policy and Forms Steering Committee Approval (if needed)	Colette Kessler: Mgr, Division Laboratory	4/5/2023
Policy and Forms Steering Committee Approval (if needed)	Gail Juleff: Project Mgr Policy	4/5/2023
Lab Chemistry Best Practice Committee	Caitlin Schein: Staff Physician	4/5/2023
Lab Chemistry Best Practice Committee	Qian Sun: Tech Dir, Clin Chemistry, Path	3/24/2023
	Colette Kessler: Mgr, Division Laboratory	3/24/2023