

Beaumont

Origination 1/5/2023
Last 1/5/2023
Approved
Effective 1/5/2023
Last Revised 1/5/2023
Next Review 1/4/2025

Document Sharon Scalise:
Contact Supv, Laboratory
Area Laboratory-
Histology
Applicability Royal Oak

Histology Special Stain - Fouchet - Royal Oak

Document Type: Procedure

I. PURPOSE AND OBJECTIVE:

Bile and lipofuchsin appear as brown-yellow pigments. The purpose of this document is to provide a procedure for the demonstration of bilirubin (bile pigment) in tissue sections. Red blood cells, when destroyed, release their hemoglobin. This hemoglobin is broken down into globin (protein) and heme (iron). When the iron is removed from the heme portion, the residue is known as biliverdin. The biliverdin is transported to the liver where it is reduced to bilirubin, where it passes into the gallbladder, then into the duodenum. In the cases of biliary obstruction or extensive liver damage, bilirubin may build up in the liver.

II. PRINCIPLE:

Bile and lipofuchsin appear as brown-yellow pigments. This stain demonstrates bilirubin (bile pigment) in tissue sections. Red blood cells, when destroyed, release their hemoglobin. This hemoglobin is broken down into globin (protein) and heme (iron). When the iron is removed from the heme portion, the residue is known as biliverdin. The biliverdin is transported to the liver where it is reduced to bilirubin, where it passes into the gallbladder, then into the duodenum. In the cases of biliary obstruction or extensive liver damage, bilirubin may build up in the liver.

III. SPECIMEN COLLECTION AND HANDLING:

A. Fixation

1. Any well-fixed tissue, except Zenker.
2. 10% neutral buffered formalin preferred.

- B. Processing
 - 1. Standard, overnight processing.
- C. Section Thickness
 - 1. Cut paraffin sections at 5 μ .
- D. Slide Drying
 - 1. 30 minutes at 60°C.
- E. Type of Slide
 - 1. Plain slides.

IV. REAGENTS:

A. 10% Ferric Chloride

Ferric chloride	10.0 gm
Distilled water	100.0 mL
OR	
29% Stock ferric chloride	3.5 mL
Distilled water	6.5 mL

Stir together with magnetic stirrer. Store at room temperature; stable for 1 month.

B. Fouchet Reagent

Trichloroacetic acid	25.0 gm
Distilled water	100.0 mL
10% ferric chloride	10.0 mL

JUST BEFORE USE, carefully add trichloroacetic acid to distilled water. Stir until dissolved. Add 10% ferric chloride. Filter before use.

C. 1% Acid Fuchsin

Acid fuchsin	1.0 gm
Distilled water	100.0 mL

Dissolve together. Store at room temperature; stable for months.

D. Saturated Aqueous Picric Acid

Use vendor pre-made saturated aqueous picric acid.

E. Van Gieson Counterstain

1% acid fuchsin	5.0 mL
Saturated aqueous picric acid	100.0 mL

Mix together. Store at room temperature; stable for 4-6 months; may be reused until weak.

V. EQUIPMENT:

- A. Balance
- B. Magnet stirrer

VI. SUPPLIES:

- A. Erlenmeyer flasks
- B. Graduated cylinders
- C. Non-metallic forceps

VII. SPECIAL SAFETY PRECAUTIONS:

- A. Trichloroacetic Acid
 - 1. Is an acid.
 - 2. Add slowly to water.
 - 3. May cause severe skin and eye burns.
 - 4. May be irritating to eyes and respiratory system.
- B. Ferric Chloride
 - 1. Is a corrosive.
 - 2. May cause skin and eye burns.
 - 3. be irritating to respiratory tract.
- C. Acid Fuchsin
 - 1. Is an irritant.
- D. Picric Acid
 - 1. Is toxic, highly reactive (4 - NFPA) and an extreme fire hazard (4 - NFPA).

VIII. QUALITY CONTROL (QC):

Section of tissue with bile.

IX. PROCEDURE:

Step	Action	Time	Notes
1	Deparaffinize and hydrate sections through graded alcohol to distilled water.		
2	Place slides in freshly filtered Fouchet reagent.	5 minutes	Make Fouchet reagent just before use and filter.
3	Rinse in distilled water, 3-5 changes.	5-10 seconds	
4	Stain in Van Gieson counterstain.	1 minute	
5	Place slides directly into 95% alcohol, 2 changes each.	5-10 seconds	Placing the slides in water may decrease the amount of counterstain.
6	Dehydrate through graded alcohols, clear	10	

	with xylene.	seconds each	
7	Coverslip.		

X. RESULTS:

- A. Bile or bilirubin - **shades of green**
- B. Collagen - **red**
- C. Cytoplasm, red blood cells - **yellow**

XI. REFERENCES:

- A. Bancroft JD, Stevens A: Theory and Practice of Histological Techniques, 3rd ed. New York, NY, Churchill Livingstone, 1990.
- B. Carson FL: Histotechnology: A Self-Instructional Text, Chicago, IL, ASCP Press, 1990. (Source of procedure.)
- C. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd edition. Columbus, Ohio, Battelle Press, 1980.
- D. Vacca L: laboratory Manual of Histochemistry. New York, NY, Raven Press, 1985.

Approval Signatures

Step Description	Approver	Date
Medical Director	Kurt Bernacki: System Med Dir, Surgical Path	1/5/2023
Policy and Forms Steering Committee (if needed)	Gail Juleff: Project Mgr Policy	1/5/2023
Policy and Forms Steering Committee (if needed)	Sharon Scalise: Supv, Laboratory	1/4/2023
	Amy Knaus: Dir, Lab Operations C	1/4/2023
	Jennifer Lehmann: Mgr Laboratory	1/3/2023
	Sharon Scalise: Supv, Laboratory	12/19/2022

Applicability

Royal Oak

COPY