

Beaumont

Origination 1/5/2023
Last 1/5/2023
Approved
Effective 1/5/2023
Last Revised 1/5/2023
Next Review 1/4/2025

Document Sharon Scalise:
Contact Supv, Laboratory
Area Laboratory-
Histology
Applicability Royal Oak

Histology Special Stain - Modified Churukian-Schenk - Royal Oak

Document Type: Procedure

I. PURPOSE AND OBJECTIVE:

The purpose of this document is to provide a procedure for the demonstration of argyrophilic granules in carcinoid tumors and neurosecretory cells and tumors. All argentaffin granules will also be stained by this manner.

II. PRINCIPLE:

Tissue components that are argyrophilic will bind silver ions from a silver nitrate solution. To reduce the silver ions to black metallic silver, a separate reducing solution must be used. Silver nitrate is the source of silver ions. Sodium sulfite and hydroquinone reduce the silver ions to a black metallic precipitate. Nuclear fast red is the counterstain.

III. SPECIMEN COLLECTION AND HANDLING:

A. Fixation

1. 10% neutral buffered formalin preferred.
2. Avoid the use of mercuric fixatives

B. Processing

1. Standard, overnight processing.

C. Section Thickness

1. Cut paraffin sections at 5 μ .

- D. Slide Drying
 - 1. 30 minutes at 60°C.
- E. Type of Slide
 - 1. Plain slides

IV. REAGENTS:

A. Stock 10% Silver Nitrate

Silver nitrate	10.0 gm
Distilled water	100.0 mL

In an acid-cleaned flask, dissolve together with the aid of a magnetic stirrer. Store in a brown glass bottle in the refrigerator(3°C); stable for 4-6 months. Discard if precipitate forms, or if solution is no longer clear.

B. Working 1% Silver Nitrate Solution

STOCK 10% silver nitrate	4.0 mL
Distilled water	36.0 mL

Mix together. May be made fresh, just before use, or may be stored in refrigerator (3°C) in brown glass bottle; stable for 2-4 months.

C. Reducing Solution

Sodium sulfite	5.0 gm
Hydroquinone	1.0 gm
Distilled water	100.0 mL

JUST BEFORE USE, dissolve together completely.

D. 5% Aluminum Sulfate

Aluminum sulfate	25.0 gm
Distilled water	500.0 mL

Dissolve together. Store at room temperature; stable for months.

E. Nuclear Fast Red

Nuclear fast red	0.1 gm
5% aluminum sulfate	100.0 mL

Dissolve together with the aid of gentle heat. Cool. Filter. Add a few crystals of thymol. Store at room temperature or in refrigerator (3°C.); stable for months.

V. EQUIPMENT:

- A. Balance
- B. 60°C oven, water bath or microwave oven
- C. Magnetic stirrer

VI. SUPPLIES:

- A. Erlenmeyer flasks
- B. Graduated cylinders

- C. Acid clean coplin jars
- D. Non-metal forceps

VII. SPECIAL SAFETY PRECAUTIONS:

- A. Silver Nitrate
 - 1. Is an oxidizer.
 - 2. Store separately from other material.
 - 3. Is poisonous and may be fatal if swallowed.
 - 4. Causes skin and eye burns.
 - 5. Is an irritant to the respiratory system.
- B. Sodium Sulfite
 - 1. Has low hazard for recommended handling.
- C. Hydroquinone
 - 1. Is irritating to eye, skin, and respiratory tract.
 - 2. May cause allergic skin reaction.
 - 3. Harmful if swallowed.
- D. Aluminum Sulfate
 - 1. Is a corrosive.
 - 2. May cause serious eye damage.
- E. Nuclear Fast Red
 - 1. Corrosive to skin.
 - 2. Irritant to eyes.
 - 3. Toxic to respiratory system.
- F. Thymol
 - 1. May cause eye burns.
 - 2. May be irritating to respiratory tract and skin.

VIII. QUALITY CONTROL(QC):

Section of tissue with argyrophilic carcinoid tumor or pancreas with islet cells.

IX. LIMITATIONS:

- A. Use acid-clean glassware, to avoid silver precipitate.
- B. Use non-metal forceps, or a silver precipitate may be formed on the slides and tissue.
- C. Use distilled water. Any metals in tap water may cause silver precipitate to be formed on the

slides and the tissue.

- D. To differentiate between argyrophil and argentaffin substances, a duplicate slide should be stained with the Fontana-Masson procedure, and the results compared.
- E. If nuclear fast red is used as a counterstain, wash slides with water after staining. Aluminum sulfate does not dissolve in alcohol. Placing slides directly from nuclear fast red into alcohol will result in a white precipitate forming on the slides, which can be removed by returning the slides to water.

X. PROCEDURE:

Step	Action	Time	Notes
1	Deparaffinize and hydrate sections through graded alcohol to distilled water.		
2	Place slides WORKING 1% silver nitrate solution.	1-2 hours	Incubate in 60°C oven. Or microwave on high for 45 seconds and allow to sit on counter for 1 minute. Save hot silver solution.
3	Rinse in distilled water, 3 changes.	5-10 seconds each	
4	Place slides in previously heated (60°C.) reducing solution.	5 minutes	
5	Wash in running tap water.	5 minutes	
6	Rinse in distilled water, 2-3 changes.	5-10 seconds	
7	Return slides to same coplin jar of 1% silver solution in 60°C oven.	10 minutes	Or heat slides in same 1% silver solution in microwave on high for 30 seconds.
8	Rinse slides in distilled water, 2-3 changes.	5-10 seconds	
9	Return slides to same coplin jar of reducing solution in 60°C oven.	1-2 hours	Or heat same reducing solution in microwave oven on high for 15 seconds. Return slides to heated solution for 1-2 minutes.
10	Rinse in distilled water, 2-3 changes each.	5-10 seconds	
11	Counterstain in nuclear fast red.	1-5 minutes	
12	Rinse in distilled water, 2 changes each.	5-10 seconds	
13	Dehydrate through graded	10	

	alcohols, clear with xylene.	seconds each	
14	Coverslip.		

XI. RESULTS:

- A. Argyrophilic granules (carcinoid tumors, neurosecretory granules) - **black**
- B. All argentaffin granules (melanin, argentaffin cells) - **black**
- C. Nuclei - **red**
- D. Background - **pale pink**

XII. REFERENCES:

- A. Bancroft JD, Stevens A: Theory and Practice of Histological Techniques, 3rd ed. New York, NY, Churchill Livingstone, 1990.
- B. Carson FL: Histotechnology: A Self-Instructional Text, Chicago, IL, ASCP Press, 1990.
- C. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd edition. Columbus, Ohio, Battelle Press, 1980.
- D. Microwave Ovens for Metallic Histologic Staining: A New Concept. Paper presented at National Society of Histotechnology Symposium. Modified by William Beaumont Hospital, Royal Oak, MI.

Approval Signatures

Step Description	Approver	Date
Medical Director	Kurt Bernacki: System Med Dir, Surgical Path	1/5/2023
Policy and Forms Steering Committee (if needed)	Gail Juleff: Project Mgr Policy	1/5/2023
Policy and Forms Steering Committee (if needed)	Sharon Scalise: Supv, Laboratory	1/4/2023
	Amy Knaus: Dir, Lab Operations C	1/4/2023
	Jennifer Lehmann: Mgr Laboratory	1/3/2023
	Sharon Scalise: Supv, Laboratory	12/19/2022

Applicability

Royal Oak

COPY