

Beaumont

Origination 1/18/2023
 Last 1/18/2023
 Approved
 Effective 1/18/2023
 Last Revised 1/18/2023
 Next Review 1/17/2025

Document Sharon Scalise:
 Contact Supv, Laboratory
 Area Laboratory-
 Histology
 Applicability Royal Oak

Histology Special Stain - Modified Gomori Trichrome - Royal Oak

Document Type: Procedure

I. PURPOSE AND OBJECTIVE:

The purpose of this document is to provide a procedure for the demonstration of connective tissue, nerves, and the myofibrils in muscle biopsies. It also demonstrates mitochondria, though they are usually not seen in normal muscles. In ragged red fibers, the accumulation of bizarre mitochondria can be seen as a thin, reddish rim, often cracked, around the inside edge of each myofibril. This stain will demonstrate cytoplasmic bodies, whose composition is unknown. Nemaline rods are only demonstrated with modified Gomori trichrome.

II. PRINCIPLE:

The reaction is based on the acid/base (+/- charges) and the hardness/softness of the tissue and the dyes. The negative charged dyes will be attracted to the positively charged tissue components. Hard dyes will be attracted to hard tissue components, while soft dyes will be attracted to soft tissue components. Chromotrope 2R will stain cytoplasm and muscle. Fast Green FCF will stain connective tissues. Phosphotungstic acid is taken up by the collagen. Fast Green FCF is subsequently bound to the phosphotungstic acid. The acetic acid changes the charges/pH of the solution. The acetic acid rinse after the stain gives a more delicate stain but will not change the ratio. Hematoxylin stains the nuclei.

III. SPECIMEN COLLECTION AND HANDLING:

A. Fixation

1. Fresh frozen tissue.

- B. Processing
 - 1. No processing.
- C. Section Thickness
 - 1. Cut frozen section muscle biopsies at 10 μ .
- D. Storage
 - 1. Store slides in refrigerator.
- E. Type of slide
 - 1. Plus slides.

IV. REAGENTS:

A. Modified Gomori Trichrome Solution

Chromotrope 2R	0.30 gm
Fast Green FCF	0.15 gm
Phosphotungstic Acid	0.40 gm
Distilled Water	50.00 ml
Acetic Acid	0.50 ml

Dissolve together. Adjust pH to 3.4 with 1 N Sodium Hydroxide. Store in refrigerator; stable for several months or until weak.

B. 0.2% Acetic Acid

Acetic Acid	60.2 ml
Distilled Water	100.0 ml

Stir together; store at room temperature; stable for several months.

C. Hematoxylin

Use hematoxylin from routine H&E set-up, Gill or Mayer preferable.

V. EQUIPMENT:

- A. Balance
- B. 60°C oven
- C. Magnet stirrer

VI. SUPPLIES:

- A. Erlenmeyer flasks
- B. Graduated cylinders
- C. Coplin jars
- D. Forceps
- E. Funnel
- F. Filter paper

VII. SPECIAL SAFETY PRECAUTIONS:

- A. Chromotrope 2R
 - 1. Is an irritant.
- B. Fast Green FCF
 - 1. Is a carcinogen and an oxidizer.
- C. Phosphotungstic Acid
 - 1. Is a corrosive.
- D. Acetic Acid
 - 1. Is an acid.

VIII. QUALITY CONTROL(QC):

Frozen section of muscle with nerve and connective tissue.

IX. PROCEDURE:

Step	Action	Time	Notes
1	Place slides in hematoxylin.	5 minutes	Incubate at room temperature.
2	Rinse in distilled water.	1 minute	
3	Stain in Modified Gomori Trichrome.	10 minutes	
4	Rinse in 0.2% Acetic Acid.		
5	Dehydrate through graded alcohols, clear with xylene.		
6	Coverslip.		

X. RESULTS:

- A. Myofibrils - **green**
(Type I muscle may appear darker than Type II)
- B. Connective tissue - **green**
- C. Cytoplasm - **green**
- D. Hypercontracted muscle fibers - **red**
- E. Mitochondria - **red**
- F. Intermyoibrillar sarcoplasmic reticulum - **red**
- G. Nuclei - **blue**

XI. REFERENCES:

- A. Carson FL: Histotechnology: A Self-Instructional Text. Chicago, IL, ASCP Press, 1990.
- B. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd edition. Columbus, Ohio, Battelle Press, 1980.

Approval Signatures

Step Description	Approver	Date
Medical Director	Kurt Bernacki: System Med Dir, Surgical Path	1/18/2023
Policy and Forms Steering Committee (if needed)	Sharon Scalise: Supv, Laboratory	1/17/2023
Policy and Forms Steering Committee (if needed)	Gail Juleff: Project Mgr Policy	1/17/2023
	Amy Knaus: Dir, Lab Operations C	1/17/2023
	Jennifer Lehmann: Mgr Laboratory	1/17/2023
	Sharon Scalise: Supv, Laboratory	1/10/2023

Applicability

Royal Oak