

Beaumont

Origination 11/4/2022
 Last 11/4/2022
 Approved
 Effective 11/4/2022
 Last Revised 11/4/2022
 Next Review 11/3/2024

Document Sharon Scalise:
 Contact Supv, Laboratory
 Area Laboratory-
 Histology
 Applicability Royal Oak

Histology Special Stain - Phosphotungstic Acid Hematoxylin (PTAH) - Royal Oak

Document Type: Procedure

I. PURPOSE AND OBJECTIVE:

The purpose of this document is to provide a procedure for the demonstration of muscle striations, fibrin, collagen, and glial fibers. Muscle striations are found in striated and cardiac muscle, and in tumors arising from such tissues. Fibrin is a specific protein formed as a result of the clotting process and is a prominent feature in renal transplant rejections and tissues that have received a physical trauma. Glial fibers are found in the brain, and in gliosarcoma, a tumor of the brain.

II. PRINCIPLE:

Tungsten is a mordant that attaches to the hematoxylin, creating a blue-colored lake. This tungsten-hematein lake will bond with certain tissue components, staining them blue. Excess phosphotungstic acid will stain fibrin and connective tissue a red-brown color. Bouin's solution is used as a post-mordant, which improves the staining.

III. SPECIMEN COLLECTION AND HANDLING:

A. Fixation

1. Any well-fixed tissue.
2. Zenker's, Bouin's or DB preferred

B. Processing

1. Standard processing.

- C. Section thickness
 - 1. Routine specimens- 5µm.
- D. Slide drying
 - 1. 60 minutes at 60°C.
- E. Type of slide
 - 1. Plain slides

IV. REAGENTS:

A. Saturated Aqueous Picric Acid

Purchased Pre-Made

B. Bouin's Fluid

Saturated aqueous picric acid	750.0 mL
Formaldehyde, 37-40%	250.0 mL
Acetic acid, concentrated	50.0 mL

Mix together, adding the formaldehyde and acetic acid slowly, as not to splash. Store at room temperature. Stable for months.

SEE CAUTION CONCERNING PICRIC ACID UNDER SPECIAL SAFETY PRECAUTIONS.

C. 5% Ferric Ammonium Sulfate

Ferric ammonium sulfate (Iron Alum)	5.0 gm
Distilled water	100.0 mL

Stir together. Store at room temperature. Stable for several months. Discard after use.

D. 5% Oxalic Acid

Oxalic acid	5.0 gm
Distilled water	100.0 mL

Dissolve together. Stable at room temperature for months. Discard after use.

E. Phosphotungstic Acid Hematoxylin (PTAH)

Hematoxylin	1.0 gm
Absolute ethanol	10.0 mL
Phosphotungstic Acid	20.0 gm
Distilled water	900.0 mL

1. Dissolve together 1 gm hematoxylin in 10 mL absolute alcohol.
2. Dissolve together 20 gm phosphotungstic acid in 900 mL hot distilled water.
3. Allow phosphotungstic acid solution to cool.
4. Mix together hematoxylin solution and phosphotungstic solution.
5. The solution ripens naturally in several months and is stable for at least one year.
6. May be stored at room temperature.
7. Discard after use.
8. If the stain is needed for immediate use, add 0.177 gm potassium permanganate to the solution. When artificially ripened solution is a reddish-brown color, it is ready for

use.

V. EQUIPMENT:

- A. Balance
- B. Magnetic stirrer

VI. SUPPLIES:

- A. Erlenmeyer flasks
- B. Graduated cylinders
- C. Acid clean coplin jars
- D. Forceps

VII. QUALITY CONTROL:

- A. To demonstrate:
 - 1. Cross-striations and connective tissue – use striated muscle.
 - 2. Fibrin – use tissue with fibrin
 - 3. Glial cells – use cerebral cortex

VIII. SPECIAL SAFETY PRECAUTIONS:

- A. Picric Acid
 - 1. Is toxic, highly reactive and an extreme fire hazard.
 - 2. Keep picric acid moist at all times.
 - 3. If a dry powder is seen around the lip of the jar, wash off with running water before opening.
 - 4. Store in an explosion proof cabinet.
- B. Formaldehyde
 - 1. Is poisonous.
 - 2. May be fatal if swallowed or cause permanent blindness.
 - 3. Possible cancer hazard.
 - 4. Irritant to eyes, skin, respiratory system.
- C. Acetic Acid
 - 1. Add drop by drop to solutions.
 - 2. May cause skin and eye burns.
- D. Ferric Ammonium Sulfate (Iron Alum)
 - 1. Is an irritant.

E. Oxalic Acid

1. Is a strong reducer.
2. Contact with other material may cause fire.
3. Store separate.
4. May cause skin and eye burns.
5. May be irritating to respiratory system.

F. Hematoxylin

1. Is incompatible with oxidizers and alkalies.
2. Store separate from these.

G. Phosphotungstic Acid

1. Is a corrosive.
2. May cause skin and eye irritation.

IX. PROCEDURE:

Step	Action	Time	Notes
1	Deparaffinize and hydrate slides through graded alcohol to distilled water.		
2	Place in Bouins solution in 60°C oven.	35-60 minutes	Steps 2-4 can be skipped if tissue was fixed in a picric acid fixative (Bouins; DB).
3	Remove from oven, and place slides in a different coplin jar containing water.		
4	Wash in running water until all the yellow color is gone.	2-10 minutes	
5	Place in Weigert hematoxylin.	10 minutes	
6	Rinse in running tap water.	1 minute	
7	Stain in Biebrich Scarlett- Acid Fuchsin.	2-3 minutes	
8	Rinse in running tap water, 2-3 changes.	5-10 seconds	
9	Place in phosphotungstic acid solution, jiggling occasionally.	1-5 minutes	
10	Rinse in distilled water, quickly, 2 changes each.	1-2 seconds	
11	Stain slides in aniline blue solution	1-10 minutes	Check the staining of the aniline blue with a microscope.

			The intensity of the red and blue should be equal. Do NOT check constantly, as this seems to make the colors muddy.
12	Rinse quickly in distilled water, 2-3 changes.	5-10 seconds	
13	Dehydrate through graded alcohols, clear with xylene.		
14	Coverslip.		

X. RESULTS:

- A. Striated muscle fibers, fibrin, glial cells - **blue**
- B. Nuclei, mitotic figures, mitochondria, keratin - **blue**
- C. Collagen, reticulin, basement membrane, cartilage - **red to red-brown**

XI. LIMITATIONS:

- A. PTAH must be either chemically oxidized or naturally ripened for the stain to work.
- B. Chemically ripened PTAH has a shorter shelf life than naturally ripened (sun, air).
- C. Over-oxidized (old) PTAH will show very pale blue staining.
 - 1. Discard over-oxidized solution.

XII. REFERENCES:

- A. Bancroft JD, Stevens A: Theory and Practice of Histological Techniques, 3rd ed, New York, NY, Churchill Livingstone, 1990.
- B. Carson FL: Histotechnology: A Self-Instructional Text. Chicago, IL, ASCP Press, 1990.
- C. Sheehan DC, Hrapchak BB: Theory and Practice of Histotechnology, 2nd edition. Columbus, Ohio, Battelle Press, 1980.
- D. Preece A: A Manual for Histologic Technicians, 3rd edition. Boston, MA, Little, Brown & Co., 1972.

Approval Signatures

Step Description	Approver	Date
Medical Director	Kurt Bernacki: System Med Dir, Surgical Path	11/4/2022

Policy and Forms Steering Committee (if needed)	Gail Juleff: Project Mgr Policy	11/4/2022
Policy and Forms Steering Committee (if needed)	Sharon Scalise: Supv, Laboratory	11/3/2022
	Amy Knaus: Dir, Lab Operations C	11/3/2022
	Jennifer Lehmann: Mgr Laboratory	11/3/2022
	Sharon Scalise: Supv, Laboratory	11/2/2022

Applicability

Royal Oak

COPY