


A. HIGH MCHC: If there is a high MCHC (>38) result on first run of sample:

1. If warming sample to 37 degrees:
 - a. Do not validate sample. Click on the Actions button and select Rerun. (All tests will be rerun if no tests are selected).



- b. After warming, run sample in manual mode on analyzer (to prevent specimen cool-down). If MCHC is acceptable, choose Run 2 from the Rerun tab.
- c. Click the Rerun hyperlink on the Results Validation screen to open the Rerun table.



- d. Select all tests by clicking on the select all checkbox, or select specific tests by using the individual check boxes.

Select all

The empty box in the column header will select all the tests in that run.

<input type="checkbox"/> RUN 1	Analyzer ID:	<input type="checkbox"/> RUN 2	Analyzer ID:		
<input type="checkbox"/>	Data Mark	Result	<input type="checkbox"/>	Data Mark	Result
<input type="checkbox"/>		2.00	<input checked="" type="checkbox"/>		3.00
<input type="checkbox"/>		2.0	<input checked="" type="checkbox"/>		3.0
<input type="checkbox"/>		3.00	<input checked="" type="checkbox"/>		3.00
<input type="checkbox"/>		1.00	<input checked="" type="checkbox"/>		3.00
<input type="checkbox"/>		2.0	<input checked="" type="checkbox"/>		3.0
<input type="checkbox"/>		3.0	<input checked="" type="checkbox"/>		3.0

- e. Add “possible cold agglutinin” comment by clicking on the comment icon in the result line of the test code.

A dark blue dialog box titled "Add New Comment for MCHC" with a close button (X) in the top right corner. Inside the dialog, there is a text input field labeled "Comment" containing the placeholder text "Enter Text". Below the input field is a blue icon of two overlapping speech bubbles followed by the text "Add Coded Comment". At the bottom right of the dialog are two buttons: "Cancel" and "Save".

- f. Validate CBC or Validate All according to the procedure.



- g. If smear does not need review, select Canceled in the pop-up box.

A dark blue dialog box titled "Select Smear Type" with a close button (X) in the top right corner. Inside the dialog, there are two radio button options: "Reviewed" and "Canceled". At the bottom of the dialog are two buttons: "Quit" and "OK" (which is highlighted in green).


2. If you perform a HCT manually:
- Replace the analyzer HCT value with the manual result by editing the analyzer results. Click on the result to be edited and backspace the result.
 - Replace the MCV and MCHC indices with the recalculated values.
 - Add internal comment: "Manual HCT. RBC indices recalculated."

d. Save and validate as appropriate.



B. If 1:5 dilution is made:

1. Do not validate sample. Click on the Actions button and select Rerun with dilution. (All tests will be rerun if no tests are selected).




2. Input the dilution factor.
3. Prepare dilution then run dilution in manual mode with the original barcode label.
4. Click the Rerun hyperlink on the Results Validation screen to open the Rerun table.




5. The dilution rerun will automatically calculate the final results. Note: The dilution rerun will not display that a dilution was performed. The only way to view if a dilution was performed is in the audit trail.
6. Check the boxes next to the values that should be released from the original and diluted run.

	RUN 1	Analyzer ID: XN905	RUN 2	Analyzer ID: XN905	RUN 3	Analyzer ID: XN905
Test Code	<input checked="" type="checkbox"/>	Data Mark	Result	<input checked="" type="checkbox"/>	Data Mark	Result
WBC	<input checked="" type="checkbox"/>	8,8,8	9.0			
RBC		@	10.35		<input checked="" type="checkbox"/>	9.85
HGB	<input checked="" type="checkbox"/>		22.3			
HCT	<input checked="" type="checkbox"/>		75.0			
MCV	<input checked="" type="checkbox"/>		75			
MCH	<input checked="" type="checkbox"/>		22			
MCHC	<input checked="" type="checkbox"/>		30			


7. Add appropriate correction comment (lipemic, icteric) as internal comment and validate as appropriate. Click on the comment icon on the left of the screen. Then click on the light blue comment icon.




8. The order comment box will appear. Click on Add New Comment.




9. Add New Comment box will appear. Select internal or report comments. Select from the list of canned comments or free text a comment. Click OK when complete.



10. Order comment appears in pop-up box. Click Close.



11. Comment now appears on the left side of the screen.




12. An Action on Smear pop-up box will appear. If slide is to be reviewed by morph bench, QUIT. If a slide does not need review, Cancel Smear.

13. Validate CBC or Validate All according to the procedure.



14. If smear does not need review, select Canceled in the pop-up box. Click OK.




15. If smear needs review, select quit for the pop-up box.

C. CBCND THAT REQUIRES WBC or PLT SLIDE REVIEW

1. If all parameters except WBC or PLT autoverify:

a. Review slide. Click the Morph tab.



b. Perform “Morphology Checks” by clicking on the three dots in the result fields and selecting appropriate responses for WBC estimate, instrument flags, RBC morphology and platelet estimate.

Test	Result
WBC ESTIMATE	...
INSTRUMENT FLAGS	...
RBC MORPHOLOGY	...
PLT ESTIMATE	...

c. If significant RBC morphology is present, you may enter appropriate grading in Diff Comment field (since morphology doesn’t go to LIS on a CBCND.)

DIFF COMMENT	...
--------------	-----

d. Click Save when complete.



e. The morph results will display on the right side of the screen.

Manual			
Diff Morph Fluid Count			
<input type="checkbox"/>	Test	Result	
<input type="checkbox"/>	WBC ESTIMATE	Confirmed ...	
<input type="checkbox"/>	INSTRUMENT FLAGS	Reviewed ...	
<input type="checkbox"/>	RBC MORPHOLOGY	Unremarkable ...	
<input type="checkbox"/>	PLT ESTIMATE	Adequate ...	
<input type="checkbox"/>	SMEAR 5	Vérifié ...	

f. Validate All.



2. If none of the CBC parameters autoverify:

a. Select the parameters to be resulted immediately before slide review.

<input type="checkbox"/>	Test	Result	Comment
<input type="checkbox"/>	WBC	6.1 ...	
<input checked="" type="checkbox"/>	RBC	5.53 ...	
<input checked="" type="checkbox"/>	HGB	15.9 ...	
<input checked="" type="checkbox"/>	HCT	49.1 ...	
<input checked="" type="checkbox"/>	MCV	89 ...	
<input checked="" type="checkbox"/>	MCH	29 ...	
<input checked="" type="checkbox"/>	MCHC	32 ...	
<input type="checkbox"/>	PLT	200 ...	
<input checked="" type="checkbox"/>	RDWCV	14 ...	

b. Click Validate selection.



c. The result fields will turn grey indicating the results crossed to the LIS.

<input type="checkbox"/>	Test	Result	Comment
<input type="checkbox"/>	WBC	6.1	...
<input type="checkbox"/>	RBC	5.53	
<input type="checkbox"/>	HGB	15.9	
<input type="checkbox"/>	HCT	49.1	
<input type="checkbox"/>	MCV	89	
<input type="checkbox"/>	MCH	29	
<input type="checkbox"/>	MCHC	32	
<input type="checkbox"/>	PLT	200	...
<input type="checkbox"/>	RDWCV	14	

d. Review slide. Click the Morph tab.

Manual

Diff | **Morph** | Fluid Count

<input type="checkbox"/>	Test	Result	Comment
<input type="checkbox"/>	SMEAR 5		...

e. Perform “Morphology Checks” by clicking on the three dots in the result fields and selecting appropriate responses for WBC estimate, instrument flags, RBC morphology and platelet estimate.

Test	Result
WBC ESTIMATE	...
INSTRUMENT FLAGS	...
RBC MORPHOLOGY	...
PLT ESTIMATE	...

f. If significant RBC morphology is present, you may enter appropriate grading in Diff Comment field (since morphology doesn’t go to LIS on a CBCND.)

DIFF COMMENT	...
--------------	-----

g. Click Save when complete.

Clear

Save

h. The morph results will display on the right side of the screen.

Manual			
Diff Morph Fluid Count			
<input type="checkbox"/>	Test	Result	
<input type="checkbox"/>	WBC ESTIMATE	Confirmed ...	
<input type="checkbox"/>	INSTRUMENT FLAGS	Reviewed ...	
<input type="checkbox"/>	RBC MORPHOLOGY	Unremarkable ...	
<input type="checkbox"/>	PLT ESTIMATE	Adequate ...	
<input type="checkbox"/>	SMEAR 5	Vérifié ...	

i. Validate All.




D. PLT < 75

1. Platelets less than 75 bill/L will automatically reflex for a fluorescent PLT (if applicable to the instruments in use).
2. The fluorescent PLT result will appear in the PLT field with a &F symbol.


&F	PLT	224
----	-----	-----

3. Note that the Op Alert states to check for a clot and review a smear.
4. Retrieve specimen and check for a clot. If no clot is found, add an internal comment. Click on the comment icon on the left of the screen. Then click on the light blue comment icon.




- a. If a clot is found, follow steps for canceling specimen both in Caresphere and the LIS. Document contact information in the LIS. See notes section of procedure for more information.


5. The order comment box will appear. Click on Add New Comment.




6. Add New Comment box will appear. Select internal or report comments. Select from the list of canned comments or free text a comment. Click OK when complete.



7. Order comment appears in pop-up box. Click Close.




8. Comment now appears on the left side of the screen.




9. The XN line will not automatically make a smear on PLTs < 75. Check previous results to determine if the patient had a previous PLT less than 75 that was verified by a smear. Be sure to scroll down to see if previous results had a smear looked at (e.g. RBC morphology, Platelet field resulted).
10. If previous results are present for a sample ID the most recent set of previous results will be displayed on the Result Validation screen.
11. Click the Prev Res hyperlink to open the Previous Results popup.




12. The list is sorted in reverse chronological order. The most recent previous result is considered as Prev Res 1.
13. If unable to determine that a smear was previously reviewed for a PLT less than 75, click on the Action button, then click Add Test.




14. Check mark the box next to Smear. Click OK.



15. Click on Smear 1 for one slide to be made. Click on Smear 2 for two slides to be made. Click OK.



16. The smear is added to the right side of the screen.



17. Retrieve specimen and place back onto the XN line so that a smear is made. Or make a manual smear.

18. Select all but the PLT and fluorescent PLT (if ran) parameters, then click on Validate Selection.

<input type="checkbox"/>	Test	Result	Comment
<input checked="" type="checkbox"/>	WBC	4.9	...
<input checked="" type="checkbox"/>	RBC	5.15	...
<input checked="" type="checkbox"/>	HGB	15.0	...
<input checked="" type="checkbox"/>	HCT	45.3	...
<input checked="" type="checkbox"/>	MCV	89	...
<input checked="" type="checkbox"/>	MCH	30	...
<input checked="" type="checkbox"/>	MCHC	34	...
<input type="checkbox"/>	PLT	250	... 1
<input checked="" type="checkbox"/>	RDWCV	13	...
<input type="checkbox"/>	NEUTRE	56.0	...
<input type="checkbox"/>	LYMPRE	33.1	...
<input type="checkbox"/>	MONORE	6.7	...
<input type="checkbox"/>	EOSIRE	2.4	...
<input type="checkbox"/>	BASORE	1.4	...
<input type="checkbox"/>	IGRE	0.4	...
<input type="checkbox"/>	NRBCRE	0.0	...
<input type="checkbox"/>	NEUTAB	2.7	...
<input type="checkbox"/>	LYMPAB	1.6	...
<input type="checkbox"/>	MONOAB	0.3	...
<input type="checkbox"/>	EOSIAB	0.1	...
<input type="checkbox"/>	BASOAB	0.1	...




19. After the smear is made and stained, access the case in Caresphere (through the To Validate portal or the Manual Review portal).
20. Perform the “Morphology Checks” by clicking on the Morph tab and entering results for WBC Estimate, Instrument Flags, RBC Morphology, and Plt Estimate.
21. Click on Validate All. This will finalize all results in Caresphere (all fields, including morphology responses will be gray) and in the LIS.




E. CLUMPED PLATELETS OR PLATELET SATELLITOSIS

1. If the XN analyzer gives any clumped platelet flags or OP Alert for Platelet Satellitosis, check sample for clots.
2. Retrieve specimen and check for a clot. If no clot is found, add an internal comment. Click on the comment icon on the left of the screen. Then click on the light blue comment icon.




- a. If a clot is found, follow steps for canceling specimen both in Caresphere and the LIS. Document contact information in the LIS. See notes section of procedure for more information.


- 3. The order comment box will appear. Click on Add New Comment.




- 4. Add New Comment box will appear. Select internal or report comments. Select from the list of canned comments or free text a comment. Click OK when complete.



- 5. Order comment appears in pop-up box. Click Close.



6. Comment now appears on the left side of the screen.



7. Select all but the PLT and fluorescent PLT (if ran) parameters, then click on Validate Selection.

<input type="checkbox"/>	Test	Result	Comment
<input checked="" type="checkbox"/>	WBC	4.9 ...	
<input checked="" type="checkbox"/>	RBC	5.15 ...	
<input checked="" type="checkbox"/>	HGB	15.6 ...	
<input checked="" type="checkbox"/>	HCT	45.3 ...	
<input checked="" type="checkbox"/>	MCV	88 ...	
<input checked="" type="checkbox"/>	MCH	30 ...	
<input checked="" type="checkbox"/>	MCHC	34 ...	
<input type="checkbox"/>	PLT	250 ... 1	
<input checked="" type="checkbox"/>	RDWCV	13 ...	
<input type="checkbox"/>	NEUTRE	56.0 ...	
<input type="checkbox"/>	LYMPRE	33.1 ...	
<input type="checkbox"/>	MONORE	6.7 ...	
<input type="checkbox"/>	EOSIRE	2.4 ...	
<input type="checkbox"/>	BASORE	1.4 ...	
<input type="checkbox"/>	IGRE	0.4 ...	
<input type="checkbox"/>	NRBCRE	0.0 ...	
<input type="checkbox"/>	NEUTAB	2.7 ...	
<input type="checkbox"/>	LYMPAB	1.6 ...	
<input type="checkbox"/>	MONOAB	0.3 ...	
<input type="checkbox"/>	EOSIAB	0.1 ...	
<input type="checkbox"/>	BASOAB	0.1 ...	

 Validate Sel

8. Review slide. If platelet clumps or platelet satellitosis are present and your platelet estimate differs from the instrument count, click on the Morph tab.

Manual

Diff | **Morph** | Fluid Count

Test	Result
<input type="checkbox"/> SMEAR 1	...

9. In the PLT estimate result field, choose the appropriate platelet estimate (increased, decreased, normal). Result all the required morphology fields appropriately.

Test	Result
WBC ESTIMATE	...
INSTRUMENT FLAGS	...
RBC MORPHOLOGY	...
PLT ESTIMATE	...

10. Click Save when complete.

11. Remove the platelet result by clicking on the platelet result field and backspacing the result.

PLT ...

12. Click on the three dots next to the platelet field to select the appropriate clumped platelet/ platelet satellitosis comment that correlates with the platelet estimate. Click Okay. Click Save.

13. Click Validate All.




F. WBC <0.4

1. All CBCWD orders with a WBC < 0.4 bill/L will not have a differential resulted.
2. The LIS order will need to be changed to a CBCND.
 - a. NOTE: If pathologist review is needed based on the smear review checklist, keep the CBCWD order and send for pathologist review.
3. Add the HE10 comment to the WBC parameter for the CBCND order in the middleware.
4. Click on the comment icon in the result line of the test code.



5. A window will open. Select Add Coded Comment to select the HE10 comment. Click OK. Click Save.



Add Coded Comment ✕

Type to search 🔍

Items Selected: 1 Unselect All

<input type="checkbox"/>	Code	Text
<input type="checkbox"/>	C	Cancel.
<input type="checkbox"/>	CLOT	Clotted.
<input type="checkbox"/>	HE01	Rechecked and Verified.
<input type="checkbox"/>	HE02	Corrected for Lipemia.
<input type="checkbox"/>	HE03	Possible cold agglutinin.
<input type="checkbox"/>	HE08	Results in question; suggest repeat if n...
<input checked="" type="checkbox"/>	HE10	Differential not performed when WBC <...
<input type="checkbox"/>	HE38	Result change from previously reported...
<input type="checkbox"/>	NCALC	not calculated
<input type="checkbox"/>	NM	#NM
<input type="checkbox"/>	PLT9	Giant Platelets Present.

Differential not performed when WBC <0.4

Cancel OK

Add New Comment for WBC ✕

Comment

Enter Text

 [Add Coded Comment](#)

Created Date: 10/10/2023 10:49:44 - Created By:
 Megan.Masakowski@beaumont.org  

Differential not performed when WBC <0.4

Cancel Save

6. Continue to result the CBCND in the middleware per procedure.

G. FIRST TIME BLASTS SEEN

1. Validate CBC or Validate Selection after addressing the flags, Op alerts, etc.



2. Review the slide under the microscope. Perform a manual differential to include the blast count. Add morphologies to the morphology tab if applicable.
3. Add a comment next to the blast percentage by clicking the blue comment icon.

4. Enter free text "First time critical blast." Click Save.
5. Add a path review in Caresphere (see step K below).
6. Validate All when complete.
7. Document all critical calls in the LIS.
8. Follow site specific workflow for sending slides and paperwork to the pathologists.

H. SECOND LEVEL REVIEW FOR BLAST FLAG

1. Validate CBC only after addressing flags, Op alerts, etc.



2. If only a scan is required:
 - a. First level tech: Perform scan and result required morphology fields (WBC estimate, Instrument flag check, RBC morphology and PLT estimate). If differential can be released without further review by pathologist and only secondary blast review is needed, Validate All.




This will send out the differential and leave only the Second Level Review field pending.


- b. Secondary Review Tech: Click the Sample Explorer tab on the top of the screen.




- c. Scan or type the order ID into the Sample ID field.

A light gray rectangular form titled "Search Criteria". It contains a label "Sample ID" above a text input field containing the value "129027651".

- d. Click Apply.



- e. The sample will appear in the search result list. Check the box next to the sample ID, click Details and Result Validation to branch to the result validation screen or click on the sample ID hyperlink (will open a new tab).



- f. The result validation screen will open.
- g. After reviewing the smear, result the Second Level Review field by clicking on the three dots next to the field.



- h. Click on Reviewed in the pop-up window. Click OK.

Coded Result	
Type to search	
Items Selected: 1 Unselect All	
Code	Text
<input type="checkbox"/>	HE10 Repeat differential not pe...
<input type="checkbox"/>	HE38 Result change from previ...
<input type="checkbox"/>	NCALC not calculated
<input type="checkbox"/>	NM #nm
<input type="checkbox"/>	PLT9 Giant Platelets Present.
<input type="checkbox"/>	PLTO PLT verified by alternate ...
<input checked="" type="checkbox"/>	R Reviewed.
<input type="checkbox"/>	TS500 Sorted by TS500
<input type="checkbox"/>	V verified.
<input type="checkbox"/>	VERIF Verified by smear.

Reviewed.

Cancel OK

- i. Click Save.



- j. Click Validate All.



3. If a differential is required:

- First level tech will perform manual differential under the Diff tab.
- Click Diff under the Manual panel in the Result Validation screen.

Manual


Diff | Morph | Fluid Count



- c. MDIFF is the default counter. Perform manual differential by selecting the appropriate key for each respective cell type. (User can select which key translates to which cell type).
1. When count limit is reached, the counting stops. No additional cells are added.
 2. The count limit default is 100 cells. If you need to count more than 100 cells edit this number before beginning to count.
 3. If the counting is stopped before the count limit is reached, you will be alerted that less than the number of cells has been reached. If accepted, the results will be averaged and rounded correctly for the total number of cells counted.
 4. In the event of a cell being misclassified, select Subtract. Remove the cell by selecting the respective key. When finished removing cell, select Add and continue enumerating cell types.



5. Additional counts can be completed by switching to Count 2. After completion, determine which results will be reported (Count 1 or Count 2).



- d. Click Save.



- e. After completion, review the Result Validation screen for any Operator Alerts related to manual differential and absolute counts.

- f. After performing the manual differential, complete the Morphology Checks.
- g. Click the Morph tab.



- h. Perform “Morphology Checks” by clicking on the three dots in the result fields and selecting appropriate responses for WBC estimate, instrument flags, RBC morphology and platelet estimate.

Test	Result
WBC ESTIMATE	...
INSTRUMENT FLAGS	...
RBC MORPHOLOGY	...
PLT ESTIMATE	...

- i. Once all results are finalized, the first level tech will Validate All. This will send out the differential and leave only the Second Level Review field pending.



- j. Secondary Review Tech: Click the Sample Explorer tab on the top of the screen.




- k. Scan or type the order ID into the Sample ID field.


Search Criteria

Sample ID
129027651

- l. Click Apply.




- m. The sample will appear in the search result list. Check the box next to the sample ID, click Details and Result Validation to branch to the result validation screen or click on the sample ID hyperlink (will open a new tab).



- n. The result validation screen will open.
- o. After reviewing the smear, result the Second Level Review field by clicking on the three dots next to the field.



- p. Click on Reviewed in the pop-up window. Click OK.



q. Click Save.



r. Click Validate All.




s. Since the second level review is internal, it will not appear in the LIS.

t. Note that if the Second Level Review is not answered, it will remain pending in Caresphere.


I. RETICULOCYTE (RET) ABN SCATTERGRAM

1. If the Retic Abn Scattergram is flagged and there are asterisks next to the RET%, RET#, IRF, and RET-He, a dilution is necessary. (If no asterisks are present, the retic parameters may be reported without further review). Refer to the workflow for proper handling of sample with this error.
2. In Caresphere under the Result Validation tab, check mark the box next to each retic parameter.

3. Click on the ACTION button and select rerun with dilution.




4. Enter dilution factor and click Save.



5. After you have prepared your dilution, run the diluted sample in the manual mode under the original barcode label. This run will cross over into Caresphere and can be seen in the Rerun tab.
6. The diluted results should automatically be calculated by Caresphere. Refer to the workflow for any additional steps to be taken before validating the results.

J. MICROORGANISMS FOUND ON PERIPHERAL SMEAR

1. In Morph tab, after reviewing the slide and resulting all of the required morphology fields (WBCEST, INSTFL check, RBCM and PLTE) click on the three dots in the Intracellular Bacteria result field.



2. Select See Below and Save.

Coded Result		
Type to search		Q
Items Selected: 1		Unselect All
Code	Text	
<input type="checkbox"/>	REMINDER	(send directly to micro for...
<input checked="" type="checkbox"/>	SB	See Below

See Below

Cancel OK

3. Validate All.



4. Send sample to microbiology following your sites specific workflow.

K. ADDING A PATHOLOGY (PATH) REVIEW IN CARESPHERE

1. Determine if the case needs to be sent for a Path Review.
2. If a path review is needed: Add a diff comment in the morph tab “Smear to be reviewed by pathologist”.

DIFF COMMENT | ...

Coded Result ✕

Type to search 🔍

Items Selected: 1 Unselect All

<input type="checkbox"/>	Code	Text
<input type="checkbox"/>	HD13	Diff performed on less th...
<input type="checkbox"/>	HE08	Results in question. Sugg...
<input type="checkbox"/>	HE10	Repeat diff not performe...
<input type="checkbox"/>	HE18	Occasional nRBC seen on ...
<input type="checkbox"/>	HE20	Echinocytes may be seen ...
<input type="checkbox"/>	HE24	Peripheral blood smear r...
<input checked="" type="checkbox"/>	HE25	Smear to be reviewed by ...
<input type="checkbox"/>	PLT6	Platelet satellitosis prese...
<input type="checkbox"/>	PLT7	Platelet satellitosis prese...
<input type="checkbox"/>	PLT8	Platelet satellitosis prese...

Smear to be reviewed by pathologist.


Cancel OK

Click Ok and Save.

3. Click the Action button and select Add Test.

The image shows a dropdown menu titled 'Action' with three options: 'Add Test', 'Rerun with Dilution', and 'Rerun'. The 'Add Test' option is highlighted in a light grey color.


4. Select Pathology Review. Click OK.



5. Path Review will display under the manual tab to the right.

Manual			
	Diff	Morph	Fluid Count
<input type="checkbox"/>	MONO ABSOLUTE	0.7	...
<input type="checkbox"/>	MONO ABSOLUTE	0.7	...
<input type="checkbox"/>	EOSI ABSOLUTE	0.1	...
<input type="checkbox"/>	BASO ABSOLUTE	0.0	...
<input type="checkbox"/>	NRBCMI ABSOLUTE	1.0	...
<input type="checkbox"/>	WBC ESTIMATE	Confirmed	...
<input type="checkbox"/>	INSTRUMENT FLAGS	Reviewed	...
<input type="checkbox"/>	RBC MORPHOLOGY	Unremarkable	...
<input type="checkbox"/>	PLT ESTIMATE	Adequate	...
<input type="checkbox"/>	Count	100	...
<input type="checkbox"/>	PATH REVIEW.		...
<input type="checkbox"/>	SMEAR 3	Verified	...
<input type="checkbox"/>	SMEAR 5	Verified	...
<input type="checkbox"/>	SECOND LEVEL REVI...		...

6. Enter Yes in the path review field by clicking on the three dots. Click Yes and OK.



7. Validate All when all results are final.



8. Follow site specific workflow for sending slides and paperwork to the pathologists.
9. Checking for a Previous Path Review in Caresphere: If you would like to know if a patient was previously sent for a Path Review:
 - a. Click the Prev Res hyperlink to open the Previous Results popup.



- b. The list is sorted in reverse chronological order. The most recent previous result is considered as Prev Res 1.

L. HOW TO REMOVE HEMCH WHEN SICKLE CELLS ARE RESULTED (IF APPLICABLE)

1. When sickle cells are called and validated in Caresphere, the order may show back up on the Caresphere list with the HEMCH field pending.

Manual		
Diff Morph Fluid Count		
Test	Result	
<input type="checkbox"/> WBC ESTIMATE	Confirmed	
<input type="checkbox"/> INSTRUMENT FLAGS	Reviewed	
<input type="checkbox"/> RBC MORPHOLOGY	See Below	
<input type="checkbox"/> PLT ESTIMATE	Adequate	
<input type="checkbox"/> SICKLE CELLS	Present	
<input type="checkbox"/> SMEAR 1	Verified	
<input type="checkbox"/> HEMCH	...	

- To cancel the HEMCH field, click on the three dots in the field and select cancel. Click Save.
- The order will autovalidate and be removed from the list.
- NOTE: If you cannot search for the sample ID afterward in Caresphere, you may need to select site R or site ALL to view the results.

M. MANUAL NRBC DISAGREEMENT FROM ANALYZER (XN SITES)

- Review Op alerts. Select the parameters that can be validated before slide review.

<input type="checkbox"/>	Test	Result	Comment	Rerun	Prev Res	Prev Com
<input type="checkbox"/>	& WBC	16.6	...		78.6	
<input checked="" type="checkbox"/>	RBC	3.22	...		4.09	
<input checked="" type="checkbox"/>	HGB	5.9	...		11.4	
<input checked="" type="checkbox"/>	HCT	28.3	...		33.7	
<input checked="" type="checkbox"/>	MCV	88	...		92	
<input checked="" type="checkbox"/>	MCH	28	...		28	
<input checked="" type="checkbox"/>	MCHC	31	...		34	
<input type="checkbox"/>	PLT	1091	...		350	
<input checked="" type="checkbox"/>	RDWCV	24	...		15	
<input type="checkbox"/>	NEUTRE	76.2	...		75.0	
<input type="checkbox"/>	LYMPRE	7.7	...		12.9	
<input type="checkbox"/>	MONORE	11.0	...		10.7	
<input type="checkbox"/>	EOSIRE	0.7	...		0.1	
<input type="checkbox"/>	BASORE	0.5	...		0.3	
<input type="checkbox"/>	IGRE	3.9	...		1.0	
<input type="checkbox"/>	NRBCRE	0.1	...		0.0	
<input type="checkbox"/>	NEUTAB	12.6	...		14.7	
<input type="checkbox"/>	LYMPAB	1.3	...		2.5	
<input type="checkbox"/>	MONOAB	1.8	...		2.1	
<input type="checkbox"/>	EOSIAB	0.1	...		0.0	
<input type="checkbox"/>	BASOAB	0.1	...		0.1	

- Validate Selection.



- This will release the CBC parameters that were checked, holding back the diff and WBC result fields.
- Review the slide and determine if a manual differential should be performed. If the number of NRBC's is in disagreement with the automated NRBC count (>10), perform a manual differential under the Diff tab.

Manual



- Complete your differential using the keypad.


Test	Key	COUNT 1		Abs 1	2
		1	%		
NEUTROPHIL	x	54	54%	...	
REACTIVE LYMPHOCYTES	c	32	32%	...	
REACTIVE LYMPHOCY...	w			...	
MONOCYTE	v	14	14%	...	
EOSINOPHIL	a		0%	...	
BASOPHIL	s		0%	...	
METAMYELOCYTE	d			...	
MYELOCYTE	f			...	
PROMYELOCYTE	g			...	
BLAST	h			...	
PLASMA CELLS	r			...	
HAIRY	t			...	
OTHER	e			...	
NRBC% M	q	10	10%		

- Click Save.
- Perform all “Morphology Checks” under the Morph tab (WBC estimate, Instrument flag check, RBC morphology and PLT estimate).


Manual



8. Return to the Result Validation tab. The WBC will need to be recalculated to correct for the manual NRBC's counted.
9. Obtain the original WBC count from the analyzer.
 - a. Select the order number from the analyzer IPU. Open the order.



- b. Select the service tab on the top and the WNR tab on the left side. The uncorrected WBC count will be listed in the TNC-N field.



10. Manually calculate the corrected WBC value using the below equation.

$$[\text{Uncorrected WBC}/(100 + \text{manual NRBC count})] \times 100$$

Example:

Uncorrected WBC: 16.599

Manual NRBC count: 10
 $[16.599 / (100 + 10)] \times 100$
 $[16.599 / 110] \times 100$
 0.1509×100
 15.09
 Corrected WBC: 15.1

11. Enter the corrected WBC into the WBC field in Caresphere. Click Save.



12. The manual diff absolutes will calculate to the right side of the screen.

Manual (Count 100)

Diff | Morph | Fluid Count

<input type="checkbox"/>	Test	Result	...	
<input type="checkbox"/>	NEUTROPHIL	54	...	
<input type="checkbox"/>	LYMPHOCYTE	32	...	
<input type="checkbox"/>	NRBC% M	10	...	
<input type="checkbox"/>	MONOCYTE	14	...	
<input type="checkbox"/>	EOSINOPHIL	0	...	
<input type="checkbox"/>	BASOPHIL	0	...	
<input type="checkbox"/>	NEUT ABSOLUTE	8.2	...	
<input type="checkbox"/>	LYMP ABSOLUTE	4.8	...	
<input type="checkbox"/>	MONO ABSOLUTE	2.1	...	
<input type="checkbox"/>	EOSI ABSOLUTE	0.0	...	
<input type="checkbox"/>	BASO ABSOLUTE	0.0	...	
<input type="checkbox"/>	NRBCMI ABSOLUTE	1.5	...	
<input type="checkbox"/>	WBC ESTIMATE	Confirmed	...	
<input type="checkbox"/>	INSTRUMENT FLAGS	Reviewed	...	

13. When all results are complete, Validate All.

- a. NOTE: If the WBC result field was released in Caresphere, the tech reviewing the slide will have to recalculate the WBC manually and make the correction in the LIS. The absolute differential results will need to be corrected as well.

N. MANUAL NRBC DISAGREEMENT FROM ANALYZER (XNL SITES)

1. Review Op alerts. Select the parameters that can be validated before slide review.

<input type="checkbox"/>	Test	Result	Comment
<input type="checkbox"/>	WBC	4.9 ...	
<input checked="" type="checkbox"/>	RBC	5.15 ...	
<input checked="" type="checkbox"/>	HGB	15.6 ...	
<input checked="" type="checkbox"/>	HCT	45.3 ...	
<input checked="" type="checkbox"/>	MCV	88 ...	
<input checked="" type="checkbox"/>	MCH	30 ...	
<input checked="" type="checkbox"/>	MCHC	34 ...	
<input checked="" type="checkbox"/>	PLT	250 ... 1	
<input checked="" type="checkbox"/>	RDWCV	13 ...	
<input type="checkbox"/>	NEUTRE	56.0 ...	
<input type="checkbox"/>	LYMPRE	33.1 ...	
<input type="checkbox"/>	MONORE	6.7 ...	
<input type="checkbox"/>	EOSIRE	2.4 ...	
<input type="checkbox"/>	BASORE	1.4 ...	
<input type="checkbox"/>	IGRE	0.4 ...	
<input type="checkbox"/>	NRBCRE	0.0 ...	
<input type="checkbox"/>	NEUTAB	2.7 ...	
<input type="checkbox"/>	LYMPAB	1.6 ...	
<input type="checkbox"/>	MONOAB	0.3 ...	
<input type="checkbox"/>	EOSIAB	0.1 ...	
<input type="checkbox"/>	BASOAB	0.1 ...	

2. Validate Selection.



3. This will release the CBC parameters that were checked, holding back the diff and WBC result fields.

4. Review the slide and determine if a manual differential should be performed. If the number of NRBC's is in disagreement with the automated NRBC count, perform a manual differential under the Diff tab.

Manual



5. Complete your differential using the keypad.

Test	Key	COUNT 1		Abs 1	2
		1	%		
NEUTROPHIL	x	54	54%		
REACTIVE LYMPHOCYTES	c	32	32%		
REACTIVE LYMPHOCY...	w				
MONOCYTE	v	14	14%		
EOSINOPHIL	a		0%		
BASOPHIL	s		0%		
METAMYELOCYTE	d				
MYELOCYTE	f				
PROMYELOCYTE	g				
BLAST	h				
PLASMA CELLS	r				
HAIRY	t				
OTHER	e				
NRBC% M	q	10	10%		

6. On the left side of the screen, select which WBC to release.

<input type="radio"/>	WBC Uncorrected	7.46
<input type="radio"/>	WBC from Analyzer	7.00
<input checked="" type="radio"/>	WBC from Manual	6.66
	% Difference in WBC Value	4.8 %

7. Click Save.

8. Perform all “Morphology Checks” under the Morph tab (WBC estimate, Instrument flag check, RBC morphology and PLT estimate).

Manual

Diff | Morph | Fluid Count

Test	Result
<input type="checkbox"/> SMEAR 1	...

9. Return to the Result Validation tab. Notice that the WBC selected from the diff tab populates in the WBC field. Validate All.

 Validate All

- a. NOTE: If the WBC result field was released in Caresphere, the tech reviewing the slide will have to recalculate the WBC manually and make the correction in the LIS.

O. OTHERS RESULTED

If others are counted in the manual differential, a comment must be added to explain what the other cells are referring to.

1. Validate CBC or Validate Selection after addressing the flags, Op alerts, etc.



2. Review the slide under the microscope. Perform a manual differential to include the others count. Add morphologies to the morphology tab if applicable.

3. Add a comment next to the others absolute result by clicking the blue comment icon.



4. Enter free text "Others are _____." Click Save.

5. Add a path review in Caresphere, if needed. (see step K above).

6. Validate All when complete.

7. Follow site specific workflow for sending slides and paperwork to the pathologists, if applicable.