

LIVER CORE BIOPSY

RA.SP.PR.GR.BX.11.r03

- I. Purpose**
To provide a procedure for the gross examination of a liver core biopsy specimen (including cores received for mass and medical liver biopsies).
- II. Principle**
To submit all tissue cores so that a diagnosis can be made microscopically by a pathologist.
- III. Equipment**
 - 1. Ruler
 - 2. Forceps
 - 3. Lens paper
- IV. Special Safety Precautions**
 - 1. PPE should be worn.
 - 2. **FORMALIN** is a known carcinogen.
- V. Supplies and Reagents**
 - 1. **10% NEUTRAL BUFFERED FORMALIN** (pH range 6.9-7.2)
 - 2. **Hematoxylin**
- VI. Quality Control**
All tissue should be submitted.
- VII. Limitations/Notes**
The following may influence the validity of test results:
 - 1. The specimen should be fixed in formalin.

LIVER CORE BIOPSY

VIII. Procedure

1. The specimen will arrive in a jar of formalin.
2. Using forceps, carefully remove all cores from the container.
3. Dictate the number, the size range (smallest length to largest length in cm.), the color, and the consistency of the cores.
4. Submit all cores by wrapping them in lens paper (preferred) or filtering in biopsy bags.
5. Ink all cores with hematoxylin.
6. If a mass is suspected, divide the largest cores into a maximum of 5 blocks. Add any extra cores to the existing 5 blocks, evenly. Serial sections X4 should be ordered on all 5 blocks (sections 1 and 4 should be stained with H&E while sections 2 and 3 should be left unstained on plus slides for possible immuno studies).
7. For all other liver biopsies (i.e. medical liver biopsies) not suspicious for mass, submit entirely with no more than two cores per block. Evenly distribute the cores, making sure the largest core is in the first block. A masson trichrome and PAS with Digestion should be ordered on the first block and an H&E on any additional blocks. Please note: the masson trichrome and PAS-D should be ordered on only one block per case, not per part.
8. The cassette or cassettes should be loaded on the appropriate processor as to allow at least 6 hours of fixation.

IX. References

1. Hruban RH, Westra, WH, Phelps, PH, & Isacson, C: Surgical Pathology Dissection An Illustrated Guide, New York, NY, Springer-Verlag Inc., 1996.
2. Lester, SC: Manual of Surgical Pathology, New York, NY, Churchill Livingstone, 2001.

X. Authorized Reviewers

1. Medical Director, Anatomic Pathology
2. Chief, Surgical Pathology

XI. Document Control

Location of Master: Master electronic file stored on the Beaumont Laboratory server under S:\AP_Grossing_Manual

Number of Controlled Copies posted for educational purposes: 0

Number of circulating Controlled Copies: 1

Location of circulating Controlled Copies: Master Grossing Manual located in Surgical Pathology

Printed copies of this document are not considered up-to-date. Please verify current version date with online document.

Document History

Printed copies of this document are not considered up-to-date. Please verify current version date with online document.

