

## ABO DISCREPANCIES

### Principle

ABO testing may give unexpected reaction patterns, some possible causes are listed below:

- ABO subgroups
- Autoagglutinins/excess plasma proteins
- Hypogammaglobulinemia
- Cold reacting allo/autoantibodies

This procedure shall serve as a guide for techs to recognize and utilize the appropriate management for ABO discrepancies that may present during routine testing.

### Clinical Significance

Forward and reverse typing need to agree so that the proper blood type is interpreted and thus a safe transfusion is assured.

### Specimen

Patient Preparation: No special preparation of the patient is required prior to specimen collection. Blood should be collected by approved techniques.

Requirements: K<sub>2</sub> EDTA pink or lavender top tube is preferred.

Minimum Volume:

- Adult: 3.0 mL whole blood
- Pediatric: 2 - K<sub>2</sub> EDTA microtainers (each with 300-500 uL) or cord blood specimen

Specimen Stability:

If stored at room temperature 15-30°C, stable for testing for 24 hours.

If stored 2-8°C, stable for testing for 72 hours.

Storage: 2-8°C for a minimum of 7 days after transfusion, or 10 days post crossmatch.

Rejection Criteria: Hemolysis. \*\*In rare occasions where sample cannot be redrawn, hemolyzed specimen may be used for testing as long as the testing personnel can accurately interpret the reactions.

### Reagents

Reagent	Storage	Stability
Isotonic Saline	15-30°C	Unopened: exp date Open: 1 month.
Biotestcell A1 and B reagent cells, BIO-RAD	2-8°C	Exp date
Anti-Human Globulin, Anti-IgG, Rabbit, BIO-RAD.	2-8°C	Exp date

### Instrumentation/Equipment

1. 10 x 75 mm tubes
2. centrifuge
3. pipettes


## Quality Control

To be performed once per day of testing. Refer to Blood Bank Quality Control Procedure for more information.

Reagent to be Tested	Positive Reactivity
A1 Reagent Cells	CorQC Reagent Antisera
B Reagent Cells	CorQC Reagent Antisera
Anti-IgG Reagent Sera	Coombs Check Cells

## Procedure

Technical or clerical errors are a common source of ABO discrepancies; repeat the ABORh with a new suspension, ensuring all tubes are labeled properly and double checking the reagents in use. Once technical or clerical error has been ruled out, it is necessary to determine the source of the discrepancy so appropriate management can take place.



*Chart is for guidance purposes, refer to following procedures for resolution*

### **Suspect discrepancy due to Rouleaux**

Excess plasma proteins can cause pseudo-agglutination, manifesting as additional reactions. Microscopic examination of the patient control will show the classic “stack of coins” appearance. **This discrepancy may appear in both the forward and the reverse typings.** The patient’s clinical diagnosis may be helpful, as abnormalities with globulin levels are often associated with conditions such as multiple myeloma or lymphoma.

1. Recentrifuge the tubes where rouleaux was observed.
2. Remove the supernatant, leaving the red cell button.
3. Replace the supernatant with equal volume saline (2 drops).
4. Resuspend the cell button gently and read for agglutination. Rouleaux will disperse when suspended in saline. True agglutination is stable in the presence of saline.
5. If agglutination persists, consider other potential causes, such as cold autoantibody or ABO subgroup.

### **Suspect discrepancy due to subgroup of A**

Subgroups of A, such as A<sub>2</sub>, may have Anti-A<sub>1</sub> present in their plasma that binds to the A<sub>1</sub> red cells ( $\leq 2+$  reactivity) used for reverse typing. A<sub>1</sub> lectin testing and reagent A<sub>2</sub> cells should be used for resolution. AB blood groups may also have subgroup of A.

**If rare subgroup showing weakened red cell reactivity (A<sub>x</sub>, A<sub>el</sub>) is suspected, notify lead.**

-Refer to UPM BB Procedure 03: Anti-A1 Lectin and Anti-A1 Testing

### **Room Temperature Incubation for hypogammaglobulinemia**

Plasma from elderly (or immune-suppressed) patients may have lower levels of antibodies, demonstrating significantly weaker reactivity in the reverse typing at immediate spin. Confirm age and diagnosis.

1. Incubate the A<sub>1</sub>, B, and patient control tubes at room temperature for 5-30 minutes.
  - a. Adding 2 extra drops of plasma may enhance reactivity
  - b. Incubating for the full 30 minutes may enhance reactivity
2. After incubation, mix well, centrifuge, and read for agglutination, using a viewer if necessary.
3. If RT incubation fails to resolve discrepancy, consider other potential causes, such as bone marrow transplantation or recent transfusion.

### **Warm saline wash for forward discrepancy due to cold autoagglutinin**

With sufficiently high titers of autoagglutinin, patient red cells may spontaneously agglutinate during centrifugation, causing a discrepancy in the forward typing. Washing the patient red cells in warmed saline may help manage interference from strong cold autoantibodies. ABO antisera should **never** be incubated at 37°C.

1. Aliquot a small amount of patient red cells for washing. A 12x75 tube may be used if desired.
2. Prepare a 3-5% cell suspension and incubate at 37°C for 15 minutes.

-If strongly reacting cold autoantibody encountered, incubate for 30-60 minutes
3. Hand wash patient red cells three times with 37°C warmed saline.
4. Decant last wash fully and prepare a 3-5% cell suspension with warmed saline for testing.
5. Perform the forward typing immediately, and read for agglutination.
6. If warm saline wash fails to resolve discrepancy, consider other potential causes, such as mixed field reactivity due to recent transfusion.

**Prewarm procedure for reverse discrepancy due to cold autoagglutinin**

1. Incubate several drops of patient plasma, one drop patient cell suspension, and one drop of reagent A<sub>1</sub> cells, and reagent B cells separately in 37°C incubator for 5 to 10 minutes.  
 -If strong antibody, incubate the full ten minutes.
2. Add 2 drops of plasma to A<sub>1</sub>, B, and patient cells without removing tubes from the incubator.
3. Incubate 30-60 minutes. If strong antibody, incubate for 60 minutes.
4. After incubation, immediately centrifuge and read for reactivity.
5. If prewarm incubation fails to resolve discrepancy, consider other potential causes, such as rouleaux or ABO subgroup.

**Reporting Results**

Determine the individual’s group based on the presence or absence of agglutination as follows:

Anti-A	Anti-B	A1 Cells	B Cells	Patient Control	Group Interpretation
NEG	NEG	POS	POS	NEG	O
POS	NEG	NEG	POS	NEG	A
NEG	POS	POS	NEG	NEG	B
POS	POS	NEG	NEG	NEG	AB


1. Enter results and interpretation into the LIS
2. Enter a comment for appropriate technique used to resolve discrepancy.

**References**

1. ARC Reference Laboratory
2. AABB. (2014). *Technical Manual* (18th ed.). Bethesda, MD: Author.
3. Bio-Rad Reagent Red Blood Cells Biotestcell A<sub>1</sub>&B insert, rev 18613/09 Aug 2014





**POLICY CREATION :**

<b>Author:</b> Kathy Maher	May 28, 2002
<b>Medical Director:</b> Douglas McGrady, MD	May 28, 2002

<b>MEDICAL DIRECTOR</b>		
DATE	NAME	SIGNATURE
	Elizabeth A. Bauer-Marsh, M.D.	

REVISION HISTORY (began tracking 2011)			
Rev	Description of Change	Author	Effective Date
1	Minor formatting changes, added Document ID to header. Added new LIS steps.	S. Schaffer	9/15/11
2	Updated reagents to BIO-RAD and changed reagent storage temp.	Kathy Turpin	12/6/13
3	Removed LIS specific reference	Vincent Strow	2/10/16
4 (Revision)	Updated Policy to cover more ABO discrepancies and maintenance thereof. Updated reverse discrepancy prewarm procedure in accordance with AABB 18 <sup>th</sup> edition technical manual. Added procedures for discrepancy maintenance in accordance with AABB 18 <sup>th</sup> edition technical manual and manufacturer's inserts.	Vincent Strow	5/4/17

**Reviewed by**

Lead	Date	Coordinator/ Manager	Date	Medical Director	Date
K. Maher	12/30/11		12/19/11	Elizabeth A. Bauer MD	10/19/11
D. Allen	1/20/14		1/20/14	Elizabeth A. Bauer MD	1/21/14
V. Strow	2/10/16		7/8/16	Elizabeth A. Bauer MD	7/11/16
V. Strow	5/8/17		5/8/17	Elizabeth A. Bauer MD	5/9/17