



Original Effective Date:	10/2016
Last Approved Review Date:	10/2016
Last Approved Revision Date:	10/2016
Responsible Department:	Laboratory
Responsible Person:	Laboratory Medical Director

Signature:

A handwritten signature in black ink, appearing to be 'C. M. M.', is written over a horizontal line.

Fetal Hemoglobin by Flow Cytometry

Purpose

This procedure provides instructions for how to send out testing for fetal hemoglobin by flow cytometry

Background

The passage of erythrocytes from a Rh positive fetus into the circulation of a Rh negative mother results in the formation of specific Rh antibodies. In subsequent pregnancies, the Rh antibodies formed in the blood serum of the Rh negative mother are readily transmissible through the placenta into the circulation of the fetus. The action of the antibodies on the Rh positive cells of the fetus may result in a disease entity recognized as iso-hemolytic disease, or erythroblastosis.

Flow Cytometry uses light that passes through or scatters and hits different points in the chamber. Based on the scattered light pattern, the cytometer program can differentiate fetal hemoglobin which has been tagged from adult hemoglobin. This test is much more precise and specific for fetal hemoglobin.

Specimen

Mother's blood

- Anticoagulated (EDTA)
- Refrigerated or ambient
- Specimens should be tested < 72 hours from collection.
- Hemolyzed specimens are not acceptable

Materials

Reagents	Supplies	Equipment
<ul style="list-style-type: none">• NA	<ul style="list-style-type: none">• Barnes-Jewish Flow Cytometry Immuno Phenotyping requisition• "Stat Testing" form• Biohazard bag	<ul style="list-style-type: none">• NA

Quality Control

NA


Fetal Hemoglobin by Flow Cytometry

BBMB 8.0

Procedure

1. Test Ordering
2. Determination of Test Flow
3. Computer Entry and Form Preparation
4. Packaging Sample
5. Arranging for Courier Service
6. Receiving and Reporting of Results

1. Test Ordering

Step	Action
1	Nursing should order Kleihauer-Betke testing on the OM side.
2	If ordered on the lab side of MCARE, order KB in the Enter/Edit tab of the Specimen tab. 

2. Determination of Test Flow



Step	Action				
1	Tech must determine if the fetal hemoglobin stain will be performed or if the testing will be sent to Barnes-Jewish hospital for flow cytometry testing. <ul style="list-style-type: none">• Note: all efforts should be made to send testing for flow as this is the better testing and technique for fetal hemoglobin determination.				
2	Refer to the chart below.				
3	<table><tr><th>Send for Flow Cytometry if any of the following...</th><th>Perform FHS if any of the following....</th></tr><tr><td><ul style="list-style-type: none">• M-F between 0600-1900• Routine• For determination of RhIg to be given after normal delivery• Requested by physician</td><td><ul style="list-style-type: none">• Weekends or holidays• After 1900 and before 0600 M-F• Doctor determines it cannot wait• It is a weekend or holiday and patient will be discharged before next regular run of flow.</td></tr></table>	Send for Flow Cytometry if any of the following...	Perform FHS if any of the following....	<ul style="list-style-type: none">• M-F between 0600-1900• Routine• For determination of RhIg to be given after normal delivery• Requested by physician	<ul style="list-style-type: none">• Weekends or holidays• After 1900 and before 0600 M-F• Doctor determines it cannot wait• It is a weekend or holiday and patient will be discharged before next regular run of flow.
Send for Flow Cytometry if any of the following...	Perform FHS if any of the following....				
<ul style="list-style-type: none">• M-F between 0600-1900• Routine• For determination of RhIg to be given after normal delivery• Requested by physician	<ul style="list-style-type: none">• Weekends or holidays• After 1900 and before 0600 M-F• Doctor determines it cannot wait• It is a weekend or holiday and patient will be discharged before next regular run of flow.				
4	If it is determined that the testing will be performed by staining, then refer to the <i>Fetal Hemoglobin Stain</i> procedure.				
5	If the testing will be sent to Barnes-Jewish hospital, continue to next section.				

3. Computer Entry and Form Preparation

Step	Action
1	Log in to MCARE

Fetal Hemoglobin by Flow Cytometry

BBMB 8.0

2	Go to the Specimen tab in the blood bank module.	
3	<p>Click on Enter Results</p> 	
4	Type in Specimen M#, patient name or BB# to retrieve testing.	
5	Select the specimen with the correct test to result	
6	<p>Enter FLOW KB in the Which Test? field, or press F9 for options.</p> 	
7	The next field will appear which is Previous BT? Enter Y if the patient has a historical blood type or N if the patient does not have a historical blood type.	
8	<p>If you entered “N”...</p> <ul style="list-style-type: none"> • The system will order an ABORH • Perform the ABORH and enter the results before proceeding. • See appropriate SOP for directions on performing and resulting ABORH. • Proceed to step 9. 	<p>If you entered “Y”....</p> <ul style="list-style-type: none"> • Proceed to step 9.
9	<p>In the PAT’S RH TYPE? field, enter the patient’s Rh type from the current ABORH or historical ABORH.</p> <ul style="list-style-type: none"> • N for negative • P for positive 	
10	Press F12 or click Save .	
11	Obtain the <i>Barnes-Jewish Flow Cytometry Immuno Phenotyping</i> requisition form and the <i>BJH Stat Testing</i> form from the forms drawer. See <i>Appendix A-BJH Flow Cytometry Immuno Phenotyping Form</i> and <i>Appendix B-BJH Stat Testing Form</i> .	
12	<p>Fill out the following at the top of the flow req.</p> <ul style="list-style-type: none"> • Patient’s first and last name • Patient’s sex • Patient’s DOB • Patient’s MRN in the space marked “Reference #” 	

Fetal Hemoglobin by Flow Cytometry

BBMB 8.0

	<p>e.g.</p>
13	On the Stat testing form, fill in <ul style="list-style-type: none"> • Patient's first and last name • Patient's DOB • Patient's MRN in the space marked "Registration #"
14	Call Flow Cytometry lab at 314-362-4628 to inform them that a sample is coming.
15	Proceed to next section.

4. Packaging Sample





Step	Action
1	Place sample in a small biohazard bag.
2	Place the requisition with the Stat testing form attached in the front pocket of the bag.
3	Seal the pocket with a Barnes-Jewish lab shipping label. See <i>Appendix C</i> .
4	Place bag in area for Lab Express pickup.

5. Arranging for Courier Service

Step	Action
1	Open Internet Explorer
2	Type www.labexp.com
3	Click on Sign In in the top right corner.


Fetal Hemoglobin by Flow Cytometry

BBMB 8.0

4	<p>The user name is memorial and the password is lab2@16. Make sure it is for “Lab”.</p> 																																																	
5	<p>Click Enter Order at the top of the screen.</p>  <p>Select Manifest Date:</p> <table><tr><th colspan="7">October 2016</th></tr><tr><th>Su</th><th>Mo</th><th>Tu</th><th>We</th><th>Th</th><th>Fr</th><th>Sa</th></tr><tr><td>25</td><td>26</td><td>27</td><td>28</td><td>29</td><td>30</td><td>1</td></tr><tr><td>2</td><td>3</td><td>4</td><td>5</td><td>6</td><td>7</td><td>8</td></tr><tr><td>9</td><td>10</td><td>11</td><td>12</td><td>13</td><td>14</td><td>15</td></tr><tr><td>16</td><td>17</td><td>18</td><td>19</td><td>20</td><td>21</td><td>22</td></tr><tr><td>23</td><td>24</td><td>25</td><td>26</td><td>27</td><td>28</td><td>29</td></tr></table>	October 2016							Su	Mo	Tu	We	Th	Fr	Sa	25	26	27	28	29	30	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29
October 2016																																																		
Su	Mo	Tu	We	Th	Fr	Sa																																												
25	26	27	28	29	30	1																																												
2	3	4	5	6	7	8																																												
9	10	11	12	13	14	15																																												
16	17	18	19	20	21	22																																												
23	24	25	26	27	28	29																																												
6	<p>Click on the arrow to the right of the Pickup Date field and choose the date.</p>																																																	
7	<p>Place cursor in the Pickup Facility field and type “Belleville Memorial” and then wait for the menu to appear.</p> <p>› Location</p> <p>count # or Search by Pickup Facility Name or Address to select an existing location. ation is not found you can enter a new location by first entering Pickup Facility name.</p> <p>Pickup Date: <input type="text"/>  </p> <p>Account #: <input type="text"/></p> <p>Notes:</p> <p>→ Pickup Facility: <input type="text" value="Search for location by typing first few letters of Location Name or Address"/></p> <p>Street Address: <input type="text"/></p>																																																	

Fetal Hemoglobin by Flow Cytometry

BBMB 8.0

8	<p>Choose the Belleville Memorial Lab option as shown below.</p> <p>Notes:</p> <p>Pickup Facility: Belleville Memorial</p> <p>Street Address:  Belleville Memorial Hospital Lab 4500 Memorial Drive Belleville, IL-62226</p> <p>Belleville Memorial Hospital Material Management 4500 Memorial Drive Belleville, IL-62226</p> <p>Belleville Memorial Hospital Pharmacy 4500 Memorial Drive Belleville, IL-62226</p> <p>Belleville Memorial Hospital Social Services 4500 Memorial Drive Belleville, IL-62226</p> <p>Belleville Memorial Hospital East Pharmacy 1404 Cross Street Shiloh, IL-62269</p> <p>Belleville Memorial Hospital Central Receiving 4500 Memorial Drive Belleville, IL-62226</p>
9	The address will populated automatically
10	In the Pickup Instructions field, type “Pick up at blood bank window”.
11	<p>In the Existing Contacts field, highlight the Blood Bank Technologist by clicking once on it.</p> <p>Existing Contacts</p> <p>Blood Bank Technologist 618-257-5093 jharris@memhosp.com</p> <p>Evan Evans 6182576605 eevans@memhosp.com</p> <p>Jan Greenwood 618-257-5101 jgreenwood@memhosp.com</p> <p>Memorial Laboratory Courier 618-257-6605 labcourier@memhosp.com</p>
12	Enter the number of specimens that you are sending in the appropriate field.
13	Check the box next to Ambient for the temperature
14	<p>Click the priority as Stat</p> <p>Number of Specimens: <input type="text" value="1"/></p> <p>Type of Specimen: <input checked="" type="checkbox"/> Ambient <input type="checkbox"/> Refrigerated <input type="checkbox"/> Frozen <input type="checkbox"/> Other</p> <p><input type="checkbox"/> Pharmacy <input type="checkbox"/> Surgical Supplies</p> <p>Type of Pickup Needed: <input type="radio"/> Regular <input checked="" type="radio"/> Stat <input type="radio"/> Special <input type="radio"/> Lockbox Sweep</p> <p><input type="radio"/> Other <input type="radio"/> Supplies <input type="radio"/> Mission Critical <input type="radio"/> Will Call</p>
15	For the delivery location, type “Barnes” in the field and wait for the menu to pop up
16	Choose the “Barnes Jewish Lab” option that matches the address on the shipping label you put on the specimen.
17	In the Pickup Instructions field, type “Deliver to Flow Cytometry Department”.
18	Click Submit Order at the bottom of the screen.
19	Lab Express should pick up specimen at blood bank window within 1 hour.

6. Receiving and Reporting of Results

Step	Action
1	The results should be called to the blood bank as well as faxed to the blood bank.
2	Once a hard copy of the results are obtained, go into the Specimen tab and then into Enter

Fetal Hemoglobin by Flow Cytometry

BBMB 8.0

	Results.																																			
3	Enter the patient's information to retrieve the KB testing screen.																																			
4	Enter the results from the Barnes report listed as Fetal RBC Pct-Flow in the FETAL RBC FLOW or FETAL RBC PCT field. (This will depend on whether the patient was Rh negative or Rh positive). See <i>Appendix D</i> for example of Barnes Report. <ul style="list-style-type: none">NOTE: Result of < 0.05% MUST be entered as "0". MCARE will not calculate correctly if you put "<" in the field.																																			
5	<p>If the patient is Rh negative, the result will look similar to this.</p> <table><tr><td>1</td><td>WHICH KB TEST?</td><td></td><td>FLOW KB</td><td></td><td></td><td></td></tr><tr><td>2</td><td>PREVIOUS BT?</td><td></td><td>YES</td><td></td><td></td><td></td></tr><tr><td>3</td><td>PAT'S RH TYPE?</td><td></td><td>NEG</td><td></td><td></td><td></td></tr><tr><td>4</td><td>FETAL RBC PCT</td><td>x</td><td>0.11</td><td></td><td>%</td><td>0.0-0.0</td></tr><tr><td>5</td><td>VIALS OF RHIG</td><td>x</td><td>1</td><td></td><td></td><td>0-1</td></tr></table>	1	WHICH KB TEST?		FLOW KB				2	PREVIOUS BT?		YES				3	PAT'S RH TYPE?		NEG				4	FETAL RBC PCT	x	0.11		%	0.0-0.0	5	VIALS OF RHIG	x	1			0-1
1	WHICH KB TEST?		FLOW KB																																	
2	PREVIOUS BT?		YES																																	
3	PAT'S RH TYPE?		NEG																																	
4	FETAL RBC PCT	x	0.11		%	0.0-0.0																														
5	VIALS OF RHIG	x	1			0-1																														
6	<p>If the patient is Rh positive, the result will look similar to this.</p> <table><tr><td>13</td><td>PAT'S RH TYPE?</td><td></td><td>POSITIVE</td><td></td><td></td><td></td></tr><tr><td>14</td><td>FETAL RBC FLOW</td><td>x</td><td>1.51</td><td></td><td></td><td></td></tr><tr><td>15</td><td>FKB QUAL INTERP</td><td>x</td><td>LARGE</td><td></td><td></td><td></td></tr></table>	13	PAT'S RH TYPE?		POSITIVE				14	FETAL RBC FLOW	x	1.51				15	FKB QUAL INTERP	x	LARGE																	
13	PAT'S RH TYPE?		POSITIVE																																	
14	FETAL RBC FLOW	x	1.51																																	
15	FKB QUAL INTERP	x	LARGE																																	
7	Press F12 or click Save .																																			

Interpretation

○ For FKB Qualitative Testing

<u>% Acid Resistant Cells</u>	<u>Results</u>
<0.05	Negative
0.06-0.80	Small
0.81-1.49	Moderate
1.50 or Greater	Large

○ For FKB Quantitative

<u>% Acid Resistant Cells</u>	<u>Vials of RhIg</u>
0.0-0.2	1
0.3-0.8	2
0.9-1.4	3
1.5-2.0	4
2.1-2.6	5
2.7-3.2	6
3.3-3.8	7
3.9-4.4	8
4.5-5.0	9
>5.1	See tech manual

Fetal Hemoglobin by Flow Cytometry

BBMB 8.0

Result Reporting

Step	Action
1	After testing has been verified, place faxed copy of report in bin at the front desk area to be scanned into the patient's chart in HIM.

References

AABB. *Standards for Blood Banks and Transfusion Services--29th Edition*. Std. 5.30.2, 5.30.3, 5.30.5. Bethesda, MD: American Association of Blood Banks; 2014
AABB. Perinatal Issues in Transfusion Practice. *Technical Manual--18th Edition*. Bethesda, MD: AABB; 2014. Pages 565-566

Related Documents

Appendix A: BJH Flow Cytometry Immuno Phenotyping Form
Appendix B: BJH Stat Testing Form
Appendix C: Barnes-Jewish Hospital Shipping Label
Appendix D: BJH Flow Cytometry Fetal Hemoglobin Report Example
Appendix E: Fetal Hemoglobin by Flow Cytometry Flowchart

Appendix A: BJH Flow Cytometry Immuno Phenotyping Form

Washington University in St. Louis SCHOOL OF MEDICINE		FLOW CYTOMETRY IMMUNO PHENOTYPING 314-362-4628 St. Louis, Missouri 63110		BARNES JEWISH Hospital	
COLLECTION INFORMATION: <input type="checkbox"/> AM <input type="checkbox"/> PM DATE: _____ TIME: _____ INITIALS: _____		PATIENT'S NAME (LAST) (FIRST) (MI) (SEX) (DATE OF BIRTH) (PATIENT'S SS #) _____			
ACCOUNT INFORMATION NAME: Memorial Hospital-Belleville ADDRESS: Laboratory-Sendouts 4500 Memorial Drive CITY: Belleville STATE: IL ZIP: 62226 PHONE: 618-257-5093 ORDERING PHYSICIAN: _____		PATIENT'S ADDRESS CITY STATE ZIP PHONE _____			
BILL TO: <input type="checkbox"/> ACCOUNT <input type="checkbox"/> PATIENT/INSURANCE <input type="checkbox"/> ALTHOUSE BLVMH		REFERENCE # _____ DIAGNOSIS _____			
SEND ADDITIONAL COPY OF REPORT TO: Blood Bank 618-257-6601 CLIENT NUMBER/PHYSICIAN NAME PHONE/FAX NUM PHYSICIAN'S ADDRESS CITY, STATE, ZIP BJH REGISTRATION # _____		PATIENT'S RELATIONSHIP TO RESPONSIBLE PARTY <input type="checkbox"/> 1-SELF <input type="checkbox"/> 2-SPOUSE <input type="checkbox"/> 3-CHILD <input type="checkbox"/> 4-OTHER NAME OF RESPONSIBLE PARTY (IF DIFFERENT FROM PATIENT) SOCIAL SECURITY (REQUIRED SSN) ADDRESS OF RESPONSIBLE PARTY APT # CITY STATE ZIP MEDICARE # STATE MEDICARE # (INCLUDE PREFIX/SUFFIX) <input type="checkbox"/> PRIMARY <input type="checkbox"/> SECONDARY MEDICARE RETIREMENT OR DISABILITY DATE INSURANCE COMPANY NAME PLAN CARRIER CODE SUBSCRIBER / MEMBER # LOCATION GROUP # INSURANCE ADDRESS PHYSICIAN'S PROVIDER # CITY STATE ZIP EMPLOYER'S NAME OR NUMBER WORKER'S COMP <input type="checkbox"/> YES <input type="checkbox"/> NO			
REGISTERED BY: _____		CLINICAL HISTORY AND DIAGNOSIS: _____ _____ _____			
SPECIMEN TYPE: <input checked="" type="checkbox"/> Peripheral Blood - Immune testing: <input type="checkbox"/> Peripheral Blood Leuk/Lym workup: <input type="checkbox"/> Peripheral Blood PNH: <input type="checkbox"/> Bone Marrow: _____ <input type="checkbox"/> Fluid: _____ <input type="checkbox"/> Tissue: _____ <input type="checkbox"/> FNA (Fine Needle Aspiration): _____ <input type="checkbox"/> Other: _____		Tube Type LAV LAV / 2 DK GRN LAV / DK GRN LAV / DK GRN / Core-Formalin FL		SAMPLE SUBMITTED: COLLECTION TIME: _____ AM _____ PM COLLECTION DATE: _____ MO _____ DAY _____ YEAR	
LABORATORY TEST <input type="checkbox"/> Lymphoma WorkUp (Lymphoproliferative disorder ex: CLL, NHL, HCL) <input type="checkbox"/> Leukemia WorkUp (Acute Leukemia ex AML, ALL, ANLL) <input type="checkbox"/> PNH Profile Includes RBC-CD59, WBC-CD59 and FLAER <input type="checkbox"/> Sezary Cell WorkUp <input type="checkbox"/> Other (Please Specify) _____ <input checked="" type="checkbox"/> Other Flow Test: Fetal Red Blood Cell Percentage <input type="checkbox"/> Lymphocyte Subpopulation 7 (CD3, CD4, CD8, CD19, CD16+56, CD2, HLA-DR) <input type="checkbox"/> Lymphocyte Subpopulation 13 (CD3, CD4, CD8, CD19, CD16+56, CD16, CD2, CD40, HLA-DR, TCR) <input type="checkbox"/> Immune Competence (CD3, CD4, CD8, CD19, CD16+56, CD4/CD8 Ratio) <input type="checkbox"/> Immune Deficiency (CD4, CD8, CD4/CD8 Ratio) <input type="checkbox"/> CD4 <input type="checkbox"/> CD45RA/CD45RO <input type="checkbox"/> Adhesion Panel (CD11a, CD11b, CD11c, CD15, CD18)					
CONTAINERS RECEIVED GRN DK GRN FL FLUID OT OTHER _____		COLLECTION TIME INITIALS _____ AM _____ PM			

Appendix B: BJH Stat Testing Form

STAT TESTING

☐

STAT TESTING/CALL BACK


☒

Patient Name: _____

DOB: _____

Registration Number: _____

Test name: **Fetal Hgb by Flow**

 Name: **Blood Bank**

Call back #: **618-257-5093**


Called Back: _____

Time/Initials: _____

To: _____

X1211-15 (6/02)

Appendix C: Barnes-Jewish Hospital Shipping Label


Barnes-Jewish Hospital Lab
Flow <u>Cytometry</u> Department
425 Euclid Suite 4701
(4 th Floor)
St. Louis, MO 63110

Fetal Hemoglobin by Flow Cytometry

BBMB 8.0

Appendix D: BJH Flow Cytometry Fetal Hemoglobin Report Example

BARNES-JEWISH
Hospital
BWH

Barnes-Jewish Hospital
One Barnes-Jewish Hospital Plaza
St. Louis, MO 63110-1003
314-362-1470
Medical Director: Dr. Charles Eby
CLIA #: 26D0438670

Memorial Hospital-Belleveille (BJH)
4500 Memorial Dr
Belleveille, IL 62226-

Patient: [REDACTED]
MRN: [REDACTED]
Patient Acct (FIN): [REDACTED]
DOB: [REDACTED]
Age/Gender: [REDACTED] Female
Location: BLVMH
Accession: 16-203-010389

Admit Date: 7/21/2016
Discharge Date:
Admitting Physician: Unknown ,Notinfile
Attending Physician: Unknown ,Notinfile
Ordering Physician: Unknown ,Notinfile

Hematology

Misc Hematology/Stains

Collected Date 7/21/2016
Collected Time 15:19 CDT

Test	Units	Reference
Fetal RBC Pct-Flow	<0.05 *	% [0.05-0.11]

Interpretive Data

i1: Fetal RBC Pct-Flow

Interpretative Data

Dosage of RhIG for the prevention of RhD alloimmunization following fetomaternal hemorrhage must be based on the individual clinical situation. The dosage recommendations provided here represent one published method and are voluntary.

Dosage of RhIG, in number of 300 mcg vials, is based on the estimated volume of fetomaternal bleed, calculated by estimating the maternal blood volume and multiplying by the percent of fetal cells in the maternal blood provided by the Kleihauer-Betke or flow cytometry test methods. The recommended dosage is the estimated volume of fetomaternal bleed, divided by 30 ml, rounded to the nearest integer, plus one. In the case of a negative screen for fetal cells, the recommended dose is one vial.

Fetomaternal bleed Recommended dose (vials)

Negative screen	1
<15 ml	1
15-44 ml	2
45-74 ml	3

As an example, if the estimated maternal blood volume is 5,000 ml and the percent of fetal cells by flow cytometry is 1.5%, the dose would be calculated as:
Estimated fetomaternal bleed = 5,000 ml X 0.015 = 75 ml
Recommend dose = 75 ml / 30 ml = 2.5, round and add one = 4 vials

These recommendations are published in the AABB Technical Manual, 17th ed, 2011, pp. 636-9. Further questions may be directed to the Laboratory Medicine Resident at 747-1320.

This test was developed and its performance characteristics determined by the Barnes Jewish Flow Cytometry laboratory. It has not been cleared or approved by the US


LEGEND: C= Critical, L=Low, H=High, *= Abnormal, ^= Corrected, F=Footnote, i=Interp Data, @=Performing Lab, OC= Order Comments

Report Request ID: 41832291

Page 1 of 2

Print Date/Time: 7/21/2016 19:00 CDT
n/a

Appendix E: Fetal Hemoglobin by Flow Cytometry Flowchart



Fetal Hemoglobin by Flow Cytometry

BBMB 8.0

PROCEDURE AND FORM CHANGE CONTROL

Title: Fetal Hemoglobin by Flow Cytometry										
Written		Validated		Path Review		Review		Effective		Reason for Revision
Date	By	Date	By	Date	By	Date	By	Date	By	
10/4/2016	JLH	10/5/16	GJM	10/14/16	ESB			10/21/16	JLH	
Revised										

Location of any copy(s) of the procedure:

Out of use:

Date: _____ By: _____ Reason: _____