

# • "Analyser" "interferences .... ?!"

- Steven Schischka
- Principal Scientist Haematology





- Delta check RBC parameters
- Sample found under centrifuge
- Staff member went through all the samples bags
- To look for the missing EDTA
- Labelled sample after finding form missing an EDTA
- Staff member was one of 3 working in the collection centre
- Why? Fear of getting in trouble



- Previously A Neg patient is not O Pos antenatal testing
- Previous collection dr's clinic
- Current collection pathology centre
- Cousin visiting from overseas and having antenatal testing done "for free" using Australian residents medicare card
- Intentional fraud



- Patient has two blood tests in one day
- One Normal FBE
- One Thalassaemic parameters
- Medicare card #1 in Melbourne
- Medicare card #2 in Sydney
- National pathology service with cumulative interstate results
- Intentional fraud

- INR result and dose phoned to patient
- Patient states this is impossible didn't have a blood test today
- Home visit pathology collector Definitely took blood of correct patient.
- Name, DOB, Address, phone number all have ticks on the form
- Pathology collector does drive by realises went to wrong house
- Home visit collector took blood off the neighbour
- Patient male not female
- Home visit INR patient too patient didn't query early visit by collector
- Failure to follow process due to familiarity





- Hospital inpatient FBE yesterday "normal" but today thalassaemic
- Collector identified patient by asking them are you "Jane Doe"
- Hospital collector took blood off "patient" after confirming ID
- "Patient" was a visitor holding the baby while mum showered
- "Patient" didn't speak English and didn't protest at having bloods taken



- Psych patient has sudden neutropenia and thrombocytopenia on clozapine
- Patient was re-bled to confirm pancytopenia FBE normal
- Entire ward bled to find the truly pancytopenic patient
- Patient not found
- In ward review the next morning the nurse admitted to removing a clot from the difficult syringe collect



- New collection system introduced
- WBIT theoretically impossible if used correctly
- WBIT identified by FBE parameters
- Two people sharing log-ins and a printer
- Post collection sign off being done in advance
- Collector(s) took wrong labels and labelled tubes away from bedside




- New collection system introduced
- WBIT theoretically impossible if used correctly
- WBIT identified by FBE parameters
- Printing devices don't have battery display and have short battery life
- When printer is recharged and turned on the Memory on the POC collection device reprinted 1 label and then 2 labels from the correct patient
- 2 other samples collected at same time were labelled correctly



- Delta check FBE on one of two trauma patients
- Unknown male and unknown female
- Unknown male has history from 4 days ago in electronic record
- Call to trauma unit patient has history, has he been re-admitted
- Trauma unit unknown male record had been able to be updated
- Patient re-identified and wrist band etc changed
- All bloods including X-match, x-rays etc re-ordered in the middle of a trauma

### The analyser lies


- Patient with fluctuating lymphocytosis
- Marked increase in target cells
- Non-lysed RBC
- Being counted as lymphocytes
- Review of scatterplots required at release




### The analyser lies



- Known CML on treatment historic basophilia
- Unusual cells on blood film degranular basophils
- Counted as lymphocytes and hiding behind 2 D scatter of lymphocytes





0.54			0.54
0.00		0.00	0.01
0.43		0.43	0.52
0.00		0.00	0.01
0.00		0.00	0.00
0.11		0.11	0.00
0.0	<<	0.0	2.4
80.0		80.0	96.4
0.0		0.0	1.1
0.0		0.0	0.1
20.0		20.0	0.0

### The computer lies

- 5 patients YYX000001 to YYX000009 have MCV's of around 62 fl,
- mildly anaemic and inverse differential, WBC scatter plots look strange
- Some previously normal and some first time in patients
- Happened just after the annual change over in lab numbers
- From YYX000001 to (YY+1)X000001
- A batch of porcine samples had been manually processed a week earlier as sample number 01-20 with no prefix
- IT system automatically inserts X00000 into the lab numbers and holds them in queue for a week
- YY was not automatically added
- Results from previous year inserted into current year due to annual cycle of lab numbers and results being held for one week







- Daily maintenance performed
- New QC loaded
- APTT, INR, FIB, TCT all incorrect
- Ddimer correct

- Routine QC reconstituted
- Ddimer QC is liquid stable
- Staff member had performed saline replacement in the middle of the night
- Saline bottle put with water bottles
- Saline used to reconstitute QC material

### **The Computer lies**



- Dr phones with complaint that abnormal blood film was not reviewed
- Lymphocytes are below normal range but film was not reviewed
- Lymphocytes are 1.5 on laboratory display with reference range 1.5 4.0
- Suspicion that results are from different laboratory
- Scientist asks Dr to fax the results in
- Results are in fact from the laboratory
  - Lymphocytes on report are 1.47 with abnormal low flag
  - Lab results uses script A to generate a cumulative display report for internal use
  - External results uses script B to generate a single display report for external use
  - Decimal places are different on the two reports from the same lab.

# Normal




Run Screen - [CBCDiff STV]	1-1-			
Open Configure Hide Display Off Arrange Print	H <u>e</u> lp			
SID P699715	Additional Platelet	Parameters	Baso	Perox
Patient ID	PDW H 67.8	%		
Aspiration Date/Time 15-Jun-05 06:41:38 PM	PCT 0.13	%	<u>-</u>	
	MPC 24.3	g/dL		
Sample Type PATIENT Rack & Position 0 0	MPM 2.39	pg		
Instrument Number BERT	Large PLT 13	x10 <sup>9</sup> cells/L		
FOR LABORATORY USE ONLY		·		
	aramatara	RBC Volume		
		KDC VOlume		
HGB 154 g/L %Hypo 0. HCT 0.47 L/L %Macro 3.				
MCV 96.4 fL %Macro 0.		RBC HC	PROVUC	PLT Scatter
MCH 31.6 pg RBC Fragments 0.			RBC V/HC	
MCHC L 328 g/L RBC Ghosts 0.				
CHCM L 329 g/L NeutX 61				
CH 31.5 pg NeutY 74		Platelet Vol		
RDW 13.9 % MNx 9.	.5			18 - 3 - 4 - 4
HDW 23.3 g/L MNV 15				A State 1
PLT 120 x10° cells/L %MN 41		- Andrewson		
MPV H 11.2 fL %PMN 57				
Cellular HGB 15		HGB Trans		
				E P P P P P P P P P P P P P P P P P P P
Routine WBC Differential Morphology F				
% # MACRO +			J	
	+	San	nple/System Flags	
Neut 54.9 * 4.86 * x10° cells/ NRBC +		N-E	RBC #LUC %LUC #LYMP %LYMP W	BC#NEUT %NEUT #MONO
Lymph 26.7 * 2.37 * x10° cells/ LARGE PLT +	Baso Rate		%MONO #EOS %EOS	
Mono H 9.2 * 0.81 * x10° cells/	HAR DURING IN			
Eos 2.7 * 0.24 * x10° cells/L				
Baso 0.7 0.06 x10° cells/L				
LUC H 5.8 * H 0.52 * x10° cells/L	Perox Rate			

### Heat affected blood



Run Screen - [CBCDiff STV]	Off Arrange Print H	elo
SI VINCENTS Playback SID P81 Patient ID Aspiration Date/Time 27-3 Sample Type PA Rack & Position 0 0 Instrument Number BEI	4878 Sep-05 10:40:32 AM FIENT	PDW     58.9     %       PCT     H     0.50       MPC     21.2     g/dL       MPM     1.47     pg       Large PLT     11     x10° cells/L
FOR LABORATORY           Roume coc           WBC         10.42         * x10° cells/L           RBC         1.0.4         * x10° cells/L           HGB         114         * g/L           HCT         0.25         * L/L           MCV         81.8         * fL           MCH         437.6         * pg           MCHC         H 460         * g/L           CHCM         331         * g/L           CH         27.0         * pg           RDW         439.9         *%           HDW         55.4         * g/L           PLT         651         * x10° cells/L           MPV         7.7         * fL	Additional Routine P           %Blast Suspect         40.4           %Hyper         7.7           %Hypo         16.4           %Macro         12.9           %Micro         23.8           RBC Fragments         0.56           RBC Ghosts         0.11           NeutX         64.3           NeutY         70.6           MNx         6.7           MNy         3.0           %MN         64.0           %PMN         25.6           Cellular HGB         82	4 9 8 6 x10 <sup>12</sup> cells/L 3 8 <b>RBC Rate</b> 0 <b>RBC Rate</b>
utine WBC Differential           %         #           VBC         10.42         x10° cd           Neut         72.2         7.52         x10° cd           mph         11.1         1.16         x10° cd           tono         H         1.2         H         1.17         x10° cd           tono         H         1.2         H         1.17         x10° cd	ells/L MACRO · ells/L HYPO · ells/L HYPER · ells/L ANISO ·	

#### Contaminated sample – glucose drip





# Patient Artefact- MPO "Deficiency"

% % g/dL



	kun Scieen - [CDCDi	11.21	vj						
<u>O</u> p				<u>A</u> rrange	<u>P</u> rint	Help	)		
	ST VINCENTS Play	/back	(				Additional	ΡĿ	atele
		SID	P70504	16					
							PDW	H	65.8
	Patie		1				PCT		0.15
	Aspiration Date/1	ime	14lun	-05 05.38	11 PM				
						-	MPC		25.0
	Sample 7	ype	PATIE	NT			LIDLI		0.15
	Rack & Pos						MPM		2.15
		IUUII	00				Lorgo DLT		11

Dun Saraan ICPCDiff STM

WBC

HGB

MCV

MCH

мснс

CH

RDW

HDW

PLT

MPV

WBC

Lymph L

Eos

Baso

KATN //

LI

Neut L 1.0

Mono H 49.3

LUC|H|46.3

3.1

0.1

0.2

0.13

0.39

0.02

0.02

1.84

00.7

L

L

H| 6.35

H 5.96 x10° cells/L

x10° cells/L

x10<sup>9</sup> cells/L

x10<sup>®</sup> cells/L

x10° cells/L

x10° cells/L

pg x10° cells L |Large PLT| |11 Instrument Number BERT FOR LABORATORY USE ONLY RBC Volum Additional Routine Parameters Róuune coc 12.88 x10<sup>9</sup> cells/L %Blast Suspect 0.4 RBC H 6.13 x10<sup>12</sup> cells/L %Hyper 0.4 176 g/L %Hypo 1.3 HCT H 0.53 L/L %Macro 0.2 RBC HC fL 86.2 %Micro 1.0 RBC V/nC PLT Scane 28.8 pg RBC Fragments 0.01 x10<sup>12</sup> cells/L 334 g/L RBC Ghosts 0.01 x10<sup>12</sup> cells/L CHCM L 326 g/L NeutX Platelet Vol 27.9 pg Neut Y 13.7 % MNx 111.7 23.9 g/L MNy| |14.6 162 x10<sup>9</sup> cells/L %MN 8.5 9.5 fL %PMN 90.7 HGB Trans Cellular HGB | 172 | g/L Routine WBC Differential Morphology Flags RBC Rate % # ATYPS +++ 12.88 x10° cells/ Sample/System Flags

Baso Rate

Perox Rate


Long Long and the

#### B-NV LI

B2.50

Perox






- Repeating samples appear to correct the error
- Faulty valve with reagent flow back
- First sample was essentially working as a reagent prime

# Patient Artefact – Dimorphic Neutrophils, Some are actually eosinophils






Patient sample results showing low SSC


Sample showing good Eosinophil differentiation

"Your analyser is performing as expected"



# **Pseudo Cold Agglutinin - Antibiotic Antibodies**





### **Pseudo** Anaemia





Syringe collect with settled sample.

Usually Emergency department – why?

Inserting access port so don't want to stab patient again.

Syringe settles whilst completing access port requirements.

Not mixed before aliquoting into collection tubes

# **Delta Check MCV**



Tha**Alfrad** 

Haematology	03-04-2019 15:21	11-03-2019 09:27	07-03-2019 09:06	05-03-2019 23:09	04-03-2019 09:16	03-03-2019 06:38	02-03-2019 03:33	01-03-2019 20:51	01-03-2019 12:12	03-05-2016 05:40
FBE										
Hb	118	99	96	85	87	82	94	90	85	127
WBC	6.08	3.13	4.17	4.87	6.81	5.93	8.50	9.20	7.43	4.34
Platelets	475	563	585	473	532	473	514	511	505	253
Hct	0.37	0.32	0.30	0.27	0.27	0.26	0.28	0.28	0.28	0.36
MCV	81	66	64	63	62	62	62	61	63	120
MCH	26.0	21.0	20.0	20.0	20.0	19.0	20.0	20.0	19.0	42.0
MCHC	316	313	314	314	319	315	332	320	307	351
RBC	4.61	4.79	4.73	4.29	4.36	4.19	4.60	4.58	4.45	3.04
RDW	40.2	23.9	22.0	21.9	21.1	21.4	21.3	21.3	21.1	11.9
Neutrophils	4.51	1.24	2.45	3.03	4.89	4.40	6.69	8.11	6.26	3.69
Lymphocytes	0.65	0.62	0.46	0.49	0.55	0.39	0.51	0.19	0.24	0.32
Monocytes	0.70	1.06	0.77	0.86	1.01	1.00	1.28	0.85	0.87	0.32
Eosinophils	0.19	0.18	0.47	0.45	0.29	0.12	0.01	0.01	0.03	0.01
Basophils	0.04	0.04	0.03	0.04	0.07	0.03	0.02	0.04	0.03	0.02

- Patient on hydroxyurea for polycythaemia rubra vera
- Patient now on venesection for polycythemia rubra vera
- Patient investigated for anaemia and given iron therapy.....

# Historic Warfarin patient presents with bleeding

<ul> <li>APTT</li> </ul>	80.0	(25 – 35)	sec
• INR	18.3	(0.8 – 1.2)	INR
• Fib	4.5	(1.5 – 4.0)	sec
• TCT	>300	(<20)	sec
<ul> <li>Anti-Xa</li> </ul>	4.0	(None or <1.0)	IU/ml



- If Anti-Xa is heparin then APTT should be >300 and INR "normal"
- If TCT is dabigatran then Anti-Xa should be 0.0
- If warfarin TCT and Anti-Xa should be normal
- Patient under cardiologist AND general practitioner. Cardiologist has changed patient to rivaroxiban and GP has changed patient to dabigatran. Patient taking both drugs fulfilled at different pharmacies.
- 2. Patient has renal failure and has been on dabigatran. Medical advice to swap from direct dabi to riva is to stop one one day and start the other the next day. Patient has not cleared dabi due to renal failure and now has riva on board too.

Unexpected pancytopenia in hematoma patient

HB	35	(115-165g/L)
WBC	1.0	(4.0-11.0)
Neutrophils	0.5	
Lymphocytes	0.3	
Monocytes	0.2	
PLTS	18	



? Bleed due to thrombocytopenia.

Call to clinician with results and clinical query.....

Clinician purposely labelled sample as FBE so lab would process it over night otherwise "they would have had to wait for microbiology in the morning"



#### Grossly abnormal urine sample over night

Blood stained urine sample is swarming with

- bacteria
- WBC
- foreign bodies



- Female patient has presented with abdominal pain
- Speaks little English
- Has been asked to provide a urine sample which staff believe she has done
- Sample is in fact fecal sample and patient has severe diarrhoea but can only explain her tummy hurts

### Natural Therapies



- Patient being treated for lead poisoning
- Chelation therapy not working
- Turns out lucky charm/crucifix that he kisses is made of lead
- Patient unable to get pregnant
- Severe lead poisoning
- Chelation therapy not working
- Turns out herbal support from "overseas" made from seaweed in estuary down stream from heavy metals factory
- Patient with psychiatric issues develops pancytopenia
- ?Therapy related cytopenia
- Munchausen syndrome self medication with Busulfan
- Patient proceeds to bone marrow transplant which fails



# The florist, the Alphabet and the trauma

- Patient being taken to theatre for operation
- · Has had pre-meds and starts laughing
- "Ha-ha" you gave me the wrong flowers this morning
- I'm Jane Smith not Joan Smith.....
- Group A blood is hung for transfusion
- Patient states "I thought I was blood group O"?
- Transfusion continues......with Incompatible blood
- Patient survives.....Just
- 3 year old child in and out of hospital with sepsis
- Multiple exploratory surgeries as commonly gut bacteria isolate
- Cyclic Pattern and suspicion develops
- Mother is found to be injecting fecal matter into the drip line.
- Ward staff split down the middle believing or not believing mother

# **Antenatal patient ? HELLP syndrome**



•	APTT	>300.0	(25 – 35)	sec
•	INR	>10	(0.8 – 1.2)	INR
•	Fib	4.5	(1.5 – 4.0)	sec
•	TCT	19	(<20)	sec
•	DDimer	0.3	(<0.25)	
•	APTT MIX	>300		
•	INR mix	>10.0		
•	Anti-Xa	0.0		IU/ml
•	PLTS	375	(150-400)	x10^9/L

- I forgot to take the blue top tube
- so I took some of the green
- and some of the purple
- and poured them into the blue top tube.

