


4 ABO(FWD/RVS)/Rh-00

Start Time: 27/11/2019 11:50 - Completion Time: 27/11/2019 11:59

ABO: ? Rh: ?

Results



Result	Anti-A	Anti-B	Anti-D	Ctrl	A1-Cells	B-Cells
Original	4+	0	4+	4+	3+	3+
Modified						

Cassette ID: 280420-00-075591-32210-3
Reagent: A1 Cells
Reagent: B Cells

Lot: 32210
Lot: 0277
Lot: 0277

Exp: 28/04/2020
Exp: 17/12/2019
Exp: 17/12/2019

Flags: indeterminate result, result expired, above/below positive reaction threshold
Accepted By: rajis

DAT Poly

Start Time: 27/11/2019 11:50 - Completion Time: 27/11/2019 12:00

Poly: POS

Results

Cassette 1



Result				Poly		
Original				4+		
Modified						

Cassette ID: 210220-22-072607-68510-0

Lot: 68510

Exp: 21/02/2020

Flags:


4 AbScr Surg Poly

Start Time: 27/11/2019 11:50 - Completion Time: 27/11/2019 12:08

ABScr: POS

Results

Cassette 1



Result	Surg 1	Surg 2	Surg 3			
Original	3+	3+	3+			
Modified						

Flags:

Presenting Complaint

Lethargy
Tiredness
Dark urine

History of Presenting Complaint

69M, background of warm autoimmune haemolytic anaemia

Presents with 3/7 history of brown/orange urine, feeling unwell and lethargic
Denies fevers
Denies infective symptoms: no cough/LUTS/diarrhoea/vomiting/constipation
No recent foreign travel

Patients historical blood group is A Pos rhesus phenotype of R1r. History of anti-C like autoantibody detected.

Screening cells and A1 cells repeated at strict 37°C. A1 cells negative, screening cells still 3+ reactions being observed.

DAT showing 4+ reactions with Poly and Anti-IgG in immediate spin, while Anti-C3d has a 1+ reaction with patient cells. After 5min Anti-C3d is observed to have a 2+ reaction with patient cells. An elution is conducted and a non-specific pan-agglutinating antibody is eluted.

To determine if the patient has any underlying clinically significant allo-antibodies potentially masked by the warm auto antibody, an auto adsorption is required.

Performing an Auto-Adsorption:

Phenocell B 3% Antigen Composition Sheet																	Batch No: 1664 197			Expiry Date: 23/12/2019							
Cell No	Reference No	Rh Phen.	RH					KEL			FY		JK		MNS			P1PK	LE		LU	Additional Typings	Cell No	Results			
			D	C	E	c	e	C ^y	K	k	Kp ^a	Fy ^a	Fy ^b	Jk ^a	Jk ^b	M	N	S	s	P1	Le ^a			Le ^b	Lu ^a	VAL	
1	4135608	R ₁ *R ₁	+	+	0	0	+	+	0	+	0	0	+	+	0	+	0	+	0	+	0	0		1			
2	2146231	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	0	+	0	+	0	+	0	+	0	2			
3	2227853	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	0	+	0	+	0	+	0	+	0	3	0	✓	
4	4398907	R ₂ R ₂	0	0	+	+	0	0	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	4	0	✓	
5	2278162	R ₂ R ₂	0	0	+	+	0	0	0	+	0	0	+	0	+	0	+	0	+	0	+	0	+	5	0	✓	
6	2172354	r'r	0	+	0	+	+	0	0	0	0	+	+	0	0	+	0	+	0	+	0	+	0	6	0	✓	
7	2138870	r'r	0	0	+	+	0	0	0	0	0	+	+	0	0	+	0	+	0	+	0	+	0	7	0	✓	
8	2259224	rr	0	0	0	+	+	0	+	0	0	+	+	0	0	+	+	0	+	0	+	0	0	8	0	✓	
9	2055924	rr	0	0	0	+	+	0	0	0	0	+	+	0	0	+	0	+	0	+	0	0	0	9	0	✓	
10	2204055	rr	0	0	0	+	+	0	0	+	0	+	+	0	0	+	0	+	0	+	0	+	0	10			
11	2117666	rr	0	0	0	+	+	0	0	+	0	+	+	0	0	+	0	+	0	+	0	+	0	11			
Auto																								Auto			
	Reference No																										
	Ab1																								Ab1	0	✓
	Ab2																								Ab2	0	✓
	Ab3																								Ab3		
		Rh Phen.	D	C	E	c	e	C ^y	K	k	Kp ^a	Fy ^a	Fy ^b	Jk ^a	Jk ^b	M	N	S	s	P1	Le ^a	Le ^b	Lu ^a		Cord		
																										A ₁	
																										A ₂	
																										B	

This product is manufactured from blood donated by voluntary donors to the Australian Red Cross Blood Service (ARCB5).

Note: Cells 6, 7 and 8 form the RhD Negative Screening Cell Subset, formulated for antenatal RhD antibody screening in the presence of Prophylactic Anti-D.

Cs: A high incidence antigen (Ref: Reid M.E., Lomas-Francis C., Olsson M.L., "The Blood Group Antigen FactsBook", 3rd Edition Pgs. 648-649).

- Once it has been confirmed that the patient has not been transfused with the last 90 days and 3 EDTA blood bank tubes have been collected the auto adsorption can begin.
- Patient cells must first be papanized. The purpose of papanizing the patient red cells is to optimise the adsorption capacity of the red cells for the auto antibody.
- Once the cells are ready the patient plasma is used in a 1:1 ratio with the papanized red cells to adsorb the auto antibody from the plasma. As this patient has a strong antibody 3 adsorption processes had to take place to ensure the auto antibody has been adsorbed.
- The adsorbed plasma is then used to set up a panel by "classic" Antihuman Globulin technique. If not enough plasma a targeted panel can be created to help exclude all possible allo antibodies against their respective antigen. Ideally if there is some plasma left it can be used for a crossmatch.
- Auto adsorption procedure has been edited and can be viewed in fast track [CD_HA_0431](#)