


FMH CHECKLIST

AlfredHealth	Alfred Pathology Services	
	Document No: CD_HA_0716	Authorised By: Gosia Gorniak
	Version No: 2.1	Issue Date: 03/09/2019

FMH CHECKLIST

Use flow-chart to assist in determining what actions are required for FMH requests

Please box to acknowledge



Pre-Delivery

- Call Sandringham maternity mid-wife to determine urgency.
- If maternity insist FMH is to be run urgently out of hours, advise them to call the on call Haematologist through switch.
- Fill out mother's group and save details on the back of the page.

Post Delivery

- Fill out section once the baby blood group has been determined.
- If sample will be tested within **70 hours of collection**, please leave paperwork in Flow, if sample will not be tested within 70 hours contact on call Haematologist.

FMH CHECKLIST

AlfredHealth	Alfred Pathology Services	
	Document No: CD_HA_0716	Authorised By: Gosia Gorniak
	Version No: 2.1	Issue Date: 03/09/2019

FMH CHECKLIST

REFER ANY ISSUES TO THE LAB HAEMATOLOGIST.

I HAVE Actioned and completed the requirements above:^{TV}..... (sign)

MOTHER

Name: Smith, Jane

UR: 7030471

Accession No: 20-020-1555

Reported ABO/Rh

O NEG

Mixed field (in ABO/Rh)

Present Not present

If Mixed Field present in mothers ABO/RH and/or there is suspicion of significant foetal blood loos:

Urgent FMH - contact lab haematologist

Antibody Screen

Positive / Negative

Antibody ID

Prophylactic-D

BABY (cord blood)

Name: Smith, Jane - Baby

UR 7030001

Accession No: 20-020-2012

CORD date & time: 20/01/2020 1400

Blood Group

O POS

Group Check

O POS

DAT

NEG

DAT Pos AND Baby's ABO group is incompatible with the mothers reverse group include cells of the baby's ABO group in the elution.

Elution result

N/A

*** If Cord Rh(D)Negative:

Perform weak D

N/A

NOTE: If the cord DAT is Positive, determination of Weak D status not possible as cells already coated with IgG. INFORM HAEMATOLOGIST AFTER ELUTING AND IDENTIFYING THE ANTIBODY.

TV 20/01/2020 (Blood Bank Scientist - sign and date)

- Ensure all sections of the form are completed.
- The cord date and time should be of the time the sample was collected.
- The form should be signed and dated by the person who has completed it

Having this at the bottom of your checklist shows that the document has been printed off fast track. H:drive should no longer be used to print out procedures or documents.