

Patient Fiona Freeman, DOB 15/02/1953, UR 555454, presents to the emergency department with fainting and fatigue. The patient's blood type is determined to be O POS and the antibody screen is positive so you decide to run a panel. Results are as below.

Cell#	Rh-hr	Rh-hr										KELL						DUFFY			KIDD			Lewis			MNS			P		Special Antigen Typing	Test Results
		D	C	E	c	e	f	Cw	V	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Jsa <sup>1</sup>	Jsb <sup>1</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Xg <sup>a</sup>	Le <sup>a</sup>	Le <sup>b</sup>	S	s	M	N	P <sub>1</sub>	Lu <sup>a</sup>	Lu <sup>b</sup>				
1	Rh-r	+	+	0	0	+	0	+	0	+	0	+	/	+	+	0	+	0	+	0	+	+	+	+	+	+	+	+	+	+	+		
2	Rh-r	+	0	+	+	0	0	0	+	+	0	+	/	+	0	+	+	0	+	0	+	+	+	+	+	+	+	+	+	+			
3	"	0	0	0	+	+	0	0	0	+	+	0	/	+	0	+	+	0	+	0	+	+	+	+	+	+	+	+	+	+			
	Patient Cells																																

Reagent Red Blood Cells  
Surgiscreen®  
Ortho Clinical Diagnostics 2010


LOT NO. 3SS993  
EXP. DATE 2021-08-17  
CCYY-MM-DD

ANTIGRAM®  
Antigen Profile  
635201091

21/2/21  
SD

FIONA FREEMAN  
15-02-1953  
SS5454  
21-105-SS55

\* represents "Not Tested" for new donors.



# Phenocell B™

Reagent Red Blood Cells  
Antibody Identification Panel

### Phenocell B 3% Antigen Composition Sheet

Cell No	Reference No	RH														KEL				FY			JK			MNS					PIPK		LE		Additional Typings	Cell No	Results	
		D	C	E	c	e	C <sup>w</sup>	K	Kp <sup>a</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	M	N	S	s	P1	Le <sup>a</sup>	Le <sup>b</sup>																		
1	2083233	R <sub>1</sub> R <sub>1</sub>	+	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0							
2	3823479	R <sub>1</sub> R <sub>1</sub>	+	+	+	0	0	+	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0					
3	1301505	R <sub>1</sub> R <sub>1</sub>	+	+	+	0	0	+	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0					
4	2008220	R <sub>2</sub> R <sub>2</sub>	+	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0			
5	5175005	R <sub>2</sub> R <sub>2</sub>	+	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0			
6	4379046	r <sup>r</sup>	+	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0			
7	5727735	r <sup>r</sup>	+	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0			
8	4086345	r <sup>r</sup>	+	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0			
9	4257497	r <sup>r</sup>	+	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0			
10	4442426	r <sup>r</sup>	+	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0			
11	2054349	r <sup>r</sup>	+	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0	+	+	0	0			
Auto																																						

Batch No: 1664 219      Expiry Date: 30/08/2021

This product is manufactured from blood donated by voluntary donors to the Australian Red Cross Lifeblood (ARCL).

Note: Cells 6, 7 and 8 form the RhD Negative Screening Cell Subset, formulated for antenatal RhD antibody screening in the presence of Prophylic Acid Anti-D.

Reference No	Cell No	Results
Ab1	1	0
Ab2	2	0
Ab3	3	0
Cord	4	0
A1	5	0
A2	6	0
B	7	0
	8	0
	9	0
	10	0
	11	0
	Auto	0